

Invacare® TDX® 4 Power Wheelchair



Yes, you can.®

Price List and Order Form Effective October 9, 2006 Revised October 1, 2007

For ease of ordering, contact Customer Service Toll Free at:
 800-333-6900 Fax: 800-678-4682
 To order Literature: 800-828-6282

www.invacare.com

Company Name _____ Account # _____

Phone # _____ P.O.# _____ Date _____

Ship To Name & Address _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

CHAIR TYPE

- TDX4 TDX4 w/MK6i™ (K0850)7,746.00
- TDX4BASE TDX4 Base Only w/MK6i.....7,006.00
- TDXSEAT TDX Seat Only740.00

Group 3 - K0850

- * TDX4 with MK6i
- * 400 lb. Weight Capacity
- * 4-Pole Heavy Duty Motors Standard
- * 1431-3 14" X 3" Tires w/Foam Filled Inserts. . . . \$250.00
- * Flip Back Adjustable Height Arm Standard
- * MK690 90 Amp Controller Standard
- * 70NHD 70° Swingaway Footrests Standard
- * 8 Amp Off Board Battery Charger Standard
- * 1515 Push Button Style Seat Positioning Strap . . Standard

FRAME

TRANSPORT TIE DOWN

- TRRO Transport Ready Option300.00
- TRBKTS Wheelchair Transport Brackets250.00

* TRANSPORT TIE-DOWN

TRRO includes four factory-installed transport brackets and a wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNA WC Vol 1. Section 19 Frontal Impact Test requirements for wheelchairs with a 168 lb. crash dummy, which corresponds to a person with a weight of 114 lbs. to 209 lbs.

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

* Transport Tie-Down option is not retrofittable to existing models and is not field serviceable.

TRRO FEATURES FOR TDX4

Features Required for TRRO	Optional Features for TRRO
♦ TDX4	♦ TrueTrack Upgrade
♦ Adjustable Rehab Seat with Adjustable Seat Tilt & Back Angle	♦ Adjustable Angle Flip UpFootplates
♦ Sectional or Straight Rear Uprights or No Push Handles	♦ Calf Strap
♦ 90° or 95° Back Angle	♦ Heel Loops w/Ankle Straps
♦ 18" - 24" ASBA Back Height	♦ Impact Guards
♦ Any Seat Width	♦ MPJ™+ Joystick
♦ Any Seat Depth	♦ All PSR / PSF+ Joysticks
♦ 400 lb. user weight limit	♦ Proportional Attendant Control
♦ Low Seat-to-Floor Height	♦ Left Hand Mounted Joystick
♦ Any Color Frame Finish	♦ "T" Handle Flexible Joystick Ext.
♦ Any Arm Type	♦ Straight Handle Flexible Joystick Ext.
♦ 70° Swingaway or 70° Tapered or Manual Center Mount Front Riggings	♦ MK6i Remote Programmer
♦ 14" x 3" Tires with Foam Filled Inserts	♦ Virtual Service Cable
♦ Please note that all TRRO chairs are set at a 10° seat angle	♦ Reinforced Back Upholstery
	♦ Group 24 Batteries and Battery Tray
	♦ Wheel Locks

FOOTNOTES

1. TDX4 comes standard with 400 lb. user weight capacity, 5 MPH top speed, heavy duty 4-pole motor package, Stability Lock. Base only includes MK690ACC Controller.

**For the most current pricing information go to www.invacare.com.
 NOTE: All specifications and dimensions are approximate.**

Invacare® TDX4

Accessories and Options

SEAT SELECTION

- ADJASBA Adj. Rehab Seat w/Adj. Seat Tilt & Back Angle.....**Standard**
- ADJRECL Manual Recline for Adjustable ASBA (E1226)744.00
- ¹ CTS Contoura™ Seating.....395.00
- OMCU Omit Contoura Seat Cushion.....Less 20.00

USER WEIGHT LIMITS

- ² U400 User Weight <400 lbs.....No Charge

TRUETRACK PACKAGE FOR TDX3

- ³ TRUTK True Track Upgrade1,899.00

WIDTH AND DEPTH

ADJUSTABLE SEAT WIDTH RANGE

- 1620W 16" to 20" Seat Width Range.....No Charge
- 2022W 20" to 22" Seat Width Range.....No Charge

⁴ ASBA SEAT WIDTH

- W16 Seat Width 16".....No Charge
- W17 Seat Width 17".....No Charge
- W18 Seat Width 18".....No Charge
- W19 Seat Width 19".....No Charge
- W20 Seat Width 20".....329.00
- W21 Seat Width 21".....329.00
- W22 Seat Width 22".....329.00

ADJUSTABLE SEAT DEPTH RANGE

- 1619D 16" to 19" Seat Depth Range.....No Charge
- 1922D 19" to 22" Seat Depth Range.....No Charge

⁵ SEAT DEPTH SETTING

- D16 Seat Depth 16".....No Charge
- D17 Seat Depth 17".....No Charge
- D18 Seat Depth 18".....No Charge
- D19 Seat Depth 19".....No Charge
- D20 Seat Depth 20".....529.00
- D21 Seat Depth 21".....529.00
- D22 Seat Depth 22".....529.00

ARMS

⁶ FIXED ARM TYPE (PAIR)

- 18F Full Length Fixed Height Conv.....No Charge
- 28F Desk Length Fixed Height Conv.....**Standard**

⁶ ADJUSTABLE ARM TYPE (PAIR)

- 19F Full Length Adj. Hgt Conv. (E0973).....230.00
- 29F Desk Length Adj. Hgt Conv. (E0973)230.00

⁶ MIXED HEIGHT ADJUSTABLE ARM TYPE (EACH)

- 19FR Full Length Adj. Hgt Right (E0973)115.00
- 29FR Desk Length Adj. Hgt Right (E0973)115.00
- 19FL Full Length Adj. Hgt Left (E0973).....115.00
- 29FL Desk Length Adj. Hgt Left (E0973).....115.00

ADJASBA SEAT BACK TYPE

- SECT Sectional Rear Uprights.....**Standard**
- SRUP Straight Rear Uprights.....No Charge
- 9720P No Push Handles (Straight Canes)82.00
- L8661 Omit Rear Uprights.....Less 20.00

SEAT BACK ANGLE DEGREE

- 90BA 90° Back Angle Setting.....**Standard**
- 95BA 95° Back Angle Setting.....No Charge

ADJASBA SEAT BACK HEIGHT

- BH12 Back Height 12".....No Charge
- BH13 Back Height 13".....No Charge
- BH14 Back Height 14".....No Charge
- BH15 Back Height 15".....No Charge
- BH16 Back Height 16".....No Charge
- BH17 Back Height 17".....No Charge
- BH18 Back Height 18".....No Charge
- BH19 Back Height 19".....No Charge
- BH20 Back Height 20".....No Charge
- BH21 Back Height 21".....No Charge
- BH22 Back Height 22".....No Charge
- BH23 Back Height 23".....No Charge
- BH24 Back Height 24".....No Charge

CONTOURA BACK HEIGHT

- TSBH25 Back Height 25".....No Charge
- TSBH295 Back Height 29.5".....No Charge

RECLINER BACK HEIGHT

- BHR18 Back Height 18".....No Charge
- BHR20 Back Height 20".....No Charge
- BHR22 Back Height 22".....No Charge
- BHR24 Back Height 24".....No Charge
- BHR26 Back Height 26".....No Charge

⁷ SEAT-TO-FLOOR HEIGHT

- LOW Seat-to-Floor Height - Low.....**Standard**
- MED Seat-to-Floor Height - Medium.....No Charge
- TALL Seat-to-Floor Height - Tall.....No Charge

WHEELS AND COLOR

WHEEL OPTION

- I430-3 14" X 3" Pneumatic Tires.....**Standard**
- I431-3 14" X 3" w/Foam Filled Inserts250.00

FENDER OPTION

- FFPKG Fender & (One-Sided) Fork Package500.00

FOOTNOTES

1. Selected seat depth on Contoura assumes 1" of back cushion compression. For example: When choosing an 18" depth the measurement from the front of the back cushion to the front of the seat pan equals 17". The Contoura back cushion is 2" thick. Please refer to page 5 for seating size offering.
2. Package includes 4 pole motors and 5 M.P.H with <400 lb. user weight limit unless TRUTK is selected.
3. Package includes TrueTrack driving technology, which enables driver to track straighter regardless of the terrain. Weight limit sets rear suspension must select MK6TT Controller.
4. Contoura™ available in even widths from 16"-22".
5. Contoura available in Seat Depths 16"-22".
6. All arms are flip-up and detachable.
7. Dimensions and specifics on page 5.

NOTE: All specifications and dimensions are approximate.

Invacare® TDX4

Accessories and Options

FRAME FINISH CONTEMPORARY

- 24P Wet Black.....No Charge
- 60P Silver Metallic.....No Charge
- 71P Silver Vein.....No Charge
- 105P Midnight Blue.....No Charge
- 115P Black Prism.....No Charge

TOO HOT

- 30P Sunny Yellow.....No Charge
- 61P Electric Red.....No Charge
- 119P Bubblegum Pink.....No Charge
- 120P Grape Madness.....No Charge

WAY COOL

- 62P Electric Blue.....No Charge
- 104P Electric Teal.....No Charge
- 121P Lolly Pop Blue.....No Charge
- 122P Cosmic Blue.....No Charge
- 125P Emerald Green.....No Charge
- 127P Grasshopper Green.....No Charge

MK6i ELECTRONICS

MK6i™ CONTROLLERS

- MK690 90 Amp Controller.....**Standard**
- MK6TT True Track Controller.....No Charge

EXPANDABLE CONTROLLER SYSTEMS

- MPJM6 MPJ+ Multiple Drive Joystick (E2377).....750.00
- PSFM6 Personal Joystick, Inductive on Front (E2399).....1,130.00
- PSRM6 Personal Joystick, Inductive on Rear (E2399).....1,130.00
- PBOD Push Button On/Off w/PB Drive Select.....No Charge
- PBSS Push Button On/Off w/Speed Select.....No Charge
- OR On/Off Mounted Right.....No Charge
- OL On/Off Mounted Left.....No Charge
- ASLRDYM6 MK6 Display Only
(No Driver Control) (E2377).....750.00

**FOR YOUR ASL DRIVE SYSTEM NEEDS
PLEASE REFER TO THE FOLLOWING:
ASL SINGLE INVOICE PRICE LIST FORM # 02-097.**

- 1500M6 Prop RIM Head Control w/Display (E2327).....2,875.00
- SMHD Small Headrest.....No Charge
- LGHD Large Headrest.....No Charge
- 1558MM6 Compact Joystick w/Display (E2373).....1,130.00

ALTERNATIVE DIGITAL CONTROL SYSTEMS

- SNPM6 Sip-N-Puff Digital Interface Box (E2325).....1,500.00
- 5018M6 Wafer Board & Digital Interface Box (E2322).....1,550.00
- 5020M6 Mini Tash Joystick & Digital Interface Box (E2321).....1,750.00

SIP-N-PUFF KIT

- PKG32666 Therafin® Sip-N-Puff Breath Tube Kit (E2326).....380.00

ADDITIONAL CONTROL CHOICES

- 1558M6 Compact Joystick - Less std Mount.....770.00
- 5018 Wafer Board.....310.00
- 5020 Mini Tash Joystick.....405.00
- PACM6 Proportional Attendant Control.....850.00
- 1552M Digital Attendant Control.....850.00

JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retract. Joystick Mnt MPJ+ (E1028).....279.00
- ARM250 Stealth Height Adjustable Joystick Mount.....210.00
- LEFT Left Handed Mounted Joystick.....No Charge
- RIMHW RIM Hardware For 1500 (E1028).....219.00
- PKG32669 Complete Bib Assembly (Tash/1558).....265.00
- GATMPJ Gatlin Midline for MPJ (E1028).....600.00
- GAT1812 Gatlin Midline for 1812 (E1028).....600.00
- RCM Rear Cane Mount for Compact Joystick.....395.00

ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer w/Pro Memory Card.....475.00
- 1813SD Pro Memory Card (USB Ready).....90.00
- 1560 "T" Handle Flexible Joystick Ext. (E2323).....77.00
- 1561 Straight Handle Flexible Joystick Ext. (E2323).....77.00
- MPJLHS MPJ+ Added Left Hand Support Platform.....49.00
- MPJRHS MPJ+ Added Right Hand Support Platform.....49.00
- 1826 Chin Cup (E2324).....48.00
- A24VPS 24V MK6i Auxiliary Power Source.....200.00
- AUX12M6 Auxiliary Module For 1 & 2.....795.00
- AUX34M6 Auxiliary Module For 3 & 4.....795.00
- MPY Monoport "Y" Cable for MPJ+ or Display.....15.00
- 2PB Dual Push Button.....147.00
- 2SW Dual Toggle.....172.00

ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Switch - Black.....65.00
- ASL304A Light Wobble Switch w/611 Mount.....320.00
- ASL202J Fiber Optic Switch.....610.00
- ASL208J Adjustable Proximity Switch.....550.00
- ASL504JM6 Remote Emergency Stop Switch.....535.00

FOOTNOTES

1. TOD-SS: Toggle On/Off-Drive Select, Speed Potentiometer Standard.
2. PSRM6 Only.
3. Drive Select through mode switch, or mode switch and Joystick Command.
4. Display Only, No Driver Control. Allows choice from order form.
5. Must also use either a Multiple Drive Joystick or MK6 Display.
6. Requires SNPM6 Sip-N-Puff Interface Box.
7. QLAM6 includes ARM250 standard.
8. For use with expandable controller systems to provide additional switch site.
9. Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options.
10. Must add MK6 Auxiliary Power Source.

NOTE: All specifications and dimensions are approximate.

Invacare® TDX4

Accessories and Options

POSTIONING STRAPS

SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Pos. Strap50.00
- 1515 Push Button Style Seat Pos. Strap**Standard**

CHEST POSITIONING STRAPS

- 7311BK Airline Buckle Chest Pos. Strap (E0960).....100.00
- 9515 Push Button Style Chest Pos. Strap (E0960).....100.00
- 7321BK Hook & Loop Chest Pos. Strap (E0960).....100.00

BATTERY

¹ BATTERY TRAY TYPE

- 24TRY 24 Group Style Tray.....**Standard**
- ² 22TRY3 22NF Style Tray for 3 Batteries.....No Charge

³ BATTERIES

- 22NFBATTERY 22NF Battery (On Chair) (E2361).....300.00
- 24BATTERY 24 Gel Battery (On Chair) (E2363).....410.00

ACCESSORIES

BACK ACCESSORIES

- L8601 Omit Back UpholsteryNo Charge
- RBP Reinforced Back Upholstery.....No Charge

MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks45.00
- SEATO2 Oxygen Holder (ASBA Seat Only) (E2208).....150.00

FRONT RIGGING OPTIONS

FRONT RIGGING

⁴ CENTER MOUNT

- ⁵ CTMT-F Center Mount Adj. Knee Angle220.00
- ⁶ CTMT-FJ Junior Center Mount Adj. Knee Angle.....220.00
- ⁷ CTMT-FL Center Mount Rigging w/Extension 16"-18.75" ..315.00

CENTER PIVOT STYLE (PAIR)

- ⁸ 70NHD 70° Swingaway Footrests**Standard**
- ⁹ 70STAPER 70° Tapered Footrests.....No Charge
- ¹⁰ AT5544 Manual ELR (E0990)255.00

¹¹ POWER LEGS MODULE CENTER PIVOT STYLE

- ¹² ELRPW Pwr ELR w/Articulation 20" Ext (Pair) (E1010)....1,895.00
- ¹³ ELRPW22 Pwr ELR w/Articulation 22" Ext (Pair) (E1010)....1,995.00

¹⁴ POWER LEGS CONTROL MODULE

- 4WSB 4 Way Sw Box for Individual Pwr Leg Control..No Charge
- S4WSB Multiple Actuator Intrfce Bx for Pwr Legs (E2311).....2,600.00
- FWT 4 Way ToggleNo Charge
- QPB 4 Push Buttons.....No Charge

For all systems with two or more actuators you must choose 4WSB at No Charge or the S4WSB at an Up Charge.

S4WSB Replaces the Four Way Switch Box and allows multiple actuator operation through the driver control.

FWT or QPB is required when choosing the 4WSB or S4WSB.

PIN STYLE HEMI (PAIR)

- PHW93 70° Swingaway Footplates.....No Charge
- PHAL4A Manual Articulating/ELR (K0053)225.00

FOOTPLATE OPTIONS (PAIR)

- ¹⁵ 1651 Composite Footplates.....No Charge
- ¹⁶ 1350 Extra Large Aluminum Footplates.....No Charge
- ¹⁷ AT5543 Adjustable Angle Flip-Up Footplates165.00
- ¹⁸ FBCTMT Footboard for Center Mount.....150.00

FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap (K0038)27.00
- 1600BK Heel Loops w/Ankle Straps (Pair) (E0951).....56.00
- LPT2 Longer Pivot Slide Tube For 1350.....83.00
- IMPACTG Impact Guards51.00
- ALPT Longer Pivot And Slide Tube.....84.00
- ALPT4 Longer Pivot And Slide Tube 4"84.00

FOOTNOTES

- 1 Due to Battery harness & tray, only batteries with terminal configuration dimension per page 4 will work with the Invacare TDX chairs.
- 2 Available when picking base only.
- 3 Sealed Gel batteries highly recommended for use in Invacare power chairs.
- 4 All front rigging options include a standard footplate with the exception of the 70NHD.
- 5 Allows 60 to 90 degree static knee angle positioning with a drop back lever for easier transferring. Extension range 12.25" - 15".
- 6 Allows 60 to 90 degree static knee angle positioning with a drop back lever for easier transferring. Extension range 9.25" - 12".
- 7 Not available with Low seat-to-floor.
- 8 With seat widths of 21" to 24" and depth of 16", must pick footplate 1350. The AT5543 is not available.
- 9 Standard with impact guards. Not available with widths 19" to 24" and 16" deep seat.
- 10 Not available with 16" wide by 16" to 18" deep or 17" wide by 16" deep and low seat-to-floor.
- 11 To operate Power Legs must select 4SWB or S4SWB.
- 12 Independent powered elevating legrests include elevating and articulating legrests with angle adjustable footplates. Extension range 13" to 20. Not available with depths of 16" or 17" with low seat-to-floor.
- 13 Independent powered elevating legrests include elevating and articulating legrests with angle adjustable footplates. Extension range 15" to 22". Not available with depths of 16" or 17" with low seat-to-floor.
- 14 These options are located in the Configurator under the MK6i Electronics Folder. The Category Header inside the folder is labeled Power Seating Accessories
- 15 1651 standard footplate on 70NHD.
- 16 Not available with AT5544 and 17" wide by 16" deep and low seat-to-floor.
- 17 Not available with AT5544 and 16" wide by 16"-18" deep or 17" wide by 16" deep and low seat-to-floor.
- 18 One piece footboard for Center Mount (CTMT-F/CTMT-FL) Riggings.

NOTE: All specifications and dimensions are approximate.

Invacare® TDX4

Accessories and Options

TDX4 Specifications	
Seat Width Range: (1" increments)	ADULT (Adjust) 16 to 24-inches
Seat Depth Range: (in 1-inch increments)	16 to 22-inches
Back Height Range:	12 to 24-inches
Back Angle Range:	85° to 105° in 5° increments
Seat-to-Floor (approximate)	
Standard:	17-inches @ 5° 16.5 @ 0°
Medium:	19-inches @ 5°
Tall:	21-inches @ 5°
Overall Width of Base:	25.5-inches (w/o joystick)
Overall Height (with 16" Back Height and 0° Seat Angle)	
Standard:	33.25-inches
Medium:	35.25-inches
Tall:	37.25-inches
Overall Length:	35.25-inches 43-inches w/o front rigging w/ 90° front rigging
Weight	
4 Pole Motor	
W/O Batteries:	186 lbs.
With Two (2) GP24 Batteries:	290 lbs.
Drive Axle:	Non-Adjustable
Drive Wheels/Tires:	(Foam Filled or Pneumatic)
Standard:	14 X 3-inches

Batteries
(not included)

GP24 batteries with terminal configuration (negative on the left and positive on the right) **MUST** be used. GP24 batteries that have the reverse terminal configuration **MUST** not be used. Terminals **MUST** have a cross hole for proper battery connection. See Owner's Manual, part number I114809. These recommendations **MUST** be followed otherwise injury and damage may occur.

Contoura Seat Sizes							
	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X
*Not available with Contoura Seat Cushion.							
Back Heights		25" or 29.5"					

Memory Card - BASIC

- Standard with ALL Expandable Control Systems (MPJ+, PSR+, PSF+, MK6i Display).
- Used to Back Up/Restore Programmed Settings for Only One Chair.
- May be given to user or kept with provider's files for safe keeping.
- Does not Contain Advanced Diagnostics/Help Screens/File Structure.
- Not compatible with SPJ+ Joysticks.

Memory Card - PROFESSIONAL

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains Advanced Diagnostics, Help Screens, Software Updates, File Storage/Retrieval.
- Used to back up multiple chairs and programming settings.
- Intended for use by qualified providers only.
- Not compatible with SPJ+ Joysticks.

Weight Limit 400 lbs
Invacare recommends all wheelchairs be fitted and programmed to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

NOTE:
All weight limits are inclusive of add-on equipment and seating systems. All specifications and dimensions are approximate.

SEATING AND POSITIONING

**FOR YOUR SEATING NEEDS
PLEASE REFER TO THE FOLLOWING:
REHAB SEATING & POSITIONING PRODUCTS
SEE PRICE LIST FORM # 07-027.**

Invacare Corporation www.invacare.com



USA
One Invacare Way
Elyria, OH USA
44035-4190
800-333-6900

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Therafin is a registered trademark of W/R Medical Electronics Co. Products designated by ASL are Manufactured by Adaptive Switch Labs, Inc. Specifications and prices are subject to change without notification.
©2007 Invacare Corporation