

Invacare® Formula™ CG Elevating Seat for Storm Series® Bases



Price List and Order Form Effective October 9, 2006 Revised October 1, 2007

For ease of ordering, contact Customer Service Toll Free at:

800-333-6900 Fax: 800-678-4682

To order Literature: 800-828-6282

www.invacare.com

Company Name _____ Account # _____

Phone # _____ P.O.# _____ Date _____

Ship To Name & Address _____

Special Note

Ordering a single power seating system on a multiple power option wheelchair base may require an affirmative determination through Advanced Determination of Medicare Coverage (ADMC). Please consult an attorney or other advisor to discuss this situation and any other specific billing situation in further detail. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program. This price list is for information purposes only and is not billing or legal advice.

POWERED SEATING BASE

- 3GAR-CG Arrow for Formula CG w/MK6i™9,325.00
- 3GTQ-CG Torque SP for Formula CG w/MK6i.....6,450.00
- 3GTQ3-CG Torque 3 for Formula CG w/MK6i.....5,600.00

Group 4 - K0884 Basic Package

- * 3GAR-CG Arrow for Formula CG with MK6i
- * 300 lb Weight Capacity
- * TrueTrack Gearless Brushless GB™ Motors.....Standard
- * 1228-3 8" X 2.25" Tires w/Foam Filled Inserts.....\$250.00
- * Flip Back Adjustable Height Arm.....Standard
- * MK6TT TrueTrack Controller.....Standard
- * 70NHD 70° Swingaway Footrests.....Standard
- * 8 Amp Off Board Battery Charger.....Standard
- * 1515 Push Button Style Seat Positioning Strap.....Standard

Seat Option

- * CGE Elevate (E2300).....\$1,895.00

Group 4- K0884 Basic Package

- * 3GTQ-CG Torque SP for Formula CG with MK6i
- * 300 lb Weight Capacity
- * Performance 4-Pole Motors.....Standard
- * 1228-3 8" X 2.25" Tires w/Foam Filled Inserts.....\$250.00
- * Flip Back Adjustable Height Arm.....Standard
- * MK690ACC Controller w/Actuator Control.....Standard
- * 70NHD 70° Swingaway Footrests.....Standard
- * 8 Amp Off Board Battery Charger.....Standard
- * 1515 Push Button Style Seat Positioning Strap.....Standard

Seat Option

- * CGE Elevate (E2300).....\$1,895.00

Group 3 - K0856 Basic Package

- * 3GTQ3-CG Torque 3 for Formula CG with MK6i
- * 300 lb Weight Capacity
- * 4-Pole Motors.....Standard
- * 1228-3 8" X 2.25" Tires w/Foam Filled Inserts.....\$250.00
- * Flip Back Adjustable Height Arm.....Standard
- * MK690ACC Controller with Actuator.....Standard
- * 70NHD 70° Swingaway Footrests.....Standard
- * 8 Amp Off Board Battery Charger.....Standard
- * 1515 Push Button Style Seat Positioning Strap.....Standard

Seat Option

- * CGE Elevate (E2300).....\$1,895.00

Invacare® Formula™ CG Elevating Seat

Accessories and Options

Base

TRANSPORT TIE DOWN

- ¹ TRBKTS Wheelchair Transport Brackets250.00

***Transport Tie Down**

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

*Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

USER WEIGHT LIMITS

ARROW

- U300 User Weight <300 lbs.....No Charge

TORQUE SP/TORQUE 3

- U300 User Weight <300 lbs.....No Charge

DRIVE TYPE OPTION

- ² TRUTK True Track Package1,899.00
- ³ PERF Performance 4 Pole Motor.....No Charge

FRAME TYPE

- ⁴ SFS Short BaseNo Charge
- LFS Long Base.....No Charge

CASTER OPTIONS

- I259 6" X 2" Semi Pneumatic CastersNo Charge
- I228-3 8" X 2.25" Semi Pneumatic Casters Lt. Grey..**Standard**
- I251-3 8" X 2" Pneumatic CastersNo Charge
- ⁵ I220-3 9" X 2.75" Pneumatic Casters.....No Charge
- ⁵ I465-3 9" X 2.75" Foam Filled Casters.....No Charge

⁶ CASTER SUSPENSION

- 9213-6 Shock Fork for 6" Casters (E1016)119.00
- ⁵ 9213-8 Shock Fork for 8" Casters (E1016)119.00

WHEEL OPTIONS

- I430-3 14" X 3" Pneumatic Tires**Standard**
- I431-3 14" X 3" Tires w/Foam Filled Inserts.....250.00
- ⁷ I440 14" X 4" Pneumatic Tires200.00
- ⁷ I441 14" X 4" Tires w/Foam Filled Inserts).....320.00

⁸ FREE WHEEL HUB

- FHUB Free Wheel Hub125.00

FRAME FINISH

CONTEMPORARY

- 24P Wet BlackNo Charge
- 60P Silver Metallic.....No Charge
- 71P Silver Vein.....No Charge
- I05P Midnight BlueNo Charge
- I15P Black Prism.....No Charge

TOO HOT

- 30P Sunny Yellow.....No Charge
- 61P Electric Red.....No Charge
- I19P Bubblegum Pink.....No Charge
- I20P Grape Madness.....No Charge

WAY COOL

- 62P Electric BlueNo Charge
- I04P Electric Teal.....No Charge
- I21P Lolly Pop Blue.....No Charge
- I22P Cosmic Blue.....No Charge
- I25P Emerald Green.....No Charge
- I27P Grasshopper GreenNo Charge

FOOTNOTES

1. Only Available with MED or TALL Seat to Floor Height.
2. TrueTrack includes Gearless Brushless GB™ motors which reaches 7 mph and has a 300 lb. weight capacity and Monroe® spring shock suspension. Must upgrade to MK6TT Controller. Available on Torque SP. **When selecting TrueTrack Motors, the top speed of the chair increases to 7 MPH and therefore no longer performs to the Group 3 code specifications.**
3. Standard on Torque SP. 6.5 MPH.
4. Not available with seat depths 19"-22". Must pick long base.
5. Not available with short base.
6. Not available with 9" casters.
7. 14" x 4" wheels increase the overall width of the base by 2.5". Does not include Storm Series design wheel or fender. Not available with (TRUTK) True Track Pkg.
8. Will add .75" to overall width. Not available with TrueTrack. Does not include Storm Series design wheel or fender.

NOTE: All specifications and dimensions are approximate.

Invacare® Formula™ CG Elevating Seat

Accessories and Options

Seat

SYSTEM TYPE

- CGE 7" Elevating Seat - up to 300 lbs. (E2300) 1,895.00

SEAT SELECTION

- ADJASBA Adj. Rehab Seat w/Adj. Seat Tilt Back Angle **Standard**
- ¹ BPO Back Posts Only - Round Canes Only No Charge
- ² CTS Contoura™ Seating 395.00
- CURVABS Curved ABS No Charge
- OMCU Omit Contura Seat Cushion Less 20.00

SEAT WIDTH RANGE

- 1620W 16" to 20" Width (Adjustable) No Charge
- 2022W 20" to 22" Width (Adjustable) No Charge

² SEAT WIDTH SETTING

- W16 16" Seat Width No Charge
- W17 17" Seat Width (not available w/Contoura) ... No Charge
- W18 18" Seat Width No Charge
- W19 19" Seat Width (not available w/Contoura) ... No Charge
- W20 20" Seat Width No Charge
- W21 21" Seat Width (not available w/Contoura) ... No Charge
- W22 22" Seat Width No Charge

SEAT DEPTH RANGE

- 1619D 16" to 19" Deep (Adjustable) No Charge
- 1922D 19" to 22" Deep (Adjustable) (Long Base) ... No Charge

SEAT DEPTH SETTING

- D16 16" Seat Depth No Charge
- ³ D17 17" Seat Depth No Charge
- ³ D18 18" Seat Depth No Charge
- D19 19" Seat Depth No Charge
- D20 20" Seat Depth No Charge
- ⁴ D21 21" Seat Depth No Charge
- ⁴ D22 22" Seat Depth No Charge

SEAT TILT

- OSA 0° **Standard**

ASBA SEAT BACK TYPE

- SECT Sectional Rear Uprights **Standard**
- SRUP Straight Rear Uprights No Charge
- 9720P No Push Handles (Straight Canes) 82.00
- L866I Omit Rear Uprights Less 20.00

BACK HEIGHT

CURVED ABS OR ADJUSTABLE ASBA

- BH16 Back Height 16" No Charge
- BH17 Back Height 17" No Charge
- BH18 Back Height 18" No Charge
- BH19 Back Height 19" No Charge
- BH20 Back Height 20" No Charge
- BH21 Back Height 21" No Charge
- BH22 Back Height 22" No Charge
- BH23 Back Height 23" No Charge
- BH24 Back Height 24" No Charge

CONTOURA

- TSBH25 Back Height 25" No Charge
- TSBH295 Back Height 29.5" No Charge

ADJUSTABLE ARM TYPE (PAIR)

- 19F Full Length Adj. Hgt. No Charge
- 29F Desk Length Adj. Hgt. No Charge

MIXED ADJUSTABLE ARM TYPE (EACH)

- 19FR Full Length Adj. Hgt. Right No Charge
- 29FR Desk Length Adj. Hgt. Right No Charge
- 19FL Full Length Adj. Hgt. Left No Charge
- 29FL Desk Length Adj. Hgt. Left No Charge

SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Pos. Strap 50.00
- 1515 Push Button Style Seat Pos. Strap **Standard**

CHEST POSITIONING STRAPS

- 7311BK Airline Buckle Chest Pos. Strap (E0960) 100.00
- 9515 Push Button Style Chest Pos. Strap (E0960) 100.00
- 7321BK Hook & Loop Chest Pos. Strap (E0960) 100.00

FOOTNOTES

1. Push handles included with BPO.
2. Contoura only available in even widths. Selected seat depth on Contoura assumes 1" of back cushion compression. For example: When choosing 18" depth the measurement from the front of the back cushion to the front of the seat pan equals 17".
The Contoura back cushion is 2" thick includes leatherette seat cushion.
3. Not available with Contoura seat cushion on 22" wide.
4. Not available with Contoura seat cushion on 16" or 18" wide.

NOTE: All specifications and dimensions are approximate.

Invacare® Formula™ CG Elevating Seat

Accessories and Options

Electronics

MK6i™ CONTROLLER

TORQUE SP/TORQUE 3

- MK690ACC MK690ACC Controller w/Actuator Control...**Standard**
- MK6TT True Track Controller.....No Charge

ARROW

- MK6TT True Track Controller.....No Charge

EXPANDABLE CONTROLLER SYSTEMS

- MPJM6 MPJ+ Multiple Drive Joystick (E2377).....750.00
- ² PSFM6 Personal Joystick, Inductive on Front (E2399) ...1,130.00
- ² PSRM6 Personal Joystick, Inductive on Rear (E2399).....1,130.00
- ³ PBOD Push Button On/Off w/PB Drive Select.....No Charge
- ^{3,4} PBSS Push Button On/Off w/Speed Select.....No Charge
- ³ OR On/Off Mounted Right.....No Charge
- ³ OL On/Off Mounted Left.....No Charge
- ⁵ ASLRDYM6 MK6 Display Only
(No Driver Control) (E2377).....750.00

**FOR YOUR ASL DRIVE SYSTEM NEEDS
PLEASE REFER TO THE FOLLOWING:
ASL SINGLE INVOICE PRICE LIST FORM # 02-097.**

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS

- 1500M6 Prop RIM Head Control w/Display (E2327)2,875.00
- SMHD Small Headrest.....No Charge
- LGHD Large Headrest.....No Charge
- 1558MM6 Compact Joystick w/Display (E2373).....1,130.00

ALTERNATIVE DIGITAL CONTROL SYSTEMS

- ⁶ SNPM6 Sip-N-Puff/Digital Interface Box (E2325).....1,500.00
- ⁶ 5018M6 Wafer Board & Digital Interface Box (E2322).....1,550.00
- ⁶ 5020M6 MiniTash Joystick & Digital Interface Box (E2321).....1,750.00

SIP N PUFF KIT

- PKG32666 Therafin® Sip-N-Puff Breath Tube Kit (E2326).....380.00

ADDITIONAL CONTROL CHOICES

- 1558M6 Compact Joystick - Less Std Mount.....770.00
- ⁷ 5018 Wafer Board310.00
- ⁷ 5020 MiniTash Joystick405.00
- PACM6 Proportional Attendant Control850.00
- ⁷ 1552M Digital Attendant Control.....850.00

JOYSTICK MOUNTING HARDWARE

- ⁸ QLAM6 Quad Link Retract. Joystick Mnt MPJ+ (E1028)....279.00
- ARM250 Stealth Height Adjustable Joystick Mount (E1028) ..210.00
- LEFT Left Handed Mounted Joystick.....No Charge
- RIMHW RIM Hardware For 1500 (E1028)219.00
- PKG32669 Complete Bib Assembly (Tash/1558)265.00
- GATMPJ6 Gatlin Midline for MPJ (E1028)600.00
- GAT1812 Gatlin Midline for 1812 (E1028)600.00
- RCM Rear Cane Mount for Compact Joystick.....395.00

ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer w/Pro Memory Card.....475.00
- 1813SD Pro Memory Card (USB Ready).....90.00
- 1560 "T" Handle Flexible Joystick Ext. (E2323)77.00
- 1561 Straight Handle Flexible Joystick Ext. (E2323).....77.00
- 1826 Chin Cup (E2324)48.00
- A24VPS 24V MK6 Auxiliary Power Source200.00
- AUX12M6 Auxiliary Module For 1 & 2795.00
- AUX34M6 Auxiliary Module For 3 & 4795.00
- ⁹ 2PB Dual Push Button147.00
- ⁹ 2SW Dual Toggle172.00

POWER SEATING ACCESSORIES

- ¹⁰ SANODE Single Actuator Operation thru Driver Control (E2310).....1,300.00

ELECTRONIC ACCESSORIES

- EGSBLK Egg Switch - Black.....65.00
- ASL304 Light Wobble Switch w/611 Mount320.00
- ¹¹ ASL202J Fiber Optic Switch.....610.00
- ¹¹ ASL208J Adjustable Proximity Switch550.00
- ¹¹ ASL504M6 Remote Emergency Stop Switch.....535.00
- ASL300A Egg Switch w/611 Mount300.00

FOOTNOTES

1. Includes Monoport "Y" Cable to provide additional switch input port.
2. TOD-SS: Toggle On/Off Drive Select, Speed, Speed Potentiometer Standard.
3. PSRM6 Only.
4. Drive Select through mode switch, or mode switch and Joystick Commands.
5. Display Only, No Driver Control. Allows choice from order form.
6. Must also use either a Multiple Drive Joystick or MK6 Display.
7. Requires SNPM6 Sip-N-Puff Interface Box.
8. QLA includes ARM250 standard.
9. Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options.
10. Required for Elevate Operation through Driver Control - Single Actuator Systems. Not compatible with SPJ+ Joysticks. The Egg Switch with the 611 Mount does not come standard when ordering a SANODE.
11. Not available with SPJ+ Joystick. Must add MK6 Auxiliary Power Source.

NOTE: All specifications and dimensions are approximate.

Invacare® Formula™ CG Elevating Seat
Accessories and Options

Miscellaneous

BATTERY BOXES

- ¹ GP24 Group 24 Battery BoxesNo Charge
- 22NF 22NF Battery Box.....No Charge

² **BATTERIES**

- ³ 22NFBATTERY 22 NF Battery (On Chair) (E2361).....300.00
- 24BATTERY 24 Gel Battery (On Chair) (E2363).....410.00

MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks45.00
- ⁴ 9205-3 Suspension Shocks (E1016).....400.00
- ⁴ 9205H-3 Suspension Shocks (Heavy Duty) (E1018)400.00
- ⁴ FDR Fender.....320.00
- ⁵ EXTP Extended Anti-Tippers.....250.00

FOOTNOTES

1. Option on Torque SP and Torque 3 (Standard on Arrow). Must select for TRBKTS.
2. Batteries will be installed by Invacare Corporation at no additional charge. Chair will be shipped complete with batteries.
3. Not available with Arrow,GB, Ranger X, or TrueTrack systems.
4. Standard with Arrow or True Track Systems.
5. Comes standard with vent tray. Adds 2.5" to overall length.
6. Not available with 19" or 20" wide and 9" casters. Standard with impact guards on the footrests.
7. Manual elevating legrest.

NOTE: All specifications and dimensions are approximate.

Front Riggings

FRONT RIGGING

CENTER PIVOT STYLE (ASB) (PAIR)

- 60HD 60° Swingaway FootrestsNo Charge
- 70NHD 70° Swingaway Footrests.....**Standard**
- ⁶ 70STAPER 70° Tapered Footrests.....No Charge
- ⁷ AT5544 Manual Elevating Legrests (E0990).....255.00

PIN STYLE (PSB) (PAIR)

- PVV93 Swingaway Footrests 13.5" - 19"No Charge
- ⁷ P904A Swingaway Elevating Legrests 16" - 19.75" (E0990).....255.00
- ⁷ PAL4A Smartleg Art. Legrests 17.25" - 21.5" (K0053).....320.00

FOOTPLATE OPTIONS (PAIR)

- 1651 Flip -Up Composite Footplates.....No Charge
- 1350 Extra Large Aluminum Footplates.....No Charge
- AT5543 Adjustable Angle Flip-Up Footplates.....165.00

FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap (E0038).....27.00
- 1600BK Heel Loops w/Ankle Straps (E0951)56.00
- LPT2 Longer Pivot Slide Tube For 1350.....83.00
- IMPACTG Impact Guards.....51.00
- ALPT Longer Pivot And Slide Tube84.00
- ALPT4 Longer Pivot And Slide Tube 4"84.00

Seating

**FOR YOUR SEATING NEEDS
PLEASE REFER TO THE FOLLOWING
REHAB SEATING & POSITIONING PRODUCTS
SEE PRICE LIST FORM # 07-027.**

Invacare® Formula™ CG Elevating Seat

Accessories and Options

Battery Boxes	
*Standard Battery Tray	
Arrow	Group 24
Torque SP/Torque 3	22NF Standard Group 24 (optional)
*Sealed Gel Cell Batteries Recommended	

Batteries

22NF batteries with terminal configuration (positive on the left and negative on the right) **MUST** be used. 22NF batteries that have the reverse terminal configuration **MUST** not be used. Terminals **MUST** have a cross hole for proper battery connection.

Contoura Seat Sizes							
	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X
*Not available with Contoura Seat Cushion							
Back Heights	25" or 29.5"						

Memory Card - BASIC

- Standard with ALL Expandable Control Systems (MPJ+, PSR+, PSF+, MK6i Display).
- Used to Back Up/Restore Programmed Settings for Only One Chair.
- May be given to user or kept with provider's files for safe keeping.
- Does not Contain Advanced Diagnostics/Help Screens/File Structure.
- Not compatible with SPJ+ Joysticks.

Memory Card - PROFESSIONAL

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains Advanced Diagnostics, Help Screens, Software Updates, File Storage/Retrieval.
- Used to back up multiple chairs and programming settings.
- Intended for use by qualified providers only.
- Not compatible with SPJ+ Joysticks.

Storm for Formula CG Elevating Seat Specifications			
	Arrow	Torque SP	Torque 3
Maximum Speed MPH			
TrueTrack	7	7 (option)	7 (option)
4-pole	N/A	N/A	4.8
4-pole (Performance)	N/A	6.5	N/A
4-pole HD	N/A	N/A	N/A
Seat-To-Floor Height	17.75"	17.75"	17.75"
Weight Limits			
<300 lb.		X	X
<300 lb. True Track	X	X	X

Formula CG Elevating Seat Features and Specifications
<ul style="list-style-type: none"> ♦ Up to 300 lb. weight capacity. ♦ Traditional Rehab Seating (ASBA). ♦ Seat-to-Floor height of 17.75" at 0° tilt. ♦ Speed inhibit feature of 20% of maximum speed while elevated.

After Market Backrests
<ul style="list-style-type: none"> ♦ The following companies also have hardware to make their backrests work with our new canes: <ul style="list-style-type: none"> ♦ AEL - planar seating system ♦ Freedom Designs® - planar seating system ♦ Canyon - planar seating system

Curved ABS or Adjustable ASBA	
Seat Widths	16" - 22"
Seat Depths	16" - 22"
Back Heights	16" - 24"

Invacare Corporation www.invacare.com



USA
 One Invacare Way
 Elyria, OH USA
 44035-4190
 800-333-6900

All rights reserved. Trademarks are identified by the symbols [™] and [®]. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Therafin is a registered trademark of W/R Medical Electronics Co. Freedom Designs is a registered trademark of Freedom Designs, Inc. Products designated by ASL are Manufactured by Adaptive Switch Labs, Inc. Specifications and prices are subject to change without notification.
 ©2007 Invacare Corporation