

Invacare® Formula™ CG Tilt Only for TDX™ SP Base



Yes, you can.®

Price List and Order Form Effective October 9, 2006 Revised October 1, 2007

For ease of ordering, contact Customer Service Toll Free at:

800-333-6900 Fax: 800-678-4682

To order Literature: 800-828-6282

www.invacare.com

Company Name _____ Account # _____

Phone # _____ P.O.# _____ Date _____

Ship To Name & Address _____

Special Note

Ordering a single power seating system on a multiple power option wheelchair base may require an affirmative determination through Advanced Determination of Medicare Coverage (ADMC). Please consult an attorney or other advisor to discuss this situation and any other specific billing situation in further detail. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program. This price list is for information purposes only and is not billing or legal advice.

POWERED SEATING BASE

TDXSP-CG TDXSP-CG for Formula CG w /MK6i™ (K0856/K0858) . . . 5,600.00

Group 3 - K0856 Basic Package

- * TDXSP-CG TDXSP for Formula CG with MK6i
- * 300 lb Weight Capacity
- * 4-Pole Motors Standard
- * 1431-3 14" X 3" Tires w/ Foam Filled Inserts . . . \$250.00
- * Flip Back Adjustable Height Arm Standard
- * MK690ACC 90 Amp Controller Standard
- * 70NHD 70° Swingaway Footrests Standard
- * 8 Amp Off Board Battery Charger Standard
- * 1515 Push Button Style Seat Positioning Strap . . Standard

Seat Option

- * CGT Tilt (E1002) \$4,500.00

Group 3 - K0858 Basic Package

- * TDXSP-CG TDXSP for Formula CG with MK6i
- * 400 lb Weight Capacity \$1,175.00
- * Heavy Duty 4-Pole Motors Standard
- * 1431-3 14" X 3" Tires w/ Foam Filled Inserts . . . \$250.00
- * Flip Back Adjustable Height Arm Standard
- * MK690ACC 90 Amp Controller Standard
- * 70NHD 70° Swingaway Footrests Standard
- * 8 Amp Off Board Battery Charger Standard
- * 1515 Push Button Style Seat Positioning Strap . . Standard

Seat Option

- * CGT Tilt (E1002) \$4,500.00

OMIT SEAT FRAME

- OSF Omit Seat Frame for Base Only No Charge

TRANSPORT TIE DOWN

- TRBKTS Wheelchair Transport Brackets 250.00

USER WEIGHT LIMITS

*Transport Tie Down

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

*Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

- ¹ U300 User Weight <300 lbs (K0856) Standard
- ² U400 User Weight <400 lbs (K0858) 1,175.00

WHEEL OPTIONS (PAIR)

- 1430-3 14" X 3" Pneumatic Tires Standard
- 1431-3 14" X 3" w/ Foam Filled Inserts 250.00

FORK OPTION

- FKPKG One Sided Fork Package 500.00

FRAME FINISH CONTEMPORARY

- 24P Wet Black No Charge
- 60P Silver Metallic No Charge
- 71P Silver Vein No Charge
- 105P Midnight Blue No Charge
- 115P Black Prism No Charge

TOO HOT

- 30P Sunny Yellow No Charge
- 61P Electric Red No Charge
- 119P Bubbegum Pink No Charge
- 120P Grape Madness No Charge

WAY COOL

- 62P Electric Blue No Charge
- 104P Electric Teal No Charge
- 121P Lolly Pop Blue No Charge
- 122P Cosmic Blue No Charge
- 125P Emerald Green No Charge
- 127P Grasshopper Green No Charge

FOOTNOTES

1. Package includes 4-pole motors and 5.8 MPH top speed with <300 lb user weight limit. Package includes Quiet Stability Lock.
2. Package includes Heavy Duty 4-pole motors and 5 MPH top speed with <400 lb user weight limit.

Note: All Specifications and dimensions are approximate.

Invacare® Formula™ CG Tilt Only

Accessories and Options

SEAT

SYSTEM TYPE

- CGT CG Tilt Only System (E1002)4,500.00
- CNVACTS Conventional Actuators.....**Standard**
- SMTACTS Smart Actuators.....500.00

SEAT-TO-FLOOR-HEIGHT

- LOW 17.5" at 0°No Charge
- MED 18.5" at 0°No Charge
- TALL 19.5" at 0°No Charge

SEAT TILT

- 0SA 0° Seat TiltNo Charge
- 5SA 5° Seat TiltNo Charge

SEAT SELECTIONS

- CTS Contoura™ Seating.....395.00
- BPO Back Posts Only/Round Canes Only.....No Charge
- CURVABS Curved ABS.....No Charge
- CNB Conventional Back.....No Charge
- OMCU Omit Contoura Seat Cushion.....Less 20.00

PUSH HANDLES

- PCNB Push Handles for Conventional.....No Charge

SEAT WIDTH RANGE

- 1216W 12" to 16" Seat Width Range (Adjustable).....No Charge
- 1620W 16" to 20" Seat Width Range (Adjustable).....No Charge
- 2024W 20" to 24" Seat Width Range (Adjustable).....No Charge

SEAT WIDTH SETTING

- W12 12" Seat Width (not available w/Contoura)No Charge
- W13 13" Seat Width (not available w/Contoura)No Charge
- W14 14" Seat Width (not available w/Contoura)No Charge
- W15 15" Seat Width (not available w/Contoura)No Charge
- W16 16" Seat WidthNo Charge
- W17 17" Seat Width (not available w/Contoura)No Charge
- W18 18" Seat WidthNo Charge
- W19 19" Seat Width (not available w/Contoura)No Charge
- W20 20" Seat WidthNo Charge
- W21 21" Seat Width (not available w/Contoura)No Charge
- W22 22" Seat WidthNo Charge
- W23 23" Seat Width (not available w/Contoura)No Charge
- W24 24" Seat WidthNo Charge

SEAT DEPTH RANGE

CURVED ABS & CONVENTIONAL

- 1218D 12" to 18" Seat Depth Range (Adjustable).....No Charge
- 1619D 16" to 19" Seat Depth Range (Adjustable)No Charge
- 1922D 19" to 22" Seat Depth Range (Adjustable)No Charge

CONTOURA

- 1619D 16" to 19" Seat Depth Range (Adjustable)No Charge
- 1922D 19" to 22" Seat Depth Range (Adjustable)No Charge

CURVED ABS, CONVENTIONAL & CONTOURA SEAT DEPTH SETTING

- D12 12" Seat DepthNo Charge
- D13 13" Seat DepthNo Charge
- D14 14" Seat DepthNo Charge
- D15 15" Seat DepthNo Charge
- D16 16" Seat DepthNo Charge
- D17 17" Seat DepthNo Charge
- D18 18" Seat DepthNo Charge
- D19 19" Seat DepthNo Charge
- D20 20" Seat DepthNo Charge
- D21 21" Seat DepthNo Charge
- D22 22" Seat DepthNo Charge

FOOTNOTES

1. Nominal measurement at front edge of seat pan to floor on 18" deep seat.
2. Selected Seat Depth on Contoura assumes 1" of back cushion compression. for example, when choosing an 18" depth the measurement from the front of the back cushion to the front of the seat pan equals 17". The Contoura back cushion is 2" thick includes leatherette seat cushion. Please refer to page 6 for seat size offerings.
3. For 12"-15" seat widths and depths; must select Conventional Back or BPO BPO standard with Push Handles.
4. Push handles will add an additional 3.5" to back-cane height.
5. Seat Widths 12"-15" and Seat Depths 12"-15" are only available with Back Post Only, (BPO) or Conventional Seating. For users under 150 lbs
6. Available only on 400 lb Package.

NOTE: All specifications and dimensions are approximate.

Invacare® Formula™ Tilt Only
Accessories and Options

SEAT

BACK HEIGHT

CURVED ABS OR CONVENTIONAL BACK

- TSBH16 Back Height 16".....No Charge
- TSBH17 Back Height 17".....No Charge
- TSBH18 Back Height 18".....No Charge
- TSBH19 Back Height 19".....No Charge
- TSBH20 Back Height 20".....No Charge
- TSBH21 Back Height 21".....No Charge
- TSBH22 Back Height 22".....No Charge
- TSBH23 Back Height 23".....No Charge
- TSBH24 Back Height 24".....No Charge

CONTOURA BACK

- TSBH25 Back Height 25".....No Charge
- TSBH295 Back Height 29.5".....No Charge

Back Height Chart	
Back Type	Back Height
CTS - Contoura	25" & 29"
BPO - Back Posts Only	16" - 24"
CURVABS - Curved ABS	16" - 24"
Conventional Back	16" - 20" on seat widths of 12" - 15" 20" - 26" on seat widths of 16" - 24"

Note: Back Heights of 21" - 26" not available with seat widths of 12" - 15"

HEADREST TYPE (OPTIONAL)

- TAR127700 Curved Headrest (E0955).....210.00
- TAR27700HW Removable & Adj Hardware (E1028).....210.00
- TAR128000 Two-Step Headrest (E0955).....210.00
- TAR128000HW Removable & Adj Hardware (E1028).....210.00
- T0022400 Large Headrest (E0955).....210.00
- TAR22400HW Removable & Adj Hardware (E1028).....210.00
- T0022600 Small Headrest (E0955).....210.00
- TAR22600HW Removable & Adj Hardware (E1028).....210.00

STEALTH HEADRESTS

- CP175 14" Multi Axis Headrest (E0955).....225.00
- CP175HW Removable & Adj Hardware (E1028).....225.00
- CP275 10" Multi Axis Headrest (E0955).....225.00
- CP275HW Removable & Adj Hardware (E1028).....225.00
- CP180 14" TWB Headrest (E0955).....225.00
- CP180HW Removable & Adj Hardware (E1028).....225.00
- CP280 10" TWB Headrest (E0955).....225.00
- CP280HW Removable & Adj Hardware (E1028).....225.00

ADJUSTABLE ARM TYPE (PAIR)

- 19F Full Length Adjustable Height.....No Charge
- 29F Desk Length Adjustable Height.....No Charge

MIXED ADJUSTABLE ARM TYPE (EACH)

- 19FR Full Length Adjustable Height Right.....No Charge
- 29FR Desk Length Adjustable Height Right.....No Charge
- 19FL Full Length Adjustable Height Left.....No Charge
- 29FL Desk Length Adjustable Height Left.....No Charge

OTTO BOCK™ ARM PADS (PAIR)

- OB20400 One Piece Channel Armpads (E2209).....230.00
- OB20500 Large Forearm Pads.....235.00
- OB20600 Medium Forearm Pads.....235.00
- OB20700 Small Forearm Pads.....235.00

OTTO BOCK HAND PADS (PAIR)

- OB23000 Large Flat Hand Pads.....100.00
- OB23300 Medium Flat Hand Pads.....100.00

LATERALS (PAIR)

3 CONVENTIONAL BACK

- TMSLM Swingaway Laterals Medium (E0956).....210.00
- TMSLMHW Removable & Adjustable Hardware (E1028).....210.00
- TMSLL Swingaway Laterals Large (E0956).....210.00
- TMSLLHW Removable & Adjustable Hardware (E1028).....210.00

3 STEALTH LATERALS FOR CONVENTIONAL BACK

- TWB-RPR Adj Swing Round Lat. Med (E0956).....320.00
- RPRHW Removable & Adjustable Hardware (E1028).....320.00
- TWB-RPRL Adj Swing Round Lat Lg (E0956).....320.00
- RPRLHW Removable & Adjustable Hardware (E1028).....320.00

4 CONTOURA SEATING

- TCSLM Swingaway Contoura Lat Med (E0956).....210.00
- TCSLMHW Removable & Adjustable Hardware (E1028).....210.00
- TCSSL Swingaway Contoura Lat Lg (E0956).....210.00
- TCSSLHW Removable & Adjustable Hardware (E1028).....210.00

4 STEALTH LATERALS FOR CONTOURA

- TWB-CPR Adj Swing Contoura Lat. Med (E0956).....320.00
- CPRHW Removable & Adjustable Hardware (E1028).....320.00
- TWB-CPRL Adj Swing Contoura Lat Lg (E0956).....320.00
- CPRLHW Removable & Adjustable Hardware (E1028).....320.00

SEAT POSITIONING STRAPS

- ⁵ 1311BK Airline Buckle Seat Pos Strap.....50.00
- 1515 Push Button Style Seat Pos Strap.....**Standard**

CHEST POSITIONING STRAPS

- 7311BK Airline Buckle Chest Pos Strap (E0960).....100.00
- 9515 Push Button Style Chest Pos Strap (E0960).....100.00
- 7321BK Hook & Loop Chest Pos Strap (E0960).....100.00

FOOTNOTES

1. Not available with seat widths 12" - 15".
2. Not compatible with quadlink retractable joystick mount.
3. For use with channel canes only.
4. For use with Contoura seating only.
5. Not available on seat widths 12" - 15".

NOTE: All specifications and dimensions are approximate.

Invacare® Formula™ Tilt Only

Accessories and Options

ELECTRONICS

MK6i CONTROLLER

- MK690ACC Controller w/ Actuator ControlStandard

STANDARD DRIVER CONTROL

- SPJAP SPJ™+ with Power Seating Switch.....Standard

ACTUATOR OPERATION THROUGH SPJ+

- SPJAPM Actuator Control through SPJ+ Joystick (E2310)1,300.00

EXPANDABLE CONTROLLER SYSTEMS

- MPJM6U MPJ™+ Multiple Drive Joystick (E2377).....750.00
- PSFM6 Personal Joystick, Inductive on Front (E2399).....1,130.00
- PSRM6 Personal Joystick, Inductive on Rear (E2399).....1,130.00
- PBOD Push Button On/Off w/ PB Drive Select.....No Charge
- PBSS Push Button On/Off w/ Speed Select.....No Charge
- OR On/Off Mounted Right.....No Charge
- OL On/Off Mounted Left.....No Charge
- ASLRDYM6 MK6 Display Only
(No Driver Control) (E2377).....750.00

**FOR YOUR ASL DRIVE SYSTEM NEEDS
PLEASE REFER TO THE FOLLOWING:
ASL SINGLE INVOICE PRICE LIST FORM # 02-097.**

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS

- I500M6 Prop RIM Head Ctrlr w/ MK6i Display (E2327)2,875.00
- SMHD Small Headrest.....No Charge
- LGHD Large Headrest.....No Charge
- I558MM6 Compact Joystick w/ Display (E2373)1,130.00

ALTERNATIVE DIGITAL CONTROL SYSTEMS

- SNPM6 Sip-N-Puff/Digital Interface Box (E2325).....1,500.00
- 5018M6 Wafer Board & Digital Interface Box (E2322)1,550.00
- 5020M6 Mini Tash Joystick & Digital Interface Box (E2321)1,750.00

SIP N PUFF KIT

- PKG32666 Therafin® Sip-N-Puff Breath Tube Kit (E2326).....380.00

ADDITIONAL CONTROL CHOICES

- I558M6 Compact Joystick - Less Std Mount770.00
- 5018 Wafer Board.....310.00
- 5020 Mini Tash Joystick.....405.00
- PACM6 Proportional Attendant Control.....850.00
- I552M Digital Attendant Control850.00

JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retract. Joystick Mnt MPJ+ (E1028)279.00
- QLA Quad Link Retract. Joystick Mnt SPJ+ (E1028).....279.00
- ARM250 Stealth Height Adjustable Joystick Mnt (E1028).....210.00
- LEFT Left Handed Mounted Joystick.....No Charge
- RIMHW RIM Hardware for I500 (E1028)219.00
- PKG32669 Complete Bib Assembly (Tash/I558).....265.00
- GATMPJM6 Gatlin Midline for MPJ (E1028).....600.00
- GAT1812 Gatlin Midline for I812 (E1028).....600.00
- RCM Rear Cane Mount for Compact Joystick.....395.00

ELECTRONIC ACCESSORIES

- I813M6 MK6i Programmer w/Pro Memory Card.....475.00
- I813SD Pro Memory Card (USB Ready).....90.00
- I560 "T" Handle Flexible Joystick Ext (E2323)77.00
- I561 Straight Handle Flexible Joystick Ext (E2323).....77.00
- I826 Chin Cup (E2324).....48.00
- A24VPS 24V MK6 Auxiliary Power Source200.00
- AUX12M6 Auxiliary Module for 1 & 2795.00
- AUX34M6 Auxiliary Module for 3 & 4795.00
- 2PB Dual Push Button.....147.00
- 2SW Dual Toggle172.00

POWER SEATING ACCESSORIES

- SANODE One Actuator Op thru Driver Control (E2310).....1,300.00
- 4SWB 4-Way Switch BoxStandard
(4WSB is standard with two or more actuators)
- S4WSB Multiple Interface Interface Box (E2311)2,600.00

For all systems with two or more actuators you must choose 4WSB at no charge or the S4WSB at an upcharge.

S4WSB Replaces the Four Way Switch Box and allows multiple actuator operation through the driver control.

FWT or QPB is required when choosing the 4WSB or S4WSB.

MULTIPLE ACTUATOR SWITCH OPTIONS

- FWT Four Way Toggle.....No Charge
- QPB Four Push ButtonsNo Charge

ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Switch - Black.....65.00
- ASL304 Light Wobble Switch w /611 Mount.....320.00
- ASL202J Fiber Optic Switch.....610.00
- ASL208J Adjustable Proximity Switch550.00
- ASL504JM6 Remote Emergency Stop Switch.....535.00
- ASL300A Egg Switch w/ 611 Mount.....300.00

FOOTNOTES

1. Allows actuator operation through the Driver Control. Includes port for the Attendant Switch.
2. TOD-SS: Toggle On/Off Drive Select, Speed Potentiometer Standard.
3. Drive Select through mode switch or mode switch and joystick Commands.
4. PSRM6 Only.
5. Display Only, No Driver Control..Allows choice from Order Form.
6. Must also use either a Multiple Drive Joystick or MK6 Display.
7. Requires SNPM6 Sip-N-Puff Interface Box.
8. QLAM6 includes ARM250 Standard.
9. Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options.
10. Required for Tilt Operation through Driver Control. Single Actuator Systems Only. Not compatible with SPJ+ Joystick. The Egg Switch with the 611 Mount does not come std. when ordering a SANODE.
11. Not available with SPJ+ Joystick. Must add MK6 Auxiliary Power Source

NOTE: All specifications and dimensions are approximate.

Invacare® Formula™ Tilt Only

Accessories and Options

MISCELLANEOUS

FRONT RIGGING OPTIONS

BATTERY TRAY TYPE

- 22TRY 22NF Style Tray**Standard**
- 24TRY 24 Group Style TrayNo Charge
- 22TRY3 22NF Style Tray for (3) BatteriesNo Charge
- 22NFTRY2BAT 22NF Style Tray for (2) BatteriesNo Charge

BATTERIES

- 22NFBATTERY 22NF Battery (On Chair) (E2361)300.00
- 24BATTERY 24 Gel Battery (On Chair) (E2363)410.00
- 22NFBATT3 22NF Battery (3) Batteries450.00

VENT TRAY

- 4 VENT Articulating Vent Tray (E1030)1,275.00

MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks45.00

FOOTNOTES

1. 22TRY Base Width 24": 24TRY and 22TRY3 Base Width 25.5".
2. Required for ventilator tray.
3. Allows for 2 Batteries in a 3 Battery Tray.
4. Not available with Contoura Seating.
5. All front riggings include a standard footplate.
6. Allows 60 degree to 90 degree static knee angle positioning with a drop back lever for easier transferring. Extension range 12.5" to 15".
7. Allows 60 degree to 90 degree static knee angle positioning with a drop back lever for easier transferring.
8. Standard with impact guards. Not available with widths 19" to 24" and 16" seat depth.
9. Not available on JR seat widths (1216W).
10. Not available with SPJ+ joystick.
11. Comes standard with a footboard, footplates.
12. Individual Footplate option for Power Center Mount. This option is only available on MPCTMT.
13. Independent powered elevating legrests including elevating and articulating legrests with angle adjustable footplates. Not available with the 16" or 17" seat depths. Extension range 13"-20". Sold as pair.
14. One piece footboard for Center Mount (CTMT-F) Riggings.

5 FRONT RIGGING

CENTER MOUNT

- 6 CTMT-F Center Mount Adj Knee Angle (12.25"-15")220.00
- 7 CTMT-FJ Center Mount Adjustable (9.25" - 12")220.00

CENTER PIVOT STYLE (PAIR)

- 70NHD 70° Swingaway Footrests**Standard**
- 70STAPER 70° Tapered FootrestsNo Charge
- AT5544 Manual Elevating Legrests (E0990)255.00

PIN STYLE-HEMI (PAIR)

- PHW93 70° Swingaway FootrestsNo Charge
- 9 PHAL4A Manual Articulating/Elevating Legrests (K0053)225.00

POWER ELEVATING LEGRESTS-ARTICULATING

- 10,11 MPCTMT Power Center Mount Foot Platform (K0108)2,995.00
- 10 MPCTMT-L Pwr Center Mount Foot Platform Long (K0108)2,995.00
- 12 IFP Individual Footplates275.00
- 10,13 ELRPW Power ELR w/ Articulation 20" Ext (E1010)1,895.00
- 10,13 ELRPW22 Power ELR w /Articulation 22" Ext (E1010)1,995.00

FOOTPLATE OPTIONS (PAIR)

- 1651 Composite Flip-Up FootplatesNo Charge
- 1350 Extra Large Aluminum FootplatesNo Charge
- AT5543 Adjustable Angle Flip-Up Footplates165.00
- 14 FBCTMT Footboard for Center Mount150.00

FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap (E0038)27.00
- 1600BK Heel Loops w/Ankle Straps (E0951)56.00
- LPT2 Longer Pivot Slide Tube for 135083.00
- IMPACTG Impact Guards51.00
- ALPT Longer Pivot and Slide Tube84.00
- ALPT4 Longer Pivot and Slide Tube 4"84.00

NOTE: All specifications and dimensions are approximate.

Motion Power Center Mount Knee to Heel			
STF Pick	STF at 0 deg	KTH w/IFP on MPCTMT	KTH w/Footboard
Low	17.5"	10.5"-14.5"	12"-13.5" MPCTMT only
Med	18.5"	11.5"-15.5"	12"-14.5" MPCTMT only
High	19.5"	11.5"-16.5"	12"-15.5" MPCTMT or MPCTMT-L
STF Pick	STF at 5 deg	KTH w/IFP on MPCTMT	KTH w/Footboard
Low	18.75"	11.5"-15.5"	12"-14" MPCTMT only
Med	19.75"	11.5"-16.5"	14"-16" MPCTMT only
High	20.75"	11.5"-17.5"	14"-17" MPCTMT or MPCTMT-L

Invacare® Formula™ Tilt Only

Accessories and Options

BATTERIES

22NF batteries with terminal configuration (positive on the left and negative on the right) **MUST** be used. 22NF batteries that have the reverse terminal configuration **MUST** not be used. Terminals **MUST** have a cross hole for proper battery connection.

GP24 batteries with terminal configuration (negative on the left and positive on the right) **MUST** be used. GP24 batteries that have the reverse terminal configuration **MUST** not be used. Terminals **MUST** have a cross hole for proper battery connection. See Owner's Manual, part number 1114809. These recommendations **MUST** be followed otherwise injury and damage may occur.

(BATTERIES ARE NOT INCLUDED)

Contoura Seat Sizes

	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X
W24	*	*	X	X	X	X	X

*Not available with Contoura Seat Cushion.

Back Heights 25" or 29.5"

Memory Card - BASIC

- Standard with All Expandable Control Systems (MPJ+, PSR+, PSF+, MK6i Display).
- Used to Back Up/Restore Programmed Settings for Only One Chair.... May be given to user or kept with provider's files for safe keeping.
- Does not Contain Advanced Diagnostics/Help Screens/File Structure.
- Not compatible with SPJ+ Joysticks.

Memory Card - PROFESSIONAL

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains Advanced Diagnostics, Help Screens, Software Updates, File Storage/Retrieval.
- Used to back up multiple chairs and programming settings.
- Intended for use by qualified providers only.
- Not compatible with SPJ+ Joysticks.

**FOR YOUR SEATING NEEDS
PLEASE REFER TO THE FOLLOWING**

**REHAB SEATING & POSITIONING PRODUCTS
SEE PRICE LIST FORM # 07-027.**

Formula_{CG} Tilt Only Features and Specifications

- ♦ Weight capacity of 300 lbs or 400 lbs (with heavy-duty 4-pole motors).
- ♦ Curved ABS or Contoura style seating.
- ♦ Adjustable seat frame.

After Market Backrests

The following companies also have hardware to make their backrests work with our new canes:

- ♦ AEL - planar seating system
- ♦ Freedom Designs® - planar seating system
- ♦ Canyon - planar seating system

Conventional Seat Sizes

Seat Widths	12" - 24"
Seat Depths	12" - 22"
Back Heights	16" - 26" 16" to 20" for 12"-15" seat widths only)

Tech Tips

Operating Powered Seating through the Driver Control
All chairs require adding or replacing one component to allow actuator operation through the driver control.

SPJ+ with PSS:

Replace joystick with SPJ+ with ACC (Actuator Control) -E2310

Single Actuator Systems with MPJ+, PSR+, PSF+ MK6i Display:

Add a SANODE (Single Actuator Node) to the system - E2310

Multiple Actuator Systems:

Replace the standard 4 Way Switch box with a Multiple Actuator Interface Box - E2311

Curved ABS

Seat Widths	16" - 24"
Seat Depths	16" - 22"
Back Heights	16" - 24"

Invacare Corporation

www.invacare.com



USA

One Invacare Way
Elyria, OH USA
44035-4190
800-333-6900

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

Therafin is a registered trademark of W/R Medical Electronics Co.

Otto Bock is a trademark of Otto Bock Orthopaedische Industries GMBH & Co.

Products designated by ASL are Manufactured by Adaptive Switch Labs, Inc.

Specifications and prices are subject to change without notification.

©2007 Invacare Corporation

Form No. 06-039 Rev 100107