

# Invacare® Storm Series® Torque™ SE Power Wheelchair



**Price List and Order Form** Effective October 9, 2006 Revised October 1, 2007

For ease of ordering, contact Customer Service Toll Free at:  
 800-333-6900 Fax: 800-678-4682  
 To order Literature: 800-828-6282

www.invacare.com

Company Name \_\_\_\_\_ Account # \_\_\_\_\_

Phone # \_\_\_\_\_ P.O.# \_\_\_\_\_ Date \_\_\_\_\_

Ship To Name & Address \_\_\_\_\_

### Special Note

**HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.**

### FRAME SELECTION

#### CHAIR TYPE

- 3GTQSE Torque SE w/MK6i™ .....5,288.00
- 3GTQSEBASE Torque SE Base Only w/MK6i.....5,048.00
- 3GSEAT Series Seat Only .....740.00

#### 3 TRANSPORT TIE-DOWN

- TRRO Transport Ready Option .....300.00
- TRBKTS Wheelchair Transport Brackets .....250.00

#### USER WEIGHT LIMITS

- U300 User Weight <300 lbs.....No Charge

#### FRAME TYPE

- SFS Short Base.....**Standard**
- LFS Long Base .....No Charge

#### 5 ADJUSTABLE SEAT WIDTH RANGE

- 1216W 12" to 16" Seat Width Range.....No Charge
- 1620W 16" to 20" Seat Width Range.....No Charge
- 2022W 20" to 22" Seat Width Range.....No Charge

#### 6 FRAME WIDTH

- W12 12" Wide.....110.00
- W13 13" Wide.....110.00
- W14 14" Wide.....110.00
- W15 15" Wide.....No Charge
- W16 16" Wide.....No Charge
- W17 17" Wide.....No Charge
- W18 18" Wide.....No Charge
- W19 19" Wide.....No Charge
- W20 20" Wide.....329.00
- W21 21" Wide.....329.00
- W22 22" Wide.....329.00

#### 4 VAN SEAT WIDTH

- VW18 18" Wide.....No Charge
- VW20 20" Wide.....329.00
- VW22 22" Wide.....329.00

#### DRIVE TYPE OPTION

- TRUTK TrueTrack upgrade.....1,899.00
- PERF Performance 4-Pole Motor.....No Charge

#### FOOTNOTES

1. Torque SE comes standard with 300 lb user weight capacity, 6.5 MPH max speed. Performance 4-pole motor package. MK6i SPJ+ joystick, MK690 Controller. Base only comes standard with MK690ACC Controller and SPJ+ w/PSS Joystick. Base only does not come standard with Seat Positioning Strap.
2. When ordering a replacement seat and changing depth to 20" or deeper, order new arms.
3. Must select Group 24 battery boxes for TRRO or TRBKTS.
4. Not available on seat depths of 19" to 22". Manual reclining Back is not available. Van Seat Depth set at the following: 18" wide set at 16" deep, 20" wide set at 18" deep, 22" wide set at 20" deep
5. Jr. seat widths require a straight rear upright and have a user weight limit of 150 lbs. Manual reclining back not available. Additional charge may apply with widths/depths less than 16".
6. Contoura™ Seating available in even widths from 16"-22".
7. TrueTrack driving technology, which enables driver to track straighter regardless of the terrain, includes Monroe® Spring Shock Suspension, <300 lb user weight limit and 7 MPH max speed.

**NOTE: All specifications and dimensions are approximate**

### \* TRANSPORT TIE DOWN

TRRO includes four factory installed transport brackets and a wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNA WC Vol I Section 19 Frontal Impact Test requirements for wheelchairs with a 168 lb crash dummy, which corresponds to a person with a weight of 114 lbs to 209 lbs.

TRBKTS includes four factory installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

\* Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

### TRRO FEATURES FOR 3G TORQUE SE

Features Required for TRRO	Optional Features for TRRO
<ul style="list-style-type: none"> <li>♦ Torque SE Base &amp; Seat</li> <li>♦ 300 lb user weight limit</li> <li>♦ Short or Long Base Frame</li> <li>♦ Adjustable Rehab Seat with Adjustable Seat Tilt &amp; Back Angle</li> <li>♦ Any Seat Width</li> <li>♦ Any Seat Depth</li> <li>♦ Low Seat-to-Floor Height</li> <li>♦ Sectional or Straight Rear Uprights</li> <li>♦ 18" - 24" ASBA Back Height</li> <li>♦ 90° or 95° Back Angle</li> <li>♦ Any Arm Type</li> <li>♦ 8" x 2.25" Semi Pneumatic Casters</li> <li>♦ Suspension Shocks</li> <li>♦ 14" x 3" Tires with Foam Filled Inserts</li> <li>♦ Rear Wheel Standard Location</li> <li>♦ Any Color Frame Finish</li> <li>♦ Black Nylon or Black Vinyl Upholstery</li> <li>♦ Center Pivot Style Bracket</li> <li>♦ 60° Swingaway, 70° Swingaway or 70° Tapered Footrests</li> <li>♦ Group 24 Battery Box</li> </ul>	<ul style="list-style-type: none"> <li>♦ TrueTrack Upgrade</li> <li>♦ Fender Option</li> <li>♦ Flip-Up Composite Footplates, Aluminum Footplates-Deluxe, Adjustable Angle Flip Up Footplates</li> <li>♦ Calf Strap</li> <li>♦ Impact Guards</li> <li>♦ All PSR / PSF Joysticks</li> <li>♦ Left Hand Mounted Joystick</li> <li>♦ "T" Handle Flexible Joystick Ext</li> <li>♦ Straight Handle Flexible Joystick Ext</li> <li>♦ MK6i Remote Programmer</li> <li>♦ SPJ™+ or MPJ™+ Joysticks</li> <li>♦ Reinforced Back Upholstery</li> <li>♦ Group 24 Batteries</li> </ul>

Please note that all TRRO chairs are set at a 10° seat angle

# Invacare® Storm Series® Torque™ SE

## Accessories and Options

### SEAT SELECTION

<input type="checkbox"/>	ADJASBA	Adj Rehab Seat w /Adj Seat Tilt & Back Angle.....	<b>Standard</b>
<input type="checkbox"/>	VAN	Van Style Captain Seat.....	100.00
<input type="checkbox"/>	ADJRECL	Manual Recline for Adj ASBA (E1226).....	744.00
<input type="checkbox"/>	CTS	Contoura Seating.....	395.00
<input type="checkbox"/>	OMCU	Omit Contoura Seat Cushion.....	Less 20.00

### ADJUSTABLE SEAT DEPTH RANGE

<input type="checkbox"/>	1218D	12" to 18" Seat Depth Range.....	No Charge
<input type="checkbox"/>	1619D	16" to 19" Seat Depth Range.....	No Charge
<input type="checkbox"/>	1922D	19" to 22" Seat Depth Range.....	No Charge

### SEAT DEPTH SETTING

<input type="checkbox"/>	D12	Seat Depth 12" (Short Base).....	110.00
<input type="checkbox"/>	D13	Seat Depth 13" (Short Base).....	110.00
<input type="checkbox"/>	D14	Seat Depth 14" (Short Base).....	110.00
<input type="checkbox"/>	D15	Seat Depth 15" (Short Base).....	No Charge
<input type="checkbox"/>	D16	Seat Depth 16" (Short Base).....	No Charge
<input type="checkbox"/>	D17	Seat Depth 17" (Short Base).....	No Charge
<input type="checkbox"/>	D18	Seat Depth 18" (Short Base).....	No Charge
<input type="checkbox"/>	D19	Seat Depth 19" (Long Base).....	No Charge
<input type="checkbox"/>	D20	Seat Depth 20" (Long Base).....	529.00
<input type="checkbox"/>	D21	Seat Depth 21" (Long Base).....	529.00
<input type="checkbox"/>	D22	Seat Depth 22" (Long Base).....	529.00

### ADJASBA SEAT TO FLOOR HEIGHT

<input type="checkbox"/>	LOW	Seat To Floor - Low.....	No Charge
<input type="checkbox"/>	MED	Seat To Floor - Med.....	No Charge
<input type="checkbox"/>	TALL	Seat To Floor - Tall.....	No Charge

### VAN SEAT STYLE

<input type="checkbox"/>	SEMIR	Semi-Recline.....	No Charge
<input type="checkbox"/>	FULLR	Full-Recline (E1226).....	469.00

### BACK SELECTION

#### ADJASBA BACK TYPE

<input type="checkbox"/>	SECT	Sectional Rear Uprights.....	No Charge
<input type="checkbox"/>	SRUP	Straight Rear Uprights.....	59.00
<input type="checkbox"/>	9720P	No Push Handles (Straight Canes).....	82.00
<input type="checkbox"/>	L866I	Omit Rear Uprights.....	Less 20.00

### RECLINER BACK HEIGHT

<input type="checkbox"/>	BHR18	Back Height 18".....	No Charge
<input type="checkbox"/>	BHR20	Back Height 20".....	No Charge
<input type="checkbox"/>	BHR22	Back Height 22".....	No Charge
<input type="checkbox"/>	BHR24	Back Height 24".....	No Charge
<input type="checkbox"/>	BHR26	Back Height 26".....	No Charge

### BACK HEIGHT

#### ADJASBA

<input type="checkbox"/>	BH12	Back Height 12".....	No Charge
<input type="checkbox"/>	BH13	Back Height 13".....	No Charge
<input type="checkbox"/>	BH14	Back Height 14".....	No Charge
<input type="checkbox"/>	BH15	Back Height 15".....	No Charge
<input type="checkbox"/>	BH16	Back Height 16".....	No Charge
<input type="checkbox"/>	BH17	Back Height 17".....	No Charge
<input type="checkbox"/>	BH18	Back Height 18".....	No Charge
<input type="checkbox"/>	BH19	Back Height 19".....	No Charge
<input type="checkbox"/>	BH20	Back Height 20".....	No Charge
<input type="checkbox"/>	BH21	Back Height 21".....	No Charge
<input type="checkbox"/>	BH22	Back Height 22".....	No Charge
<input type="checkbox"/>	BH23	Back Height 23".....	No Charge
<input type="checkbox"/>	BH24	Back Height 24".....	No Charge

#### CONTOURA

<input type="checkbox"/>	TSBH25	Back Height 25".....	No Charge
<input type="checkbox"/>	TSBH295	Back Height 29.5".....	No Charge

### BACK ANGLE DEGREE

<input type="checkbox"/>	85BA	85°.....	No Charge
<input type="checkbox"/>	90BA	90°.....	No Charge
<input type="checkbox"/>	95BA	95°.....	No Charge
<input type="checkbox"/>	100BA	100°.....	No Charge
<input type="checkbox"/>	105BA	105°.....	No Charge

### BACK ACCESSORIES

<input type="checkbox"/>	L860I	Omit Back Upholstery.....	No Charge
<input type="checkbox"/>	RBP	Reinforced Back Upholstery.....	No Charge

### ARM SELECTION

#### ARM TYPE

<input type="checkbox"/>	18F	Full Length Fixed Hgt Conv Flip Up.....	No Charge
<input type="checkbox"/>	28F	Desk Length Fixed Hgt Conv Flip Up.....	<b>Standard</b>
<input type="checkbox"/>	19F	Full Length Adj Hgt Conv Flip Up (E0973).....	230.00
<input type="checkbox"/>	29F	Desk Length Adj Hgt Conv Flip Up (E0973).....	230.00

### MK6i ELECTRONICS

#### MK6i CONTROLLERS

<input type="checkbox"/>	MK690	90 Amp Controller.....	<b>Standard</b>
<input type="checkbox"/>	MK6TT	TrueTrack Controller.....	No Charge

#### STANDARD DRIVER CONTROL

<input type="checkbox"/>	SPJM6	SPJ+ Single Drive Joystick.....	<b>Standard</b>
--------------------------	-------	---------------------------------	-----------------

#### BASE ONLY SPJ+ OPTIONS

<input type="checkbox"/>	SPJAP	SPJ+ with Powered Seating Switch.....	<b>Standard</b>
<input type="checkbox"/>	SPJAPM	Single Actuator Operation through SPJ+ (E2310).....	1,300.00

### EXPANDABLE CONTROLLER SYSTEMS

<input type="checkbox"/>	MPJM6U	MPJ+ Multiple Drive Joystick (E2377).....	750.00
<input type="checkbox"/>	PSFM6	Personal Joystick, Inductive on Front (E2399).....	1,130.00
<input type="checkbox"/>	PSRM6	Personal Joystick, Inductive on Rear (E2399).....	1,130.00
<input type="checkbox"/>	PBOD	Push Button On/Off w/PB Drive Select.....	No Charge
<input type="checkbox"/>	PBSS	Push Button On/Off w/Speed Select.....	No Charge
<input type="checkbox"/>	OR	On/Off Mounted Right.....	No Charge
<input type="checkbox"/>	OL	On/Off Mounted Left.....	No Charge
<input type="checkbox"/>	ASLRDYM6	MK6 Display Only (No Driver Control) (E2377).....	750.00

#### FOR YOUR ASL DRIVE SYSTEM NEEDS

#### PLEASE REFER TO THE FOLLOWING:

#### ASL SINGLE INVOICE PRICE LIST FORM # 02-097.

<input type="checkbox"/>	1500M6	Prop RIM Head Control w/ Display (E2327).....	2,875.00
<input type="checkbox"/>	SMHD	Small Headrest.....	No Charge
<input type="checkbox"/>	LGHD	Large Headrest.....	No Charge
<input type="checkbox"/>	1558MM6	Compact Joystick w/ Display (E2373).....	1,130.00

### FOOTNOTES

- Selected seat depth on Contoura assumes 1" of back cushion compression. For example: when choosing an 18" depth the measurement from the front of the back cushion to the front of the seat pan equals 17". The Contoura back cushion is 2" thick. Please see page 5 for seat size offerings. Contoura seating available in depths of 16" to 22".
- Must select 12"-18" depth range (1218D) when ordering the 12"-16" seat width range (1216W).
- Manual reclining back not available.
- Must use long base when selecting 19" - 22" depth.
- Based on 5° tilt and 16" seat depth on ADJASBA add .25" for additional 2" of seat depth. See page 5 for nominal Seat-To-Floor Heights.
- Used only with manual reclining seat selection.
- All arms are flip-up and detachable.
- Not available w /seat widths 12" - 15" must use adjustable height.
- Base Only Systems for Multiple Actuator Powered Seating Systems will require MPJ+, PSR+, PSF+, or MK6i Display & Alternative Driver Control. Includes Egg Switch
- Allows Actuator Operation through Driver Control.
- TOD-SS: Toggle On/Off Drive Select, Speed Potentiometer Standard.
- PSRM6 Only.
- Drive select through mode switch, or mode switch and Joystick command.
- Display only, allows choice from order form.

**NOTE: All specifications and dimensions are approximate.**

# Invacare® Storm Series® Torque™ SE

## Accessories and Options

### ALTERNATIVE DIGITAL CONTROL SYSTEMS

- SNPM6 Sip-N-Puff Digital Interface Box (E2325) .....1,500.00
- 5018M6 Wafer Board & Digital Interface Box (E2322).....1,550.00
- 5020M6 Mini Tash Joystick & Digital Interface Box (E2321) .1,750.00

### SIP-N-PUFF KIT

- PKG32666 Therafin® Sip-N-Puff Breath Tube Kit (E2326)....380.00

### ADDITIONAL CONTROL CHOICES

- 1558M6 Compact Joystick - Less Std Mount.....770.00
- 5018 Wafer Board .....310.00
- 5020 Mini Tash Joystick .....405.00
- PACM6 Proportional Attendant Control .....850.00
- 1552M Digital Attendant Control .....850.00

### JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retract. Joystick Mnt MPJ+ (E1028) .....279.00
- QLA Quad Link Retract. Joystick Mnt SPJ+ (E1028).....279.00
- QLM6 Quad Link Retract. Jystck Mnt MPJ+ Van (E1028)...230.00
- QL Quad Link Retract. Jystck Mnt SPJ+ Van (E1028)....230.00
- ARM250 Stealth Height Adjustable Joystick Mount .....210.00
- LEFT Left Handed Mounted Joystick.....No Charge
- RIMHW RIM Hardware for 1500 (E1028) .....219.00
- PKG32669 Complete Bib Assembly (Tash/1558) .....265.00
- GATMPJ6 Gatlin Midline for MPJ+ (E1028).....600.00
- GAT1812 Gatlin Midline for 1812 (E1028) .....600.00
- RCM Rear Cane Mount for Compact Joystick.....395.00

### ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer w/ Pro Memory Card .....475.00
- 1813SD Pro Memory Card (USB Ready) .....90.00
- 1560 "T" Handle Flexible Joystick Ext (E2323) .....77.00
- 1561 Straight Handle Flexible Joystick Ext (E2323) .....77.00
- 1826 Chin Cup (E2324).....48.00
- A24VPS 24V MK6 Auxiliary Power Source .....200.00
- AUX12M6 Auxiliary Module for 1 & 2 .....795.00
- AUX34M6 Auxiliary Module for 3 & 4 .....795.00
- MPY Monoport "Y" Cable for MPJ+ or Display .....15.00
- 2PB Dual Push Button.....147.00
- 2SW Dual Toggle .....172.00

### ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Switch Black.....65.00
- ASL304A Light Wobble Switch w /611 Mount.....320.00
- ASL202J Fiber Optic Switch.....610.00
- ASL208J Adjustable Proximity Switch .....550.00
- ASL504JM6 Remote Emergency Stop Switch.....535.00

### FRAME AND UPHOLSTERY COLOR

#### FRAME FINISH CONTEMPORARY

- 24P Wet Black.....No Charge
- 60P Silver Metallic .....No Charge
- 71P Silver Vein.....No Charge
- 105P Midnight Blue.....No Charge
- 115P Black Prism.....No Charge

#### TOO HOT

- 30P Sunny Yellow .....No Charge
- 61P Electric Red.....No Charge
- 119P Bubblegum Pink .....No Charge
- 120P Grape Madness .....No Charge

#### WAY COOL

- 62P Electric Blue.....No Charge
- 104P Electric Teal .....No Charge
- 121P Lolly Pop Blue .....No Charge
- 122P Cosmic Blue.....No Charge
- 125P Emerald Green.....No Charge
- 127P Grasshopper Green.....No Charge

### UPHOLSTERY TYPE

- U240 Black Nylon Upholstery .....**Standard**
- U67 Black Vinyl.....No Charge
- GC Grey Cloth .....No Charge
- GV Grey Vinyl .....No Charge
- TV Tan Vinyl .....No Charge

### CASTER AND WHEEL SELECTION

#### CASTER OPTIONS

- 1259 6" X 2" Semi-Pneumatic Casters.....No Charge
- 1228-3 8" X 2 1/4" Semi Pneu. Casters Lt Grey ....**Standard**
- 1251-3 8" X 2" Pneumatic Casters .....No Charge
- 1220-3 9" X 2 3/4" Pneumatic Casters .....No Charge
- 1465-3 9" X 2 3/4" Foam Filled Casters .....No Charge

#### 9 CASTER SUSPENSION

- 9213-6 Shock Fork For 6" Casters (E1016).....119.00
- 9213-8 Shock Fork For 8" Casters (E1016).....119.00

#### WHEEL OPTION

- 1430-3 14" X 3" Pneumatic Tires .....No Charge
- 1431-3 14" X 3" Tires w/Foam Filled Inserts .....250.00
- 1440 14" X 4" Pneumatic Tires.....200.00
- 1441 14" X 4" Tires w/Foam Filled Inserts .....320.00

#### 13 FREE WHEEL HUB

- FHUB Free Wheel Hub.....125.00

#### 14 REAR WHEEL LOCATION

- RWF2 Rear Wheel Forward 2".....No Charge
- RWF1 Rear Wheel Forward 1" .....No Charge
- RWS Rear Wheel Standard .....**Standard**

#### 15 FENDER OPTION

- FDR Fender .....320.00

### FOOTNOTES

1. Must also use either a Multiple Drive Joystick or MK6i Display.
2. Requires SNPM6 Sip-N-Puff Interface Box.
3. QLAM6 and QLA include ARM250 standard.
4. Only available with Van Seating.
5. For use with with expandable controller systems to provide additional switch sites.
6. Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options.
7. Not available with SPJ+ Joystick. Must add MK6 Auxiliary Power Source.
8. Used only with Van Style Seat Selection.
9. Not available with TRRO or TRBKTS.
10. Not available with short base. Not available with 9" casters.
- 11 14" x 4" wheels increase overall width of the base by 2.5". Not available with TrueTrack (TRUTK).
12. Does not include Storm Series designed mag wheel or fender.
13. Free wheel hub increases overall width .75. Not needed or available with TRUTK TrueTrack package.
14. This location changes the center-of-gravity and performance. Forward wheel position will make turning and climbing easier, but more prone to "weeely".
15. Standard with TrueTrack (TRUTK) option.

**NOTE: All specifications and dimensions are approximate.**

# Invacare® Storm Series® Torque™ SE

## Accessories and Options

### POSITIONING STRAPS

#### SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Positioning Strap.....50.00
- 1515 Auto Style Buckle Seat Positioning Strap .....**Standard**

#### CHEST POSITIONING STRAPS

- 7311BK Airline Buckle Chest Positioning Strap (E0960)..100.00
- 9515 Auto Style Chest Positioning Strap (E0960).....100.00
- 7321BK Hook & Loop Chest Positioning Strap (E0960)..100.00

### MISC ACCESSORIES

#### MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks .....45.00
- VAN02 02 Holder - Van (E2208).....150.00
- SEATO2 Oxygen Holder (ASBA Seat Only) .....150.00
- 9205-3 Suspension Shocks (E1016) .....400.00
- 9205H-3 Suspension Shocks Heavy Duty .....400.00

### BATTERY SELECTION

#### BATTERY BOXES

- 22NF 22NF Battery Box .....No Charge
- GP24 Group 24 Battery Boxes .....No Charge

#### BATTERIES

- 22NBATTERY 22NF Battery (on chair) (E2361) .....300.00
- 24BATTERY 24Gel Battery (on chair) (E2363).....410.00

### HANGER BRACKET SELECTION

#### HANGER BRACKET TYPE

- ASB Center Pivot Style Bracket .....No Charge
- PSB Pin Style Bracket .....No Charge
- NHB No Hanger Bracket .....No Charge

### FRONT RIGGING OPTIONS

#### FRONT RIGGING

##### CENTER PIVOT

- 60HD 60° Swingaway Footrests.....No Charge
- 70NHD 70° Swingaway Footrests.....**Standard**
- 70STAPER 70° Tapered Footrests.....No Charge
- AT5544 Elevating Legrests (E0990) .....255.00
- 90SWING 90° Swingaway Footrests.....No Charge

##### PIN STYLE

- PAS4A Smartleg® Art. Legrests 17 1/4" - 21 1/2" (K0053).....320.00
- PVW93 Swingaway Footrests. 13 1/2" - 19".....No Charge
- P904A Swingaway Elevating Legrests 16"-19 3/4" (E0990)..255.00
- PAL4A Smartleg Art. Legrests 17 1/4" - 21 1/2" (K0053).....320.00

#### POWER LEGS

- 6.7 ELRPW Power Elevating Legrests w/Articulation (E1010).....1,895.00

#### POWER LEGS CONTROL MODULE

- 9 4WSB Four Way Switch Box.....No Charge
- 10 S4WSB Multiple Actuator Interface Box (E2311) .....2,600.00
  - FWT 4 Way Toggle.....No Charge
  - QPB 4 Push Buttons.....No Charge

**For all systems with two or more actuators you must choose 4WSB at no charge or the S4WSB at an upcharge.**

**S4WSB Replaces the Four Way Switch Box and allows multiple actuator operation through the driver control.**

**FWT or QPB is required when choosing the 4WSB or S4WSB.**

### FOOTPLATE OPTIONS

- 1651 Flip-Up Composite Footplates.....No Charge
- 1350 Extra Large Aluminum Footplates .....No Charge
- 1450 Aluminum Footplates-Deluxe .....No Charge
- AT5043 Flip Up Footplates.....No Charge
- S1263 One Piece Wooden Footboard.....No Charge
- AT5543 Adjustable Angle Flip Up Footplates.....165.00

### FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap (K0038) .....27.00
- 1600BK Heel Loops with Ankle Straps (pair) (E0951).....56.00
- IMPACTG Impact Guards .....51.00
- A9342 "H" Heel Strap (K0039) .....51.00
- LPT2 Longer Pivot Slide Tube For 1350.....83.00
- ALPT Longer Pivot And Slide Tube.....84.00
- ALPT3 3" Extension For 90° Footrest .....31.00
- ALPT4 Longer Pivot And Slide Tube 4" .....84.00

### FOOTNOTES

1. 1311BK not available in widths 12" - 15" wide.
2. Not available with TRUTK package selected.
3. Batteries will be installed by Invacare Corporation at no additional charge. Chair will be shipped complete w/batteries.
4. Standard with impact guards.
5. For seat widths greater than or equal to 16".
6. Independent powered legrests include elevating and articulating legrests with angle adjustable footplates. Extension range of 13" - 20".
7. To operate power legs must choose: 4WSB - 4 Way Switch Box or S4WSB Multiple Actuator Interface Box. Must include with pick either FWT- 4-Way Toggle or QPB - 4 Push Buttons.
8. These options are located in the Configurator under the MK6i Electronics Folder. The Category Header inside the folder is labeled Power Seating Accessories.
9. Required with all multiple actuator systems. Provides a D-9 port for the 4-way toggle or 4 push buttons.
10. Replaces the Four Way Switch Box and allows multiple actuator operation through the driver control. Not compatible with SPJ+ joysticks.

**NOTE: All specifications and dimensions are approximate.**

### SEATING AND POSITIONING

**FOR YOUR SEATING NEEDS  
PLEASE REFER TO THE FOLLOWING**

**REHAB SEATING & POSITIONING PRODUCTS  
SEE PRICE LIST FORM # 07-027.**

# Invacare® Storm Series® Torque™ SE

## Specifications

### Chair Features and Options

- ♦ Independent rear wheel suspension with new rubber element design
- ♦ 6.5 mph standard speed range
- ♦ Programmable performance to meet individuals needs
- ♦ 17.5" seat-to-floor Standard  
19.75" and 21.5" Optional
- ♦ Seat location system allows for balanced weight distribution
- ♦ Continuous tilt seat 0° - 10°
- ♦ Dynamic anti tippers
- ♦ Rigid seat pan
- ♦ Convenient rear battery access
- ♦ 22NF batteries (not supplied) in a single 22NF box  
Single battery connector
- ♦ Arms are flip back/removable
- ♦ Adjustable depth rigging attachment 2" (1" increments)
- ♦ Automatic dual mode battery charger (Lester Brand)
- ♦ Compact size: 24.75" overall base width w/14" x 3" wheels
- ♦ Weight limit - 300 lbs

### Weight Limit 300 Pounds (Standard)

Invacare recommends all wheelchairs be fitted and programmed to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

**NOTE:** All weight limits are inclusive of add-on equipment and seat systems

### Footrest Extension Range

	AT5043	AT5543	1651
<b>60HD</b>	N/A	12.25"-16.50"	13.50"-17.50"
<b>70MFX</b>	6.75"-12.50"	4.50"-10"	N/A
<b>70NHD</b>	N/A	10.75"-15 75"	11 75"-17"
<b>70STAPER</b>	N/A	12 50"-16 50"	13.25"-17.25"
<b>90SWING</b>	4"-9"	5"-8"	N/A

Measured from Seat Pan to Back of Footplate

### Memory Card - BASIC

- Standard with ALL Expandable Control Systems (MPJ+, PSR+, PSF+, MK6i Display).
- Used to Back Up/Restore Programmed Settings for Only One Chair.
- May be given to user or kept with provider's files for safe keeping.
- Does not Contain Advanced Diagnostics/Help Screens/File Structure.
- Not compatible with SPJ+ Joysticks.

### Memory Card - PROFESSIONAL

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains Advanced Diagnostics, Help Screens, Software Updates, File Storage/Retrieval.
- Used to back up multiple chairs and programming settings.
- Intended for use by qualified providers only.
- Not compatible with SPJ+ Joysticks.

### Torque SE Specifications

* Overall Width	24.75" 27.25"	w/14"x3" wheel w/14"x4" wheel
**Overall Length	29.5" 32.5"	Short Base Long Base
***Seat-to-Floor	17.50" Standard 19.75" Optional	Low Medium
Battery	22NF deep cycle; gel/sealed (not supplied)	

Weight w/o batteries 149 lb. (40 - 45) lb. each battery)

Shipping Weight 266 lbs (18" X 18" seat w/Group 24 batteries).

\*Measurement based on wheel base. Seat width is width picked plus 7.5" for the armrest and joystick.

\*\*Measured from headtube to anti-tipper.

\*\*\*Based on 5° tilt and 16" seat depth; add .25" for additional 2" of seat depth

**Note: All specifications and dimensions are approximate.**

### Manual Recline Specifications

Seat Widths	14"-22" (1" increments)		
Seat Depths	16"-22" (1" increments) (Adjustable 16"-19", Fixed 20"-22")		
Back Heights	24" Standard, 18.50", 20", 22", 24", 26"		
Weight Limits	300 lbs with Standard Cylinders		
Headrest	10" Height		
Recline Range	90° to 170° Approximate		
Overall Back Height	Low Seat	Med. Seat	High Seat
{w/o extension}			
Standard Base	41.50"-	43.50"-	45.25"-
	{Measured with 5° Tilt and 10" for Extension}		

### Contoura Seat Sizes

	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X
W24	*	*	X	X	X	X	X

\*Not available with Contoura Seat Cushion.

Back Heights 25" or 29.5"

Invacare Corporation [www.invacare.com](http://www.invacare.com)



**USA**  
One Invacare Way  
Elyria, OH USA  
44035-4190  
800-333-6900

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Therafin is a registered trademark of W/R Medical Electronics Co. Monroe is a registered trademark of Tenneco Automotive, Inc. Products designated by ASL are Manufactured by Adaptive Switch Labs, Inc. Specifications and prices are subject to change without notification.  
©2007 Invacare Corporation

Form No. 06-067 Rev 100107