

Invacare® Formula™ CG Recline for TDX® SP Base



Yes, you can.®

Price List and Order Form Effective March 1, 2007 Revised October 1, 2007

For ease of ordering, contact Customer Service Toll Free at:
 800-333-6900 Fax: 800-678-4682
 To order Literature: 800-828-6282

www.invacare.com

Company Name _____ Account # _____

Phone # _____ P.O.# _____ Date _____

Ship To Name & Address _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

POWERED SEAT BASE

TDXSP-CG TDXSP for Formula CG w/MK6i™5,600.00

Group 3 - K0856

Basic Packages

- * TDXSP-CG w/Formula CG with MK6i
- * 300 lb Weight Capacity
- * 4-Pole Motors.....Standard
- * 1431-3 14" X 3" Tires w/Foam Filled Inserts.....\$250.00
- * Flip Back Adjustable Height Arm.....Standard
- * MK690 90 Amp Controller.....Standard
- * 70NHD 70° Swingaway Footrests.....Standard
- * 8 Amp Off Board Battery Charger.....Standard
- * 1515 Push Button Style Seat Positioning Strap.....Standard

Seat Option

- * CGR Recline (E1004).....\$5,355.00

Group 3 - K0858

Basic Packages

- * TDXSP-CG w/Formula CG with MK6i
- * 400 lb Weight Capacity
- * 4-Pole Motors.....Standard
- * 1431-3 14" X 3" Tires w/Foam Filled Inserts.....\$250.00
- * Flip Back Adjustable Height Arm.....Standard
- * MK690 90 Amp Controller.....Standard
- * 70NHD 70° Swingaway Footrests.....Standard
- * 8 Amp Off Board Battery Charger.....Standard
- * 1515 Push Button Style Seat Positioning Strap.....Standard

Seat Option

- * CGR Recline (E1004).....\$5,355.00

Base

OMIT SEAT FRAME

OSF Omit Seat Frame for Base OnlyNo Charge

TRANSPORT TIE-DOWN

TRBKTS Wheelchair Transport Brackets250.00

*Transport Tie Down

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

*Transport Tie-Down option is not retrofittable to existing models and is not field serviceable.

USER WEIGHT LIMITS

- U300 User Weight <300 lbs. (K0856).....Standard
- U400 User Weight <400 lbs. (K0858).....1,175.00

WHEEL OPTIONS (PAIR)

- 1430-3 14" X 3" Pneumatic Tires.....Standard
- 1431-3 14" X 3" Tires w/Foam Filled Inserts250.00

FORK OPTION

FKPKG One-Sided Fork Package500.00

FRAME FINISH

CONTEMPORARY

- 24P Wet Black.....No Charge
- 60P Silver Metallic.....No Charge
- 71P Silver Vein.....No Charge
- 105P Midnight Blue.....No Charge
- 115P Black Prism.....No Charge

TOO HOT

- 30P Sunny Yellow.....No Charge
- 61P Electric Red.....No Charge
- 119P Bubblegum Pink.....No Charge
- 120P Grape Madness.....No Charge

WAY COOL

- 62P Electric Blue.....No Charge
- 104P Electric Teal.....No Charge
- 121P Lolly Pop Blue.....No Charge
- 122P Cosmic Blue.....No Charge
- 125P Emerald Green.....No Charge
- 127P Grasshopper Green.....No Charge

FOOTNOTES

1. Package includes 4-pole motors and 5.8 MPH top speed with <300 lbs User Weight limit. Package includes Quiet Stability Lock.
2. Package includes Heavy-Duty 4-pole motors and 5 MPH top speed with <400 lb User Weight limit. Package includes Quiet Stability Lock.

NOTE: All specifications and dimensions are approximate.

Invacare® Formula™ CG Recline Only

Accessories and Options

Seat

SYSTEM TYPE

- CGR CG Recline Only System (E1004)5,355.00
- CNVACT Conventional Actuator.....**Standard**
- SMTACT Smart Actuator.....500.00

¹ SEAT-TO-FLOOR HEIGHT

- LOW 17.5" at 0°.....No Charge
- MED 18.5" at 0°.....No Charge
- TALL 19.5" at 0°.....No Charge

SEAT TILT

- OSA 0° Seat Tilt.....No Charge
- 5SA 5° Seat Tilt.....No Charge

SEAT SELECTION

- ² CTS Contoura™ Seating.....395.00
- CNB Conventional Back.....No Charge
- OMCU Omit Contoura Seat CushionLess 20.00

³ PUSH HANDLES

- PCNB Push Handles for Conventional.....No Charge

SEAT WIDTH RANGE

- 1620W 16" to 20" Seat Width Range (Adjustable).....No Charge
- 2024W 20" to 24" Seat Width Range (Adjustable).....No Charge

SEAT WIDTH SETTING

- W16 16" Seat Width.....No Charge
- W17 17" Seat Width (not available w/Contoura).....No Charge
- W18 18" Seat Width.....No Charge
- W19 19" Seat Width (not available w/Contoura).....No Charge
- W20 20" Seat Width.....No Charge
- W21 21" Seat Width (not available w/Contoura).....No Charge
- W22 22" Seat Width.....No Charge
- ⁴ W23 23" Seat Width (not available w/Contoura).....No Charge
- ⁴ W24 24" Seat Width.....No Charge

SEAT DEPTH RANGE

CONVENTIONAL & CONTOURA

- 1619D 16" to 19" Seat Depth Range (Adjustable).....No Charge
- 1922D 19" to 22" Seat Depth Range (Adjustable).....No Charge

CONVENTIONAL & CONTOURA

SEAT DEPTH SETTING

- D16 16" Seat Depth.....No Charge
- D17 17" Seat Depth.....No Charge
- D18 18" Seat Depth.....No Charge
- D19 19" Seat Depth.....No Charge
- D20 20" Seat Depth.....No Charge
- D21 21" Seat Depth.....No Charge
- D22 22" Seat Depth.....No Charge

BACK HEIGHT

CONVENTIONAL BACK

- TSBH20 Back Height 20".....No Charge
- TSBH21 Back Height 21".....No Charge
- TSBH22 Back Height 22".....No Charge
- TSBH23 Back Height 23".....No Charge
- TSBH24 Back Height 24".....No Charge
- TSBH25 Back Height 25".....No Charge
- TSBH26 Back Height 26".....No Charge

CONTOURA

- TSBH25 Back Height 25".....No Charge
- TSBH295 Back Height 29.5".....No Charge

HEADREST TYPE (OPTIONAL)

- TAR127700 Curved Headrest (E0955).....210.00
- TAR27700HW Removable & Adjustable Hardware (E1028).....210.00
- TAR128000 Two-Step Headrest (E0955).....210.00
- TAR128000HW Removable & Adjustable Hardware (E1028).....210.00
- T0022400 Large Headrest (E0955).....210.00
- TAR22400HW Removable & Adjustable Hardware (E1028).....210.00

HEADREST TYPE (CONTINUED)

- T0022600 Small Headrest (E0955).....210.00
- TAR22600HW Removable & Adjustable Hardware (E1028).....210.00

STEALTH HEADRESTS

- CP175 14" Multi Axis Headrest (E0955).....225.00
- CP175HW Removable & Adjustable Hardware (E1028).....225.00
- CP275 10" Multi Axis Headrest (E0955).....225.00
- CP275HW Removable & Adjustable Hardware (E1028).....225.00
- CP180 14" TWB Headrest (E0955).....225.00
- CP180HW Removable & Adjustable Hardware (E1028).....225.00
- CP280 10" TWB Headrest (E0955).....225.00
- CP280HW Removable & Adjustable Hardware (E1028).....225.00

ADJUSTABLE ARM TYPE (PAIR)

- RA19FR Reclining Full Length Adj. Hgt. Right.....No Charge
- RA19FL Reclining Full Length Adj. Hgt. Left.....No Charge
- RA29FR Reclining Desk Length Adj. Hgt. Right.....No Charge
- RA29FL Reclining Desk Length Adj. Hgt. Left.....No Charge

⁵ OTTO BOCK™ ARM PADS (PAIR)

- OB20400 One Piece Channel Armpads (E2209).....230.00
- OB20500 Large Forearm Pads.....235.00
- OB20600 Medium Forearm Pads.....235.00
- OB20700 Small Forearm Pads.....235.00

OTTO BOCK HAND PADS (PAIR)

- OB23000 Large Flat Hand Pads.....100.00
- OB23300 Medium Flat Hand Pads.....100.00

LATERALS (PAIR)

CONVENTIONAL BACK

- TMSLM Swingaway Laterals Medium (E0956).....210.00
- TMSLMHW Removable & Adjustable Hardware (E1028).....210.00
- TMSLL Swingaway Laterals Large (E0956).....210.00
- TMSLLHW Removable & Adjustable Hardware (E1028).....210.00

STEALTH LATERALS FOR CONVENTIONAL BACK

- TWB-RPR Adj. Swing Round Lat. Med (E0956).....320.00
- RPRHW Removable & Adjustable Hardware (E1028).....320.00
- TWB-RPRL Adj. Swing Round Lat. Lg (E0956).....320.00
- RPRLHW Removable & Adjustable Hardware (E1028).....320.00

CONTOURA SEATING

- TCSLM Swingaway Contoura Lat. Med (E0956).....210.00
- TCSLMHW Removable & Adjustable Hardware (E1028).....210.00
- TCSSL Swingaway Contoura Lat. Lg (E0956).....210.00
- TCSSLHW Removable & Adjustable Hardware (E1028).....210.00

STEALTH LATERALS FOR CONTOURA

- TWB-CPR Adj. Swing Contoura Lat. Med (E0956).....320.00
- CPRHW Removable & Adjustable Hardware (E1028).....320.00
- TWB-CPRL Adj. Swing Contoura Lat. Lg (E0956).....320.00
- CPRLHW Removable & Adjustable Hardware (E1028).....320.00

SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Pos. Strap.....50.00
- 1515 Push Button Style Seat Pos. Strap.....**Standard**

CHEST POSITIONING STRAPS

- 7311BK Airline Buckle Chest Pos. Strap (E0960).....100.00
- 9515 Push Button Style Chest Pos. Strap (E0960).....100.00
- 7321BK Hook & Loop Chest Pos. Strap (E0960).....100.00

FOOTNOTES

1. Nominal measurement at front edge of seat pan to floor on 18" deep seat.
2. Selected Seat Depth on Contoura assumes 1" of back cushion compression. For example: When choosing an 18" depth the measurement from the front of the back cushion to the front of the seat pan equals 17". The Contoura back cushion is 2" thick includes leatherette seat cushion. Please refer to page 6 for seat size offerings.
3. Push handles will add an additional 3.5" to back-cane height.
4. Available only on 400 lb. package.
5. Not compatible with quadlink retractable joystick mount.

NOTE: All specifications and dimensions are approximate.

Invacare® Formula™ CG Recline Only
Accessories and Options

Electronics

MK6i™ CONTROLLER

- MK690ACC Controller w/Actuator ControlStandard

STANDARD DRIVER CONTROL

- SPJAP SPJ™+ with Power Seating Switch.....Standard
- SPJAPM Actuator Control through SPJ+ Joystick (E2310) ..1,300.00

EXPANDABLE CONTROLLER SYSTEMS

- MPJM6U MPJ™+ Multiple Drive Joystick (E2377).....750.00
- PSM6 Personal Joystick, Inductive on Front (E2399) ...1,130.00
- PSRM6 Personal Joystick, Inductive on Rear (E2399).....1,130.00
 - PBOD Push Button On/Off w/PB Drive Select.....No Charge
 - PBSS Push Button On/Off w/Speed Select.....No Charge
 - OR On/Off Mounted Right.....No Charge
 - OL On/Off Mounted Left.....No Charge
- ASLRDYM6 MK6 Display Only
(No Driver Control) (E2377).....750.00

**FOR YOUR ASL DRIVE SYSTEM NEEDS
PLEASE REFER TO THE FOLLOWING:
ASL SINGLE INVOICE PRICE LIST FORM # 02-097.**

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS

- I500M6 RIM Head Control w/MK6i Display (E2327)2,875.00
 - SMHD Small Headrest.....No Charge
 - LGHD Large Headrest.....No Charge
- I558MM6 Compact Joystick w/Display (E2373).....1,130.00

ALTERNATIVE DIGITAL CONTROL SYSTEMS

- SNPM6 Sip-N-Puff/Digital Interface Box (E2325).....1,500.00
- 5018M6 Wafer Board & Digital Interface Box (E2322)1,550.00
- 5020M6 Mini Tash Joystick & Digital Interface Box (E2321)..1,750.00

SIP N PUFF KIT

- PKG32666 Therafin® Sip-N-Puff Breath Tube Kit (E2326).....380.00

ADDITIONAL CONTROL CHOICES

- I558M6 Compact Joystick - Less Std Mount770.00
- 5018 Wafer Board.....310.00
- 5020 Mini Tash Joystick.....405.00
- PACM6 Proportional Attendant Control.....850.00
- I552M Digital Attendant Control850.00

JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retract. Joystick Mnt MPJ+ (E1028)279.00
- QLA Quad Link Retract. Joystick Mnt SPJ+ (E1028).....279.00
- ARM250 Stealth Height Adjustable Joystick Mnt (E1028).....210.00
- LEFT Left Handed Mounted JoystickNo Charge
- RIMHW RIM Hardware For I500 (E1028).....219.00
- PKG32669 Complete Bib Assembly (Tash/I558).....265.00
- GATMPJM6 Gatlin Midline for MPJ (E1028).....600.00
- GAT1812 Gatlin Midline for I812 (E1028).....600.00
- RCM Rear Cane Mount for Compact Joystick.....395.00

ELECTRONIC ACCESSORIES

- I813M6 MK6i Programmer w/Pro Memory Card.....475.00
- I813SD Pro Memory Card (USB Ready).....90.00
- I560 “T” Handle Flexible Joystick Ext. (E2323).....77.00
- I561 Straight Handle Flexible Joystick Ext. (E2323).....77.00
- I826 Chin Cup (E2324).....48.00
- A24VPS 24V MK6 Auxilliary Power Source200.00
- AUX12M6 Auxillary Module For 1 & 2795.00
- AUX34M6 Auxillary Module For 3 & 4795.00
- 2PB Dual Push Button.....147.00
- 2SW Dual Toggle172.00

POWER SEATING ACCESSORIES

- SANODE One Actuator Op thru Driver Control (E2310)...1,300.00
- 4WSB Four Way Switch Box.....No Charge
- S4WSB Multiple Interface Box (E2311).....2,600.00
 - FWT Four Way Toggle.....No Charge
 - QPB Four Push ButtonsNo Charge

**For all systems with two or more actuators you must choose
4WSB at No Charge or the S4WSB at an Up Charge.**

**S4WSB Replaces the Four Way Switch Box and allows
multiple actuator operation through the driver control.**

FWT or QPB is required when choosing the 4WSB or S4WSB.

ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Reset-Black.....65.00
- ASL304 Light Wobble Switch w/611 Mount.....320.00
- ASL202J Fiber Optic Switch.....610.00
- ASL208J Adjustable Proximity Switch550.00
- ASL504JM6 Remote Emergency Stop Switch.....535.00
- ASL300A Egg Switch w/611 Mount300.00

FOOTNOTES

1. Allows Actuator Operation through the Driver Control. include Port for Attendant Switch.
2. TOD-SS: Toggle On/Off Drive Select, Speed Potentiometer Standard.
3. Drive Select through mode switch, or mode switch and Joystick Commands.
4. PSRM6 Only.
5. Display Only, No Driver Control. Allows choice from Order Form.
6. Must also use either a Multiple Drive Joystick or MK6i Display.
7. Requires SNPM6 Sip-N-Puff Interface Box.
8. QLAM6 and QLA include ARM250 standard.
9. Options to replace Mono Port “Y” cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options.
10. Required for Recline Operation through Driver Control - Single Actuator Systems. Not compatible with SPJ+ Joysticks. The Egg Switch with the 611 Mount does not come standard when ordering a SANODE.
11. Not available with SPJ+ Joystick. Must add MK6 Auxilliary Power Source.

NOTE: All specifications and dimensions are approximate.

Invacare® Formula™ CG Recline Only
Accessories and Options

Miscellaneous

BATTERY TRAY TYPE

- 22TRY 22NF Style Tray**Standard**
- 24TRY 24 Group Style Tray.....No Charge
- 22TRY3 22NF Style Tray for (3) BatteriesNo Charge
- 22NFTRY3BAT 22NF Style Tray for (2) BatteriesNo Charge

BATTERIES

- 22NFBATTERY 22NF Battery (On Chair) (E2361).....300.00
- 24 BATTERY 24 Gel Battery (On Chair) (E2363).....410.00
- 24NFBATT3 22NF Battery (3) Batteries450.00

VENT TRAY

- 4 VENT Articulating Vent Tray (E1030)1,275.00

MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks45.00

FOOTNOTES

1. 22TRY Base Width 24": 24TRY and 22TRY3 Base Width 25.5".
2. Required for ventilator tray.
3. Allows for 2 batteries in a 3 battery tray.
4. Not available with Contoura Seating.
5. All front riggings include a standard footplate.
6. Allows 60 degree to 90 degree static knee angle positioning with a drop back lever for easier transferring. Extension range 12 1/2" to 15". Only available on Low Seat-to-Floor.
7. Allows 60 degree to 90 degree static knee angle positioning with a drop back lever for easier transferring. Extension range 9 1/4" to 12". Available with Low Seat-to-Floor height only.
8. Standard with impact guards. Not available with widths 19" to 24" and 16" seat depth. Standard with composite footplates only.
9. See page 6 for MPCTMT and MPCTMFL extension range chart. Center mount availability chart on page 6.
10. Not available with SPJ+ joystick.
11. Comes standard with a single footboard.
12. Individual Footplate option for Power Center Mount. This option is only available on MPCTMT.
13. Independant powered elevating legrests including elevating and articulating legrests with angle adjustable footplates. Not available with the 16" or 17" seat depths. Extension range 13"-20". Sold as pair.
14. One piece footboard for Center Mount (CTMT-F) Front Rigging.

NOTE: All specifications and dimensions are approximate.

Front Rigging

HANGER BRACKET TYPE

- ASB Center Pivot Style Bracket.....No Charge
- PSB Pin Style Bracket.....No Charge

5 FRONT RIGGING

CENTER MOUNT (PAIR)

- 6 CTMT-F Center Mount Adjustable Knee Angle220.00
- 7 CTMT-FJ Jr. Center Mount Adjustable Knee to Heel.....220.00

CENTER PIVOT STYLE (PAIR)

- 70NHD 70° Swingaway Footrests.....**Standard**
- 8 70STAPER 70° Tapered Footrests.....No Charge
- AT5544 Manual Elevating Legrests (E0990).....255.00

PIN STYLE- HEMI (PAIR)

- PHW93 70°Swingaway Footrests.....No Charge
- PHAL4A Manual Art./Elev Legrests (E0990).....225.00

9 POWER ELEVATING LEGRESTS - ARTICULATING

- 10,11 MPCTMT Power Center Mount Foot Platform (K0108).....2,995.00
- 10 MPCTMT-L Pwr Center Mount Foot Platform Long (K0108).....2,995.00
- 12 IFP Individual Footplates.....275.00
- 10,13 ELRPW Power ELR w/Articulation 20" Ext. (E1010)1,895.00
- 10,13 ELRPW22 Power ELR w/Articulation 22" Ext. (E1010)1,995.00

FOOTPLATE OPTIONS (PAIR)

- 1651 Composite FootplatesNo Charge
- 1350 Extra Large Aluminum Footplates.....No Charge
- AT5543 Adjustable Angle Flip-Up Footplates.....165.00
- 14 FBCTMT Footboard for Center Mount150.00

FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap (E0038).....27.00
- 1600BK Heel Loops w/Ankle Straps (E0951)56.00
- LPT2 Longer Pivot Slide Tube For 1350.....83.00
- IMPACTG Impact Guards.....51.00
- ALPT Longer Pivot And Slide Tube84.00
- ALPT4 Longer Pivot And Slide Tube 4"84.00

Motion Power Center Mount Knee to Heel			
STF Pick	STF at 0 deg	KTH w/IFP on MPCTMT	KTH w/Footboard
Low	17.5"	10.5"-14.5"	2"-13.5" MPCTMT only
Med	18.5"	11.5"-15.5"	12"-14.5" MPCTMT only
High	19.5"	11.5"-16.5"	12"-15.5" MPCTMT or MPCTMT-L
STF Pick	STF at 5 deg	KTH w/IFP on MPCTMT	KTH w/Footboard
Low	18.75"	11.5"-15.5"	12"-14" MPCTMT only
Med	19.75"	11.5"-16.5"	14"-16" MPCTMT only
High	20.75"	11.5"-17.5"	14"-17" MPCTMT or MPCTMT-L

Invacare® Formula™ CG Recline Only
Accessories and Options

Seating

**FOR YOUR SEATING NEEDS
PLEASE REFER TO THE FOLLOWING**

**REHAB SEATING & POSITIONING PRODUCTS
SEE PRICE LIST FORM # 07-027.**

Batteries (Not Included)

22NF batteries with terminal configuration (positive on the left and negative on the right) **MUST** be used. 22NF batteries that have the reverse terminal configuration **MUST** not be used. Terminals **MUST** have a cross hole for proper battery connection.

GP24 batteries with terminal configuration (negative on the left and positive on the right) **MUST** be used. GP24 batteries that have the reverse terminal configuration **MUST** not be used. Terminals **MUST** have a cross hole for proper battery connection. See Owner's Manual, part number 1114809. These recommendations **MUST** be followed otherwise injury and damage may occur.

Contoura Seat Sizes

	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X
W24	*	*	X	X	X	X	X

*Not available with Contoura Seat Cushion

Back Heights	25" or 29.5"
--------------	--------------

Formula_{CG} Recline Only Features and Specifications

- ♦ Weight capacity of 300 lb. or 400 lb. (with heavy-duty 4-pole motors).
- ♦ Curved ABS or Contoura style seating.
- ♦ Adjustable seat frame.

After Market Backrests

- ♦ The following companies also have hardware to make their backrests work with our new canes:
 - ♦ AEL - planar seating system
 - ♦ Freedom Designs® - planar seating system
 - ♦ Canyon - planar seating system

Conventional Seat Sizes

Seat Widths	16" - 24"
Seat Depths	16" - 22"
Back Heights	20" - 26"

Tech Tips

Operating Powered Seating through the Driver's Control
All chairs require adding or replacing one component to allow actuator operation through the driver control.

SPJ+ with PSS:

Replace joystick with SPJ+ with ACC (Actuator Control) -E2310

Single Actuator Systems with MPJ+, PSR+, PSF+ MK6i Display:

Add a SANODE (Single Actuator Node) to the system - E2310

Multiple Actuator Systems:

Replace the standard 4 Way Switch box with a Multiple Actuator Interface Box - E2311

Memory Card - BASIC

- Standard with ALL Expandable Control Systems (MPJ+, PSR+, PSF+, MK6i Display).
- Used to Back Up/Restore Programmed Settings for Only One Chair.
- May be given to user or kept with provider's files for safe keeping.
- Does not Contain Advanced Diagnostics/Help Screens/File Structure.
- Not compatible with SPJ+ Joysticks.

Memory Card - PROFESSIONAL

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains Advanced Diagnostics, Help Screens, Software Updates, File Storage/Retrieval.
- Used to back up multiple chairs and programming settings.
- Intended for use by qualified providers only.
- Not compatible with SPJ+ Joysticks.

Invacare Corporation www.invacare.com



USA
One Invacare Way
Elyria, OH USA
44035-4190
800-333-6900

All rights reserved. Trademarks are identified by the symbols TM and [®]. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.
Therafin is a registered trademark of W/R Medical Electronics Co.
Otto Bock is a trademark of Otto Bock Orthopaedische Industrien GMBH & Co.
Products designated by ASL are Manufactured by Adaptive Switch Labs, Inc.
Specifications and prices are subject to change without notification.
©2007 Invacare Corporation