

INVACARE®
BAR600IVC
BARIATRIC BED



The **Invacare® BAR600IVC Bariatric Bed** is a heavy-duty full-electric bed frame designed for bariatric individuals. It is capable of supporting users who weigh up to 600 lb. The bed frame extends to the edges of the bed deck for better support than most other split-spring bariatric beds. Using the hand pendant, the user or caregiver is able to change the positioning of the head and foot sections, as well as the bed height. Despite being substantial in size, bed setup within the home is easy – no special tools are needed.

Invacare beds, mattresses and rails, excluding therapeutic support surfaces and bariatric beds, when used as a system and sold new and unused after August 1, 2007, meet the voluntary FDA Guidelines for bed rail entrapment. Visit www.invacare.com/bedsafety for additional safety and entrapment information.

Features

Heavy-duty frame design ensures added strength for client support

Heavy-duty bed sleep surface is 17% larger than conventional beds to ensure comfortable positioning

Split-spring design with universal bed ends for ease in setup

No special tools necessary for in-home setup

Easy-release motor mounts allow quick motor removal

Robotically welded, steel frame construction

Pendant allows for dynamic client positioning

Quiet, smooth operation

Use with Invacare bariatric mattress (BARMATT42) and bed rails (BAR6640IVC):

- BARMATT42 meets fire safety regulations 16CFR1633 and 16CFR1632

600 lb. Client weight capacity

HCPCS Code: E0301



1. Hand pendant

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2021

Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.
Form No. 00-258
Rev. 03/21
210354



Model #	Dimensions	Sleep Surface	Weight Capacity	Shipping Weight	Carton Dimensions	Regulation Compliance	Limited Warranty
BAR600IVC	88" L 14"-22" H 42" W	80" L 42" W	Safe Working Load: 700 lb.* Max. User Weight: 600 lb.	N/A	N/A	UL962 CSA Certified	5 Years Welds/ Frame/Structure 1 Years Mechanical/Electrical
BAR5000IVC Headspring				114 lb.	42" L 10" W 43" D		5 Years Welds/ Frame/Structure 2 Years Mechanical/Electrical
BAR5490IVC Footspring				97 lb.	48" L 9" W 43" D		5 Years Welds/ Frame/Structure 2 Years Mechanical/Electrical
BAR5301IVC Bed End				55 lb.	46" L 7" W 22" D		

*Includes user, mattress and all accessories.

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.