

INVACARE®

Personal Care Products



Yes, you can.®

Invacare offers products that improve patient's lives and increase their freedom, because personal care is just that, it's personal. This is reflected in products designed with unique patient needs in mind and centered on creating safe and comfortable environments. From shower chairs to crutches to bariatric products, Invacare is known for innovation, quality and reliability. The Personal Care line continues to expand and enhance its offerings, because proper patient care can make all the difference.

Invacare® I•Class™ Brand

Personal Care Products

Invacare's best products. They provide the highest level of comfort, performance, features and value.



Latex Free Logo

Invacare's Personal Care Products represented in this catalog are manufactured by Invacare and/or their 3rd party suppliers with latex free components.



Tool-Less Assembly Logo

Certain Personal Care Products do not require tools to assemble before use or to disassemble.



Retail Packaging

Certain Personal Care Products are offered in attractive, consumer-oriented packaging. They can be displayed on pre-designed plan-o-grams to help build retail business and increase consumer awareness. Look for the Retail Packaging symbol within this catalog to easily identify these products.



Invacare® Bariatric Solutions Logo

Invacare's Bariatric Solutions icon represents the company's commitment to providing bariatric solutions for providers and clinicians, as well as their patients. Invacare achieves this through offering a breadth of the highest quality bariatric products, as well as providing solutions by educating providers and clinicians on Invacare products and the bariatric market as a whole. Invacare's products and programs improve the quality of life for patients. An uplifting and active presence is symbolized through the bariatric icon, a presence that Invacare wishes to convey to all bariatric clients through these products and services.

Table of Contents

Ambulatory

Walkers	3 - 5
Walker Accessories	5 - 7
Invacare® WalkLite™ Walkers	8
Invacare® Rollite® Rollators	9
Invacare® Value-Line Rollators	10 - 11
Canes	12 - 13
Crutches and accessories	14 - 15
Knee Walker	15

Bath Safety

Invacare® Rio™ Bath Lift	17
Transfer Benches and Shower Chairs	18 - 20
Grab Bars	21
Tub Bar, Bath Mat, Bath Board and Hand Showers	22
Raised Toilet Seats and Toilet Safety Frame	23

Commodes

Standard Commodes	25 - 26
Folding Commode	26
Drop-Arm Commode	27
Bariatric and Heavy-Duty Commode	28
Commode Accessories	28
Invacare® Mariner™ Rehab Shower Commode Chairs	29

Merchandising Program

With fluctuating reimbursement challenges, a focus on retail sales is an invaluable way to enhance your business model. The Invacare Merchandising Program offers attention-grabbing tools to draw prospects in and maximize profitability.


From specially-designed plan-o-grams to packaging and signage, the Merchandising Program makes it easy to turn your store into a showroom. Exclusively designed to give health care providers the tools and support necessary to build consumer sales, this program maximizes the visual impact, emphasis and appeal of your retail section and lets customers know what Invacare products are available to them. Use these tools to target the expansive baby boomer generation and increase your sales and profit margins.

- Retail Product Packaging
- Plan-o-grams
- Point-of-Sale Materials

Plan-o-grams provide the health care provider with pre-designed product layouts in 4', 8' and 12' sections. Plan-o-grams provide recommended product location, quantity and hook size, as well as product displays and signage.

For additional information and support on building consumer sales, contact your Territory Business Manager or

- 1.) go to www.invacare.com
- 2.) Search plan-o-gram using the search function on left
- 3.) Click on 'plan-o-gram' icon
- 4.) Browse merchandising display options



12' x 8' Plan-o-gram



Invacare Bath Safety Trifold



Ceiling Dangler




Wobbler



Poster



Trifold Holder w/ Header



Window Cling



INVACARE®

Ambulatory Products

- Invacare® Walkers
- Invacare® Rollators
- Invacare® Canes
- Invacare® Crutches
- Invacare® Knee Walker

Folding Mechanisms


Our folding mechanisms make it easy to open and close the walkers even from a seated position.

Walkers fold easily for transport or storage.



Single-Release Button

- One button to press for folding
- Button provides both visual and audible "locked" cues



Invacare® Dual Blue-Release™ Button

- Dual Blue-Release buttons provide both visual and audible "locked" cues



Dual-Release Paddle

- Ergonomic paddles "pop up" when walker is opened
- Nylon paddles are resistant to cracking or shattering

Invacare contemporary walkers feature wide, deep frames with a large number of height adjustments and a lower side brace for added stability. They are light and easy to lift and maneuver. Walkers include PVC handgrips that provide comfort and long-lasting wear and anti-rattle silencers for quiet operation.



Invacare® Single-Release Walkers

Model no. 6281-A - Single-Release Paddle Adult (shown)
Model no. 6281-JR - Single-Release Paddle Junior



Invacare® Blue-Release™ Walkers

Model no. 6240-A - Dual Blue-Release Adult (shown)
Model no. 6240-JR - Dual Blue-Release Junior



Invacare® I-Class™ Paddle Walkers

Model no. 6291-A - Dual-Release Paddle Adult (shown)
Model no. 6291-1 - Dual-Release Paddle Adult
Model no. 6291-JR - Dual-Release Paddle Junior

Invacare® Wheeled Walkers

Model no. 6240-5F - Dual Blue-Release Adult Walker
Model no. 6240-JR5F - Dual Blue-Release Junior Walker

Invacare® Wheeled Walkers

Model no. 6291-3F - Dual-Release Paddle Adult Walker
Model no. 6291-JR3F - Dual-Release Paddle Junior Walker
Model no. 6291-5F - Dual-Release Paddle Adult Walker
Model no. 6291-JR5F - Dual-Release Paddle Junior Walker



Invacare® Heavy-Duty Paddle Walker

Model no. 6291-HDA - Heavy-Duty Dual-Release Paddle Adult

- Compatible with 6375 bariatric 5" fixed wheels and leg extensions



Invacare® Bariatric Dual Release Walker

Model no. 6441-A

- Compatible with Invacare model no. 6372 bariatric 5" fixed wheels and leg extensions

Specifications

Model no.		Width inside hand-grips	Width inside base legs	Depth opened at base	Depth folded	Height adjustment min.- max.	*Approx. client height min. - max.	Product weight	Weight capacity	Units/ carton	Limited warranty
Single Release Walker											
6281-A	Single-Release Adult	17"	21"	17.5"	4.5"	30.375" – 37.375"	5' 3" – 6' 4"	6 lb. 2 oz.	300 lb.	4	5 years
6281-JR	Single-Release Junior	17"	20.5"	16.5"	4.5"	25.375" – 32.375"	4' 4" – 5' 7"	5 lb. 14 oz.	300 lb.	4	5 years
P421B	Two Step Walker	16.5"	21.5"	19"	5"	31" – 36"	5' 3" – 6' 4"	5 lb. 8 oz.	286 lb.	5	5 years
Invacare® Dual Blue-Release™ Walker											
6240-A	Dual Blue-Release Adult	17"	21"	16"	4"	32" – 39"	5' 4" – 6' 6"	5 lb. 7 oz.	300 lb.	4	5 years
6240-JR	Dual Blue-Release Junior	17"	20"	16"	4"	25" – 32"	4' 4" – 5' 7"	5 lb. 5 oz.	300 lb.	4	5 years
6240-5F	Wheeled Dual Blue-Release Adult	17"	21"	18.5"	4"	32" – 39"	5' 4" – 6' 6"	6 lb. 6 oz.	300 lb.	4	5 years
6240-JR5F	Wheeled Dual Blue-Release Junior	17"	20"	18"	4"	25" – 32"	4' 4" – 5' 7"	6 lb. 12 oz.	300 lb.	4	5 years
Dual-Release Paddle Walker											
6291-A (4); 6291-1	Dual-Release Paddle Adult	17"	21.5"	17.5"	4"	30.375" – 37.375"	5' 3" – 6' 4"	5 lb. 5 oz.	300 lb.	4	5 years
6291-JR	Dual-Release Paddle Junior	17"	20.5"	16.5"	4"	25.375" – 32.375"	4' 4" – 5' 7"	5 lb. 1 oz.	300 lb.	4	5 years
6291-3F	Wheeled Dual-Release Paddle Adult	17"	21.5"	17.5"	4.5"	32" – 38"	5' 4" – 6' 4"	5 lb. 5 oz.	300 lb.	4	5 years
6291-JR3F	Wheeled Dual-Release Paddle Junior	17"	20.25"	16.5"	4.5"	27" – 33"	4' 6" – 5' 6"	5 lb. 1 oz.	300 lb.	4	5 years
6291-5F**	Wheeled Dual-Release Paddle Adult	17"	21.5"	17.5"	4.5"	33" – 39"	5' 6" – 6' 6"	5 lb. 6 oz.	300 lb.	4	5 years
6291-JR5F**	Wheeled Dual-Release Paddle Junior	17"	20.25"	16.5"	4.5"	28" – 34"	4' 4" – 5' 6"	5 lb. 6 oz.	300 lb.	4	5 years
6291-HDA	Heavy-Duty Dual-Release Paddle Adult	20"	24"	17.5"	4"	30.375" – 37.375"	5' 3" – 6' 4"	6 lb. 11 oz.	500 lb.	2	5 years
6441-A	Bariatric Dual-Release	20"	23.5"	16"	4.25"	31.25" – 41.25"	5' 3" – 6' 10"	7 lb.	700 lb.	2	5 years
Hemi-Walker											
6252	Hemi Walker Adult	n/a	n/a	n/a	3"	29.5" – 34.5"	5' 1" – 5' 10"	2 lb. 14 oz.	250 lb.	2	5 years

*All client height ranges are approximate. Each individual should always consult with his/her physician or therapist to determine proper adjustment and usage.

"F" designation in Model no. 3" or 5" fixed wheels with glide tips. ** With 13" rear leg extensions.



Invacare® Adult Hemi Walker/Sidestepper

Model no. 6252 - Hemi Walker

- Wider base provides more stability than a quad cane
- 14.75" x 17.75" width at high height setting
- 12.25" x 16.5" width at low height setting



New

Invacare® Two Step Walker

Model no. P421B

- Combination single release folding walker and sit to stand assistance device
- Can also be used as a toilet safety frame
- Adjustable height settings
- Accessories currently not available

Walker Accessories



Walker Basket

Model no. 6098

(Compatible with all walkers except model nos. 6291-HDA, 6252, 6300-JRA, 6300-ATA)

- Mounts easily with snap hooks
- Rust-resistant, gray vinyl coating
- Basket can be installed on the inside or outside of the walker
- 5 lb. weight capacity
- One per carton
- One year limited warranty



Walker Tray

Model no. 6007

(Compatible with all walkers except model nos. 6291-HDA, 6252, 6300-JRA, 6300-ATA, 6240, 6281 walkers)

- Provides convenient assistance for a variety of everyday activities
- Molded with details for cups and pen holders
- Folds out of the way when not in use
- 5 lb. weight capacity
- One per carton
- One year limited warranty



Walker Pouch

Model no. 6015-4

(Compatible with all walkers except model nos. 6291-HDA, 6252, 6300-JRA, 6300-ATA)

- Made of forest green/deep sea blue denier nylon
- Inside and outside pockets
- Easily attaches over front crossbar
- Lightweight and water-repellent
- 5 lb. weight capacity
- Four per carton
- One year limited warranty



Platform Attachment

Model no. 6027

(Compatible with all walkers except model nos. 6291-HDA, 6441-A, 6252 or 6300-JRA, 6300-ATA)

- Thick, vinyl-covered pad with heavy-duty Velcro® strap
- 300 lb. weight capacity
- Five year limited warranty



Walker Tray

Model no. 6002

(Only to be used on model nos. 6240-A, 6240-JR, 6240-5F, 6240-JR5F walkers)



Walker Basket

Model no. 6096

(For model nos. 6240-A, 6240-JR, 6240-5F, 6240-JR5F walkers)

- Snap hooks for easy attachment
- Rust resistant vinyl coating cleans easily
- One year limited warranty

Invacare® walker accessories enhance walkers for safety, maneuverability and durability. Wheel rims are made of durable, easy-to-clean polypropylene; tires are made of non-marring thermoplastic rubber; wheels and bearings are corrosion protected.



Invacare® Court-Side Glide™ Client-Walker Height Range Table

Walker Model No.	Approx. Client Height Range w/Court-Side Glides*
6291, 6240, 6281 w/ 5" fixed wheels	5' 6" – 6' 4"
6291, 6240, 6281 w/ 3" fixed wheels	5' 4" – 6' 4"
6291JR, 6240JR, 6281JR w/ 5" fixed wheels	4' 8" – 5' 6"
6291JR, 6240JR, 6281JR w/ 3" fixed wheels	4' 6" – 5' 6"
6300-JRA Adult to Junior WalkLite	4'11" – 5' 6"
6300-ATA Tall Adult WalkLite	5' 4" – 6' 0"

* All client height ranges are approximate. Do not use on 6291HDA

Model no. 6290: Sold as a carton; 6 pairs per carton. Each pair complete with aluminum leg extensions, housings and tennis balls.

Swivel Wheel Attachments



5" Swivel Wheels with 13" Cambered Invacare® SureGlide™ Brakes

Model no. 6268



5" Swivel Wheel Attachments with Rear Glide Tips

Model no. 6267



3" Swivel Wheels with 13" Cambered Invacare® SureGlide™ Brakes

Model no. 6272



3" Swivel Wheel Attachments with Rear Glide Tips

Model no. 6273

Fixed-Wheel Attachments



3" Single Fixed-Wheel Attachments with Rear Glide Tips

Model no. 6270**



5" Single Fixed Wheel Attachments with Rear Glide Tips

Model no. 6271**

- Use 13" rear leg extensions for maximum adjustability



Tall 3" Single Fixed-Wheel Attachments with Tall Rear Legs and Glide Tips

Model no. 6278 with 13" rear leg extensions

- Adjusts to accommodate persons up to 6' 6" tall



Bariatric Wheel Kit

Model no. 6372 for use with 6441-A. 6375 for 6291-HDA not shown

- Dual 5" fixed wheels for added stability
- Large number of height adjustments

**Also available in retail packaging with a minimum purchase of four sets. Order as Model nos. 6270-4 and 6271-4.

Rear Brake Attachments

Brake attachments operate without interfering with regular walking motion.



10" Rear Wheel Brakes

Model no. 6264

- Use with 3" or 5" single fixed-wheel attachments

13" Rear Wheel Brakes

Model no. 6265 (not shown)

- Use with 3" or 5" swivel wheel attachments or tall 3" single fixed-wheel attachments

Leg Extensions



Tall Leg Extensions

Model no. 6275

- Adjusts to accommodate persons up to 6' 6" tall



13" Rear Leg Extensions with Glide Tips

Model no. 6266

- Use with swivel or fixed wheels



Glide Tips

Model no. 6280

- Packaged 24 per carton
- Durable construction for long wear

Walker Accessory Specifications*

*Walker accessories are not to be used with Invacare® WalkLite™ Walker.

Height adjustment

Adult walkers

6291-1, 6291-3F, 6291-5F,
6291-A, 6281-A, 6240-A,
6240-5F

Junior walkers

6291-JR, 6281-JR, 6291-JR3F,
6291-JR5F, 6240-JR, 6240-JR5F

Model no.		min. – max.	min. – max.	Limited warranty
Fixed-wheel				
6270	3" single fixed-wheels w/rear glide tips	31" – 37"	26" – 32"	5 years
6271	5" single fixed-wheels w/rear glide tips	33" – 39"	28" – 33"	5 years
6277	tall 3" single fixed-wheels w/glide tips (not shown)	31" – 37"	26" – 32"	5 years
6278	tall 3" single fixed-wheels w/tall rear leg ext. w/glide tips	33" – 40"	31" – 35"	5 years
6372*	double 5" fixed-wheels for 6441-A bariatric walker	33.5" – 43.5"	–	5 years
6375*	double 5" fixed-wheels for 6291-HDA	33.5" – 43.5"	–	5 years
Swivel wheel				
6267	5" swivel wheels w/rear glide tips	37" – 40"	32" – 35"	5 years
6268	5" swivel wheels w/13" cambered Invacare® SureGlide™ brakes	37" – 40"	32" – 35"	5 years
6269	3" swivel wheels w/13" rear leg extensions w/glide tips (not shown)	34" – 40"	29" – 35"	5 years
6272	3" swivel wheels w/13" cambered Invacare SureGlide brakes	34" – 40"	29" – 35"	5 years
6273	3" swivel wheels w/rear glide tips	34" – 37"	29" – 32"	5 years
Leg extensions				
6266	13" rear leg extensions w/glide tips	33" – 40"	29" – 35"	5 years
6275	tall front and rear leg extensions	33" – 40"	28" – 35"	5 years
Rear brake attachments				
6264	10" rear wheel brakes	33" – 37"	28" – 32"	5 years
6265**	13" rear wheel brakes	33" – 40"	29" – 35"	5 years

**For use with 3" or 5" swivel wheels or tall 3" fixed wheels.

The Invacare® WalkLite™ Walker revolutionized the world of walkers by addressing the needs of providers, clinicians and consumers at an affordable price. The line of walkers covers a wide range of consumer heights from Junior to Tall Adult. The WalkLite Walker continues to feature all the benefits of the original WalkLite Walker – a comfortable flip-up seat, durable front 5" fixed wheels, rear spring-loaded Invacare® SureGlide™ brakes, a lightweight aluminum frame with flexible backrest and ergonomic dual-paddle folding mechanisms.



Model no. 6300-JRA
(Junior to Adult)



Model no. 6300-ATA
(Adult to Tall)



Folds into a compact unit for storage. Built-in seat handle allows for easy transport.

Features

- Lightweight aluminum frame
- Flip-up seat with built-in handle
- Ergonomic dual-paddle folding mechanisms
- Folds into a compact unit for storage and transport
- Flexible backrest for comfort
- Rear spring-loaded Invacare® SureGlide™ brakes
- Adjustable seat-to-floor height

Basket and Tray Accessories

- Basket (model no. 65110 (shown) for 6300-ATA; 65109 for 6300-JRA)
 - 10 lb. weight capacity
 - Position basket underneath seat to secure valuables
- Tray (model no. 65408)
 - For food and beverage placement
 - Tray fits on both 6300-JRA and 6300-ATA



Flip-up seat allows user to step inside the frame.



Rear spring-loaded SureGlide brakes provide additional security when seated.**

** The rear SureGlide brakes engagement may vary depending upon the weight distributed on the seat and the positioning of the user on the seat.

Specifications

	6300-JRA	6300-ATA
Depth opened at base	23"	23"
Depth folded	7"	7"
Width inside handgrips	18"	18"
Width inside base of legs	20"	24.5"
Width overall	27"	27"
Height adjustment min. - max.	29.25" – 36.5"	32" – 39"
Approx. client height* min. - max.	4' 11" – 6' 0"	5' 5" – 6' 5"
Seat-to-floor height min. - max.	16.625" – 23.625"	19" – 25.75"
Seat width x depth	17" x 9.5"	17" x 9.5"
Wheels	5" fixed	5" fixed
Product weight	11 lb.	11 lb.
Weight capacity	300 lb.	300 lb.
Units/carton	1	1
Limited warranty	5 years	5 years

*Client height is approximate. Each individual should always consult with his/her physician or therapist to determine proper adjustment.

The Invacare® Rollite® Rollator line is unlike any other rollator line in the industry. Invacare has taken the Rollite Rollator one step further, offering the sizes, wheel types and colors that consumers want. The Rollite line offers adult and junior sizes in different colors. For improved mobility over rougher terrain, large 8" wheel models are available in adult and tall sizes (Midnight Blue color). Along with additional features, the Rollite line provides the comfort, durability and quality that consumers have come to expect from Invacare products.



Color: Electric Red
Model nos. 65100R Rollite Adult,
65100R-JR Rollite Junior (shown)



Color: Midnight Blue
Model nos. 65100 Rollite Adult (shown)
65100-JR Rollite Junior



Color: Midnight Blue
Model nos. 68100 Rollite Adult w/8" wheels,
68100-TA Rollite Tall Adult w/8" wheels (shown).

Features

- Lightweight aluminum frame
- Flip-up seat with built-in handle
- Ergonomic dual-paddle folding mechanism
- Folds into a compact unit for storage and transport
- Flexible backrest for comfort
- Ergonomic handbrakes lock the rear wheels for security
- Adjustable seat-to-floor height

Basket and Tray Accessories

- Basket (model no. 65110) For use on Rollite models 65100, 65100R and 68100-TA
- Tray (model no. 65408)
 - For food and beverage placement
 - For use on all Rollite models
- Basket (model no. 65109) For use on Rollite models 65100-JR, 65100R-JR and 68100
 - 10 lb. weight capacity
 - Position basket underneath seat to secure valuables



Folds into a compact unit for storage. Built-in seat handle allows for easy transport.

Specifications

	65100R Rollite Adult (Electric Red)	65100R-JR Rollite Junior (Electric Red)	65100 Rollite Adult (Midnight Blue)	65100-JR Rollite Junior (Midnight Blue)	68100 Rollite Adult w/8" wheels (Midnight Blue)	68100-TA Rollite Tall Adult w/8" wheels (Midnight Blue)
Depth opened at base	24"	24"	24"	24"	26"	26"
Depth folded	7"	7"	7"	7"	7"	7"
Width inside handgrips	18.25"	18.25"	18.25"	18.25"	18.25"	18.25"
Width inside base of legs	22"	22"	22"	22"	22"	22"
Width overall	27"	27"	27"	27"	27"	27"
Height adjustment min. - max.	32" – 36"	29" – 33"	32" – 36"	29" – 33"	32" – 36"	35" – 39"
*Approx. client height min. - max.	5' 4" – 6'	4' 11" – 5' 7"	5' 4" – 6'	4' 11" – 5' 7"	5' 4" – 6'	5' 11" – 6' 6"
Seat-to-floor height min. - max.	19" – 23"	16.5" – 20.5"	19" – 23"	16.5" – 20.5"	19.5" – 23.5"	22" – 26"
Seat width x depth	17" x 9.5"	17" x 9.5"	17" x 9.5"	17" x 9.5"	17" x 9.5"	17" x 9.5"
Wheels	5"	5"	5"	5"	8"	8"
Product weight	13.5 lb.	13.5 lb.	13.5 lb.	13.5 lb.	13.5 lb.	16.5 lb.
Weight capacity	300 lb.	300 lb.	300 lb.	300 lb.	300 lb.	300 lb.
Units/carton	1	1	1	1	1	1
Limited warranty						
brakes	5 years	5 years	5 years	5 years	5 years	5 years
other components**	1 year	1 year	1 year	1 year	1 year	1 year

*Client height is approximation. Each individual should always consult with his/her physician or therapist to determine proper adjustment.

**No tires/wheels.

Invacare® Value-Line Rollators will meet all your needs with a focus on style and affordability. These rollators feature easy-to-operate ergonomic hand brakes and easy folding for storage and transport. They offer great maneuverability and stability, and are height adjustable for a custom fit. Basket is standard on all rollators.



65350GR Junior Rollator

Model no. 65350GR

- Easy-to-operate, ergonomic hand brakes lock the rear wheels for security
- Flip-up seat
- Lightweight aluminum frame
- Basket hooks easily on front seat brace
- Removable foam padded straight backrest enhances comfort and can be easily removed



65851R / 65851B Adult Rollator

Model no. 65851R Red, 65851B Blue

- Easy-to-operate, ergonomic hand brakes lock the rear wheels for security
- Flip-up seat
- Lightweight aluminum frame
- Basket hooks easily on front seat brace
- Removable foam-padded straight backrest enhances comfort and can be easily removed



65650 / 65650R Rollator

Model no. 65650 Blue, 65650R Red

- Easy-to-operate, ergonomic hand brakes lock the rear wheels for security
- Flip-up padded seat with zipper compartment for storage of valuables
- Removable foam padded curved backrest enhances comfort and can be easily removed
- Lightweight aluminum frame
- Rollator folds easily with basket attached
- Universal basket fits underneath seat or hooks easily on front seat brace



New

Invacare® FR300 Rollator

Model no. 67100

- Tray/Basket included
- Cross Brace
- Streamlined design
- Curved backrest
- Curb climber



65510 Rollator

Model no. 65510

- Easy-to-operate, ergonomic hand brakes
- Rollator folds easily with tote bag attached
- Indoor/outdoor 7½" x 2" solid rubber wheels
- Basket standard, attaches easily to the frame
- Tray standard, attaches easily to the basket
- Tote bag standard, attaches easily to the frame



65550 Rollator

Model no. 65550

- Easy-to-operate, ergonomic hand brakes lock the rear wheels for security
- Sturdy seat
- Rollator folds easily with basket attached
- Basket and tray standard

Baskets are standard on model nos. 65410, 65350GR, 65450, 65650, 65650R, 65851R and 66500. Tray is standard on model nos. 65410 and 65450.



66550 Bariatric Rollator

Model no. 66550

- 500 lb. weight capacity
- Wide, deep frame
- Flip-up padded seat
- Removable padded backrest
- Folds for ease of storage and transport
- Large 8" wheels for improved mobility over rougher terrain
- Height adjustable loop handbrakes with foam handgrips
- Durable steel construction
- Pouch included

Specifications

	65510	65350GR	65550	65650/65650R	65851R/65851B	66550	67100
Depth opened	24"	24.5"	27"	28.25"	27"	29.5"	26.3"
Width opened	26"	26"	23.75"	24.25"	28.5"	27"	24.4"
Adj. handle height min.-max.	31.25" – 37"	29.5" – 34"	28" – 38.5"	32.6875" – 37.3125"	34" – 38"	35.75" – 38"	30" - 38.2"
Client approx. height min.-max.	5' 3" – 6'	5' – 5' 6"	5' 3" – 6' 5"	5' 5" – 6' 2"	5' 6" – 6' 3"	5' 8" – 6' 5"	
Distance between handles	14"	13.5"	16.5"	18"	16"	22"	16"
Wheels	7.5" x 2"	6"	7.5"	7.5"	6"	7.5"	8"
Seat width x depth	N/A	14" x 13"	15.25" x 6.375"	13.75" x 12.75"	14" x 13"	18" x 14"	15.5" x 8"
Seat-to-floor height	15" – 21"	19.5"	24"	23.25"	22"	22"	23.6"
Product weight	17 lb.	13 lb.	25 lb.	19 lb.	13.4 lb.	26 lb.	17.5 lb.
Weight capacity	295 lb.	250 lb.	285 lb.	300 lb.	300 lb.	500 lb.	285 lb.
Limited warranty							
brakes	5 years	5 years	5 years	5 years	5 years	5 years	5 years
other components*	1 year	1 year	1 year	1 year	1 year	1 year	1 year
Units/carton	1	1	1	1	1	1	1

*No tires/wheels.

Invacare offers a stylish and unique line of canes. With a wide assortment of handgrip styles, you can fit virtually any of your clients with ease. The quad canes come in medium and large base varieties, have a 250 lb. weight capacity and allow for proper fit with right- or left-hand use. Invacare Bariatric Canes are built with durable steel construction that provides maximum stability for individuals weighing up to 700 lb. All of the canes are height adjustable for a proper fit to the consumer and are designed with a twist locking ring providing additional security during height adjustment.



The Invacare® Grip

- The familiar triangular-shaped grip. The flat top helps reduce the stress a user's hand may experience through everyday cane use.



Cane Tip and Collar

- Locking adjustment collar secures cane at desired length. Tip features reinforced metal insert for longer wear.



Invacare® Standard Crook Canes

Model no. 8916 – Black

- Invacare Grip



Invacare® Offset Canes

Model no. 40918-6 – Black

Model no. 8917 – Silver

- Invacare Grip with Strap

Cane Specifications

Model no.	Style	Color	Grip type	Adjustment	*Client height (approx.)	Weight	Weight capacity	Units/ carton	Limited warranty
8916	crook	black	Invacare	30" – 39"	4' 6" – 6' 6"	15 oz.	250 lb.	6	3 years
40918-6	offset	black	Invacare	29" – 38"	4' 4" – 6' 4"	15 oz.	250 lb.	6	3 years
8917	offset	silver	Invacare	29" – 38"	4' 4" – 6' 4"	15 oz.	250 lb.	6	3 years
8940-6	offset	black		29" – 38"	4' 10" – 6' 4"	3 lb.	700 lb.	6	3 years

*All client height ranges are approximate. Each individual should always consult with his/her physician or therapist to determine proper adjustment and usage.

Quad Cane Specifications

Model no.	Base	Color	Grip type	Adjustment	*Client height (approx.)	Footprint	Weight	Weight capacity	Units/ carton	Limited warranty
3920-2	medium	silver-vein finish	Invacare	29.25" – 38.25"		10.25" x 6"	1 lb. 13 oz.	250 lb.	2	3 years
7827	large	silver-vein finish	Invacare	29.25" – 38.5"		13" x 7.75"	2 lb. 10 oz.	250 lb.	2	3 years
7831-2	small	black		29.5" – 38.5"	4' 11" – 6' 5"	8" x 6"	3 lb.	700 lb.	2	3 years
7835-2	large	black		29.5" – 38.5"	4' 11" – 6' 5"	12" x 8.5"	3 lb.	700 lb.	2	3 years

Orthopedic Grip Cane Specifications

This unique orthopedic grip is ergonomically shaped to spread pressure evenly across the palm.

Model no.	Handle style	Color	Adjustment	height* (approx.)	Weight	Weight capacity	Units/ carton	Limited warranty
8929	orthopedic (right)	black	31" – 40"	4' 6" – 6' 6"	15 oz.	250 lb.	6	3 years

*Client height is approximation. Each individual should always consult with his/her physician or therapist to determine proper adjustment.



Medium Base Quad Cane

Model no. 3920-2 - Invacare Grip

Large Base Quad Cane

Model no. 7827 - Invacare Grip

Invacare® Orthopedic Grip Cane

Model no. 8929 - Black (Right Hand)

- This unique orthopedic grip is ergonomically shaped to spread pressure evenly across the palm



Bariatric Quad Canes

Model no. 7831-2 small base, 7835-2 large base

- 700 lb. weight capacity
- Right- or left-hand use
- Height adjustable
- Comfortable PVC hand grip
- Locking ring for additional stability
- Durable steel construction

Bariatric Offset Cane

Model no. 8940-6

- 700 lb. weight capacity
- Height adjustable
- Comfortable PVC hand grip
- Locking ring for additional stability
- Durable steel construction

Small base quad cane model no. 7831-2

Large base quad cane model no. 7835-2

Offset base quad cane model no. 8940-6



Cane Tips

Model nos. 40850-4, 40851-4, 40852-4, 40853-4

- Each pair of cane tips is retail packaged

Specifications

Model no.	Units/ carton	Fits cane diameter
40850-4	8	1"
40851-4	8	.875"
40852-4	8	.75"
40853-4	8	.625"

**350 lb.
Weight
Capacity**



Invacare® Quick-Change Crutches

Model no. 8115-T – Tall Adult, 8115-A – Adult, 8115-J – Junior

Invacare® quick-change crutches offer a solid one-piece crutch frame in lightweight anodized aluminum with push-button height adjustment and contoured handgrips.

Invacare® Quick-Adjust Crutches

Model no. 8120-T – Tall Adult, 8120-A – Adult, 8120-J – Junior

The one-piece handgrip assembly eliminates the need for wing nuts and other assembly pieces. Plus, the push-button height adjustment and patient height guide allow the crutch to be fitted to the proper user height. This patented model comes completely assembled for the consumer's convenience.

New

Invacare® Forearm Crutches

Model no. 8153-T – Tall, 8153-A – Adult, 8153-J – Junior

Forearm crutches feature 1" solid composite cuffs, cushioned vinyl handgrips and extra large crutch tips.

Quick-Change/Quick-Adjust Crutch Specifications

Model no.	Size	*Approximate client height	Adjustment range	Weight per pair	Weight capacity	Units/carton	Limited warranty
8115-J	junior	4' 6" – 5' 2"	36" – 44"	3 lb. 4 oz.	350 lb.	8 pairs	3 years
8115-A	adult	5' 2" – 5' 10"	44" – 52"	3 lb. 8 oz.	350 lb.	8 pairs	3 years
8115-T	tall adult	5' 10" – 6' 6"	52" – 60"	3 lb. 10 oz.	350 lb.	8 pairs	3 years
8120-J	junior	4' 6" – 5' 2"	36" – 44"	3 lb. 7 oz.	250 lb.	8 pairs	3 years
8120-A	adult	5' 2" – 5' 10"	44" – 52"	3 lb. 11 oz.	250 lb.	8 pairs	3 years
8120-T	tall adult	5' 10" – 6' 6"	52" – 60"	4 lb.	250 lb.	8 pairs	3 years

Forearm Crutch Specifications

Model no.	Size	*Approximate client height	Height adj. to center of handle	Top of cuff to center of handle	Weight per pair	Weight capacity	Units/carton	Limited warranty
8153-J	junior	4' 5" – 5' 2"	23" – 32"	10.5" – 13.5"	4 lb. 8 oz.	300 lb.	6 pairs	3 years
8153-A	adult	5' 2" – 5' 11"	29" – 38"	10.5" – 13.5"	4 lb. 8 oz.	300 lb.	6 pairs	3 years
8153-T	tall adult	5' 11" – 6' 8"	32" – 41"	11" – 14"	4 lb. 8 oz.	300 lb.	6 pairs	3 years

*All client height ranges are approximate. Each individual should always consult with his/her physician or therapist to determine proper adjustment and usage.

Invacare® Crutch Accessories for Model No. 8115 and 8120 Crutches



Model nos. 6134, 6135

Rubber Tips

Model no. 6134
Rubber Tips – Large 1.75" O.D. Base
Model no. 6135
Rubber Tips – Extra Large 2" O.D. Base

- Help grip the surface and absorb shock
- Reinforced with metal inserts for longer-lasting wear



Model no. 6131

Arm Pads

Model no. 6131 – Arm Pads

- Add comfort to armrests



Model nos. 6125*, 6126** (not shown)

Handgrips

Model nos. 6125, 6126 – Solid Handgrips

- Padded to keep hands and wrists comfortable
- Contoured grips that limit rotation

Model no.	Only available for Model nos.
* 6125	8120-J, 8120-A and 8120-T crutches
**6126	8115-J, 8115-A and 8115-T crutches

Bariatric Crutches

Strong and comfortable, Invacare's bariatric crutches range in size from adult to tall adult. Large, double push buttons allow for height adjustments in 1" increments. Crutch frame offers durable steel construction.



Invacare® Bariatric Crutches

Model nos. 8130-A, 8130-T

- 1000 lb. weight capacity
- Push-button height adjustment
- Patient height guide on crutch for easy adjustment
- Contoured underarm design
- Arm pads (model no. 7131), handgrips (model no. 7126) and 2" crutch tips (model no. 6135) also available as accessories

Invacare® Bariatric Crutch Accessories For Model No. 8130 Crutches



Rubber Tips

Model no. 6135

- Extra large, 2" outside diameter (O.D.) base



Arm Pads

Model no. 7131

- Contoured underarm design



Handgrips

Model no. 7126

- Solid handgrips

Bariatric Crutch Specifications

Model no.	Size	*Approximate client height	Adjustment range	Weight per crutch	Weight capacity	Units/carton	Limited warranty
8130-A	adult	5' 2" – 5' 10"	44" – 52"	4 lb. 6 oz.	1000 lb.	1 pair	3 years
8130-T	tall adult	5' 10" – 6' 6"	52" – 60"	5 lb.	1000 lb.	1 pair	3 years

* Client height is approximate. Each individual should always consult with his/her physician or therapist to determine proper adjustment.



Model no. 65960



The Invacare® Knee Walker is a great alternative for those who are weary of the discomfort related to crutches. It features a comfortable knee pad for resting the wounded limb, and it is easily maneuverable with its 5" front swivel wheels. Its easy-to-operate ergonomic hand brakes keep the consumer in control. In addition, the Knee Walker folds easily for storage or transport, and it comes with a basket for the consumer's convenience.

Specifications

	65960
Depth opened	29"
Width opened	18"
Adj. handle height min.-max.	31" – 39"
Client approx. height min.-max.	5' 2" – 5' 10"
Wheels	6"
Seat width x depth	7.5" x 14"
Seat-to-floor height	17" – 20.5"
Product weight	16.5 lb.
Weight capacity	300 lb.
Limited warranty	1 year
Units/carton	1



INVACARE®

Bath Safety Products

- Invacare® Rio™ Bath Lift
- Invacare® Transfer Benches
- Invacare® Shower Chairs
- Invacare® Grab Bars
- Invacare® Bath Safety Accessories
- Invacare® Raised Toilet Seats and Toilet Safety Frames



1. Remote Control, Battery Operated
2. Textured Seating
3. Self-Release Suction Cups

New

Invacare® Rio™ Bath Lift

With the rise of new and better technology, there are more innovative and attractive products that are available to provide safety in the bathroom. Bath lifts that gently and safely lower and raise a person into a bathtub in particular are growing in popularity. Bath lifts are attractive because they also provide a cost-effective alternative to expensive and labor-intensive bathroom renovations. Installing simple bath safety accessories such as bath lifts will aid in the safety and security in the bathroom and promote a sense of empowerment for seniors.

Specifications

	H605
Seat width without side flaps	14.6"
Seat width with side flaps	28.3"
Seat depth	17.3"
Seat height max.	17.3"
Seat height min.	2.36"
Width of backrest	13.8"–15.7"
Height of backrest	24.2"
Total depth	26.8"
Total height with backrest	40.6"
Base area	26.4" L x 14.6" W
Total weight with battery	21 lb.
Weight, seat	3.5 lb.
Weight, backrest	12.5 lb.
Weight, battery	3.5 lb.
Weight capacity	264 lb.
Color	White

Invacare® Transfer Benches help make transfers in and out of the bathtub safer and more comfortable. The 1" anodized aluminum frame is rust-resistant and features an arm rail for support and non-marring rubber tips.

Invacare® I-Fit™ and Invacare® CareGuard™ Shower Chairs offer stability, comfort, and tool-less assembly in seconds. Invacare® shower chairs are designed for individuals who require seating support when bathing or showering and provide comfort and safety in the bathtub. Available with a composite or vinyl padded seat, shower chairs fit in most standard tubs and showers.



Model no. 98071

Invacare® CareGuard™ Tool-less Transfer Bench

Model no. 98071

- Compact design
- Three-piece blow-molded seat
- Textured seat with drain holes
- No tools required to assemble
- Hand hole in seat back for easy transport
- Suction cups
- Half-inch increments and adjustments



Model no. 9670C



Model no. 9670



Invacare® I•Class™ Heavy-Duty Transfer Benches

Model nos. 9670 (partially assembled), 9670U (unassembled not shown), 9670C (partially assembled with commode opening)

- Sturdy, durable, one-piece blow-molded seat
- Textured seat with drain holes
- Ergonomically shaped, reversible backrest snaps in and out easily without tools
- Built-in soap dish and hand-held shower holder



Model no. 96-2

Model no. 91-2

Invacare® CareGuard™ Tool-less Shower Chairs

Model no. 91-2 and 96-2

- Tool-less assembly
- Deep wide seating surface
- Contoured seat with ergonomic backrest
- Textured seat with drain holes
- Hand hole built into seat back for easy transport
- Aluminum legs with half-inch snap button adjustments to help with stability
- Wide-based crutch tips to help increase stability



Model no. 9872

Invacare® Vinyl Padded Shower Chair

Model no. 9872 (assembled)

- Vinyl seat and back are padded with 1.25" open-cell foam for comfort
- Nine height adjustments in half-inch increments
- Non-marring rubber tips



Model no. 9781

Model no. 9780

Invacare® I-Fit™ Shower Chairs

Model no. 9781 (with back), Model no. 9780 (without back)

- Tool-less assembly
- Deep and wide seating surface
- Contoured seat with drain holes and ergonomic backrest
- Plastic legs with twist height adjustments
- Numbering system built into legs to help keep chair level



Model no. 9871

Invacare® Vinyl Padded Transfer Benches

Model nos. 9871 (partially assembled), 9871-1 (partially assembled)

- Vinyl seat and back are padded with 1.25" open-cell foam for comfort

The Invacare® Bariatric Shower Chairs provide a wide and deep seating surface for those who need extra stability in the shower. These 700 lb. weight capacity shower chairs come with contoured seat and backrest and adjustable height legs. The durable aluminum frame with special cross brace design provides maximum stability for larger individuals.



Model no. 9675-2



Model no. 9785-2

Invacare® Bariatric Transfer Bench

Model no. 9675-2

- 700 lb. weight capacity
- Wide, deep seating surface
- Contoured seat and backrest
- Adjustable height legs
- Durable aluminum frame with dual cross brace for added stability
- Backrest easily snaps into assembled seat frame with no tools required

Invacare® Bariatric Shower Chair

Model no. 9785-2

- 700 lb. weight capacity
- Wide, deep seating surface
- Contoured seat and backrest
- Adjustable height legs
- Durable aluminum frame with dual cross brace for added stability
- Backrest easily snaps into assembled seat frame with no tools required

Specifications Transfer Benches

Model no.	Seat height min. - max.	Seat depth	Seat width	Overall height min. - max.	Clearance beneath frame min. - max.	Overall width at base min. - max.	Overall depth at base min. - max.	Height adjustments	Weight	Weight capacity	Units/carton	Limited warranty
98071	18" – 22"	16.25"	27"	30" – 34"	14.75" – 18.75"	29"	15¼"	9**	10 lb.	315 lb.	2	3 years
9670, 9670U, 9670C	18" – 22.5"	18"	33.25"	31.25" – 35.75"	14.5" – 19"	30.5" – 31"	16.25" – 17"	9*	12 lb.	400 lb.	2, 1, 2	3 years
9871, 9871-1	17.5" – 22"	16"	23"	28.5" – 33"	14.5" – 19"	30.5" – 31"	16.25" – 17"	9*	13 lb.	300 lb.	2, 1	3 years
9675-2	16.25" – 20.25"	16"	27"	28" – 32.25"	14" – 18.25"	29.5" – 30.5"	21" – 23"		16 lb.	700 lb.	2	3 years

* in 1" increments ** in half inch increments

Specifications Shower Chairs

Model no.	Seat height min. - max.	Seat depth	Seat width	Back height from seat	Overall width at base min. - max.	Overall depth at base min. - max.	Weight	Weight capacity	Units/carton	Limited warranty
91-2	16.125" – 22"	12.75"	21"	n/a	17.25"	9.5" – 18.875"	5 lb. 14 oz.	315 lb.	2	3 years
96-2	16.1225" – 22"	12.75"	21"	17.25"	17.25"	13.75" – 18.875"	6 lb. 6 oz.	315 lb.	2	3 years
9781	15" – 21"	18"	20"	14.25"	16.25"	17.25"	10 lb.	400 lb.	2	3 years
9780	15" – 21"	18"	20"	n/a	16.25"	17.25"	8 lb. 4 oz.	400 lb.	2	3 years
9872	17.5" – 22"	16"	15.5"	14.75"	19.25" – 19.75"	16.25" – 17"	9 lb.	250 lb.	2	3 years
9785-2	16.25" – 20"	16.25"	17"	13.5"	18.25" – 18.75"	21" – 22"	11 lb.	700 lb.	2	3 years

91-2 and 96-2 – half inch increments on leg extensions

Invacare® Folding Shower Chair With Back

Model no. 9981 (with back - assembled)

Invacare® folding shower chairs offer bathing comfort and support and fold easily for storage or travel. Features include a contoured blow-molded seat and optional backrest, 1" anodized aluminum frame and fixed-height legs.

- Wide contoured seat and back
- Ergonomically positioned drain holes allow consistent drainage
- Hand-held shower holders are molded in for convenience
- Non-marring rubber tips



Model no. 9981



Specifications

Model no.	Seat height	Seat depth	Seat width	Back height from seat	Overall width at base	Overall depth at base	Folded length	Folded width	Folded depth	Weight	Weight capacity	Units/carton	Limited warranty
9981	17.5"	17"	19.75"	12"	16.75"	17.75"	25.5"	19.75"	7.5"	6 lb. 8 oz.	250 lb.	2	3 years

Invacare® Mobile Shower Chair

Model no. 6358 (partially assembled)

The Invacare® mobile shower chair is designed for safe transport to and from the shower. Features include a flexible nylon backrest, anodized, non-corrosive 1" aluminum tubing and 5" swivel casters with rear wheel safety locks. Once pail is removed, chair easily rolls over most toilets.



Model no. 6358

Specifications

Model no.	Seat height	Seat-to-armrest height	Width between arms	Height overall	Width overall	Length overall	Casters	Weight	Weight capacity	Units/carton	Limited warranty
6358	21"	8"	18"	38"	22.375"	24"	5"	11 lb.	250 lb.	2	3 years

New

Invacare® Suction Cup Grab Bar

Model no. 218-3, 222-3

Comes in two sizes, 17.5" and 21.5". Other features include visual adhesion control, tool-less assembly, a non-slip surface that offers a secure grip and a bar that mounts in a variety of positions.



Specifications

Model no.	Gripping surface length	Gripping surface diameter	Wall clearance	Overall length	Suction cup diameter
218-3	6.5"	1"	1"– 1.6"	17.5"	3.7"
222-3	10.5"	1.4"	1"– 1.6"	21.5"	3.7"



Knurled Grab Bars (Residential)

Model nos. 512-3, 516-3, 518-3, 524-3, 532-3

These knurled grab bars are constructed of heavy-gauge steel for added stability with a knurled, chrome-plated finish that is rust-resistant. They mount in a variety of angles to a wall's structural supports, with rubber gaskets to protect the wall surface.



Firmgrip Grab Bars (Residential)

Model nos. 312-3, 316-3, 318-3, 324-3, 332-3

These grab bars include a soft-touch, Firmgrip white powder coating surface, durable steel construction with 3.75" wall clearance and a 1" O.D. gripping area.

Grab Bar Specifications *

Model no.	Length	Gripping diameter	Clearance from wall	Product description	Weight	Weight capacity	Units/ carton	Limited warranty
312-3, 512-3	12"	1"	3.75"	Firmgrip/knurled	1 lb.	250 lb.	3	3 years
316-3, 516-3	16"	1"	3.75"	Firmgrip/knurled	1 lb. 4 oz.	250 lb.	3	3 years
318-3, 518-3	18"	1"	3.75"	Firmgrip/knurled	1 lb. 8 oz.	250 lb.	3	3 years
324-3, 524-3	24"	1"	3.75"	Firmgrip/knurled	1 lb. 12 oz.	250 lb.	3	3 years
332-3, 532-3	32"	1"	3.75"	Firmgrip/knurled	2 lb.	250 lb.	3	3 years

* All models are not ANSI or ADA compliant.



Invacare® I•Class™ Tub Bar

Model no. 710-4

The Invacare I•Class tub bar provides two gripping surfaces. The upper gripping position helps consumers enter the tub without turning the wrist or body. The lower gripping position helps consumers to lower into or rise from a seated position in the tub.




New

Invacare® Bath Board

Model no. H112-5

The lightweight Invacare Bath Board not only features an easy to clean surface and a stable handle for secure transfer, it also looks great.



New

Invacare® Bath Mat

Model no. H190

Designed to provide added safety and security, it is held firmly in place with suction cups and features a non-slip surface for added peace of mind.

Specifications

Model no.	Height	Length	Width	Adjustment range	Weight	Weight capacity	Units/ carton	Limited warranty
710-4	16"	8.25"	—	1.25" – 6"	5 lb.	250 lb.	4	3 years
H112-5	6.25"	27"	13"	15" – 24.75"	3.5 lb.	330 lb.	5	3 years
H190	.4"	32"	16"	—	3.75 lb.	286 lb.	6	1 year

Invacare® Hand Showers

Invacare hand showers are lightweight and feature a wall bracket for mounting on the shower wall and a white finish.

Deluxe Hand Shower

Model no. 828-4



- Adjustable head provides a choice of fine or firm spray, pulsating massage spray or a combination spray
- Comes standard with 7' chrome shielded hose
- On/Off switch

Standard Hand Shower

Model no. 826-4



- Adjustable head provides a choice of fine, firm or pulsating spray
- Comes standard with 7' nylon reinforced hose

Specifications

Model no.	Units/carton	Limited warranty
828-4	4	1 year

Specifications

Model no.	Units/carton	Limited warranty
826-4	4	1 year

Invacare® Toilet Safety products make the bathroom a safer environment for individuals with special needs. Choose from products with a variety of features including a weight capacity of 300 lb. and adjustable arms on the 1302RTS.



New

Invacare® Clamp-On Raised Toilet Seats

Model no. 1300RTS

- Recommended for use with toilet safety frame
- Attaches without tools and fits on most round bowl configurations
- Rear seat wings for more secure attachment
- Contoured seat for internal rotation of hips
- Ample depth of front recess makes hygiene performance more convenient



Model no. 1302RTS

New

Invacare® Raised Toilet Seat

Model no. 1301RTS without arms
Model no. 1302RTS with arms (shown)



- Comfortable, contoured seat is made of easy-to-clean, durable molded composite
- Clamp-on locking mechanism tightens easily with one hand
- Fits on most round bowl configurations
- New adjustable locking arms (1302RTS)

Invacare® Toilet Safety Frame (not shown)

Model no. 1392KD



- Arms fold back for cleaning and are easily removed for transfers
- Armrests provide secure hand support
- Frame is made of anodized, non-corrosive 1" aluminum tubing
- Aluminum bracket is easy to attach

Toilet Seat Specifications

Model no.	Seat height	Seat width	Seat depth	Arm height above seat	Overall width with arms	Product weight	Weight capacity	Units/carton	Limited warranty
1300RTS	5"	14.75"	15.5"	—	—	3 lb. 6 oz.	300 lb.	4	3 years
1301RTS	5"	15.8"	16.5"	—	—	4.5 lb.	300 lb.	2	3 years
1302RTS	5"	15.8"	16.5"	10"	19.5"	5.5 lb.	300 lb.	2	3 years

Toilet Safety Frame

Model no.	Arm height from floor	Distance between arms	Weight	Weight capacity	Units/carton	Limited warranty
1392KD	25.75" – 30" min. - max.	17.5" – 20.5" min. - max.	2 lb. 5 oz.	250 lb.	2	3 years



INVACARE®

Commode

Products

- Invacare® Commodes
- Invacare® Mariner™ Rehab Shower Commode Chair

The Invacare® All-in-One Commode offers consumers the comfort and stability they need. The frame is lightweight. With its 16" wide seat, consumers gain additional comfort and support. The Invacare All-in-One Commode can be used bedside, or can act as a toilet safety frame or raised toilet seat.



Model no. 9630-1



Large 16" wide snap-on seat for added comfort and support.

Invacare® I-Class™ All-In-One Steel Commode

Model no. 9630-1 (shown, gray-coated steel, one per carton), 9630-4 (gray-coated steel, four per carton)

- Can be used bedside, or with backrest removed can act as a toilet safety frame or raised toilet seat
- Large 16" wide snap-on seat with lid for added comfort and support
- Comes complete with commode bucket, bucket cover and splash shield
- Wing nuts allow for removal of backrest when needed



Model no. 9650-4



Large 16" wide snap-on seat for added comfort and support.

Invacare® All-In-One Aluminum Commode

Model no. 9650-4

- Lightweight aluminum frame is rust resistant
- Can be used bedside, or with backrest removed can act as a toilet safety frame or raised toilet seat
- Large 16" wide snap-on seat with lid for added comfort and support
- Comes complete with commode bucket, bucket cover and splash shield
- Wing nuts allow for removal of backrest when needed

Specifications

	Base width overall	Width between arms	Seat surface width	Seat opening width	Base depth overall	Seat surface depth	Seat opening depth	Height overall min.-max.	Seat-to-floor height min.-max.	Seat-to-arm height	Product weight	Weight capacity	Units/ carton	Limited warranty
9630-1, 9630-4	25"	19"	16"	8.5"	19"	14"	10"	23" – 29"	16" – 22"	7.5"	12 lb.	300 lb.	4	5 years
9650-4	25"	19"	16"	8.5"	19"	14"	10"	23" – 29"	16" – 22"	7.5"	8 lb.	300 lb.	4	5 years



Model no. 9610-4

Invacare® Safeguard Commode

Model no. 9610-4 (gray-coated steel, four per carton)

- Can be used as a raised toilet seat or a toilet safety frame
- Equipped with a splash shield
- Snap-on seat

The Invacare folding commode offers all the great features that you expect. It folds up in seconds for storage or to transport. It comes with armrests, angled legs for stability, non-marring rubber tips, anti-rattle collars that reduce noise, snap-on removable back, pail, splash shield and lid.



Model no. 9830-4



Specifications

	Base width overall	Width between arms	Seat surface width	Seat opening width	Base depth overall	Seat surface depth	Seat opening depth	Height overall min.-max.	Seat-to-floor height min.-max.	Seat-to-arm height	Product weight	Weight capacity	Units/ carton	Limited warranty
9610-4	25"	19"	16"	8.5"	19"	14"	10"	23" – 29"	16" – 22"	7.5"	12 lb.	300 lb.	4	5 years
9830-4	21.5"	18"	13.25"	8.5"	18.5"	14.5"	9.75"	25" – 31"	16.5" – 21.5"	8"	13 lb.	300 lb.	4	5 years



Model no. 9669

Invacare® Drop-Arm Commode

Model no. 9669 (one per carton)

- Can be used over the toilet with optional splash shield
- Accommodates transfers in restricted areas
- Drop-arm lever operates simply: lift to lock or push to release
- Gray powder coated, 1" welded steel tubing
- Snap-on seat and lid

The Invacare® Bariatric Drop-Arm Commode is built with durable steel construction and is designed with a wide and deep seating surface for individuals weighing up to 1000 lb. The commode comes with easy to operate drop arms that accommodate transfers in and out of the commode. The drop-arm commode also has an extra-wide snap-on seat and backrest and also has height adjustable legs.

**New
Elongated Seat**



Model no. 6599

Invacare® Bariatric Drop-Arm Commode

Model no. 6599

- 1000 lb. weight capacity
- Easy-to-operate drop-arm mechanism
- Extra-wide snap-on seat
- Removable padded backrest
- Adjustable height legs
- Durable steel frame design
- Comes complete with pail and lid

Specifications

	Base width overall	Width between arms	Seat surface width	Seat opening width	Base depth overall	Seat surface depth	Seat opening depth	Height overall min.-max.	Seat-to-floor height min.-max.	Seat-to-arm height	Product weight	Weight capacity	Units/carton	Limited warranty
9669	20.25"	17.5"	14"	8.25"	20.5"	15"	9.5"	27" – 31.5"	17" – 21"	8"	22 lb.	250 lb.	1	5 years
6599	28"	25"	23.5"	8.5"	20.5"	18"	12"	27" – 31"	19.25" – 23.25"	7.75"	21 lb.	1,000 lb.	1	5 years

New



Model no. 6500-BHD



Invacare® Bariatric Commode

Model no. 6500-BHD

- 650 lb. weight capacity
- Extra-wide, snap-on seat
- Rectangular pan can be removed from the front or rear of the commode
- Adjustable seat-to-floor height

Invacare® Bariatric and Heavy-Duty Commodes combine comfort with strength and durability. They each feature armrests, angled legs for stability, non-marring rubber tips and anti-rattle collars that reduce noise.



Model no. 6497



Invacare® Heavy-Duty Commode

Model no. 6497

- 450 lb. weight capacity
- Designed for bedside use
- Rectangular pan and pan holder
- Armrests for added comfort

Specifications

	Base width overall	Width between arms	Seat surface width	Seat opening width	Base depth overall	Seat surface depth	Seat opening depth	Height overall min.-max.	Seat-to-floor height min.-max.	Seat-to-arm height	Product weight	Weight capacity	Units/carton	Limited warranty
6500-BHD	36"	30"	27.75"	9"	21"	18.5"	12"	25" – 29"	17" – 21"	8.25"	20 lb.	650 lb.	2	5 years
6497	30"	21.75"	15.5"	8.25"	19"	15.25"	12"	25" – 31"	16.25" – 22"	9"	16 lb.	450 lb.	2	5 years

Invacare® Pan and Holder

Part no. 13245X000 (pan)
Part no. 1027004 (holder)

- 8-quart pan capacity



Specifications

Part no.	Length overall	Length inside	Width overall	Width inside
13245X000	16.75"	14"	10"	9"

Invacare® Pail with Lid

Model no. 6317 (12 per carton)

- 12-quart pail capacity



Invacare® Splash Shield

Model no. 6319 (12 per carton)



Invacare® Mariner™ Rehab Shower Commode Chairs offer aluminum frames and stainless steel hardware that are rust-resistant, making them ideal for use in the shower. These chairs fit over most standard and elongated bowls for over-the-toilet use as well. They feature flip-back padded arms that lock down for safety, swing-away front riggings with tool-less adjustable height footrests and a four-position padded seat. The Invacare Mariner chair offers various seat widths and wheel options for a custom fit along with additional features and benefits.



Features

- Aluminum frame and stainless steel hardware
- Fits over most standard and elongated toilet bowls
- Padded four-position seat
- Flip-back padded arms that lock down for safety
- Swing-away front riggings
- Tool-less adjustable height footrests
- Folds into a compact unit for storage and transport
- Removable pan and pan holder
- Safety strap included
- 23" treaded urethane tires or 5" casters

Model nos. 6895 - 23" treaded urethane tires with 18.5" seat (shown).
 6795 - 23" treaded urethane tires with 16" seat
 6891 - Four 5" locking casters with 18.25" seat



Flip-back padded arms that lock down for safety.



Padded four-position seat.



Swing-away front riggings with tool-less adjustable height footrests.

Specifications	6895	6795	6891 w/5" casters
Overall chair width (handrim to handrim)	26.5"	24.5"	22.5" (no handrim)
Overall back width	18"	16"	18"
Overall chair height (floor to top of push handle)	39"	39"	39"
Overall depth w/o front riggings	32"	32"	29"
Seat width	18.25"	16"	18.25"
Seat depth	18.25"	16"	18.25"
Width between arms	18.5"	16.5"	18.5"
Folded width	12"	12"	8"
Seat-to-floor height	22.75"	22.5"	22.75"
Seat-to-armrest height	8"	8"	8"
Back height (top of seat to top of upholstery)	15.5"	15"	15.5"
Footrest extension min. - max. (from top of seat)	13" – 18"	13" – 18"	13" – 18"
Product weight	41 lb.	39 lb.	32 lb.
Weight capacity	300 lb.	300 lb.	300 lb.
Units/carton	1 per carton	1 per carton	1 per carton
Limited warranty	3 years	3 years	3 years
Clearance for over-the-toilet use			
Depth (front cross brace to edge of rear tire)	26.5"	26.5"	23.5"
Depth (rear wheel axle bolt to edge of tire)	11.5"	11.5"	–
Height (floor to seat base)	20.25"	20.5"	20.75"
Inside width (between side frames)	16"	14"	16"

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

© 2011

Invacare Corporation
All rights reserved.
Trademarks are identified
by the symbols ™ and ®.
All trademarks are owned
by or licensed to Invacare
Corporation unless otherwise
noted. Velcro is a registered
trademark of Velcro
Industries B.V.
Specifications are subject
to change without
notification. Patents:
5,529,425 D480,
666, D480, 995 and D483,
702.5, 529,425 D480,666
D480,995 and D483,702.
Other Patents pending.

Form No. 00-267
Rev. 09/11
110837

