

Top End® XLT Jr. Handcycle

Price List and Order Form

Effective January 1, 2006 Revised July 1, 2014

For accuracy of ordering, please fax order to Top End at 727-522-1007 or e-mail to TopEndOrders@invacare.com

For questions call 800-532-8677 or 727-522-8677.

To order Literature: 800-828-6282

www.topendwheelchair.com



Provider Name _____ Date _____ Top End Job Number _____
 Account # _____ P.O. # _____ Phone # _____ Fax # _____
 Address _____ City _____ State _____ Zip Code _____
Provider Primary Contact: _____ **Email:** _____ **Phone No** _____
 Athlete Name _____ Email _____ Phone No _____
 Address _____ City _____ State _____ Zip Code _____
Ship order to: Provider Customer **Order will ship to Provider unless otherwise noted.**

CHAIR TYPE

XLTJR Top End Exceleator XLT Jr. Handcycle.....\$2395.00

Seat Width: 14"

Frame Color: Grasshopper Green

OTHER THAN STANDARD FRAME COLOR:

P3.....\$100.00

CONTEMPORARY COLORS

Wet Black White Snow Silver Metallic

Stealth Black

WAY COOL COLORS


Royal Blue Midnight Blue

TOO HOT COLORS

Tangerine Red Raspberry Red Grape Madness
 Safety Yellow Orange Bubblegum Pink

SAFETY AND COMFORT OPTIONS:

- | | | |
|--------------------------|--|----------|
| <input type="checkbox"/> | Seat and Back Cushions | Standard |
| <input type="checkbox"/> | Chain Guard..... | Standard |
| <input type="checkbox"/> | Safety Flag..... | Standard |
| <input type="checkbox"/> | Chest restraint: hook and loop..... | Standard |
| <input type="checkbox"/> | Quad twist shifter option | 50.00 |
| <input type="checkbox"/> | Horizontal foam covered handpedals..... | 75.00 |
| <input type="checkbox"/> | Mirror..... | 30.00 |
| <input type="checkbox"/> | Water Bottle & Cage | 35.00 |
| <input type="checkbox"/> | Helmet: choose size..... | 75.00 |
| | <input type="checkbox"/> Small/ Medium <input type="checkbox"/> Medium/Large | |
| <input type="checkbox"/> | Front Safety Light | 25.00 |
| <input type="checkbox"/> | Rear Flashing Light..... | 25.00 |
| <input type="checkbox"/> | Crutch Holder & Strap | 125.00 |
| <input type="checkbox"/> | Leg Guard Attachment | 75.00 |
| <input type="checkbox"/> | Bike Rack for Trailer Hitch..... | 450.00 |
| | 1 1/4" (Standard) <input type="checkbox"/> 2" (No Charge) | |
| <input type="checkbox"/> | Freight* XLT 20 | \$ _____ |



Includes:

- Aluminum Frame/Chrome-moly Fork
- Adjustable Seat
- Adjustable Back Angle
- Adjustable Crank Height
- Chain Guard, cushion, safety flag
- Chain Tension Idler
- Adjustable Footrest System
- Shimano® 7 Speed Hub w/Reversing Drum Brake
- 155 mm Crankset w/one Chainring
- Parking Brake
- 20" Cruiser Wheels and Tires
- Quick Release Axles
- 15° Camber
- Road Crown Compensator
- Reflectors
- Chest Restraint

* Indicates an item where the standard discount does not apply or there are additional shipping charges or both. Please request a quote.

Top End® Excelerator XLT™ Jr. Handcycle Specifications

Frame – 2.5 x .065" ovalized 6061 aluminum aircraft tubing.

Frame Design – Single main tube, recumbent position.

Fork Assembly – 4130 Chrome-moly for strength/stiffness. Raked angle for cornering stability.

Wheelbase – Overall length 64".

Weight – Overall weight 30 – 35 lb.

Turning Radius – 10 feet for 360° circle.

Static Tip Angle - 39° for 250 lbs. subject.

Transfer Seat Height – Seat height front w/o cushion is 12" and rear seat height is 10" w/o cushion.

Back Style/Height-Standard narrow back is 7" wide. Back height is fixed at 16.5".

Upholstery – Nylon w/hoop and loop adjustable tension seat sling with seat and back cushion. Hook and loop chest/seat restraint included.

Crankset – Sugino 155mm w/chainring, adjustable height. Includes chain tension idler.

Chainguard- Full length shock resistant composite construction

Steering – Pivot steering w/raked fork for controlled turning and stability. Road crown compensator is standard.

Standard Propulsion System - -150mm Crankset w/one chainring, adjustable height crank, Shimano 7 speed front hub w/internal hub brake and 7 speed twist shifter.

Braking System – Shimano internal drum brake (coaster brake) and Chain, activated by pedaling backwards.

Gear Shifter- Shimano Tourney revoshifter indexed twist grip with gear display. For customers with limited hand function the Quad twist shifter option, a 3 pronged adaptation that is secured around the twist shifter is recommended.

Handpedals – Oval Ergonomic handpedals. Optional horizontal style is available.

Options for impaired hand function–The quad twist shifter adaptation option will allow the users to easily change gears.

Parking Brake – Caliper brake with locking lever.

Footrest – Multi-adjustable for and aft w/straps to securely cradle calves and feet. Leg guard attachments are recommended for persons who have no leg movement. The leg guard option will protect the legs from getting into the chain.

Performance Wheels – Rear hubs are precision heat treated, silver anodized w/quick release stainless axles. Quick release axles are easiest for transporting the handcycle. 20" x 1.75" treaded Tire,
20 x 1.5" Tube 40 psi,
20" Spoke Wheels

Camber - 15° for stability.

Recommended Options – Helmet, Cushion, and Water Bottle and Cage(as this is custom to fit on this bike).

Warranty – 3 year on frame to original owner / limited parts warranty.

Assembly – Shipped completely assembled except for rear wheels. Shipping weight is 50 lb.

Options

- Cordless Computer/Speedometer
- Water Bottle and Cage
- Helmet
- Quad twist shift adaptation
- Safety Lights
- Mirror
- Leg Guard Attachment
- Crutch Holder
- Bike Rack



Invacare Corporation www.invacare.com

TOP END
Sports and Recreation
Products
4501 63rd Circle North
Pinellas Park, FL USA
33781-5914

All rights reserved. Trademarks are identified by the symbols TM and [®]. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Shimano is a registered trademark of Shimano, Inc. Specifications and prices are subject to change without notification. ©2014 Invacare Corporation