

# Invacare® 3G Arrow® Power Wheelchair w/TrueTrack



## Price List and Order Form

Effective May 18, 2009

For ease of ordering, contact Customer Service Toll Free at:

800-333-6900 Fax: 800-678-4682

To order Literature: 800-828-6282

[www.invacare.com](http://www.invacare.com)

Company Name \_\_\_\_\_ Account # \_\_\_\_\_

Phone \_\_\_\_\_ P.O. # \_\_\_\_\_ Date \_\_\_\_\_

Ship to Name & Address \_\_\_\_\_

### Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

### CHAIR TYPE

- 3GAR Arrow Base & Seat w/MK6i™ (K0868/K0870).. 9,775.00 <sup>1</sup>
- 3GARBASE Arrow Base Only w/MK6i..... 9,020.00 <sup>1</sup>
- 3GSEAT Seat Only ..... 755.00 <sup>2</sup>

### Group 4 – K0868 Basic Package

- 3GAR Arrow with MK6i
- 300 lb Weight Capacity
- 1431-3 14" X 3" Tires w/Foam Filled Inserts..... No Charge
- TrueTrack Gearless Brushless GB™ Motors.....Standard
- 1228-3 8" X 2.25" Semi Pneumatic Casters.....Standard
- 28F Desk Length Fixed Height Arm.....Standard
- MK6TT TrueTrack Controller.....Standard
- 70NHD 70° Swingaway Footrests.....Standard
- 8 Amp Off Board Battery Charger.....Standard
- 1515 Push Button Style Seat Positioning Strap.....Standard

### Group 4 – K0870 Basic Package

- 3GAR Arrow with MK6i
- 400 lb Weight Capacity
- 1431-3 14" X 3" Tires w/Foam Filled Inserts.....No Charge
- TrueTrack Gearless Brushless GB™ Motors..... Standard
- 1228-3 8" X 2.25" Semi Pneumatic Casters..... Standard
- 28F Desk Length Fixed Height Arm..... Standard
- MK6TT TrueTrack Controller..... Standard
- 70NHD 70° Swingaway Footrests..... Standard
- 8 Amp Off Board Battery Charger..... Standard
- 1515 Push Button Style Seat Positioning Strap..... Standard

### TRANSPORT TIE DOWN

- TRRO Transport Ready Option..... 300.00
- TRBKTS Wheelchair Transport Brackets..... 250.00

#### TRANSPORT TIE-DOWN

TRRO includes four factory-installed transport brackets and a wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNA WC Vol 1. Section 19 Frontal Impact Test requirements for wheelchairs with a 168 lb crash dummy, which corresponds to a person with a weight of 114 lbs to 209 lbs.

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry.

Invacare cannot and does not recommend any wheelchair transportation systems.

\* Transport Tie-Down option is not retrofittable to existing models and is not field serviceable.

#### TRRO FEATURES FOR 3G ARROW

Features Required for TRRO	Optional Features for TRRO
<ul style="list-style-type: none"> <li>• Storm Series Arrow Base &amp; Seat</li> <li>• 300 lb or 400 lb user weight limit</li> <li>• Short or Long Base Frame</li> <li>• Adjustable Rehab Seat with Adjustable Seat Tilt &amp; Back Angle</li> <li>• Low Seat-to-Floor Height</li> <li>• Sectional or Straight Rear Uprights</li> <li>• 18"- 24" ASBA Back Angle</li> <li>• 8" x 2.25" Semi Pneumatic Casters</li> <li>• Reinforced Back Uph.</li> <li>• 14" x 3" Tires w/Foam Filled Inserts</li> <li>• Rear Wheel Standard Location</li> <li>• Black Nylon or Black Vinyl Uph.</li> <li>• Center Pivot Style Bracket</li> <li>• 60° Swingaway, 70° Swingaway or 70° Tapered Footrests</li> </ul>	<ul style="list-style-type: none"> <li>• Fender</li> <li>• Flip-Up Composite Footplates, Aluminum Footplates-Deluxe or Adjustable Angle Flip Up Footplates</li> <li>• Calf Strap</li> <li>• Impact Guards</li> <li>• All PSR/PSF Joysticks</li> <li>• Left Hand Mounted Joystick</li> <li>• "T" Handle Flexible Joystick Ext.</li> <li>• Straight Handle Flexible Joystick Ext.</li> <li>• MK6i Remote Programmer</li> <li>• Group 24 Batteries</li> </ul>
<p><b>Please note that all TRRO chairs are set at 10° seat angle</b></p>	

### FOOTNOTES

<sup>1</sup> Arrow comes standard with 300 lb or 400 lb user weight capacity, 7 MPH top speed, TrueTrack driving technology w/GB™ motors, MK6i MPJ™+ joystick, MK6TT controller, Base only does not come standard with Seat Positioning Strap.

<sup>2</sup> When ordering a replacement seat and changing depth to 20" or deeper, must also order new arms.

**NOTE: All specifications and dimensions are approximate.**  
For the most current pricing information go to [www.invacare.com](http://www.invacare.com).

# Invacare® 3G Arrow® Power Wheelchair w/TrueTrack Accessories and Options

## USER WEIGHT LIMITS

- U300 User Weight <300 lb (K0868)..... No Charge
- U400 User Weight <400 lb (K0870)..... No Charge

## FRAME TYPE

- SFS Short Frame ..... No Charge <sup>3</sup>
- LFS Long Frame ..... No Charge

## SEAT SELECTION

### ADJUSTABLE SEAT WIDTH RANGE

- 1216W 12" to 16" Seat Width Range ..... No Charge
- 1620W 16" to 20" Seat Width Range ..... No Charge
- 2024W 20" to 24" Seat Width Range ..... No Charge

### FRAME WIDTH

- W12 12" Wide ..... 110.00 <sup>4</sup>
- W13 13" Wide ..... 110.00 <sup>4</sup>
- W14 14" Wide ..... 110.00 <sup>4</sup>
- W15 15" Wide ..... No Charge <sup>4</sup>
- W16 16" Wide ..... No Charge
- W17 17" Wide ..... No Charge
- W18 18" Wide ..... No Charge
- W19 19" Wide ..... No Charge
- W20 20" Wide ..... No Charge
- W21 21" Wide ..... 329.00
- W22 22" Wide ..... 329.00
- W23 23" Wide ..... 329.00 <sup>5</sup>
- W24 24" Wide ..... 469.00 <sup>5</sup>

### VAN SEAT WIDTH <sup>6</sup>

- VW18 18" Wide ..... No Charge
- VW20 20" Wide ..... 329.00
- VW22 22" Wide ..... 329.00

### SEAT SELECTION

- ADJASBA Adj. Rehab Seat w/Adj. Seat Tilt & Back Angle ..... **Standard** <sup>7</sup>
- VAN Van Style Captain Seat ..... 100.00
- ADJRECL Manual Recline for Adjustable ASBA (E1226) .. 744.00
- CURV Curved ..... No Charge
- CTB Contoura Seating (E2620/E2621) ..... 625.00
- OMCU Omit Contoura Seat Cushion..... Less 20.00

### ADJUSTABLE SEAT DEPTH RANGE

- 1218D 12" to 18" Seat Depth Range..... No Charge <sup>8</sup>
- 1619D 16" to 19" Seat Depth Range..... No Charge
- 1922D 19" to 22" Seat Depth Range..... No Charge

## FOOTNOTES

- <sup>3</sup> Not available on seat depths of 1922D (19" to 22" seat depth range).
- <sup>4</sup> Jr. Seat widths/depths require a straight rear upright and a user weight limit of 150 lbs. Manual Reclining back not available. Additional charges apply to widths/depths less than 16". Jr. seat widths have Jr. arms with an adjustment range of 6"-10" inches in height.
- <sup>5</sup> The 23" & 24" seat widths are not available with low seat-to-floor. Must pick Heavy Duty Pkg to order 23" & 24" seat widths.
- <sup>6</sup> Manual Reclining Back not available. Van Seat Depth set at the following:  
18 wide set @ 16" deep  
20 wide set @ 18" deep  
22 wide set @ 20" deep
- <sup>7</sup> Full Length adjustable arms are standard on Van Style Seat.
- <sup>8</sup> Must select 12"-18" depth range(1218D) when ordering the 12"-16" seat width range (1216W). 12" & 13" widths are not available with 20"-24" back heights.

## SEAT DEPTH SETTING <sup>9</sup>

- D12 Seat Depth 12" (Short Base)..... 110.00 <sup>4</sup>
- D13 Seat Depth 13" (Short Base)..... 110.00 <sup>4</sup>
- D14 Seat Depth 14" (Short Base)..... 110.00 <sup>4</sup>
- D15 Seat Depth 15" (Short Base)..... No Charge <sup>4</sup>
- D16 Seat Depth 16" (Short Base)..... No Charge
- D17 Seat Depth 17" (Short Base)..... No Charge
- D18 Seat Depth 18" (Short Base)..... No Charge
- D19 Seat Depth 19" (Long Base)..... No Charge
- D20 Seat Depth 20" (Long Base)..... No Charge
- D21 Seat Depth 21" (Long Base)..... 529.00
- D22 Seat Depth 22" (Long Base)..... 529.00

## SEAT-TO-FLOOR HEIGHT <sup>10</sup>

- LOW Seat-To-Floor – 17.75"..... No Charge
- MED Seat-To-Floor – 20" ..... No Charge
- TALL Seat-To-Floor – 21.25"..... No Charge

## VAN SEAT STYLE <sup>11</sup>

- SEMIR Semi-Recline..... No Charge
- FULLR Full-Recline (E1126) ..... 469.00

## ADJASBA BACK TYPE

- SECT Sectional Rear Uprights..... No Charge
- SRUP Straight Rear Uprights ..... No Charge
- 9720P No Push Handles (Straight Canes) ..... 82.00
- L8661 Omit Rear Uprights..... Less 20.00

## RECLINER BACK HEIGHT <sup>12</sup>

- BHR18 Back Height 18"..... No Charge
- BHR20 Back Height 20"..... No Charge
- BHR22 Back Height 22"..... No Charge
- BHR24 Back Height 24"..... No Charge
- BHR26 Back Height 26"..... No Charge

## BACK SELECTION

### ADJASBA BACK HEIGHT

- BH12 Back Height 12"..... No Charge
- BH13 Back Height 13"..... No Charge
- BH14 Back Height 14"..... No Charge
- BH15 Back Height 15"..... No Charge
- BH16 Back Height 16"..... No Charge
- BH17 Back Height 17"..... No Charge
- BH18 Back Height 18"..... No Charge
- BH19 Back Height 19"..... No Charge
- BH20 Back Height 20"..... No Charge
- BH21 Back Height 21"..... No Charge
- BH22 Back Height 22"..... No Charge
- BH23 Back Height 23"..... No Charge
- BH24 Back Height 24"..... No Charge

### CONTOURA

- TSBH25 Back Height 25"..... No Charge
- TSBH28 Back Height 28"..... No Charge

## FOOTNOTES

- <sup>9</sup> Contoura Seating available in depths 16" - 22".
- <sup>10</sup> Based on 5° tilt and seat depth range of 16-19" : add .25" for seat depth range of 19-22"
- <sup>11</sup> Used only with Van style seat selection.
- <sup>12</sup> Used only with manual reclining seat selection.

# Invacare® 3G Arrow® Power Wheelchair w/TrueTrack Accessories and Options

## BACK ANGLE DEGREE

<input type="checkbox"/> 85BA	85°	No Charge
<input type="checkbox"/> 90BA	90°	No Charge
<input type="checkbox"/> 95BA	95°	No Charge
<input type="checkbox"/> 100BA	100°	No Charge
<input type="checkbox"/> 105BA	105°	No Charge

## BACK ACCESSORIES

<input type="checkbox"/> L8601	Omit Back Upholstery	No Charge
<input type="checkbox"/> RBP	Reinforced Back Upholstery	No Charge

## ARM SELECTION

### ARM TYPE <sup>13</sup>

<input type="checkbox"/> 18F	Full Length Fixed Hgt Conv Flip Up	No Charge
<input type="checkbox"/> 28F	Desk Length Fixed Hgt Conv Flip Up	<b>Standard</b>
<input type="checkbox"/> 19F	Full Length Adj Hgt Conv Flip Up (E0973)	230.00
<input type="checkbox"/> 29F	Desk Length Adj Hgt Conv Flip Up (E0973)	230.00

### CANTILEVER ARM TYPE (PAIR) <sup>14</sup>

<input type="checkbox"/> CLVF	Full Length Arm	\$350.00
<input type="checkbox"/> CLVD	Desk Length Arm	\$350.00

### CANTILEVER MIXED ARM TYPE (EACH)

<input type="checkbox"/> CLVFR	Full Length Arm Right	\$175.00
<input type="checkbox"/> CLVDR	Desk Length Arm Right	\$175.00
<input type="checkbox"/> CLVFL	Full Length Arm Left	\$175.00
<input type="checkbox"/> CLVDL	Desk Length Arm Left	\$175.00

### ARMPAD OPTIONS (CANTILEVER ARMS)

<input type="checkbox"/> WF	Waterfall Armpads	<b>Standard</b>
<input type="checkbox"/> RCL	Reclining Armpads	No Charge

## MK6i ELECTRONICS

### MK6i CONTROLLER

<input type="checkbox"/> MK6TT	TrueTrack Controller	<b>Standard</b>
--------------------------------	----------------------	-----------------

### EXPANDABLE CONTROLLER SYSTEMS

<input type="checkbox"/> PWH	Harness req. for Expandable Electronics (E2313)	350.00	<sup>15</sup>
<input type="checkbox"/> MPJM6	MPJ+ Multiple Drive Joystick (E2377)	750.00	
<input type="checkbox"/> PSFM6	Personal Joystick, Inductive on Front	1130.00	<sup>16</sup>
<input type="checkbox"/> PSRM6	Personal Joystick, Inductive on Rear	1130.00	<sup>17</sup>
<input type="checkbox"/> PBOD	Push Button On/Off w/PB Drive Select	No Charge	<sup>17</sup>
<input type="checkbox"/> PBSS	Push Button On/Off w/Speed Select	No Charge	<sup>18 19</sup>
<input type="checkbox"/> OR	On/Off Mounted Right	No Charge	<sup>18</sup>
<input type="checkbox"/> OL	On/Off Mounted Left	No Charge	<sup>18</sup>
<input type="checkbox"/> ASLRDYM6	MK6 Display Only (No Driver Control) (E2377)	750.00	<sup>18</sup>

**FOR YOUR ASL DRIVE SYSTEM NEEDS  
PLEASE REFER TO THE FOLLOWING:  
ASL SINGLE INVOICE PRICE LIST FORM # 02-097.**

<input type="checkbox"/> 1500M6	Prop RIM Head Control w/MK6i Display (E2327)	2875.00
<input type="checkbox"/> SMHD	Small Headrest	No Charge
<input type="checkbox"/> LGHD	Large Headrest	No Charge
<input type="checkbox"/> 1558MM6	Compact Joystick w/Display (E2373)	1130.00

## FOOTNOTES

- <sup>13</sup> All arms are flip up and detachable. Seat widths 12-15" not available with fixed height arms must use adjustable height arm style only.
- <sup>14</sup> Cantilever Arms only available with 300 lb. rated systems. Not available with BPO, Otto Bock or 12"-15" seat width. Must use hard back. See chart on last page for arm height adjustment range.
- <sup>15</sup> Harness is required to interface joystick or MK6i display with expandable controllers. Option is located under Electronic Accessories in configurator.
- <sup>16</sup> TOD-SS: Toggle On/Off-Drive Select, Speed Potentiometer Standard.
- <sup>17</sup> PSRM6 Only.
- <sup>18</sup> Display only, allows choice from order form.

## ALTERNATIVE DIGITAL CONTROL SYSTEMS

<input type="checkbox"/> SNPM6	Sip-N-Puff Digital Interface Box (E2325)	1,500.00	<sup>19</sup>
<input type="checkbox"/> 5018M6	Wafer Board & Digital Interface Box (E2322)	1,550.00	<sup>20</sup>
<input type="checkbox"/> 5020M6	Mini Tash Joystick & Digital Interface Box (E2321)	1,750.00	<sup>20</sup>

## SIP-N-PUFF KIT

<input type="checkbox"/> PKG32666	Therafin® Sip-N-Puff Breath Tube Kit (E2326)	380.00
-----------------------------------	--	--------

## ADDITIONAL CONTROL CHOICES

<input type="checkbox"/> 1558M6	Compact Joystick - Less Std Mount	770.00	
<input type="checkbox"/> 5018	Wafer Board	310.00	<sup>20</sup>
<input type="checkbox"/> 5020	Mini Tash Joystick	405.00	<sup>21</sup>
<input type="checkbox"/> PACM6	Prop. Attendant Control w/quick release mount	850.00	<sup>21</sup>
<input type="checkbox"/> 1552M	Digital Attendant Control	850.00	<sup>21</sup>

## JOYSTICK MOUNTING HARDWARE

<input type="checkbox"/> QLAM6	Quad Link Retract. Joystick Mnt MPJ+ (E1028)	279.00	<sup>21</sup>
<input type="checkbox"/> QLM6	Quad Link Retract. Joystick Mnt MPJ+ Van (E1028)	230.00	<sup>22</sup>
<input type="checkbox"/> ARM250	Stealth Height Adjustable Joystick Mount (E1028)	210.00	
<input type="checkbox"/> LEFT	Left Handed Mounted Joystick	No Charge	
<input type="checkbox"/> RIMHW	RIM Hardware For 1500 (E1028)	219.00	
<input type="checkbox"/> PKG32669	Complete Bib Assembly (Tash/1558)	265.00	
<input type="checkbox"/> GATMPJ6	Gatlin Midline for MPJ (E1028)	600.00	
<input type="checkbox"/> GAT1812	Gatlin Midline for 1812 (E1028)	600.00	
<input type="checkbox"/> RCM	Rear Cane Mount for Compact Joystick	395.00	

## ELECTRONIC ACCESSORIES

<input type="checkbox"/> 1813M6	MK6i Programmer w/Pro Memory Card	475.00	
<input type="checkbox"/> 1813SD	Pro Memory Card (USB Ready)	90.00	
<input type="checkbox"/> 1560	"T" Handle Flexible Joystick Ext. (E2323)	77.00	
<input type="checkbox"/> 1561	Straight Handle Flexible Joystick Ext. (E2323)	77.00	
<input type="checkbox"/> 1826	Chin Cup (E2324)	48.00	
<input type="checkbox"/> A24VPS	24V MK6 Auxiliary Power Source	200.00	
<input type="checkbox"/> AUX12M6	Auxiliary Module For 1 & 2	795.00	
<input type="checkbox"/> AUX34M6	Auxiliary Module For 3 & 4	795.00	
<input type="checkbox"/> MPY	Monoport "Y" Cable for MPJ+ or Display	15.00	<sup>22</sup>
<input type="checkbox"/> 2PB	Dual Push Button	147.00	<sup>23</sup>
<input type="checkbox"/> 2SW	Dual Toggle	172.00	<sup>23</sup>
<input type="checkbox"/> PBMT	Push Button Mounting Kit	252.00	<sup>24</sup>
<input type="checkbox"/> IR	Infrared Module with Wireless Mouse Control	995.00	<sup>24</sup>

## ASL OPTIONAL ACCESSORIES

<input type="checkbox"/> EGSBLK	Egg Switch - Black	156.00	
<input type="checkbox"/> ASL304A	Light Wobble Switch w/611 Mount	432.00	<sup>25</sup>
<input type="checkbox"/> ASL202J	Fiber Optic Switch	610.00	<sup>25</sup>
<input type="checkbox"/> ASL208J	Adjustable Proximity Switch	550.00	<sup>25</sup>
<input type="checkbox"/> ASL504M6	Remote Emergency Stop Switch	1,185.00	
<input type="checkbox"/> ASL611	Adjustable Switch Mounting Bracket	252.00	

## FOOTNOTES

- <sup>19</sup> Must also use either a Multiple Drive Joystick or MK6 Display.
- <sup>20</sup> Requires SNPM6 Sip-N-Puff Interface Box.
- <sup>21</sup> QLA includes Arm250 standard.
- <sup>22</sup> For use with MPJ+ or MK6i Display to provide additional switch site.
- <sup>23</sup> Option to replace the Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ joystick options.
- <sup>24</sup> Capable of controlling up to six remote Infra-Red controlled devices and also Computer Mouse operation through the driver control. Includes stereo switch (with hook/loop mount) for mouse clicks
- <sup>25</sup> Not available with SPJ+ Joystick. Must add MK6 Auxiliary Power Source

**NOTE: All specifications and dimensions are approximate.  
For the most current pricing information go to [www.invacare.com](http://www.invacare.com).**

# Invacare® 3G Arrow® Power Wheelchair w/TrueTrack Accessories and Options

## FRAME AND UPHOLSTERY COLOR

### FRAME FINISH

#### CONTEMPORARY

- 24P Wet Black ..... No Charge
- 60P Silver Metallic ..... No Charge
- 71P Silver Vein ..... No Charge
- 105P Midnight Blue ..... No Charge
- 115P Black Prism ..... No Charge

#### TOO HOT

- 30P Sunny Yellow ..... No Charge
- 61P Electric Red ..... No Charge
- 119P Bubblegum Pink ..... No Charge
- 120P Grape Madness ..... No Charge

#### WAY COOL

- 62P Electric Blue ..... No Charge
- 104P Electric Teal ..... No Charge
- 121P Lolly Pop Blue ..... No Charge
- 122P Cosmic Blue ..... No Charge
- 125P Emerald Green ..... No Charge
- 127P Grasshopper Green ..... No Charge

### UPHOLSTERY TYPE

- U240 Black Nylon Upholstery ..... **Standard**
- U67 Black Vinyl ..... No Charge <sup>26</sup>
- GC Grey Cloth ..... No Charge <sup>26</sup>
- GV Grey Vinyl ..... No Charge <sup>26</sup>
- TV Tan Vinyl ..... No Charge <sup>26</sup>

## CASTER AND WHEEL SELECTION

### CASTER OPTIONS

- 1259 6" X 2" Semi Pneumatic Casters ..... No Charge <sup>27</sup>
- 1228-3 8" X 2.25" Semi Pneumatic Casters Lt Grey ..... **Standard**
- 1251-3 8" X 2" Pneumatic Casters ..... No Charge
- 1220-3 9" X 2.75" Pneumatic Casters ..... No Charge <sup>28</sup>
- 1465-3 9" X 2.75" Foam Filled Casters ..... No Charge <sup>28</sup>

### CASTER SUSPENSION <sup>29</sup>

- 9213-6 Shock Fork For 6" Casters (E1016) ..... 119.00
- 9213-8 Shock Fork For 8" Casters (E1016) ..... 119.00

### WHEEL OPTION

- 1430-3 14" X 3" Pneumatic Tires ..... **Standard**
- 1431-3 14" X 3" Tires w/Foam Filled Inserts ..... No Charge

### REAR WHEEL LOCATION

- RWS Rear Wheel Standard ..... **Standard**

### FENDER OPTION

- L8603 Omit Fender ..... Less 50.00

### FOOTNOTES

- <sup>26</sup> Used only with Van style seat selection.
- <sup>27</sup> Only available with U300 lb weight capacity.
- <sup>28</sup> Not available with short base. . Not available with Transport Tiedown options
- <sup>29</sup> Not available with 9" casters.

## POSITIONING STRAPS

### SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Positioning Strap ..... 50.00 <sup>30</sup>
- 1515 Push Button Style Seat Positioning Strap ..... **Standard**

### CHEST POSITIONING STRAPS

- 7311BK Airline Buckle Chest Positioning Strap (E0960) ..... 100.00 <sup>30</sup>
- 9515 Push Btn Style Chest Positioning Strap (E0960) ..... 100.00
- 7321BK Hook & Loop Chest Pos. Strap (E0960) ..... 100.00

## MISC ACCESSORIES

### MISCELLANEOUS ACCESSORIES

- VANO2 O2 Holder - Van (E2208) ..... 150.00
- SEATO2 Oxygen Holder (ADJASBA Seat Only) (E2208) ..... 150.00
- 1303 Wheel Locks ..... 45.00 <sup>31</sup>
- 9205H-3 Heavy Duty Suspension Shock ..... 400.00
- ADD Adductors Small (E0956) ..... 200.00
- ADDHW Hardware for ADD (E1028) ..... 400.00
- ADD L Adductors Large (E0956) ..... 200.00
- ADDHWL Hardware for ADD L (E1028) ..... 400.00

## BATTERY SELECTION

### BATTERY BOXES

- GP24 Group 24 Battery Boxes ..... **Standard**

### BATTERIES

- 24BATTERY 24 Gel Battery (On Chair) (E2363) ..... 410.00

### BATTERY CHARGER OPTION

- 110CHARGER 110 Volt Battery Charger ..... No Charge
- LESTR Lester 8 AMP Battery Charger (E2367) ..... 320.00
- L8900 Omit Battery Charger ..... Less 80.00

### FOOTNOTES

- <sup>30</sup> Not available with 12"-15" seat widths.
- <sup>31</sup> Not available with pneumatic tires.

**NOTE: All specifications and dimensions are approximate.**  
*For the most current pricing information go to [www.invacare.com](http://www.invacare.com).*

# Invacare® 3G Arrow® Power Wheelchair w/TrueTrack Accessories and Options

## HANGER BRACKET TYPE

### HANGER BRACKET TYPE

- ASB Center Pivot Style Bracket ..... No Charge
- PSB Pin Style Bracket ..... No Charge
- NHB No Hanger Bracket..... No Charge

## FRONT RIGGING OPTIONS

### FRONT RIGGING

#### CENTER PIVOT

- 60HD 60° Swingaway Footrests 13.5"-17.5" ..... No Charge
- 70NHD 70° Swingaway Footrests 13.25"-17" ..... No Charge
- 90SWING 90° Swingaway Footrests ..... No Charge
- 70STAPER 70° Tapered Footrests 13"-17" ..... No Charge <sup>32</sup>
- AT5544 Elevating Legrests (E0990)..... 295.00

#### PIN STYLE

- PAS4A Smartleg® Art. Legrests 17.25"- 21.5" (K0053)..... 320.00
- PW93 Swingaway Footrests 13.5" - 19" ..... No Charge
- P904A Swingaway Elevating Legrests 16" - 19.75" (E0990) ... 303.00
- PAL4A Smartleg Art. Legrests 17.25"- 21.5" (K0053)..... 320.00 <sup>33</sup>

### POWER LEGS <sup>34</sup>

- ELRPW Power ELR w/Articulation 20" Ext (E1010).....1,895.00

### POWER LEGS CONTROL MODULE <sup>35</sup>

- 4WSB Four Way Switch Box ..... No Charge
- FWT 4 Way Toggle ..... No Charge
- QPB 4 Push Buttons ..... No Charge
- S4WSB Multiple Actuator Interface Box (E2311) ..... 2,600.00
- FWT 4 Way Toggle ..... No Charge
- QPB 4 Push Buttons ..... No Charge



FWT



QPB

**For systems with two or more actuators, NOT operating through driver control, choose 4WSB with either 4-way toggle or 4-push button.**

**For systems with 2 or more actuators OPERATING THROUGH DRIVER CONTROL, choose S4WSB. Uses 4-way toggle or 4-push button as attendant control.**

## FOOTPLATE OPTION

- 1651 Flip-Up Composite Footplates ..... No Charge
- 1350 Extra Large Aluminum Footplates ..... No Charge
- 1450 Aluminum Footplates-Deluxe ..... No Charge
- AT5543 Adjustable Angle Flip Up Footplates (K0040) .... 165.00
- AT5043 Flip Up Footplates ..... No Charge
- S1263 One Piece Wooden Footboard..... No Charge

## FOOTREST/LEGREST ACCESSORIES

- LPT2 Longer Pivot Slide Tube For 1350 ..... 83.00
- 1337 Calf Strap (K0038) ..... 27.00
- IMPACTG Impact Guards..... 51.00
- ALPT Longer Pivot and Slide Tube ..... 84.00
- A9342 "H" Heel Strap (K0039) ..... 51.00

## FOOTNOTES

<sup>32</sup> Standard with Impact guards.

<sup>33</sup> For seat widths greater than or equal to 16".

<sup>34</sup> Must pick 4WSB or S4WSB to operate power legs. Independent powered elevating articulating legrests with adjustable angle footplates. Extension range is 13"-20". Sold as pair.

<sup>35</sup> These options are located in the Configurator under the MK6i Electronics Folder. The Category Header inside the folder is labeled Power Seating Accessories.

# Invacare® 3G Arrow® Power Wheelchair w/TrueTrack Accessories and Options

Arrow Specifications	
Maximum Speed	7 MPH
Seat Widths	12" - 24" 1" increments
Seat Depths	12" - 19" 1" increments
	20" - 22" 1" increments
Overall Width	25"
Overall Length	29.5"
	32.5"
Seat-to-Floor	Standard Base
	17.5" Standard
	19.75" Optional
	21" Optional
Weight Limits	300 lb
	400 lb
Measurement based on wheel base. Seat width is width picked plus 7.5" for armrest and joystick	
Measured from headtube to anti-tipper	
Based on 5° Tilt & 16" Seat Depth. Add .25" for additional 2" of Seat Depth. Field adjustable for 10° of total pre-tilt.	
<b>Note: All specifications and dominions are approximate</b>	

Contoura Seat Sizes							
	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X
W24	*	*	X	X	X	X	X
*Not available with Contoura Seat Cushion							
Back Heights	25" or 28"						

Memory Card - BASIC
<ul style="list-style-type: none"> <li>Standard with ALL Expandable Control Systems (MPJ+, PSR+, PSF+, MK6i Display).</li> <li>Used to Back Up/Restore Programmed Settings for Only One Chair.</li> <li>May be given to user or kept with provider's files for safe keeping.</li> <li>Does not Contain Advanced Diagnostics/Help Screens/File Structure.</li> <li>Not compatible with SPJ+ Joysticks.</li> </ul>

Memory Card - PROFESSIONAL
<ul style="list-style-type: none"> <li>Standard with ALL MK6i Programmers.</li> <li>Available also with USB Card Reader on Order Form.</li> <li>Contains Advanced Diagnostics, Help Screens, Software Upgrades, File Storage/Retrieval.</li> <li>Used to back up multiple chairs and programming settings.</li> <li>Intended for use by qualified providers only.</li> <li>Not compatible with SPJ+ Joysticks.</li> </ul>

**FOR YOUR SEATING NEEDS  
PLEASE REFER TO THE FOLLOWING  
REHAB SEATING & POSITIONING PRODUCTS  
See Price List Form 07-027**

Chair Features and Options
<ul style="list-style-type: none"> <li>Independent Suspension for RWD for Clinical and Functional Comfort</li> <li>Full Electronic Communication &amp; Switch Option Capability</li> <li>Programmable Performance to Meet Individuals Needs</li> <li>Four Drives 1, 2, 3 &amp; 4 – Four Chairs in One</li> <li>5" Ground Clearance (RWD)</li> <li>Solid Seat Pan</li> <li>Convenient Rear Battery Access</li> <li>Automatic Self Connecting Battery Hookup</li> <li>Adjustable Depth Rigging Attachment 2" {1" increments}</li> <li>Automatic Dual Mode Battery Charge (Lester Brand)</li> <li>Multiple Colors to Choose From</li> <li>Option of Adjustable Tilt and Back Angle Seat</li> </ul>

Manual Recline Specifications			
Seat Widths	14" – 24" 1" increments		
Seat Depths	16" – 24" 1" increments {Adjustable 16" – 19", Fixed 20" – 22"}		
Back Heights	24" Standard, 18.5", 20", 22", 24", 26"		
Weight Limits	300 lb or 400 lb		
Headrest	10" Height		
Recline Range	90° to 170° Approximate		
Overall Back Height (w/o extension)	<u>Low Seat</u>	<u>Med Seat</u>	<u>High Seat</u>
Standard Base	41.5"	43.5"	45.25"
Heavy Duty Base		42.75"	44.62"
(Measured with 5° Tilt and 10" for Extension)			

Footrest Extension Range			
	AT5043	AT5543	1651
60HD	N/A	12.25" – 16.5"	13.5" – 17.5"
70MFX	6.75" – 12.5"	4.5" – 10"	N/A
70NHD	N/A	10.75" – 15.75"	11.75" – 17"
70STAPER	N/A	12.5" – 16.5"	13.25" – 17.25"
Measured From Seat Pan to Back of Footplate			

Weight Limit up to 400 Pounds
Invacare recommends all wheelchairs be fitted and programmed to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

**NOTE:**  
All weight limits are inclusive of add-on equipment and seating systems

Cantilever Arm Height Adjustment			
Back Type	Height	Adjustment	Adjustment with laterals
CURV	20" or greater	11" to 16"	11" to 14"
CURV	19"	11" to 15"	11" to 13"
CURV	18"	11" to 14"	11" to 12"
CURV	17"	11" to 13"	no adjustment
CURV	16"	11" to 12"	not available
Contoura	25" & 28"	12.5" to 14.5"	no adjustment

**Invacare Corporation** [www.invacare.com](http://www.invacare.com)



**USA**  
One Invacare Way  
Elyria, OH USA  
44035-4190  
800-333-6900

All rights reserved. Trademarks are identified by the symbols <sup>TM</sup> and <sup>®</sup>. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Therafin is a registered trademark of W/R Medical Electronics Co. Products designated by ASL are manufactured by Adaptive Switch Labs, Inc. Specifications and prices are subject to change without notification.  
©2008 Invacare Corporation