



Yes, you can.®

### Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email\* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

\*Saving and emailing the eForms is available with Acrobat Professional software

#### If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

**Please note:** The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

#### If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit [www.adobe.com](http://www.adobe.com).

# Invacare® Nutron® R51™ LXP Power Wheelchair

## Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



Yes, you can.®

### TRANSACTION TYPE:

- Order  Quote

Date \_\_\_\_\_ Account # \_\_\_\_\_ P.O. # \_\_\_\_\_  
Company Name \_\_\_\_\_ Phone # \_\_\_\_\_  
Address \_\_\_\_\_ City \_\_\_\_\_ State \_\_\_\_\_ Zip Code \_\_\_\_\_ Country \_\_\_\_\_  
Contact Name \_\_\_\_\_ Comments \_\_\_\_\_

### Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

## FRAME SELECTION

### POWER CHAIR TYPE

- R51LXP Nutron R51LXP with MK5™ NX™ Electronics ..... \$5,195.00

### FRAME WIDTH

- W16 16" Seat Width ..... Standard  
 W18 18" Seat Width ..... No Charge  
 W20 20" Seat Width ..... \$329.00

### SEAT DEPTH

- D16 16" Seat Depth ..... Standard  
 D17 17" Seat Depth ..... No Charge  
 D18 18" Seat Depth ..... No Charge

### SEAT-TO-FLOOR HEIGHT

- STF1734 17 3/4" Seat-to-Floor ..... No Charge  
 STF1934 19 3/4" Seat-to-Floor ..... Standard

## FRAME FINISH SELECTIONS

### BASIC COLORS

Clear

- 24P Wet Black ..... No Charge  
 71P Silver Vein ..... No Charge  
 125P Emerald Green ..... No Charge

### CANDY COLORS

- 61P Electric Red ..... No Charge  
 62P Electric Blue ..... No Charge

### UPHOLSTERY

- U240 Black Nylon Upholstery ..... No Charge

## BACK AND SEAT SELECTIONS

### BACK TYPE

- AHAB Adjustable Height/Angle Back ..... No Charge  
 FAHAB Fold Down Adj. Ht. Angle Back ..... \$85.00

### BACK HEIGHT

- BH16 16" Back Height ..... No Charge  
 BH17 17" Back Height ..... Standard  
 BH18 18" Back Height ..... No Charge  
 BH19 19" Back Height ..... No Charge

### SEAT POSITIONING STRAPS

Clear

- 1515 Push Button Style Seat Pos. Strap ..... Standard  
 1311BK Airline Buckle Seat Pos. Strap ..... \$50.00

### SEAT ACCESSORIES

- 9120 Solid Seat Insert ..... \$98.00  
 L8651 Omit Seat Upholstery ..... Less \$10.00

## ARM SELECTIONS

### ARM TYPE (PAIR)

Arms are flip-up and detachable.

- 28F Flip Up Desk Length Fixed Ht Conv ..... No Charge  
 29F Flip Up Desk Length Adj. Ht Conv ..... \$230.00 E0973  
 18F Flip Up Full Length Fixed Ht Conv ..... No Charge  
 19F Flip Up Full Length Adj. Ht Conv ..... \$230.00 E0973

### ARM ACCESSORIES

- L8401 Omit Arms ..... Less \$20.00  
 L8405 Omit Arm Pads ..... Less \$10.00

**NOTE: All specifications and dimensions are approximate.**  
*For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).*

# Invacare® Nutron® R51™ LXP Power Wheelchair

## CASTER AND REAR WHEEL OPTIONS

### CASTERS

- 1251BK 8" x 2" Pneumatic ..... No Charge
- 1228 8" x 1 3/4" Semi Pneumatic Lt Grey ..... Standard
- 1259 6" x 2 " Semi Pneumatic ..... No Charge

### REAR WHEEL TYPE

- 1476 12 1/2" x 2 1/4" Pneumatic Tires ..... Standard
- 1477 12 1/2" x 2 1/4" Flat Free Tires ..... \$150.00

## BATTERY OPTIONS

### BATTERY CHARGER OPTIONS

- 110CHARGER 110 Volt Battery Charger ..... No Charge
- LESTR Lester 8-Amp Battery Charger ..... \$320.00
- L8900 Omit Charger ..... Less \$80.00

### BATTERIES

- Clear
- 22NFMKG-2 MK Battery Gel (2 each) ..... \$300.00

## ELECTRONICS

### JOYSTICK MOUNTING HARDWARE

- QL Quad Link Retractable Joystick Mount ..... \$230.00 E1028
- LEFT Left Hand Mounted Joystick ..... No Charge

### ELECTRONIC ACCESSORIES

- Clear
- 1560 "T" Handle Flexible Joystick Ext ..... \$59.00
- 1561 Straight Handle Flexible Joystick Ext ..... \$39.00
- 1813M6 MK6i Programmer ..... \$350.00

## MISCELLANEOUS ACCESSORIES

### MISCELLANEOUS ACCESSORIES

- 1304 Wheel Locks ..... \$45.00
- 1209 Tire Pump ..... \$35.00
- NO2 Oxygen Holder ..... \$150.00

## FRONT RIGGING OPTIONS

### FRONT RIGGINGS

- PW93 Swingaway Footrests 13 1/2" - 19" ..... \$165.00
- PW93ST Swingaway Tapered Footrests 14 3/4" - 19 3/4" ..... \$165.00
- P904A Swingaway Elev. Legrests 16" - 19 3/4" ..... \$255.00 E0990
- PAL4A Smartleg Art. Legrests 17 1/4" - 21 1/2" ..... \$320.00 K0053

### FOOTPLATE OPTIONS

- 1650 Aluminum Footplates ..... No Charge
- 1350 Extra Large Aluminum Footplates ..... \$115.00  
**Not available on 16" wide frame.**
- AT5543 Adjustable Angle Flip-Up Footplates ..... \$165.00
- 1450 Aluminum Footplates Deluxe ..... No Charge

### FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap ..... \$36.00  
**2" wide hook and loop closure.**
- HLOOP Heel Loop ..... \$45.00 E0951
- 9342 "H" Strap ..... \$56.00
- IMPACTG Impact Guards ..... \$51.00

#### Seating and Positioning

For your seating needs please refer to the following:  
Rehab Seating and Positioning Products  
see Price List Form 07-027

**NOTE: All specifications and dimensions are approximate.**  
**For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

# Invacare® Nutron® R51™ LXP Power Wheelchair

R51 LXP Specifications	
<b>Length*</b>	40 3/4"
<b>Width</b>	20" - 26"
<b>Incline Capability</b>	9°
<b>Ground Clearance</b>	3 1/2"
<b>Weight Capacity</b>	300 lbs.
<b>Motor</b>	650 watt
<b>Maximum Speed</b>	6 mph
<b>Turning Radius**</b>	
Front (footrest in place)	33"
Rear (anti tippers in place)	12"
<b>Heights</b>	
17 3/4" seat-to-floor	35" - 38"
19 3/4" seat-to-floor	37" - 40"
Back angle	-5° - 20°
<b>Weight</b>	
With 22NF batteries	170 lbs.
Without batteries	105 lbs.
<b>Battery Options</b>	22NF
<b>Wheel and caster size</b>	
Front casters (standard)	8" x 2"
Front casters (optional)	6" x 2"
Drive wheels	12 1/2" x 2 1/4"
<b>Warranty</b>	
Base frame	60 months
Seat frame	36 months
Electronics	12 months
Motor/gearbox	12 months
* Measured with PW93 front riggings and anti tippers with a 16" seat depth. ** Measured with 16" W x 16" D seat.	

Comments

**NOTE: All specifications and dimensions are approximate.**  
*For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).*



**Invacare Corporation**  
**USA**  
 One Invacare Way  
 Elyria, OH USA  
 44035-4190  
 800-333-6900

[www.invacare.com](http://www.invacare.com)

©2013 Invacare Corporation  
 All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.