

# Top End® XLT PRO Handcycle

## Price List and Order Form

Effective January 1, 2006 Revised July 1, 2014

For accuracy of ordering, please fax order to Top End at 727-522-1007 or e-mail to TopEndOrders@invacare.com

For questions call 800-532-8677 or 727-522-8677.

To order Literature: 800-828-6282

www.topendwheelchair.com



Provider Name \_\_\_\_\_ Date \_\_\_\_\_ Top End Job Number \_\_\_\_\_  
 Account # \_\_\_\_\_ P.O. # \_\_\_\_\_ Phone # \_\_\_\_\_ Fax # \_\_\_\_\_  
 Address \_\_\_\_\_ City \_\_\_\_\_ State \_\_\_\_\_ Zip Code \_\_\_\_\_  
**Provider Primary Contact:** \_\_\_\_\_ **Email:** \_\_\_\_\_ **Phone No** \_\_\_\_\_  
 Athlete Name \_\_\_\_\_ Email \_\_\_\_\_ Phone No \_\_\_\_\_  
 Address \_\_\_\_\_ City \_\_\_\_\_ State \_\_\_\_\_ Zip Code \_\_\_\_\_  
**Ship order to:**  Provider  Customer **Order will ship to Provider unless otherwise noted.**

### Chair Type:

- XLTPRO** Shimano Rapid fire 27 speed external gearing ..... 3895.00  
 11-34 cassette operated by hands-on shifter w/caliper brake mounted on right handpedal

### STANDARD SEAT WIDTH AND FRAME COLOR: (Select One)

- 15" Tangerine Red Frame
- 17" Safety Yellow Frame

### OTHER THAN STANDARD FRAME COLOR:

- P3.....\$100.00

### CONTEMPORARY COLORS

- Wet Black  White Snow  Silver Metallic
- Stealth Black

### WAY COOL COLORS

- Royal Blue  Midnight Blue  Grasshopper Green


### TOO HOT COLORS

- Raspberry Red  Grape Madness
- Orange  Bubblegum Pink

### SAFETY AND COMFORT OPTIONS:

- |                          |                                                                                                |                 |
|--------------------------|------------------------------------------------------------------------------------------------|-----------------|
|                          | Seat and Back Cushions.....                                                                    | <b>Standard</b> |
|                          | Chain Guard.....                                                                               | <b>Standard</b> |
|                          | Safety Flag.....                                                                               | <b>Standard</b> |
|                          | Chest restraint: hook and loop.....                                                            | <b>Standard</b> |
| <input type="checkbox"/> | <b>XLT84</b> Crank width adaptor set (adds 1.5" per set).....                                  | 20.00           |
| <input type="checkbox"/> | <b>XLT12</b> Mirror, leg rest mount.....                                                       | 30.00           |
| <input type="checkbox"/> | <b>XLT63</b> Mirror, handle mount.....                                                         | 35.00           |
| <input type="checkbox"/> | <b>XLT13</b> Water bottle & cage.....                                                          | 35.00           |
| <input type="checkbox"/> | <b>XLT14</b> Cordless computer.....                                                            | 125.00          |
| <input type="checkbox"/> | <b>XLT15</b> Helmet: choose size.....                                                          | 75.00           |
|                          | <input type="checkbox"/> Small/ Medium <input type="checkbox"/> Medium/Large                   |                 |
| <input type="checkbox"/> | <b>XLT6</b> Front safety light.....                                                            | 25.00           |
| <input type="checkbox"/> | <b>XLT7</b> Rear flashing light.....                                                           | 25.00           |
| <input type="checkbox"/> | <b>XLT17</b> Crutch holder & strap.....                                                        | 125.00          |
| <input type="checkbox"/> | <b>XLT9</b> Leg guard attachment.....                                                          | 75.00           |
| <input type="checkbox"/> | <b>XLT19</b> Tow bar for wheelchair.....                                                       | 250.00          |
| <input type="checkbox"/> | <b>XLT20*</b> Bike rack for trailer hitch.....                                                 | 450.00          |
|                          | <input type="checkbox"/> 1 1/4" (std) <input type="checkbox"/> 2"                              |                 |
| <input type="checkbox"/> | <b>XLT28</b> Backpack hydration system.....                                                    | 100.00          |
| <input type="checkbox"/> | <b>TER126</b> Click strap large (pad length~12").....                                          | 86.00           |
| <input type="checkbox"/> | <b>TER125</b> Click strap med (pad length~10").....                                            | 80.00           |
| <input type="checkbox"/> | <b>TER124</b> Click strap small (pad length~8").....                                           | 75.00           |
| <input type="checkbox"/> | <b>XLT 33</b> Alignment Gauge.....                                                             | 75.00           |
| <input type="checkbox"/> | <b>RACE22*</b> Top End T-Shirt <input type="checkbox"/> L <input type="checkbox"/> XL.....     | 10.00           |
| <input type="checkbox"/> | <b>RACE13*</b> Top End Bike Jersey <input type="checkbox"/> L <input type="checkbox"/> XL..... | 70.00           |
| <input type="checkbox"/> | <b>Freight*</b> XLT 20.....                                                                    | \$ _____        |

\* Indicates an item where the standard discount does not apply or there are additional shipping charges or both. Please request a quote.



### Includes:

- Aluminum frame/chrome-moly fork
- Adjustable seat and back angle
- Narrow back style
- Adjustable crank height
- Chain guard, cushions and safety flag
- Adjustable footrest system
- Ovalized ergonomic hand pedals
- Parking brake
- 26" High Performance wheels with high pressure clincher tires
- Quick release axles
- 15° camber
- Road crown compensator
- Reflectors/Chest restraint

**Top End® XLT PRO™  
Handcycle Specifications**

**Frame** - 2.5" x .065 ovalized 6061 aluminum aircraft tubing.

**Frame Design** - Single main tube, recumbent position.

**Fork Assembly** - 4130 Chrome-moly for strength/stiffness. Raked angle for cornering stability.

**Wheelbase** - Overall length 72".

**Weight** - Overall weight 30 - 35 lbs.

**Turning Radius** - 18 feet for 360° circle, 12 feet w/one backing maneuver.

**Static Tip Angle** - 39° for 250 lbs. subject.

**Transfer Seat Height** - Seat height Front 13", Seat Height Rear 9".

**Back Angle** - Adjustable angle back.

**Back Style/Height** -Standard narrow back is ideal for competitive cyclists. Seat width minus 7". Back Height is fixed at 20".

**Upholstery** - Nylon w/hoop and loop adjustable tension seat sling with seat and back cushion. Hook and loop chest/seat restraint included.

**Crank** - Adjustable height.

**Parking Brake:** mounted on the right side of the handcycle. This brake is available if the primary brake fails and can be moved for easy access.

**Chainguard-** Black plastic.

**Steering** - Pivot steering w/raked fork for controlled turning and stability. Road crown compensator is standard.

**27 Speed Propulsion and Braking System** – 180mm V Crankset x 20" wide w/three chainrings, features 27 speeds, 48-39-24 chainrings and 11-34 cassette. The 11/34 cassette is operated by the Hands-on Rapid-Fire unit which is mounted on the right handpedal. It operates both the 11/34 cassette and the hands-on caliper brake. The manual derailleur operates the top chain rings utilizing a manual system which eliminates cable and their maintenance.

**Not recommended for customers with limited hand function because they typically cannot operate this option safely.**

**Handpedals** - Shape and angle is ergonomically designed for maximum power and are a lightweight ovalized aluminum shape with frictionless bearings.

**Footrest** - Multi-adjustable fore and aft w/straps to securely cradle calves and feet. Leg guard attachments are recommended for persons who have no leg movement. The leg guard will protect the legs from getting into the chain.

**Performance Wheels** - Rear hubs are precision heat treated, silver anodized w/quick release stainless axles:

26" High Performance Wheels include 26" x 1" High Performance Clincher Tire/Tube 100 psi

**Camber** – 15° for stability.

**Warranty** - 3 year on frame to original owner / limited parts warranty.

**Assembly** - Shipped completely assembled except for rear wheels. Shipping weight is 50 lbs.

Options:

Cordless Computer/Speedometer

Water Bottle and Cage

Helmet

Safety Lights, Mirror

Leg Guard Attachment

Click straps

Crutch Holder

Bike Rack

Cushion

Tow bar for Wheelchair

Backpack Hydration System

Alignment gauge



**Invacare Corporation**

**TOP END**

Sports and Recreation Products  
4501 63<sup>rd</sup> Circle North  
Pinellas Park, FL USA  
33781-5914  
800-532-8677

**www.invacare.com**

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

Shimano is a trademark of Shimano, Inc.

Specifications and prices are subject to change without notification.

©2014 Invacare Corporation