

Invacare® Pronto® M51™ Power Wheelchair



Yes, you can.®

Seating • POWER • CUSTOM  
REHAB  
Sports  
Pronto® M51™



*Anything else is a compromise.™*



Yes, you can.®

Versatility comes standard with the Invacare® Pronto® M51™ power wheelchair with SureStep®. The M51 boasts true center-wheel-drive performance for exceptional maneuverability and intuitive driving. New seat-to-floor heights allow for ideal seating and positioning selections. The new captain's seat (optional) offers a plush charcoal grey vinyl semi-recline seat and adjustable height arms, providing optimal consumer comfort. The M51-R also offers a standard rehab (Adjustable Seat and Back Angle – ASBA) seat. The M51 provides a smooth ride over transitions and thresholds up to two-inches while maintaining stability. The M51 design offers sleek styling and comfort while optimizing maneuverability and accessibility.



### Pronto M51 with Van Seat

- SureStep technology allows for smooth driving over thresholds and transitions
- Stylish design with choice of Candy Red or Deep Blue shroud colors
- Compact size for the tightest overall turning radius
- True center-wheel-drive system for intuitive driving
- Six wheels on the ground provide optimum stability
- Fully programmable MK5™ NX™ electronics for maximum function and precise driver control
- Adjustable joystick mount
- Light grey medium back van seat with flip-up, adjustable width arms
- Charcoal grey semi-recline seat (optional) with standard flip-up, adjustable width arms
- Rehab (ASBA) seat available on M51-R
- Powerful 350-watt, high-torque motors for reliable performance
- Regenerative braking and automatic motor locks for dependable stopping
- 3-amp on-board charger for convenience
- 300 lb. weight capacity

### Center-Wheel Drive



Center-wheel-drive offers the tightest turning radius, as well as improved tracking and ease of driving, over all other Invacare configurations.

### SureStep®



This patented automatic system ensures that everyday obstacles and thresholds are handled elegantly for a smooth ride.

### Optional Rehab Seating Available (ASBA) on M51-R



Customize the M51-R for consumers with moderate seating and positioning needs.

### Specifications

#### Overall dimensions

Length (with footboard folded)	34"
Length (with footboard down)	38" - 40"
Width (excluding joystick)	24"
Height	35" - 37"

#### Seat dimensions

Width	16", 18", 20"
Depth	16" - 18"

#### Weight

With batteries	198 lb.
Without batteries	151 lb.

#### Wheel and caster sizes

Front casters	6" x 2"
Drive wheels (flat-free)	10.5" x 3.5"
Rear casters	6" x 2"

#### Maximum speed

4.25 mph

#### Incline capability

9°

#### Weight capacity

300 lb.

#### Turning radius

Front (with footboard down)	19.5"
Rear	19.5"

#### Ground clearance

2.6"

#### Limited warranties

Base frame	60 months
Seat frame	60 months
Electronics	12 months
Motor/gearbox	12 months



"Invacare says 'Yes, you can.®' and I couldn't have said it better myself."