

# Invacare® IVC™ Tracer EX2 Wheelchair



Yes, you can.®

## Price List and Order Form

Effective July 1, 2008, Revised July 22, 2008

For ease of ordering, Contact Customer Service Toll Free at:

800-333-6900 Fax: 800-678-4682

To order Literature: 800-828-6282

www.invacare.com

Company Name \_\_\_\_\_ Account # \_\_\_\_\_

Phone \_\_\_\_\_ P.O. # \_\_\_\_\_ Date \_\_\_\_\_

Ship to Name & Address \_\_\_\_\_

### Special Note

**HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.**

### STOCK WHEELCHAIR CONFIGURATIONS

	Suggested List Price
<input type="checkbox"/> TREX26R (9153637773) <b>K0001/K0002</b> 16" x 16" Frame with Removable Fixed Height, Desk Length Arms .....	\$395.00
<input type="checkbox"/> TREX26RF (9153645183) <b>K0001/K0002</b> 16" x 16" Frame with Removable Fixed Height, Full Length Arms.....	395.00
<input type="checkbox"/> TREX28FF (9153637776) <b>K0001/K0002</b> 18" x 16" Frame with Permanent Arms & Fixed Footrests.....	410.00
<input type="checkbox"/> TREX28P (9153637782) <b>K0001/K0002</b> 18" x 16" Frame with Permanent Arms.....	395.00
<input type="checkbox"/> TREX28R (9153637774) <b>K0001/K0002</b> 18" x 16" Frame with Removable Fixed Height, Desk Length Arms.....	395.00
<input type="checkbox"/> TREX28RF (9153637784) <b>K0001/K0002</b> 18" x 16" Frame with Removable Fixed Height, Full Length Arms.....	395.00
<input type="checkbox"/> TREX20P (9153645182) <b>K0001/K0002/E2201</b> 20" x 16" Frame with Permanent Arms.....	395.00
<input type="checkbox"/> TREX20R (9153637775) <b>K0001/K0002/E2201</b> 20" x 16" Frame with Removable Fixed Height, Desk Length Arms.....	395.00
<input type="checkbox"/> TREX20RF (9153645181) <b>K0001/K0002/E2201</b> 20" x 16" Frame with Removable Fixed Height, Full Length Arms.....	395.00

### STOCK FRONT RIGGING OPTIONS

	Suggested List Price
<input type="checkbox"/> T93HE (9153637918) Economy Hemi Footrests; Composite Footplates without Heel Loops .....	\$65.00
<input type="checkbox"/> T93HC (9153637786) <b>E0951</b> Hemi Footrests; Composite Footplates with Heel Loops.....	93.00
<input type="checkbox"/> T93HA (9153637801) <b>E0951</b> Hemi Footrests; Aluminum Footplates with Heel Loops.....	93.00
<input type="checkbox"/> T94HE (9153637919) <b>E0990</b> Economy Hemi Elevating Legrests; Composite Footplates & Calf Pads.....	110.00
<input type="checkbox"/> T94HC (9153637788) <b>E0990</b> Hemi Elevating Legrests; Composite Footplates & Padded Calf Pads.....	180.00
<input type="checkbox"/> T94HA (9153637804) <b>E0990</b> Hemi Elevating Legrests; Aluminum Footplates & Padded Calf Pads.....	180.00
<input type="checkbox"/> AHL4A (9153619911) <b>K0053</b> Hemi Smartleg Articulating Legrests; Aluminum Footplates & Padded Calf Pads.....	320.00

# Invacare® IVC™ Tracer EX2 Wheelchair

## ACCESSORIES and OPTIONS

### STOCK OPTIONS & ACCESSORIES

**Suggested List Price**

#### ARMS

- 19ST (9153645340) **E0973**  
Full Length Adjustable Height Arms (10" – 14" Adjustable Range) ..... \$150.00
- 29ST (9153645339) **E0973**  
Desk Length Adjustable Height Arms (10" – 14" Adjustable Range)..... 150.00

#### ANTI-TIPPERS

- 9758ST (9153645343) **E0971**  
Rear Anti-Tippers for TREX28FF..... 82.00
- <sup>1</sup>  1360 (9153629710) **E0971**  
Adjustable Anti-Tippers (Does not fit fixed frame model)..... 82.00

#### MISCELLANEOUS

- 1329BKST (9153645341) **E0961**  
8" Wheel Lock Extensions..... 39.00
- <sup>2</sup>  1335ST (9153645342)  
Non-Folding Device..... 40.00
- 1514 (9153602056)  
Overhead Anti-Theft Device..... 46.00
- 1515ST (9153645344) **E0978**  
Push Button Style Seat Positioning Strap..... 50.00
- 1496 (9153632659) **E2208/K0105**  
IV/Oxygen Cylinder Holder Package Assembly..... 158.00
- 1497 (9153632671) **E2208**  
Oxygen Cylinder Holder Package Assembly..... 151.00
- <sup>1</sup>  1820ST (9153645345)  
Non-Removable Device for Front Riggings..... 20.00

#### Footnotes

1. Not available on Fixed Frame model.
2. Not available on 20" widths.

#### All Stock Tracer EX2 Wheelchairs Include: Standard

- \* Dual Axle Positions (17 ½" and 19 ½" Seat-to-Floor)
- \* 16" Back Height
- \* 8" x 1" Solid Rubber Casters
- \* 24" Composite Wheels with Urethane Tires
- \* Composite Handrims
- \* Midnight Blue Vinyl Upholstery (U550)
- \* Padded Armpads
- \* Wheel Locks

#### Specifications

Seat Widths	16", 18", or 20"
Overall Width	Permanent Arm – Seat Width plus 7 ¼" Removable Arm – Seat Width plus 8 ½"
Seat Depth	16"
Overall Depth	31 7/8" (18x16 without Front Riggings)
Back Height	16"
Seat-to-floor	17 ½" or 19 ½" (21" on Fixed Frame)
Weight Limit	250 lbs.
Chair Weight	36 lbs. (18x16 without Front Riggings)
Shipping Weight	45 ½ lbs.

For the most current pricing information go to [www.invacare.com](http://www.invacare.com)



**Invacare Corporation**  
USA  
One Invacare Way  
Elyria, OH USA  
44035-4190  
800-333-6900

[www.invacare.com](http://www.invacare.com)  
All rights reserved. Trademarks are identified by the symbols ™ and ®.  
All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.  
Specifications and prices are subject to change without notification.  
©2008 Invacare Corporation