



Yes, you can.®

Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.

Invacare® Pronto® M91™ Van Seat with SureStep™

Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



Yes, you can.®

TRANSACTION TYPE:

- Order Quote

Date _____ Account # _____ P.O. # _____

Company Name _____ Phone # _____

Address _____ City _____ State _____ Zip Code _____ Country _____

Contact Name _____ Comments _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

FRAME SELECTION

POWER CHAIR TYPE

- M91 Pronto M91 w/SureStep™ with MK₅™ NX75 Controller SPJ
+Joystick(push button Style seat positioning strap included) \$5,495.00 K0849
- M91SEAT Pronto M91 Seat Only..... \$740.00

SEAT SELECTION

- VAN Van Style Seat Standard

VAN SEAT (Standard Dimensions: 18"W X 16"D)

- VW16 16" Wide Van Seat No Charge
- VD16 16" Depth No Charge
- VD17 17" Depth No Charge
- VD18 18" Depth No Charge
- VW18 18" Wide Van Seat Standard
- VD16 16" Depth Standard
- VD17 17" Depth No Charge
- VD18 18" Depth No Charge
- VW20 20" Wide Van Seat \$329.00
- VD18 18" Depth No Charge
- VD19 19" Depth No Charge
- VD20 20" Depth \$529.00
- VW22 22" Wide Van Seat \$529.00
- VD18 18" Depth No Charge
- VD19 19" Depth No Charge
- VD20 20" Depth \$529.00

VAN UPHOLSTERY COLORS

- GV Gray Vinyl Upholstery No Charge

MOTOR PACKAGES

- U300 User Weight Up to 300 lbs. (STDPKG) Standard K0849
- U400 User Weight 301 - 400 lbs. (HDPKG) \$400.00 K0851

SEAT SELECTION

- SEMIR Semi-Recline Van Seat Standard
For Semi-Recline Van Seat with Solid Base see M91 Rehab Version price list.

ARM SELECTIONS

VAN ARM TYPE (PAIR)

Arms are flip-up and adjustable width.

- VAH Full Length, Adjustable Height \$230.00 E0973
- VFH Full Length Fixed Height Standard
- VAHD Desk Length, Adjustable Height \$230.00 E0973
- VFHD Desk Length, Fixed Height No Charge

POSITIONING STRAPS

SEAT POSITIONING STRAPS

Clear

- 1515 Push Button Style Seat Pos. Strap Standard
- 1311BK Airline Buckle Seat Pos. Strap No Charge

DRIVE WHEEL TYPE

DRIVE WHEEL TYPE

- DW143 14" x 3" Pneumatic Tires Standard
- DW143FF 14" x 3" Flat Free Tires \$250.00

SHROUD COLORS

- RD Candy Red Pearl No Charge
- BL Deep Blue Pearl No Charge
- MB Magic Black No Charge

ACCESSORIES

MISCELLANEOUS ACCESSORIES

- VANO2 Oxygen Holder (Van Seat Only) \$150.00 E2208
- HFO2-ML6 Home Fill O2 Holder-ML6 Only \$150.00 E2208
- 1209 Tire Pump \$35.00
- LESTR Lester 8 Amp Battery Charger \$320.00 E2367

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare Pronto M91 Van Seat with SureStep

ELECTRONICS

ELECTRONIC ACCESSORIES

1813M6 MK6i Programmer \$350.00

JOYSTICK MOUNTING HARDWARE

QL Quad Link Retractable Joystick Mount SPJ+ \$230.00 E1028

LEFT Left Hand Mounted Joystick No Charge

BATTERY OPTIONS

BATTERIES

22NFINSTAL 22NF Battery (On-Chair)(1 pair) \$300.00 E2361

FRONT RIGGING OPTIONS

FOOTBOARD

90FB Footboard Standard

L90FB Omit Footboard Less \$40.00

FOOTREST/LEGREST (HEMI)

Clear

PHW93 Swingaway Footrests 13 1/2" - 19" \$165.00

PH904A Swingaway ELR 16" - 19 3/4" \$260.00 E0990

PHAL4A Smartleg® Art. Legrests 17 1/4" - 21 1/2" \$225.00 K0053

FOOTPLATE OPTIONS

1650 Aluminum Footplates No Charge

1350 Extra Large Aluminum Footplates \$115.00

Not available on 14" or 16" wide frame.

AT5543 Adjustable Angle Flip-Up Footplates \$165.00

1450 Aluminum Footplates Deluxe No Charge

Not available on 14" wide frame.

FOOTREST/LEGREST ACCESSORIES

Clear

LPT2 Longer Pivot Slide Tube For 1350 \$83.00

ALPT Longer Pivot And Slide Tube \$84.00

ALPT4 Longer Pivot And Slide Tube 4" \$84.00

A1337 Calf Strap \$27.00 K0038

IMPACTG Impact Guards \$51.00

HLOOP Heel Loops \$45.00 E0951

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.**

Invacare Pronto M91 Van Seat with SureStep

Specifications	
Dimensions	
Lengths (with Footboard folded)	39"
Lengths (with Footboard down)	40" - 42"
Widths (excluding Joystick)	25.9"
Seat dimensions	
	Van Seat
	Semi-Recline
Widths	18" - 22"
Depths	37" - 40"
Overall Height	37" - 40"
Back Height	18"
with headrest	25"
Back Angle	80° - 115°
Seat to Floor Heights	
Range	
Van	21" - 23"
Weight	
With Batteries	273 lbs.
Without Batteries	199 lbs.

Specifications	
Wheel and Caster Sizes	
Front Casters	6" x 2"
Drive Wheels	14" x 3"
Rear Casters	6" x 2"
Maximum Speed	
Standard	6 mph
Heavy-Duty 400 lbs Weight Capacity	5 mph
Incline Capability	
	9 degrees
Weight Capacity	
	300 lbs. (400 lbs. option)
Turning Radius	
Front (with footboard down)	19 1/2"
Rear	21 1/2"
Ground Clearance	3 1/4"
Warranties	
Base frame	60 months
Seat Frame	36 months
Electronics	12 months
Motor/gearbox	12 months

Comments

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.



Invacare Corporation
USA
 One Invacare Way
 Elyria, OH USA
 44035-4190
 800-333-6900

www.invacare.com

©2013 Invacare Corporation
 All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.