



Yes, you can.®

Need a lightweight folding wheelchair with simple adjustability and reliable performance at an economical price? The Invacare® Patriot™ is the model of choice. This classic manual rehab wheelchair offers valuable features including: a lightweight, aluminum frame; adjustable seat-to-floor heights; durable center mount front riggings; and adjustable angle head tubes all in a compact folding design.



Features

- Lightweight, aluminum frame
- Adjustable seat-to-floor height and rear axle positions
- Adjustable angle head tubes
- Durable center mount front riggings
- Extra smooth precision head tube bearings
- Compact for transport (11.75" wide folded with wheels)
- Adjustable-height back available in 1" increments, 15"-19"
- Multiple arm, footplate and tire options
- 21 color options
- Quick-release axle option



Simple, triangular axle bracket offers 8 different axle positions – 4 axle positions in front of back canes and 4 axle positions behind back canes.



Precision head tube bearings allow casters to turn smoothly.

Specifications

Overall		
Width open	7.50" plus seat width*	
Width closed	11.75"**	
Depth (w/rigging)	43" or 45"	
Seat widths	14", 16", 17", 18", 20"	
Seat depths	16", 18"	
Knee-to-heel**	Min.	Max.
(range adjusts in 1" increments)		
70 degree tapered	14.75"	18.75"
70 degree	12.50"	17.50"
60 degree	13.75"	17.50"
Elevating	14.25"	20.25"
Arm Types		
Conventional dual-point		
Desk or Full Length/Fixed or Adjustable Height		
Space-Saver dual-point		
Desk or Full Length/Fixed or Adjustable Height		
Back heights	15"-20"	
Weight capacity	250 lb.	
Weight***	29 lb.	
Limited warranty	Lifetime on frame and crossbraces	

Front Seat-to-Floor

Rear Wheel Size	Seat-to-Floor	Caster	Front Caster		Rear Wheel	
			Mounting Position	Mounting Position	Mounting Position	Mounting Position
20"	15.50"	6"	top	middle	upper	
	16.50"	6"	†	refer to service manual		
	17.50"	6"	bottom	bottom	lower	
	17.50"	8"	top	bottom	lower	
	18.50"	8"	†	refer to service manual		
	18.50"	8"	†	refer to service manual		
22"	15.50"	6"	†	refer to service manual		
	16.50"	6"	†	refer to service manual		
	17.50"	6"	bottom	bottom	lower	
	17.50"	8"	top	bottom	lower	
	18.50"	8"	†	refer to service manual		
	19.50"	8"	†	refer to service manual		
24"	16.50"	6"	†	refer to service manual		
	17.50"	6"	bottom	middle	upper	
	17.50"	8"	top	middle	upper	
	18.50"	8"	†	refer to service manual		
	19.50"	8"	bottom	bottom	lower	
	19.50"	8"	bottom	bottom	lower	

* Measurements taken from outside of handrim to outside of handrim with 0 degrees camber

** Measurement taken from top of seat rail to back of footrest

*** 18" x 16" with 8" x 1 1/4" composite urethane casters, 24" composite urethane rear wheels. 70° swingaway footrest, flip-up composite footplates, desk length fixed height space-saver arm, adjustable back with 10° bend, aluminum anodized handrim, black sail cloth adjustable tension seat upholstery, auto style seat positioning strap and pull-to-lock wheellock.

† Not available as a factory setting

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

© 2009 Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and ©. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

Specifications are subject to change without notification. Accessories and seating products shown in the photographs are not necessarily included in the base price of wheelchair. Refer to current price list located at www.invacare.com.

Form No. 03-051
Rev. 04/09



Yes, you can.®

HomeCare
Bed • T.S.S. • Personal Care • Manual Wheelchair

Patriot™ Wheelchair

