



Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.

Invacare® P9000™ XDT

Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



Yes, you can.®

TRANSACTION TYPE:

Order Quote

Date _____ Account # _____ P.O. # _____

Company Name _____ Phone # _____

Address _____ City _____ State _____ Zip Code _____ Country _____

Contact Name _____ Comments _____

STOCK MODELS

P9000XDT1816 18"W x 16"D \$5,035.00 K0822

P9000XDT1818 18"W x 18"D \$5,035.00 K0822

OPTIONS AND ACCESSORIES

UIBH-2 UI Battery with handle 2 Each \$232.00

P900XDT does NOT come standard with batteries. Must select UIBH-2 if you would like batteries.

P9000XDT Features

- Choice of 18"W x 16"D or 18"W x 18"D Black Nylon Sling Seat
- MK₅™ Electronics with CSPJ Joystick
- Flip Up, Adjustable Height, Desk Length Arms
- Swingaway Footrests
- Drive Wheel - Flat Free Tires
- Casters - Semi-Pneumatic
- Wheel Locks
- Auto Style Positioning Strap
- 8 Amp Off-Board Charger
- Wet Black Frame Color

Seating and Positioning

For your seating needs please refer to the following:
Rehab Seating and Positioning Products
see Price List Form 07-027

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® P9000™ XDT

P9000XDT Specifications - specifications are approximate		
Dimensions	English	Metric
Length*	47"	(119.38 cm)
Width (excluding joystick)	24 1/2"	(62.23 cm)
Height (back at 17")	38"	(96.52 cm)
Seat-to-Floor	19 1/2"	(49.53 cm)
Seat Dimensions		
Width	18"	(45.72 cm)
Depths	16", 18"	(40.64 cm, 45.72 cm)
Back Height	17", 18", 19"	(43.18 cm, 45.72 cm, 48.26 cm)
Seat Angle	3°	3°
Weight Capacity	250 lbs.	(113.4 kg)
Wheel and Caster Sizes		
Front Casters	8" x 1 3/4"	(20.32 cm x 4.45 cm)
Drive Wheels	12 1/2" x 2 1/4"	(31.75 cm x 5.72 cm)
Maximum Speed	4.0 mph	(6.4 kph)
Motor	350 watts	(350 watts)
Incline Capability	9° grade	9° grade
Turning Radius		
Front (footrests in place)	36"	(86.36 cm)
Rear (anti-tippers in place)	16"	(33.02 cm)
Weight		
with UI Batteries	142 lbs.	(64.4 kg)
w/o Batteries & Front Riggings	82 lbs.	(37.2 kg)
Ground Clearance	3 1/4"	(8.26 cm)
Warranties	See Owner's Manual	
* Measured with front riggings and anti-tippers with 16" seat depth.		

Comments

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare Corporation

www.invacare.com

USA

One Invacare Way
 Elyria, OH USA
 44035-4190
 800-333-6900

©2010 Invacare Corporation
 All rights reserved. Trademarks are identified by the symbols ™ and ®.
 All trademarks are owned by or licensed to Invacare Corporation unless
 otherwise noted. Therafin is a registered trademark of W/R Medical
 Electronics Co. Specifications and prices are subject to change without
 notification.

