

Invacare® P9000™ XDT Power Wheelchair

REHAB
POWER
CUSTOM
SPORTS
SEATING

P9000 XDT™



Yes, you can.®



Anything else is a compromise.™

Invacare® P9000™ XDT Power Wheelchair

The Invacare® P9000™ XDT is a traditional rear-wheel drive chair that easily maneuvers and negotiates obstacles both inside and outside the home. It offers all the features consumers want at an affordable price. And it easily disassembles and folds for easy transport.



Yes, you can.®



Features

- Rear-wheel drive and compact size provide optimum maneuverability and accessibility
- Invacare® MK₅™ NX electronics featuring CSPJ joystick
- Folds to narrow 13" for transport
- Adjustable-height back provides comfort and support
- Powerful 350-watt, high-torque motors ensure dependable performance
- Telescoping front riggings allow optimal foot placement
- Wheel locks for added safety
- 250 lb. user weight capacity
- Available in two seat sizes: 18" W x 16" D, 18" W x 18" D

Specifications (All weight/ measurements are approximate.)

Dimensions	
Length	47"
Width (excluding joystick)	24.5"
Height	38"
Seat-to-floor height	19.5"
Seat Dimensions	
Width	18"
Depth	16", 18"
Back height	17", 18", 19"
Seat angle	3°
Weight capacity	250 lb.
Wheel & Caster Sizes	
Front/Rear casters	8" x 1.75"
Drive wheels	12.5" x 2.25"
Maximum Speed	
	4.0 mph
Incline Capability	
	9°
Turning Radius	
Front (footrest in place)	36"
Rear (anti-tippers in place)	16"
Weight	
with UI batteries	142 lb.
without batteries and front riggings	82 lb.
Ground Clearance	
	3.25"
Limited Warranties	
Base frame	36 months
Seat frame	36 months
Electronics	12 months
Motor/Gearbox	12 months

Measured with front riggings and anti-tippers with 16" seat depth.

Invacare® P9000™ XDT Power Wheelchair



Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Invacare, the Blue Medallion Design and Yes, you can, are registered trademarks and P9000 and MK₅ are trademarks of Invacare Corporation.

© 2006
Invacare Corporation
Form No. 03-067
Rev. 10/06



"Invacare says 'Yes, you can.®' and I couldn't have said it better myself."

Arnold Palmer

