

INVACARE®
P9000® XDT
POWER WHEELCHAIR



The **Invacare® P9000™ XDT Power Wheelchair** is a traditional rear-wheel drive chair that maneuvers easily and negotiates obstacles both inside and outside the home. It offers all the features consumers want at an affordable price and it easily disassembles and folds for transport.

Features

- Rear-wheel drive and compact size provide optimum maneuverability and accessibility
- Invacare® MK5™ NX™ Electronics featuring the Invacare® SPJ™ + Joystick
- Folds to narrow 13" for transport
- Adjustable-height back provides comfort and support
- Powerful 350-watt, high-torque motors ensure dependable performance
- Telescoping front riggings allow optimal foot placement
- Wheel locks for added safety
- Available in two seat sizes: 18" W x 16" D, 18" W x 18" D
- 8" casters

HCPCS Code: K0822



1. Disassembles for transport
2. Folds for transport



Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2012
Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and
®. All trademarks are
owned by or licensed
to Invacare Corporation
unless otherwise noted.
Specifications are subject
to change without
notification.
Form No. 03-067
Rev. 09/12
120352

Model #	Overall Dimensions	Seat Dimensions	Seat-to-Floor Height	Drive Wheels	Ground Clearance	Turning Radius	Incline
P9000XDT1816, P9000XDT1818	47.6" L* 25.4" W 35.2" H	18" W 16", 18" D	17.75"-19.75"	12.5"	3.25"	35.6"	9°
	Total Weight	Weight Capacity	Motor	Maximum Speed	Battery	Battery Range	Limited Warranty
	142 lb. w/ U1 batteries	250 lb.	2-Pole	4 mph	U1, 22NF	Up To 12-14 Miles w/ U1 batteries	3 Years Base Frame 1 Year Seat Frame 1 Year Electronics/ Motor/Gearbox

* Including front riggings

Do not operate on roads, streets or highways.

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

