

Stationary Oxygen • Aerosol  
RESPIRATORY  
Ambulatory Oxygen

# Platinum™ XL



Yes, you can.®





Yes, you can.®



The Invacare® Platinum™ XL Concentrator features an upgraded compressor that extends the time before preventative maintenance is required by 40% vs. the Platinum 5\*. Platinum XL offers a limited warranty – 5 years with no hour limit on the compressor. The unit features quiet operation, a top handle and bottom recesses for easy handling and a side filter access door. It is also compatible with the HomeFill® Oxygen System, and is available with SensO<sub>2</sub>™ (oxygen purity sensor).

**Invacare Corporation**  
www.invacare.com

**USA**

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

**Canada**

570 Matheson Blvd. E.  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

©2008

Invacare Corporation.  
All rights reserved.  
Trademarks are identified by the symbols ™ and ®.  
All Trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.  
Specifications are subject to change without notification.

Form No. 03-085  
Rev. 02/08

## Invacare® Platinum™ XL Oxygen Concentrator

Model no. IRC5LX 5-liter concentrator  
Model no. IRC5LX02 5-liter concentrator w/ SensO<sub>2</sub>®

### Features

- Upgraded compressor – time before rebuild extended by 40%
- Limited warranty – Full five years with no hour limit on the compressor



HomeFill® connection.



Filter access and adapter storage.

### Key Features

Top handle	For ease of mobility
Quick access door	Provides easy access to the inlet filter, humidifier adapter and diagnostics display
Servicing parts compatibility	Compatible with Platinum 5 and 10-Liter concentrators
Alarms and safety systems	Battery-free power-loss alarm; Sieve-GARD™ monitors sieve performance; SensO <sub>2</sub> ® oxygen monitor; Compressor 35 psi pressure-relief valve, low flow alarm
Pressure-based system	No timing adjustments needed
Self-diagnostic electronics	Minimizes troubleshooting time
Simplified P.M. schedule	Replace P.M. components as needed
Limited warranty	5 years with no hour limit on the compressor
Humidifier adapter storage	Located inside access door
Quiet operation	Spring mounted compressor = reduced noise and low vibration
One-piece hood	More durable than multiple piece hoods
Superior life cycle costs	The most efficient concentrator on the market *
HomeFill® compatibility	Reduce service costs when connected to HomeFill

\* Based on data gathered from independent, third party repair centers (as of 8/26/04)

### Specifications

Electrical requirements	120 VAC +/- 10% (132 VAC/108 VAC) 60Hz
Rated current input	4.0 A
Sound level	50 dB average
Oxygen output concentration levels	95.6% to 87% at all flow rates
Maximum outlet pressure	5 +/- 0.5 psi.
Flow range	1/2 - 5 L/min.
Average power consumption	390W
Pressure relief mechanism operational	35 psi +/- 3.5 psi
Filters	Cabinet, outlet HEPA, compressor inlet
Dimensions	
Width	18 3/8"
Height	26 3/8"
Depth	14 3/8"
Weight	52 lb.
Shipping weight	58 lb.
Limited warranty	5 years (no hour limit on compressor)
Approval	ETL