

Stationary Oxygen • Aerosol
RESPIRATORY
Ambulatory Oxygen

Platinum™ XL



Yes, you can.®





Yes, you can.®



The Invacare® Platinum™ XL Concentrator features an upgraded compressor that extends the time before preventative maintenance is required by 40% vs. the Platinum 5*. Platinum XL offers a limited warranty – 5 years with no hour limit on the compressor. The unit features quiet operation, a top handle and bottom recesses for easy handling and a side filter access door. It is also compatible with the HomeFill® Oxygen System, and is available with SensO₂™ (oxygen purity sensor).

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2008

Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ®.
All Trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.
Specifications are subject to change without notification.

Form No. 03-085
Rev. 02/08

Invacare® Platinum™ XL Oxygen Concentrator

Model no. IRC5LX 5-liter concentrator
Model no. IRC5LX02 5-liter concentrator w/ SensO₂®

Features

- Upgraded compressor – time before rebuild extended by 40%
- Limited warranty – Full five years with no hour limit on the compressor



HomeFill® connection.



Filter access and adapter storage.

Key Features

Top handle	For ease of mobility
Quick access door	Provides easy access to the inlet filter, humidifier adapter and diagnostics display
Servicing parts compatibility	Compatible with Platinum 5 and 10-Liter concentrators
Alarms and safety systems	Battery-free power-loss alarm; Sieve-GARD™ monitors sieve performance; SensO ₂ ® oxygen monitor; Compressor 35 psi pressure-relief valve, low flow alarm
Pressure-based system	No timing adjustments needed
Self-diagnostic electronics	Minimizes troubleshooting time
Simplified P.M. schedule	Replace P.M. components as needed
Limited warranty	5 years with no hour limit on the compressor
Humidifier adapter storage	Located inside access door
Quiet operation	Spring mounted compressor = reduced noise and low vibration
One-piece hood	More durable than multiple piece hoods
Superior life cycle costs	The most efficient concentrator on the market *
HomeFill® compatibility	Reduce service costs when connected to HomeFill

* Based on data gathered from independent, third party repair centers (as of 8/26/04)

Specifications

Electrical requirements	120 VAC +/- 10% (132 VAC/108 VAC) 60Hz
Rated current input	4.0 A
Sound level	50 dB average
Oxygen output concentration levels	95.6% to 87% at all flow rates
Maximum outlet pressure	5 +/- 0.5 psi.
Flow range	1/2 - 5 L/min.
Average power consumption	390W
Pressure relief mechanism operational	35 psi +/- 3.5 psi
Filters	Cabinet, outlet HEPA, compressor inlet
Dimensions	
Width	18 3/8"
Height	26 3/8"
Depth	14 3/8"
Weight	52 lb.
Shipping weight	58 lb.
Limited warranty	5 years (no hour limit on compressor)
Approval	ETL