

Top End® TRANSFORMER™ All Sport

Price List and Order Form Effective January 1, 2006 Revised July 1, 2007

For accuracy of ordering, fax your order to Top End at 727-522-1007.

For questions call 1 800-532-8677 or 727-522-8677.

To order Literature: 800-828-6282

www.invacare.com



Provider Name _____ Date _____ Top End Job Number _____
 Account # _____ P.O. # _____ Phone # _____ Fax # _____
 Address _____ City _____ State _____ Zip Code _____
 Customer Name _____ Phone No _____
 Address _____ City _____ State _____ Zip Code _____
Ship order to: Provider Customer **Order will ship to Provider unless otherwise noted.**

Chair Type
 TRAF Transformer All Sport.....\$2,395.00

Seat width: 12" 13" 14" 15" 16" 17" 18" 19" 20"
Seat depth: 12" 13" 14" 15" 16" 17" 18" 19" 20"
Adjustable frame height: choose tall or short
 Tall: Front seat height 18-21" set at _____
 Rear seat height 15-21" set at _____
 Short: Front seat height 16-18" set at _____
 Rear seat height 15-18" set at _____

Back type / angle: (choose one)
 Folding back / Adjustable angle
 Fixed back / 90° Angle

Back height: (choose one)
 8-11" 10-14" 12-16" 14-18"

Wing/bumper: (choose one)
 Full wing (available with camber of 12°, 15°, 18° or 20°)
 Bumper
 TRA190 Both wing and bumper100.00

Back upholstery
 Fold over **Standard**
 TRA49 Adjustable tension (w/straps)100.00

Camber system
 Single camber (0 -20) with camber tube **Standard**
 0° 3° 6° 9° 12° 15° 18° 20°
 TRA64 Extra camber tube150.00
 0° 3° 6° 9° 12° 15° 18° 20°
 TRA66 Adjustable camber system (choose 2)250.00
 Primary angle
 0° 3° 6° 9° 12°
 Secondary angle
 0° 3° 6° 9° 12°
 TRA46 Additional camber inserts (choose 2)100.00
 0° 3° 6° 9° 12°
 TRA48 Quick release lever30.00

Frame Color
UPHOLSTERY Black std UPHI Camouflage \$200.00
CONTEMPORARY COLORS
 Wet Black Black Prism White Snow
 Silver Metallic Silver Vein Copper Penny
WAY COOL COLORS
 Electric Blue Midnight Blue Cosmic Blue
 Electric Teal Emerald Green Grasshopper Green
TOO HOT COLORS
 Tangerine Red Electric Red Grape Madness
 Safety Yellow Sunkissed Orange Fruit Fusion

FOOTNOTES
 1. Rear seat height is dependent on rear wheel size. 15" RSH requires 22" or 24" rear wheels.
 2. If no color is chosen, black will be shipped.
 3. Crossed Spokes not available with Spineray wheels.
 4. Titanium Pushrims not available w/700C or 22" wheels.
 5. Needed for lower leg lengths of less than 15" and 12" (tall and short frames respectively).
 6. Kenda tires are available in 24", 25" and 26" sizes.
 7. Scissor style low-mount not available with camber over 6 degrees. Other restrictions may apply depending on configuration.
 8. Folding back only. Folding sideguard rigid w/slot not avail. with hi-mount wheellocks.

High performance wheels:.....Standard
 24", 25", 26" and 700 C wheels are double wide rims, high flange hub equipped with high pressure clincher tires and quick release axles.

Rear wheel size (choose one)
 22" (25-501) 24" (25-540)
 25" (25-559) 26" (25-590) 700 C

TRA67 Spineray® Spox wheels upgrade.....695.00
 22" (25-520) 24" (25-540) 25" (25-559)
 26" (25-590) 700C
 2 Spoke color Black Yellow

TRA68 6 Kenda® tires Red Yellow75.00
 TRA6 3 Crossed spokes upgrade125.00
 TRA21 2 Spoke guards.....75.00
 Black Clear

TRA27 Quad release axles.....125.00
 TRA76 Titanium axles.....150.00

Pushrim type:
 Aluminum **Standard**
 TRA8 Coated pushrims100.00
 TRA75 4 Titanium pushrim400.00

Pushrim mount:
 Close mount Standard Mount

Front wheels:
 3" Solid urethane w/ matching mini-fork **Standard**
 TRA77 Titanium fork stems100.00
 TRA15 Lighted front wheels and anti-tip100.00

Footrest:
 Adjustable height/angle one piece footplate **Standard**
 TRA111 Adj. ht/angle one-piece footplate w/toe clips250.00
 TRA191 5 Adaptive footplate spacers95.00
 2" 4" 6"
 Calf strapNo charge

Rear anti-tip
 Single swivel **Standard**
 TRA193 Double swivel100.00

Options:
 TRA14 Cushion: soft.....75.00
 2" 3" 4"
 TRA33 Cushion: hard125.00
 2" 3" 4"
 TRA194 Solid seat50.00
 TRA18 Removable sideguards.....100.00
 TRA51 8 Folding sideguards: Rigid with slot175.00
 TRA121 8 Fold-in style sideguards.....125.00
 TRA35 Frame protector (rear)20.00
 TRA37 Frame protector pair (front)40.00
 TRA94 Alignment gauge.....75.00
 TRA16 Click straps (set of 2 pair: large and small).....100.00
 TRA17 Positioning straps - **Hook and Loop**40.00
 TRA19 Positioning straps - **Auto Style Buckle**50.00
 TRA2 8 Wheel locks - **Push to Lock Hi-Mount**.....110.00
 TRA3 7 Wheel locks - **Scissor Style Low-Mount**.....250.00
 Race22 Top End T-Shirt, Blue L XL10.00
 TRA120 Fencing arm150.00

Total _____

All prices subject to change w/o notice. All prices in U.S. Dollars. All orders F.O.B. Pinellas Park, Florida, USA. For the most current pricing information refer to www.invacare.com.

NOTE: All specifications and dimensions are approximate.

Top End® TRANSFORMER™ All Sport

TRANSFORMER PRODUCT PARAMETERS	
Adjustable 4 Wheel Multi-Use Sports Chair	
Weight Limit:	250 lbs.
Material:	Aluminum Frame/Chrome-moly Wing or bumper.
Seat Widths:	12" - 20"
Seat Depths:	12" - 20"
Seat Upholstery:	Adjustable Tension
Back Upholstery:	Fold Over
Rear Wheels:	22-700c High Performance
Tires:	High Performance Clincher
Front Wheels:	3" Solid
Adjustable Rear Seat Height:	Tall or Short
Adjustable Front Seat Height:	Tall or Short
Back:	Folding Back Standard
Fixed Back:	No Charge Option
Back Height:	Adjustable
Back Angle:	Fixed Back is 90° to Floor Folding Back is Adjustable in 5° Increments
Center of Gravity:	Adjustable
Camber:	Camber Bar 2000 0° - 20°
Single Swivel Anti-Tip:	Standard
Adjustable Height/ Angle Footplate:	Standard
Removable Full Wing or Bumper:	Standard
Calf Strap:	Standard
Paint:	All Basic Top End Colors
Frame Length is determined by Seat Depth plus 9" for Folding Back	
Frame Length is determined by Seat Depth plus 7" for Fixed Back	
Footrest Width is determined by Seat Width Minus 5"	

TRANSFORMER OPTIONS	
*Removable Sideguards	*Spokeguard
*Folding Sideguards	*Spinergy Wheels
*Double Swivel Anti-Tippers	*Titanium Axles
*Clamp-on Raised Footrest	*Titanium Handrims
*Footplate Spacers	*Titanium Fork Stems
*Cushions	*Coated Handrims
*Solid Seat	*Frame Protectors
*Adjustable Tension Back Upholstery	*Extra Camber Inserts

NOTE: All specifications and dimensions are approximate.

Invacare Corporation www.invacare.com



TOP END
Sports and Recreation Products
4501 63rd Circle North
Pinellas Park, FL USA
33781-5914

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Spinergy is a registered trademark of Spinergy Corp. Kenda is a registered trademark of Kenda Rubber Industrial Co. Ltd. Specifications and prices are subject to change without notification. ©2007 Invacare Corporation

Form No. 03-111 Rev 0707
P2003_111H