

**CareGuard™ by Invacare®**  
**K2000 Standard**  
**K3000 Lightweight**  
**K4000 High Strength Lightweight**  
**K7000 Heavy Duty**  
**Wheelchairs**



*Yes, you can.*

# Parts Catalog Usage Guide

The information contained in this document is subject to change without notice. This catalog contains service parts for Form No. 03-037. If you have any questions regarding this catalog, please contact one of our Customer Service Representatives.

**For Customer Service Support Call:**

**1-800-333-6900**

**For Faxed In Orders Call:**

**1-800-678-4682 or**

**1-440-366-0709**

**For Technical Service Call:**

**1-800-832-4707**

**Literature Can Be Found at:**

**www.invacare.com or**

**1-800-828-6272**

## SYMBOL DEFINITION

Throughout this catalog you will notice letters to draw your attention to specific footnotes which are located at the bottom of the page or below specified charts.

Examples:

A,B,C etc. = Special Comments

NOTE: = General information pertaining to that specific page

NLA = No Longer Available

## HARDWARE MEASUREMENTS

The measurements for the hardware will be listed as follows:

Screws, Bolts Diameter X Thread Count/Pitch X Length

Washers Inside Diameter X Outside Diameter X Thickness

Spacers Inside Diameter X Outside Diameter X Length

For your convenience, we have supplied an area for you to list your Invacare Representatives.

Name	Phone Numbers

©2019 Invacare Corporation, Elyria, Ohio. Invacare is a registered trademark and Yes, you can, CareGuard, K2000, K3000, K4000 and K7000 are trademarks of the Invacare Corporation.


# TABLE of CONTENTS

<b>Description</b>	<b>Page Number</b>
Accessories .....	15
Arms - Desk Length & Full Length .....	10
Casters - Front, Forks and Mounting Hardware .....	11
Chair Components - Models K2000, K3000, K4000, K7000 .....	3
Crossbraces with Hardware .....	6
Frames - Side with Hardware .....	5
Front Riggings - Footrests & Legrests .....	14
Upholstery - Back .....	8
Upholstery - Seat .....	9
Wheel Locks .....	13
Wheels - Rear Wheels .....	12

# Chair Components

Model - K2000, K3000, K4000, K7000

**Chair Shown  
is for Reference Only**



# Chair Components

Model - K2000, K3000, K4000, K7000

CareGuard™ Chair Components
Back Upholstery
Seat Upholstery
Frames with Hardware
Crossbraces with Hardware
Arms - Desk Length & Full Length
Forks & Caster with Hardware
Wheel Locks
Rear Wheels
Footrest & Legrest with Footplate
Accessories - Telescopic IV Rod / IV / O2 Holder / O2 Holder

# Side Frames with Hardware

## Model K2000

Side Frame w/ Hardware	w/ Removable Arm	w/ Permanent Arm
Right	1123171-NLA See Note(s) A	1123173-NLA See Note(s) A
Left	1123172-NLA See Note(s) A	1123174-NLA See Note(s) A

NOTE: A - Complete frame assembly includes handgrip, caster journal bearings, all tube plugs/caps and the lower frame insert tube with attaching hardware.

## Model K3000

Side Frame w/ Hardware	Part Number
Right	1123191-NLA See Note(s) A
Left	1123192-NLA See Note(s) A

NOTE: A - Complete side frame assembly includes handgrip, caster journal bearings, plastic h-blocks, rivets & washers and all tube plugs/caps.

## Model K4000

Side Frame w/ Hardware	Part Number
Right	1123204-NLA See Note(s) A
Left	1123205-NLA See Note(s) A

NOTE: A - Complete side frame assembly includes caster journal bearings, plastic h-blocks, rivets & washers, all tube plugs/caps and back cane assembly with handgrip, attaching screw and knob.

## Model K7000

Side Frame w/ Hardware	Part Number
Right	1123210-NLA See Note(s) A
Left	1123211-NLA See Note(s) A

NOTE: A - Complete side frame assembly includes handgrip, caster journal bearings, all tube plugs/caps and the lower frame insert tube with attaching hardware.

# Crossbraces with Hardware

## Crossbraces with Hardware - Model K2000

Crossbrace Assembly Includes:	16" Wide	18" Wide
Right & Left X-brace w/tube plugs (1) Pivot Bolt Screw & Washer (1) Guide Tubes (2) Plastic Seat Guides (4)	1123175-NLA  See Note(s) A	1123176-NLA  See Note(s) A
NOTE: A - Assemblies for both removable and permanent arm chairs include hardware.		

## Crossbraces with Hardware - Model K3000

Crossbrace Assembly Includes:	16" Wide	18" Wide	20" W
Right & Left X-brace w/tube plugs (1) Pivot Bolt Screw & Washers (1) Pivot Link Assemblies (4) Pivot Link attaching Screws & Washers (2) Saddle Block and Plastic Liner (2) Attaching Screws & Washers (2)	1123193-NLA	1123194-NLA	1123195-NLA

## Crossbraces with Hardware - Model K4000

Crossbrace Assembly Includes:	16" Wide	18" Wide
Right & Left X-brace w/tube plugs (1) Pivot Bolt Screw & Washers (1) Pivot Link Assemblies (4) Pivot Link attaching Screws & Washers (2) Saddle Block and Plastic Liner (2) Attaching Screws & Washers (2)	1123193-NLA	1123194-NLA

# Crossbraces with Hardware

## Crossbraces with Hardware - Model K7000

Crossbrace Assembly Includes:	22" Wide	24" Wide
Right & Left X-brace w/tube plugs (1) Pivot Bolt Screw & Washers (1) Guide Tubes (4)	1123212-NLA	1123213-NLA

## Seat Guide - Model K2000 Only

Part Number	Description
1123170-NLA	NLA - Guide, Seat



# Back Upholstery

## Vinyl - Model K2000

Chair Width	16" Height	Hardware Quantity
16" W	1123160-NLA See Note(s) A	Screws & Washers (8)
18" W	1123161-NLA See Note(s) A	Screws & Washers (8)

NOTE: A - Includes metal straps, screws and washers.

## Nylon - Model K3000

Chair Width	16" Height	Hardware Quantity
16" W	1123182-NLA See Note(s) A	Screws & Washers (8)
18" W	1123183-NLA See Note(s) A	Screws & Washers (8)
20" W	1123184-NLA See Note(s) A	Screws & Washers (8)

NOTE: A - Includes metal straps, screws and washers.

## Nylon, Wrap Around - Model K4000

Chair Width	16" Height	Hardware Quantity
16" W	1123196-NLA See Note(s) A	Screws & Washers (2)
18" W	1123197-NLA See Note(s) A	Screws & Washers (2)

NOTE: A - Includes screws and washers.

## Vinyl - Model K7000

Chair Width	16" Height	Hardware Quantity
22" W	1123208-NLA See Note(s) A	Screws & Washers (8)
24" W	1123209-NLA See Note(s) A	Screws & Washers (8)

NOTE: A - Includes metal straps, screws and washers.

# Seat Upholstery

## Vinyl - Model K2000

Chair Width	16" Depth	Hardware Quantity
16" W	1123158-NLA See Note(s) A	Screws & Washers (8)
18" W	1123159-NLA See Note(s) A	Screws & Washers (8)

NOTE: A - Includes metal straps, screws and washers.

## Nylon - Model K3000

Chair Width	16" Depth	Hardware Quantity
16" W	1123179-NLA See Note(s) A	Screws & Washers (8)
18" W	1123180-NLA See Note(s) A	Screws & Washers (8)
20" W	1123181-NLA See Note(s) A	Screws & Washers (8)

NOTE: A - Includes metal straps, screws and washers.

## Nylon - Model K4000

Chair Width	16" Depth	Hardware Quantity
16" W	1123179-NLA See Note(s) A	Screws & Washers (8)
18" W	1123180-NLA See Note(s) A	Screws & Washers (8)

NOTE: A - Includes metal straps, screws and washers.

## Vinyl - Model K7000

Chair Width	18" Depth	Hardware Quantity
22" W	1123206-NLA See Note(s) A	Screws & Washers (10)
24" W	1123207-NLA See Note(s) A	Screws & Washers (10)

NOTE: A - Includes metal straps, screws and washers.

# Arms - Desk Length & Full Length Assemblies

## Arm Assembly - Model K2000 & K7000

Arm Style	Arm Assembly - Right	Arm Assembly - Left	Arm Assembly - Pair
Desk Length	1123163-NLA See Note(s) A	1123164-NLA See Note(s) A	1123162-NLA See Note(s) B
Full Length	1123166-NLA See Note(s) A	1123167-NLA See Note(s) A	1123165-NLA See Note(s) B

NOTE: A - Armrest assembly includes hardware, armpad and guard.  
B - Includes Right and Left assemblies.

## Arm Assembly - Model K3000

Arm Style	Arm Assembly - Right	Arm Assembly - Left	Arm Assembly - Pair
Desk Length	1123186-NLA See Note(s) A	1123187-NLA See Note(s) A	1123185-NLA See Note(s) B
Full Length	1123189-NLA See Note(s) A	1123190-NLA See Note(s) A	1123188-NLA See Note(s) B

NOTE: A - Armrest assembly includes hardware, armpad and guard.  
B - Includes right and left assemblies.

## Arm Assembly - Model K4000

Arm Style	Arm Assembly - Right	Arm Assembly - Left	Arm Assembly - Pair
Desk Length	1123199-NLA See Note(s) A	1123200-NLA See Note(s) A	1123198-NLA See Note(s) B
Full Length	1123202-NLA See Note(s) C, A	1123203-NLA See Note(s) A	1123201-NLA See Note(s) B

NOTE: A - Armrest assembly includes armpad, guard and hardware (including rear plastic socket, pivot screw, and nut).  
B - Includes right and left armrest assemblies.  
C - No longer available, use part number 1123201 as replacement

## Armrest Pad w/Hardware - Model K2000, K3000, K4000, K7000

Arm Style	Armrest Pad w/Hardware	Quantity Required
Desk Length	1134674 See Note(s) A	2
Full Length	1134675-NLA See Note(s) A	2

NOTE: A - Armpad assembly includes upholstered armpad, base and hardware.

# Casters, Forks and Mounting Hardware

## Model K2000, K3000, K4000

Part Number	Fork and Stem Assembly includes:
1123177-NLA  See Note(s) A	Fork & Stem Assembly (1) 8" x 1" Caster Assembly (1) Fork Stem Locknut (1) Fork Stem Washer (1)
NOTE: A - 8x1" caster includes tire, bearings, hub spacer, axle bolt and locknut.	

## Model K7000

Part Number	Fork and Stem Assembly includes:
1123214-NLA  See Note(s) A	Fork & Stem Assembly (1) 8 x 2" Caster Assembly (1) Fork Stem Locknut (1) Fork Stem Washer (1)
NOTE: A - 8x2" caster includes tire, bearings, hub spacer, axle bolt and locknut.	

**Rear Wheels**  
**with Handrim and Hardware**  
**Models - K2000, K3000, K4000**

Part Number	Rear Wheel Assembly includes:
1123178-NLA	24" Rear wheel with tire, bearing, hub spacer, composite handrim and hardware which includes axle bolt, axle, spacer, locknut and hubcap.

**Models - K7000**

Part Number	Rear Wheel Assembly includes:
1123215-NLA	24" Rear wheel with tire, bearing, hub spacer, chrome handrim and hardware which includes axle bolt, axle, spacer, locknut and hubcap.

# Wheel Locks

## Wheel Locks - Conventional, Removable Arm Models

### Model - K2000, K3000, K7000

Wheel Lock - Pair	Right	Left
1123293	1101452	1101453
See Note(s) A, B	See Note(s) A	See Note(s) A
NOTE: A - Black painted finish. B - Includes right and left wheel lock assemblies.		

## Wheel Locks - Space Saver, Permanent Arm Models

### Model - K2000, K4000

Wheel Locks - Pair	Right	Left
1123294-NLA	1111924	1111925
See Note(s) B, A	See Note(s) A	See Note(s) A
NOTE: A - Black painted finish. B - Includes right and left wheel lock assemblies.		

# Footrest & Legrest with Footplate

## Footrest & Legrest with Composite Footplate

### Models K2000, K3000, K4000

Front Riggings	Right	Left	Pair
Footrest	1123147-NLA	1123148-NLA	1123146-NLA
Legrest	1123153-NLA See Note(s) A	1123154-NLA See Note(s) A	1123152-NLA See Note(s) A

NOTE: A - Includes padded calf pads.

## Footrest & Legrest with Aluminum Footplate

### Model K7000

Front Riggings	Right	Left	Pair
Footrest	1123150-NLA	1123151-NLA	1123149-NLA
Legrest	1123156-NLA See Note(s) A	1123157-NLA See Note(s) A	1123155-NLA See Note(s) A

NOTE: A - Includes padded calf pads.  
NOTE: Aluminum recommended for K7000 Heavy Duty Chairs.

## Accessories

Part Number	Description
1101760-NLA See Note(s) A	Telescopic IV Rod
1070899 See Note(s) C, B	IV / O2 Holder
1070901 See Note(s) C, B	O2 Holder
1487 See Note(s) E	IV / O2 Holder
1488 See Note(s) D	O2 Holder
1122138-NLA	Owners Manual
511110-NLA	Seat Positioning Straps
623070-NLA	Rear Anti-Tippers - K2000/K7000
623074-NLA	Rear Anti-Tippers - K3000/K4000
<p>NOTE: A - Telescopic IV Rod is used with all CareGuard models.            B - I.V./O2 Holder and O2 Holder models 1496 &amp; 1497 are used on K2000 and K7000 only.            C - I.V./O2 Holder or O2 Holder cannot be used when rear anti-tippers are used on K7000 model chairs.            D - Oxygen Holder is used on models K3000 and K4000, however, it cannot be used when rear anti-tippers are in place.            E - I.V. / Oxygen Holder is used on models K3000 and K4000, however, it cannot be used when rear anti-tippers are in place.</p>	



**Invacare Corporation**

[www.invacare.com](http://www.invacare.com)



**Yes, you can.**

One Invacare Way  
Elyria, Oh USA  
44036-2125  
800-333-6900

Technical Support  
800-832-4707  
Florida Technical Support  
407-321-5630

Customer Service  
800-333-6900

© 2019 Invacare Corporation. All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by <sup>TM</sup> and ©. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted. (Revised 01/30/19)

Form No. 03-037