

INVACARE®

At'm

Take Along Wheelchair



Yes, you can.®

Wherever you're going, the Invacare® At'm Take Along Wheelchair goes, too. It's the power chair lightweight and compact enough to fit in the trunk of a car; simple enough to assemble in 60 seconds; and powerful enough for a full day of adventure.

How easy is it to set up an At'm power chair? If you can set up a lawn chair, you're ready. No mobility lifts, no tools – just three easy pieces to assemble. Open the seat, snap it on the lightweight base, add the battery and go! No wonder we call the At'm power chair the Take Along Chair.



Invacare® At'm Take Along Chair

The At'm chair goes wherever you go. It folds together like a lawn chair, and sets up in seconds. How easy is it to drive? Flip the switch and go.

Available Colors



Candy Red Pearl



Deep Blue Pearl

Specifications

Dimensions

Overall width	25.5"
Overall length	35.5"
Seat style	Sling style
Seat depth	16"
Seat width	18"
Seat back height	21"

Tires (flat-free std.)

Front tire size	1" x 6"
Rear tire size	2.25" x 12.5"

Weight

Base section	34 lb.
Seat and arms	14 lb.
Battery pack	26 lb.
Total weight	74 lb.
Weight Limit	250 lb.

Mechanical/Performance

Max. speed	3.5 mph
Incline capability	6°
Ground clearance	2"
Charger	Off-board
Battery size	17 Ah
Battery range	10-12 miles *

Limited Warranty

Frame	36 months
Electronics	12 months

*Range varies with battery conditions surface, terrain and operator weight. The Invacare® At'm power chair is not intended for street or highway use. MINI is a registered trademark of Bayerische Motoren Werke Aktiengesellschaft Corporation and is not affiliated with Invacare Corporation. All weight/measurements are approximate.

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

© 2010

Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ®.
All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

Form No. 04-159
Rev 08/10
100502

