

The best just got bigger...introducing the new Platinum™ 10 concentrator. The Platinum 10 delivers the reliability, performance and value advantage you have come to expect from the Platinum line of concentrators and provides an attractive alternative to liquid systems. The Platinum 10 is available at up to 10 LPM for high flow patients, and has less than 10 parts different from the Platinum™ XL 5 Liter concentrators.



Invacare® Platinum™ 10 Concentrator

Model no. IRC10LX 10-Liter Concentrator
Model no. IRC10LXO2 10-Liter Concentrator w/
SensO₂ (shown above)

FEATURES:

- Delivers continuous oxygen up to 10 L/min
- Alternative to liquid systems
- Same platform as the Platinum XL you know and trust
- Common maintenance and spare parts with the Platinum 5 and XL
- 3 year full warranty
- Superior life cycle costs



Filter access and adapter storage.

ADDED FEATURES

Top handle	For ease of mobility
One-piece hood	More durable than multiple piece hoods
Lifetime 4-way valve	Maintenance free – never needs servicing
Quick access door	The only concentrator with access to BOTH the compressor inlet and outlet HEPA® filters
Humidifier adapter storage	Located in access door
Service parts compatibility	Compatible with existing Platinum 5 and Platinum XL concentrators
Alarms and safety systems	Battery-free power-loss alarm; Sieve-GARD™ monitors sieve performance; SenseO ₂ oxygen monitor; Compressor 35 psi pressure-relief valve; Low-flow alarm
Pressure-based system	Superior to timing-based systems
HomeFill™ II compatibility	Same connection as on XL
Simplified P.M. schedule	Replace P.M. components as needed

SPECIFICATIONS

Electrical Requirements	120 VAC +/- 10% (132 VAC/108 VAC) 60 Hz
Rated current input	4.9 A
Sound level	58 dB Average
Oxygen output concentration levels	94% to 87% at all flow rates
Maximum outlet pressure	5+/- 0.5 psi
Flow range	2-10 L/min
Average power consumption	585 W
Pressure relief mechanism operational	35 psi +/-3.5 psi
Filters	Cabinet, outlet HEPA, compressor inlet
Dimensions	
Width	18 3/8"
Height	26 3/8"
Depth	14 3/8"
Weight	54 lb.
Shipping weight	60 lb.
Limited warranty	3 years



Yes, you can.™

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Invacare and the blue medallion design are registered trademarks and Yes, you can., HomeFill, Platinum and Sieve-GARD are trademarks of Invacare Corporation. HEPA is a registered trademark of HEPA Corporation. Specifications are subject to change without notification.

© 2004

Invacare Corporation
Form No. 04-169



"Invacare says 'Yes, you can.' and I couldn't have said it better myself."



Invacare® Platinum™ 10 Concentrator



Yes, you can.™