



Yes, you can.®

### Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email\* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

\*Saving and emailing the eForms is available with Acrobat Professional software

#### If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

**Please note:** The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

#### If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit [www.adobe.com](http://www.adobe.com).

# Invacare® Pronto® M91™ with SureStep® (Rehab Version)

## Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



Yes, you can.®

### TRANSACTION TYPE:

- Order       Quote

Date \_\_\_\_\_ Account # \_\_\_\_\_ P.O. # \_\_\_\_\_

Company Name \_\_\_\_\_ Phone # \_\_\_\_\_

Address \_\_\_\_\_ City \_\_\_\_\_ State \_\_\_\_\_ Zip Code \_\_\_\_\_ Country \_\_\_\_\_

Contact Name \_\_\_\_\_ Comments \_\_\_\_\_

### Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

## FRAME SELECTION

### POWER CHAIR TYPE

- M91R M91 with SureStep with MK6i™ electronics .. \$5,600.00 K0848  
**Push button style seat positioning strap included.**
- M91-C M91 Complete with SureStep Base Only, with MK6i electronics ..... \$5,255.00
- M91-M M91 Motion Concepts with SureStep Base Only, with MK6i electronics ..... \$5,255.00
- M91SEAT Pronto M91 Seat Only ..... \$740.00

### TRANSPORT TIE DOWN

- TRRO Transport Ready Option ..... \$300.00
- TRBKTS Wheelchair Transport Brackets ..... \$250.00

#### TRRO Features for M91 with SureStep

- Pronto M91R with MK6i
- 300 lbs. or 400 lbs. User Weight Capacity
- Adjustable Rehab Seat with Adj Seat and Tilt Back Angle
- Any Seat Width
- Any Seat Depth
- Sectional Rear Uprights
- 18" - 20" Back Height
- Any Arm Type
- Black Nylon Upholstery
- 14" x 3" Flat Free Tires
- Any Shroud Color
- Footboard
- MPJ+™
- All PSR/PSF Joysticks
- Left Handed Mounted Joystick Option
- "T" Handle Flexible Joystick Ext.
- Straight Handle Flexible Joystick Extension
- MK6i Remote Programmer
- Reinforced Back Upholstery
- 22NF Batteries
- Tire Pump
- Lester 8 amp Battery Charger

#### Transport Tie Down

TRRO includes four factory-installed transport brackets and a wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNAWC Vol 1. Section 19.5.3 Frontal Impact Test requirements for wheelchairs with a 168 lbs. crash dummy, which corresponds to a person with a weight of 114 lbs. to 209 lbs.

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

\*Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

## SEAT SELECTION

- ADJASBA Seat with Adj. Seat Tilt and Back Angle ..... Standard
- CTB Contoura™ Seating ..... \$625.00  
**Selected Seat Depth on Contoura assumes 1" of back cushion compression. For example: When choosing 18" depth the measurement from the front of the back cushion to the front of the seat pan equals 17". The Contoura back cushion is 2" thick.**
- OMCU Omit Contoura Seat Cushion ..... Less \$30.00
- SEMIRB Van Seat - Semi Recline with Solid Base ..... \$70.00

## MOTOR PACKAGES

- U300 User Weight Up to 300lbs. (STDPKG) ..... Standard K0849  
**Standard motor package includes 4-pole motors with a weight limit of 300 lbs.**
- U400 User Weight 301 - 400 lbs. (HDPKG) ..... \$400.00 K0851

**NOTE: All specifications and dimensions are approximate.**  
**For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

# Invacare® Pronto® M91™ with SureStep® (Rehab Version)

## ADJUSTABLE ASBA SEAT SELECTIONS

### ADJUSTABLE SEAT WIDTH RANGE

- 1620W 16" to 20" Seat Width Range ..... No Charge
- 2024W 20" to 24" Seat Width Range ..... No Charge

### SEAT WIDTH

- W16 16" Seat Width ..... No Charge
- W17 17" Seat Width ..... No Charge
- W18 18" Seat Width ..... No Charge
- W19 19" Seat Width ..... No Charge
- W20 20" Seat Width ..... No Charge
- W21 21" Seat Width ..... \$329.00
- W22 22" Seat Width ..... \$329.00
- W23 23" Seat Width ..... \$329.00  
**Only available with U400 motor package.**
- W24 24" Seat Width ..... \$469.00  
**Only available with U400 motor package.**

### ADJUSTABLE SEAT DEPTH RANGE

- 1619D 16" to 19" Seat Depth Range ..... No Charge
- 1922D 19" to 22" Seat Depth Range ..... No Charge

### SEAT DEPTH

- D16 16" Seat Depth ..... Standard
- D17 17" Seat Depth ..... No Charge
- D18 18" Seat Depth ..... No Charge
- D19 19" Seat Depth ..... No Charge
- D20 20" Seat Depth ..... \$529.00
- D21 21" Seat Depth ..... \$529.00
- D22 22" Seat Depth ..... \$529.00

### UPHOLSTERY

- U240 Black Nylon Upholstery ..... No Charge

## BACK SELECTIONS

### SEAT BACK HEIGHT

- BH16 16" Back Height ..... No Charge
- BH17 17" Back Height ..... No Charge
- BH18 18" Back Height ..... No Charge
- BH19 19" Back Height ..... No Charge
- BH20 20" Back Height ..... No Charge

### CONTOURA

- TSBH25 25" Back Height ..... No Charge
- TSBH28 28" Back Height ..... No Charge

### SEAT BACK TYPE

- SECT Sectional Rear Uprights ..... No Charge

## ARM SELECTIONS

### ARM TYPE (PAIR)

- 28F Desk Length Fixed Ht Conv ..... Standard
- 29F Desk Length Adj. Ht Conv ..... \$230.00 E0973
- 18F Full Length Fixed Ht Conv ..... No Charge
- 19F Full Length Adj. Ht Conv ..... \$230.00 E0973

### MIXED ARMS (EACH)

- 19FR Right Arm, Full Length, Adj. Height ..... \$115.00 E0973
- 29FR Right Arm, Desk Length, Adj. Height ..... \$115.00 E0973
- 19FL Left Arm, Full Length, Adj. Height ..... \$115.00 E0973
- 29FL Left Arm, Desk Length, Adj. Height ..... \$115.00 E0973

## VAN SEAT - WITH SOLID BASE

### VAN SEAT WIDTH AND DEPTH

- VW16 16" Wide Van Seat ..... No Charge
  - VD16 16" Depth ..... No Charge
  - VD17 17" Depth ..... No Charge
  - VD18 18" Depth ..... No Charge
- VW18 18" Wide Van Seat ..... Standard
  - VD16 16" Depth ..... No Charge
  - VD17 17" Depth ..... No Charge
  - VD18 18" Depth ..... No Charge
- VW20 20" Wide Van Seat ..... \$329.00
  - VD18 18" Depth ..... No Charge
  - VD19 19" Depth ..... No Charge
  - VD20 20" Depth ..... \$529.00
- VW22 22" Wide Van Seat ..... \$529.00
  - VD18 18" Depth ..... No Charge
  - VD19 19" Depth ..... No Charge
  - VD20 20" Depth ..... \$529.00

### VAN UPHOLSTERY COLORS

- GV Gray Vinyl Upholstery ..... No Charge

### VAN ARM TYPE (PAIR)

- VAH Full Length, Adj. Height ..... \$230.00 E0973
- VFH Full Length, Fixed Height ..... Standard
- VAHD Desk Length, Adj. Height ..... \$230.00 E0973
- VFHD Desk Length, Fixed Height ..... No Charge

**NOTE: All specifications and dimensions are approximate.**  
**For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

# Invacare® Pronto® M91™ with SureStep® (Rehab Version)

## POSITIONING STRAPS AND WHEEL OPTIONS

### SEAT POSITIONING STRAPS

Clear

- 1515 Push Button Style Seat Pos. Strap ..... Standard  
 1311BK Airline Buckle Seat Pos. Strap ..... No Charge

### CHEST POSITIONING STRAPS

Clear

- 7311BK Airline Buckle Chest Pos. Strap ..... \$100.00 E0960  
 9515 Push Button Style Chest Pos. Strap ..... \$100.00 E0960  
 7321BK Hook & Loop Chest Pos. Strap ..... \$100.00 E0960

### SHROUD COLORS

- BL Deep Blue Pearl ..... No Charge  
 RD Candy Red Pearl ..... No Charge  
 MB Magic Black ..... No Charge

### DRIVE WHEEL TYPE

- DW143 14" x 3" Pneumatic Tires ..... Standard  
 DW143FF 14" x 3" Flat Free Tires ..... \$149.00

## ACCESSORIES

### BACK ACCESSORIES

- L8601 Omit Back Upholstery ..... No Charge  
 RBP Reinforced Back Upholstery ..... No Charge

### ARM ACCESSORIES

Clear

- L8401 Omit Arms ..... Less \$20.00  
 L8405 Omit Arm Pads ..... Less \$10.00

### MISCELLANEOUS ACCESSORIES

- SEAT02 Oxygen Holder (ASBA Seat Only) ..... \$150.00 E2208  
 HFO2-ML6 HomeFill™ O2 Holder-ML6 Only ..... \$150.00 E2208  
 1209 Tire Pump ..... \$35.00  
 LESTR Lester 8 Amp Battery Charger ..... \$320.00 E2367  
 VANO2 Oxygen Holder (Van Seat Only) ..... \$150.00 E2208  
 CTE Single Mode Charger ..... \$205.00 E2366

## BATTERY OPTIONS

### INSTALLATION OF BATTERIES

- 22NFINSTAL 22NF Battery(On-Chair)(1 pair) ..... \$300.00 E2361

## MK6i ELECTRONICS

### ELECTRONIC CONTROLLERS

- MK690 MK6i 90 Amp Controller ..... Standard

### STANDARD DRIVER CONTROL

- SPJM6 SPJ™+ Single Drive Joystick ..... Standard

### BASE ONLY SPJ+ OPTIONS

Base Only Systems for Multiple Actuator Powered Seating Systems will require MPJ+, PSRM6, PSFM6, or MK6i Display and Alternative Driver Control.

- SPJAP SPJ+ with Powered Seating Switch ..... Standard  
 SPJAPM Single Actuator Operation thru SPJ+ ..... \$1,300.00 E2310  
Allows Actuator Operation through Driver Control. For Single Actuator Powered Seating Systems.

## EXPANDABLE CONTROLLER SYSTEMS

- PWH Harness req. for Expandable Electronics ..... \$350.00 E2313  
Harness is required to interface joystick or MK6i display with expandable controllers. Option is located under Electronic Accessories in configurator.
- CMPJM6U MPJ+ Multiple Drive Joystick with Color Display ... \$750.00 E2377
- PSFM6 Personal Joystick, Inductive on Front ..... \$1,130.00 E2377  
TOD-SS: Toggle On/Off Drive Select, Speed Potentiometer Standard.
- PSRM6 Personal Joystick, Inductive on Rear ..... \$1,130.00 E2377  
TOD-SS: Toggle On/Off Drive Select, Speed Potentiometer Standard.
- PBOD Push Button On/Off with PB Drive Select . No Charge  
Only available with PSRM6 option.
- PBSS Push Button On/Off with Speed Select ..... No Charge  
Only available with PSRM6 option.  
Drive Select through mode switch or mode switch and Joystick Command.
- OR On/Off Mounted Right ..... No Charge  
Only available with PSRM6 option.
- OL On/Off Mounted Left ..... No Charge  
Only available with PSRM6 option.
- ASLRDYM6 MK6 Display Only (No Driver Control) ..... \$750.00 E2377  
Display Only, No Driver Control. Allows choice from order form.
- 1500M6 Proportional RIM Head Control ..... \$2,875.00 E2327
- SMHD Small Headrest ..... No Charge
- LGHD Large Headrest ..... No Charge
- 1558MM6 Compact Joystick with Display ..... \$1,130.00 E2377

**For your ASL Drive System Needs**  
Please refer to the following form:  
ASL Single Invoice Price List Form # 02-097

## ALTERNATE DIGITAL CONTROL SYSTEMS

Must also use either a Multiple Drive Joystick or MK6i Display.

- SNPM6 Sip-N-Puff Interface Box ..... \$1,500.00 E2325
- 5020M6 Mini Tash Joystick and D/I Box ..... \$1,750.00 E2321

## SIP N PUFF KIT

- PKG32666 Therafin® Sip-N-Puff Breath Tube Kit ..... \$380.00 E2326

## ADDITIONAL CONTROL CHOICES

- 1558M6 Compact Joystick - Less Std Mount ..... \$770.00
- 5020 Mini Tash Joystick ..... \$405.00  
Requires SNPM6 Sip-N-Puff Interface Box.
- PACM6 Proportional Attendant Control ..... \$850.00  
Requires SNPM6 Sip-N-Puff Interface Box.
- 1552M Digital Attendant Control ..... \$850.00

**NOTE: All specifications and dimensions are approximate.**  
For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).

# Invacare® Pronto® M91™ with SureStep® (Rehab Version)

## MK6i ELECTRONICS continued

### JOYSTICK MOUNTING HARDWARE

- QL Quad Link Retract Joystick Mnt SPJ+ ..... \$230.00 E1028  
**QL and QLM6. ARM250 does not come standard with Quad Links on Pronto Series.**
- QLM6 Quad Link Retract Joystick Mnt MPJ+ ..... \$279.00 E1028  
**QL and QLM6. ARM250 does not come standard with Quad Links on Pronto Series.**
- ARM250 Stealth Height Adjustable Joystick Mount ..... \$210.00
- RIMHW Cane Mounted Head Control Hardware ..... \$219.00 E1028
- GATMPJ6 Gatlin Midline for MPJ+ ..... \$600.00 E1028
- GAT1812 Gatlin Midline for Compact Joystick ..... \$600.00 E1028
- PKG32669 Complete Bib Assembly ..... \$265.00
- RCM Rear Cane Mount for Compact Joystick ..... \$395.00
- LEFT Left Handed Mounted Joystick ..... No Charge

### ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer with Pro Memory Card ..... \$475.00
- 1813SD Pro Memory Card (USB Ready) ..... \$90.00
- 1560 "T" Handle Flexible Joystick Ext. .... \$77.00 E2323
- 1561 Straight Handle Flexible Joystick Ext. .... \$77.00 E2323
- 1826 Chin Cup ..... \$48.00 E2324
- A24VPS 24V MK6 Auxiliary Power Source ..... \$200.00
- AUX12M6 Auxiliary Module For 1 and 2 ..... \$795.00
- AUX34M6 Auxiliary Module For 3 and 4 ..... \$795.00
- MPY Monoport "Y" Cable ..... \$15.00
- RFM Wireless (RF) Mouse Only Module ..... \$750.00  
**RFM includes Transmitter, USB Radio Frequency Mouse Dongle (Receiver), and Stereo Switch with hook and loop mount for mouse clicks. Not compatible with SPJ+.**
- IR Infrared Module with Wireless RF Mouse Control ..... \$995.00  
**Capable of operating up to six remote Infra-Red controlled devices. Also allows wireless (RF) Mouse operation through the driver control. Compatible with RFM module for control of two wireless mice from one chair. Includes two button stereo switch (with hook/loop mount) for mouse clicks. Not compatible with SPJ+ Joysticks.**

### ASL OPTIONAL ACCESSORIES

- EG5BLK Egg Switch - Black ..... \$156.00
- ASL304A Light Wobble Switch with 611 Mount ..... \$180.00
- ASL202J Fiber Optic Switch ..... \$732.00
- ASL208J Adjustable Proximity Switch ..... \$660.00
- ASL504M6 Remote Emergency Stop Switch ..... \$1,185.00

## FRONT RIGGING OPTIONS

### FOOTBOARD

- 90FB Footboard ..... Standard
- L90FB Omit Footboard ..... Less \$40.00

### FOOTREST/LEGREST (HEMI)

- PHW93 Swingaway Footrests 13 1/2" - 19" ..... \$165.00
- PH904A Swingaway ELR 16" - 19 3/4" ..... \$260.00 E0990
- PHAL4A Smartleg® Art. Legrests 17 1/4" - 21 1/2" ..... \$225.00 K0053

### FOOTPLATE OPTIONS

- 1650 Aluminum Footplates ..... No Charge
- 1350 Extra Large Aluminum Footplates ..... \$115.00
- AT5543 Adjustable Angle Flip-Up Footplates ..... \$165.00
- 1450 Aluminum Footplates Deluxe ..... No Charge

### FOOTREST/LEGREST ACCESSORIES

- Clear
- LPT2 Longer Pivot Slide Tube For 1350 ..... \$83.00
- ALPT Longer Pivot And Slide Tube ..... \$84.00
- ALPT4 Longer Pivot And Slide Tube 4" ..... \$84.00
- A1337 Calf Strap ..... \$27.00 K0038
- IMPACTG Impact Guards ..... \$51.00
- HLOOP Heel Loops ..... \$45.00 E0951

#### Memory Card - PROFESSIONAL

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains Advanced Diagnostics, Help Screens, Software Upgrades, File Storage/Retrieval
- Used to back up multiple chairs and programming settings.
- Intended for use by qualified providers only.
- Not compatible with SPJ+ Joysticks.

#### Memory Card - BASIC

- Standard with ALL Expandable Control Systems (MPJ+, PSR+, PSF+, MK6i Display).
- Used to Back Up/Restore Programmed Settings for Only One Chair.
- May be given to user or kept with provider's files for safe keeping.
- Does not Contain Advanced Diagnostics/Help Screens/File Structure.
- Not compatible with SPJ+ Joysticks.

#### Seating and Positioning

For your seating needs please refer to the following:  
Rehab Seating and Positioning Products  
see Price List Form 07-027

**NOTE: All specifications and dimensions are approximate.**  
**For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

# Invacare® Pronto® M91™ with SureStep® (Rehab Version)

Specifications	
<b>Dimensions</b>	
Lengths (with Footboard folded)	39"
Lengths (with Footboard down)	40" - 42"
Widths (excluding Joystick)	25.9"
Height	37" - 40"
<b>Seat dimensions</b>	
<b>ADJASBA Seat</b>	
Widths	16" - 24"
Depths	16" - 22"
Back Height	16" - 20"
Seat to Floor Heights	
<b>ADJASBA</b>	18 1/2" - 21"
<b>Weight</b>	
With Batteries	200 lbs.
Without Batteries	150 lbs.
<b>Weight of Components Disassembled</b>	
Battery and Cover Assembly (2 per unit)	25 lbs.
Battery tray, charger, controller shrouds	21 lbs.
Side frame, motor gearbox, side shrouds (2 per unit)	39 lbs. (ea.)
Footboard and mounting tube	5 lbs.

Specifications	
<b>Wheel and Caster Sizes</b>	
Front Casters	6" x 2"
Drive Wheels	14" x 3"
Rear Casters	6" x 2"
<b>Maximum Speed</b>	
Standard	6 mph
Heavy Duty	5 mph
<b>Incline Capability</b>	
	9°
<b>Weight Capacity</b>	
	300 lbs. (400 lb. option)
<b>Turning Radius</b>	
Front (with footboard down)	19 1/2"
Rear	21 1/2"
Ground Clearance	3 1/4"
<b>Warranties</b>	
Base frame	60 months
Seat Frame	36 months
Electronics	12 months
Motor/gearbox	12 months

Comments

**NOTE: All specifications and dimensions are approximate.**  
*For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).*

**Invacare Corporation**

[www.invacare.com](http://www.invacare.com)

**USA**

One Invacare Way  
 Elyria, OH USA  
 44035-4190  
 800-333-6900

©2010 Invacare Corporation  
 All rights reserved. Trademarks are identified by the symbols ™ and ®.  
 All trademarks are owned by or licensed to Invacare Corporation unless  
 otherwise noted. Products designated by ASL are Manufactured by  
 Adaptive Switch Labs, Inc. Therafin is a registered trademark of W/R  
 Medical Electronics Co.

