

Invacare eForms

For additional information, visit www.invacare.com/eForms

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order Literature: 800-828-6282

Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.

Please consider the environment before printing this page with your order form.
For more information on features and benefits of electronic order forms, please visit www.invacare.com/eForms

Invacare® Pronto® M91™ with SureStep® (Rehab Version)

Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



TRANSACTION TYPE:

- Order Quote

Date _____ Account # _____ P.O. # _____

Company Name _____ Phone # _____

Address _____ City _____ State _____ Zip Code _____ Country _____

Contact Name _____ Comments _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

FRAME SELECTION

POWER CHAIR TYPE

- M91R M91 with SureStep with MK6i™ electronics .. \$5,600.00 K0848
Push button style seat positioning strap included.
- M91-C M91 Complete with SureStep Base Only, with MK6i electronics \$5,255.00
- M91-M M91 Motion Concepts with SureStep Base Only, with MK6i electronics \$5,255.00
- M91SEAT Pronto M91 Seat Only \$740.00

TRANSPORT TIE DOWN

- TRRO Transport Ready Option \$300.00
- TRBKTS Wheelchair Transport Brackets \$250.00

TRRO Features for M91 with SureStep

- Pronto M91R with MK6i
- 300 lbs. or 400 lbs. User Weight Capacity
- Adjustable Rehab Seat with Adj Seat and Tilt Back Angle
- Any Seat Width
- Any Seat Depth
- Sectional Rear Uprights
- 18" - 20" Back Height
- Any Arm Type
- Black Nylon Upholstery
- 14" x 3" Flat Free Tires
- Any Shroud Color
- Footboard
- MPJ+™
- All PSR/PSF Joysticks
- Left Handed Mounted Joystick Option
- "T" Handle Flexible Joystick Ext.
- Straight Handle Flexible Joystick Extension
- MK6i Remote Programmer
- Reinforced Back Upholstery
- 22NF Batteries
- Tire Pump
- Lester 8 amp Battery Charger

Transport Tie Down

TRRO includes four factory-installed transport brackets and a wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNAWC Vol 1. Section 19.5.3 Frontal Impact Test requirements for wheelchairs with a 168 lbs. crash dummy, which corresponds to a person with a weight of 114 lbs. to 209 lbs.

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

*Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

SEAT SELECTION

- ADJASBA Seat with Adj. Seat Tilt and Back Angle Standard
- CTS Contoura™ Seating \$395.00
Selected Seat Depth on Contoura assumes 1" of back cushion compression. For example: When choosing 18" depth the measurement from the front of the back cushion to the front of the seat pan equals 17". The Contoura back cushion is 2" thick.
- OMCU Omit Contoura Seat Cushion Less \$30.00
- SEMIRB Van Seat - Semi Recline with Solid Base \$70.00

MOTOR PACKAGES

- U300 User Weight Up to 300lbs. (STDPKG) Standard K0849
Standard motor package includes 4-pole motors with a weight limit of 300 lbs.
- U400 User Weight 301 - 400 lbs. (HDPKG) \$400.00 K0851

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.**

Invacare® Pronto® M91™ with SureStep® (Rehab Version)

ADJUSTABLE ASBA SEAT SELECTIONS

ADJUSTABLE SEAT WIDTH RANGE

- 1620W 16" to 20" Seat Width Range No Charge
- 2024W 20" to 24" Seat Width Range No Charge

SEAT WIDTH

- W16 16" Seat Width No Charge
- W17 17" Seat Width No Charge
- W18 18" Seat Width No Charge
- W19 19" Seat Width No Charge
- W20 20" Seat Width No Charge
- W21 21" Seat Width \$329.00
- W22 22" Seat Width \$329.00
- W23 23" Seat Width \$329.00
Only available with U400 motor package.
- W24 24" Seat Width \$469.00
Only available with U400 motor package.

ADJUSTABLE SEAT DEPTH RANGE

- 1619D 16" to 19" Seat Depth Range No Charge
- 1922D 19" to 22" Seat Depth Range No Charge

SEAT DEPTH

- D16 16" Seat Depth Standard
- D17 17" Seat Depth No Charge
- D18 18" Seat Depth No Charge
- D19 19" Seat Depth No Charge
- D20 20" Seat Depth \$529.00
- D21 21" Seat Depth \$529.00
- D22 22" Seat Depth \$529.00

UPHOLSTERY

- U240 Black Nylon Upholstery No Charge

BACK SELECTIONS

SEAT BACK HEIGHT

- BH16 16" Back Height No Charge
- BH17 17" Back Height No Charge
- BH18 18" Back Height No Charge
- BH19 19" Back Height No Charge
- BH20 20" Back Height No Charge

CONTOURA

- TSBH25 25" Back Height No Charge
- TSBH28 28" Back Height No Charge

SEAT BACK TYPE

- SECT Sectional Rear Uprights No Charge

ARM SELECTIONS

ARM TYPE (PAIR)

- 28F Desk Length Fixed Ht Conv Standard
- 29F Desk Length Adj. Ht Conv \$230.00 E0973
- 18F Full Length Fixed Ht Conv No Charge
- 19F Full Length Adj. Ht Conv \$230.00 E0973

MIXED ARMS (EACH)

- 19FR Right Arm, Full Length, Adj. Height \$115.00 E0973
- 29FR Right Arm, Desk Length, Adj. Height \$115.00 E0973
- 19FL Left Arm, Full Length, Adj. Height \$115.00 E0973
- 29FL Left Arm, Desk Length, Adj. Height \$115.00 E0973

VAN SEAT - WITH SOLID BASE

VAN SEAT WIDTH AND DEPTH

- VW16 16" Wide Van Seat No Charge
 - VD16 16" Depth No Charge
 - VD17 17" Depth No Charge
 - VD18 18" Depth No Charge
- VW18 18" Wide Van Seat Standard
 - VD16 16" Depth No Charge
 - VD17 17" Depth No Charge
 - VD18 18" Depth No Charge
- VW20 20" Wide Van Seat \$329.00
 - VD18 18" Depth No Charge
 - VD19 19" Depth No Charge
 - VD20 20" Depth \$529.00
- VW22 22" Wide Van Seat \$529.00
 - VD18 18" Depth No Charge
 - VD19 19" Depth No Charge
 - VD20 20" Depth \$529.00

VAN UPHOLSTERY COLORS

- GV Gray Vinyl Upholstery No Charge

VAN ARM TYPE (PAIR)

- VAH Full Length, Adj. Height \$230.00 E0973
- VFH Full Length, Fixed Height Standard
- VAHD Desk Length, Adj. Height \$230.00 E0973
- VFHD Desk Length, Fixed Height No Charge

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Pronto® M91™ with SureStep® (Rehab Version)

POSITIONING STRAPS AND WHEEL OPTIONS

SEAT POSITIONING STRAPS

Clear

- 1515 Push Button Style Seat Pos. Strap Standard
 1311BK Airline Buckle Seat Pos. Strap No Charge

CHEST POSITIONING STRAPS

Clear

- 7311BK Airline Buckle Chest Pos. Strap \$100.00 E0960
 9515 Push Button Style Chest Pos. Strap \$100.00 E0960
 7321BK Hook & Loop Chest Pos. Strap \$100.00 E0960

SHROUD COLORS

- BL Deep Blue Pearl No Charge
 RD Candy Red Pearl No Charge
 MB Magic Black No Charge

DRIVE WHEEL TYPE

- DW143 14" x 3" Pneumatic Tires Standard
 DW143FF 14" x 3" Flat Free Tires \$149.00

ACCESSORIES

BACK ACCESSORIES

- L8601 Omit Back Upholstery No Charge
 RBP Reinforced Back Upholstery No Charge

ARM ACCESSORIES

Clear

- L8401 Omit Arms Less \$20.00
 L8405 Omit Arm Pads Less \$10.00

MISCELLANEOUS ACCESSORIES

- SEAT02 Oxygen Holder (ASBA Seat Only) \$150.00 E2208
 HFO2-ML6 HomeFill™ O2 Holder-ML6 Only \$150.00 E2208
 1209 Tire Pump \$35.00
 LESTR Lester 8 Amp Battery Charger \$320.00 E2367
 VANO2 Oxygen Holder (Van Seat Only) \$150.00 E2208
 CTE Single Mode Charger \$205.00 E2366

BATTERY OPTIONS

INSTALLATION OF BATTERIES

- 22NFINSTAL 22NF Battery(On-Chair)(1 pair) \$300.00 E2361

MK6i ELECTRONICS

ELECTRONIC CONTROLLERS

- MK690 MK6i 90 Amp Controller Standard

STANDARD DRIVER CONTROL

- SPJM6 SPJ™+ Single Drive Joystick Standard

BASE ONLY SPJ+ OPTIONS

Base Only Systems for Multiple Actuator Powered Seating Systems will require MPJ+, PSRM6, PSFM6, or MK6i Display and Alternative Driver Control.

- SPJAP SPJ+ with Powered Seating Switch Standard
 SPJAPM Single Actuator Operation thru SPJ+ \$1,300.00 E2310
Allows Actuator Operation through Driver Control. For Single Actuator Powered Seating Systems.

EXPANDABLE CONTROLLER SYSTEMS

- PWH Harness req. for Expandable Electronics \$350.00 E2313
Harness is required to interface joystick or MK6i display with expandable controllers. Option is located under Electronic Accessories in configurator.
- CMPJM6U MPJ+ Multiple Drive Joystick with Color Display ... \$750.00 E2377
- PSFM6 Personal Joystick, Inductive on Front \$1,130.00 E2377
TOD-SS: Toggle On/Off Drive Select, Speed Potentiometer Standard.
- PSRM6 Personal Joystick, Inductive on Rear \$1,130.00 E2377
TOD-SS: Toggle On/Off Drive Select, Speed Potentiometer Standard.
- PBOD Push Button On/Off with PB Drive Select . No Charge
Only available with PSRM6 option.
- PBSS Push Button On/Off with Speed Select No Charge
Only available with PSRM6 option.
Drive Select through mode switch or mode switch and Joystick Command.
- OR On/Off Mounted Right No Charge
Only available with PSRM6 option.
- OL On/Off Mounted Left No Charge
Only available with PSRM6 option.
- ASLRDYM6 MK6 Display Only (No Driver Control) \$750.00 E2377
Display Only, No Driver Control. Allows choice from order form.
- 1500M6 Proportional RIM Head Control \$2,875.00 E2327
- SMHD Small Headrest No Charge
- LGHD Large Headrest No Charge
- 1558MM6 Compact Joystick with Display \$1,130.00 E2377

For your ASL Drive System Needs
Please refer to the following form:
ASL Single Invoice Price List Form # 02-097

ALTERNATE DIGITAL CONTROL SYSTEMS

Must also use either a Multiple Drive Joystick or MK6i Display.

- SNPM6 Sip-N-Puff Interface Box \$1,500.00 E2325
- 5018M6 Wafer Board and D/I Box \$1,550.00 E2322
- 5020M6 Mini Tash Joystick and D/I Box \$1,750.00 E2321

SIP N PUFF KIT

- PKG32666 Therafin® Sip-N-Puff Breath Tube Kit \$380.00 E2326

ADDITIONAL CONTROL CHOICES

- 1558M6 Compact Joystick - Less Std Mount \$770.00
- 5018 Wafer Board \$310.00
Requires SNPM6 Sip-N-Puff Interface Box.
- 5020 Mini Tash Joystick \$405.00
Requires SNPM6 Sip-N-Puff Interface Box.
- PACM6 Proportional Attendant Control \$850.00
Requires SNPM6 Sip-N-Puff Interface Box.
- 1552M Digital Attendant Control \$850.00

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Pronto® M91™ with SureStep® (Rehab Version)

MK6i ELECTRONICS continued

JOYSTICK MOUNTING HARDWARE

- QL Quad Link Retract Joystick Mnt SPJ+ \$230.00 E1028
QL and QLM6. ARM250 does not come standard with Quad Links on Pronto Series.
- QLM6 Quad Link Retract Joystick Mnt MPJ+ \$279.00 E1028
QL and QLM6. ARM250 does not come standard with Quad Links on Pronto Series.
- ARM250 Stealth Height Adjustable Joystick Mount \$210.00
- RIMHW Cane Mounted Head Control Hardware \$219.00 E1028
- GATMPJ6 Gatlin Midline for MPJ+ \$600.00 E1028
- GAT1812 Gatlin Midline for Compact Joystick \$600.00 E1028
- PKG32669 Complete Bib Assembly \$265.00
- RCM Rear Cane Mount for Compact Joystick \$395.00
- LEFT Left Handed Mounted Joystick No Charge

ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer with Pro Memory Card \$475.00
- 1813SD Pro Memory Card (USB Ready) \$90.00
- 1560 "T" Handle Flexible Joystick Ext. \$77.00 E2323
- 1561 Straight Handle Flexible Joystick Ext. \$77.00 E2323
- 1826 Chin Cup \$48.00 E2324
- A24VPS 24V MK6 Auxiliary Power Source \$200.00
- AUX12M6 Auxiliary Module For 1 and 2 \$795.00
- AUX34M6 Auxiliary Module For 3 and 4 \$795.00
- MPY Monoport "Y" Cable \$15.00

ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Switch - Black \$65.00
- ASL304A Light Wobble Switch with 611 Mount \$320.00
- ASL202J Fiber Optic Switch \$610.00
- ASL208J Adjustable Proximity Switch \$550.00
- ASL504M6 Remote Emergency Stop Switch \$1,185.00

FRONT RIGGING OPTIONS

FOOTBOARD

- 90FB Footboard Standard
- L90FB Omit Footboard Less \$40.00

FOOTREST/LEGREST (HEMI)

- PHW93 Swingaway Footrests 13 1/2" - 19" \$165.00
- PH904A Swingaway ELR 16" - 19 3/4" \$165.00 E0990
- PHAL4A Smartleg® Art. Legrests 17 1/4" - 21 1/2" \$225.00 K0053

FOOTPLATE OPTIONS

- 1650 Aluminum Footplates No Charge
- 1350 Extra Large Aluminum Footplates \$115.00
- AT5543 Adjustable Angle Flip-Up Footplates \$165.00
- 1450 Aluminum Footplates Deluxe No Charge

FOOTREST/LEGREST ACCESSORIES

- Clear
- LPT2 Longer Pivot Slide Tube For 1350 \$83.00
- ALPT Longer Pivot And Slide Tube \$84.00
- ALPT4 Longer Pivot And Slide Tube 4" \$84.00
- A1337 Calf Strap \$27.00 K0038
- IMPACTG Impact Guards \$51.00
- 1600BK Heel Loops with Ankle Straps \$56.00 E0951

5 STAR PREMIUM SERVICE PLAN

- 5STAR1 Invacare 5 Star Premium Service Plan \$179.00
5STAR1 Guarantees predictable service costs for one year on Invacare Consumer Power Wheelchairs. Visit www.invacare.com/partner for more information.

Memory Card - PROFESSIONAL

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains Advanced Diagnostics, Help Screens, Software Upgrades, File Storage/Retrieval
- Used to back up multiple chairs and programming settings.
- Intended for use by qualified providers only.
- Not compatible with SPJ+ Joysticks.

Memory Card - BASIC

- Standard with ALL Expandable Control Systems (MPJ+, PSR+, PSF+, MK6i Display).
- Used to Back Up/Restore Programmed Settings for Only One Chair.
- May be given to user or kept with provider's files for safe keeping.
- Does not Contain Advanced Diagnostics/Help Screens/File Structure.
- Not compatible with SPJ+ Joysticks.

Seating and Positioning

For your seating needs please refer to the following:
Rehab Seating and Positioning Products
see Price List Form 07-027

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Pronto® M91™ with SureStep® (Rehab Version)

Specifications	
Dimensions	
Lengths (with Footboard folded)	39"
Lengths (with Footboard down)	40" - 42"
Widths (excluding Joystick)	25.9"
Height	37" - 40"
Seat dimensions	
ADJASBA Seat	
Widths	16" - 24"
Depths	16" - 22"
Back Height	16" - 20"
Seat to Floor Heights	
ADJASBA	18 1/2" - 21"
Weight	
With Batteries	200 lbs.
Without Batteries	150 lbs.
Weight of Components Disassembled	
Battery and Cover Assembly (2 per unit)	25 lbs.
Battery tray, charger, controller shrouds	21 lbs.
Side frame, motor gearbox, side shrouds (2 per unit)	39 lbs. (ea.)
Footboard and mounting tube	5 lbs.

Specifications	
Wheel and Caster Sizes	
Front Casters	6" x 2"
Drive Wheels	14" x 3"
Rear Casters	6" x 2"
Maximum Speed	
Standard	6 mph
Heavy Duty	5 mph
Incline Capability	
	9°
Weight Capacity	
	300 lbs. (400 lb. option)
Turning Radius	
Front (with footboard down)	19 1/2"
Rear	21 1/2"
Ground Clearance	3 1/4"
Warranties	
Base frame	60 months
Seat Frame	36 months
Electronics	12 months
Motor/gearbox	12 months

Comments

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare Corporation

www.invacare.com

USA

One Invacare Way
 Elyria, OH USA
 44035-4190
 800-333-6900

©2010 Invacare Corporation
 All rights reserved. Trademarks are identified by the symbols ™ and ®.
 All trademarks are owned by or licensed to Invacare Corporation unless
 otherwise noted. Products designated by ASL are Manufactured by
 Adaptive Switch Labs, Inc. Therafin is a registered trademark of W/R
 Medical Electronics Co.

