



## Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email\* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

\*Saving and emailing the eForms is available with Acrobat Professional software

### If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

**Please note:** The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

### If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit [www.adobe.com](http://www.adobe.com).

# Invacare® Pronto® M51™ with SureStep® (Rehab Version)

## Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



Yes, you can.®

### TRANSACTION TYPE:

Order  Quote

Date \_\_\_\_\_ Account # \_\_\_\_\_ P.O. # \_\_\_\_\_

Company Name \_\_\_\_\_ Phone # \_\_\_\_\_

Address \_\_\_\_\_ City \_\_\_\_\_ State \_\_\_\_\_ Zip Code \_\_\_\_\_ Country \_\_\_\_\_

Contact Name \_\_\_\_\_ Comments \_\_\_\_\_

### Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

## PRECONFIGURED STOCK SELECTIONS

### STOCK MODEL WITH SEMI-RECLINE VAN SEAT WITH SOLID SEAT PAN

|                          |               |                         |            |       |
|--------------------------|---------------|-------------------------|------------|-------|
| <input type="checkbox"/> | M51PRSOLID16R | 16"W x 16"D, Red .....  | \$5,986.00 | K0822 |
| <input type="checkbox"/> | M51PRSOLID16B | 16"W x 16"D, Blue ..... | \$5,986.00 | K0822 |
| <input type="checkbox"/> | M51PRSOLID18R | 18"W x 18"D, Red .....  | \$5,986.00 | K0822 |
| <input type="checkbox"/> | M51PRSOLID18B | 18"W x 18"D, Blue ..... | \$5,986.00 | K0822 |
| <input type="checkbox"/> | M51PRSOLID20R | 20"W x 18"D, Red .....  | \$5,986.00 | K0822 |
| <input type="checkbox"/> | M51PRSOLID20B | 20"W x 18"D, Blue ..... | \$5,986.00 | K0822 |

#### \*Total MSRP (Base Model plus Pre-Configured Features)

- Save time and money with stock models available for quick shipment!
- Pronto M51 Power Base w/ Solid Seat Pan... (MSRP) \$5454.00 K0822
  - Semi-Recline Van Seat w/ Solid Seat Pan, High Back, Light Gray Vinyl ..... (MSRP) \$70.00
- Note: All stocked M51 (and M41) seats swivel, have a center post seat mount, and are interchangeable.*
- U1 Batteries Installed ..... (MSRP) \$232.00 E2365
  - Full Length, Adjustable Height Arms .... (MSRP) \$230.00 E0973
  - Flat Free Tires
  - Flip-Up Footboard
  - MK5 NX controller with SPJ+ Joystick
  - Push-Button Style Seat Positioning Strap

### FRONT RIGGING OPTIONS

|  |                |                                     |                |
|--|----------------|-------------------------------------|----------------|
| <input type="checkbox"/>   | FRAKH          | Front Rigging Attachment Kit .....  | \$55.00        |
| <b>FRAKH required for use with front riggings on stock models.</b> |                |                                     |                |
| <input type="checkbox"/>   | PHW93/34P      | Swing Away Footrests .....          | \$165.00       |
| <input type="checkbox"/>   | PH904A/U67/34P | Swing-Away Elevating Legrests ..... | \$260.00 E0990 |

### ABSOLUTE™ CUSHION General Purpose Cushion

|                          |      |                              |          |       |
|--------------------------|------|------------------------------|----------|-------|
| <input type="checkbox"/> | EC66 | Absolute Cushion 16x16 ..... | \$104.00 | E2601 |
| <input type="checkbox"/> | EC88 | Absolute Cushion 18x18 ..... | \$104.00 | E2601 |
| <input type="checkbox"/> | EC08 | Absolute Cushion 20x18 ..... | \$104.00 | E2601 |

### COMFORT-MATE EXTRA CUSHION Skin Protection Cushion

|                          |        |                                |          |       |
|--------------------------|--------|--------------------------------|----------|-------|
| <input type="checkbox"/> | CMEX66 | Comfort Mate Extra 16x16 ..... | \$208.00 | E2603 |
| <input type="checkbox"/> | CMEX88 | Comfort Mate Extra 18x18 ..... | \$208.00 | E2603 |
| <input type="checkbox"/> | CMEX08 | Comfort Mate Extra 20x18 ..... | \$208.00 | E2603 |

### SWINGAWAY JOYSTICK KITS FOR STOCKED MODELS

|                          |         |  |          |       |
|--------------------------|---------|--|----------|-------|
| <input type="checkbox"/> | QLRIGHT | Quad Link Joystick Mount Kits, Right Side ...  | \$230.00 | E1028 |
| <input type="checkbox"/> | QLLEFT  | Quad Link Joystick Mount Kits, Left Side ..... | \$230.00 | E1028 |

**NOTE: All specifications and dimensions are approximate.**  
For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).

# Invacare® Pronto® M51™ with SureStep® (Rehab Version)

## FRAME SELECTION

### POWER CHAIR TYPE-CONFIGURABLE

- M51PR M51 with SureStep with MK5™ NX w/SPJ+ ..... \$5,454.00 K0822
- M50SEAT Pronto M51 Seat Only ..... \$740.00

### TRANSPORT TIE DOWN

- TRRO Transport Ready Option ..... \$300.00  
**Not available with SEMIRB seat.**
- TRBKTS Wheelchair Transport Brackets ..... \$250.00

#### Transport Tie Down

TRRO includes four factory-installed transport brackets and a wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNAWC Vol 1, Section 19.5.3 Frontal Impact Test requirements for wheelchairs with a 168 lbs. crash dummy, which corresponds to a person with a weight of 114 lbs. to 209 lbs.

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

\*Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

## SEAT SELECTIONS

### SEAT SELECTION

*Note: Both ADJASBA and SEMIRB seat styles on model M51PR have a 4-post seat mount for adjustable seat-to-floor height, adjustable seat slope, and can be easily retrofitted with a power tilt kit.*

- ADJASBA Seat with Adj. Seat Tilt and Back Angle ..... Standard
- SEMIRB Van Seat - Semi Recline with Solid Base ..... \$100.00

## ADJUSTABLE ASBA SEAT SELECTIONS

### ADJUSTABLE SEAT WIDTH RANGE

- 1620W 16" to 20" Seat Width Range ..... No Charge

### SEAT WIDTH

- W16 16" Seat Width ..... No Charge
- W17 17" Seat Width ..... No Charge
- W18 18" Seat Width ..... No Charge
- W19 19" Seat Width ..... No Charge
- W20 20" Seat Width ..... No Charge

### ADJUSTABLE SEAT DEPTH RANGE

- 1619D 16" to 19" Seat Depth Range ..... No Charge
- 1920D 19" to 20" Seat Depth Range ..... No Charge

## SEAT DEPTH

- D16 16" Seat Depth ..... Standard
- D17 17" Seat Depth ..... No Charge
- D18 18" Seat Depth ..... No Charge
- D19 19" Seat Depth ..... No Charge
- D20 20" Seat Depth ..... \$529.00

## UPHOLSTERY COLORS

- U240 Black Nylon Upholstery ..... No Charge
- L8601 Omit Back Upholstery ..... Less \$10.00
- RBP Reinforced Back Upholstery ..... No Charge

## SEAT BACK TYPE

- SECT Sectional Rear Uprights ..... No Charge

## SEAT BACK HEIGHT

- BH16 16" Back Height ..... No Charge  
**Not available with TRRO.**
- BH17 17" Back Height ..... No Charge  
**Not available with TRRO.**
- BH18 18" Back Height ..... No Charge
- BH19 19" Back Height ..... No Charge
- BH20 20" Back Height ..... No Charge

## ARM TYPE (PAIR)

**Arm Types are flip-up and detachable.**

- 28F Flip Up Desk Length Fixed Ht Conv ..... No Charge
- 29F Flip Up Desk Length Adj. Ht Conv ..... \$230.00 E0973
- 18F Flip Up Full Length Fixed Ht Conv ..... No Charge
- 19F Flip Up Full Length Adj. Ht Conv ..... \$230.00 E0973

## MIXED ARMS (EACH)

**The below arm selections are flip-up and detachable.**

- 19FR Right Arm, Full Length, Adj. Height ..... \$115.00 E0973
- 29FR Right Arm, Desk Length, Adj. Height ..... \$115.00 E0973
- 19FL Left Arm, Full Length, Adj. Height ..... \$115.00 E0973
- 29FL Left Arm, Desk Length, Adj. Height ..... \$115.00 E0973

## ARM ACCESSORIES

- Clear
- L8401 Omit Arms ..... Less \$20.00
- L8405 Omit Arm Pads ..... Less \$10.00

**NOTE: All specifications and dimensions are approximate.**  
**For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

# Invacare® Pronto® M51™ with SureStep® (Rehab Version)

## VAN SEAT WITH SOLID SEAT PAN

### VAN UPHOLSTERY COLORS

- GV Gray Vinyl Upholstery ..... No Charge

### VAN SEAT WIDTH AND DEPTH

- VW16 16" Wide Van Seat ..... No Charge
- VD16 16" Depth ..... No Charge
  - VD17 17" Depth ..... No Charge
  - VD18 18" Depth ..... No Charge
- VW18 18" Wide Van Seat ..... Standard
- VD16 16" Depth ..... Standard
  - VD17 17" Depth ..... No Charge
  - VD18 18" Depth ..... No Charge
- VW20 20" Wide Van Seat ..... \$329.00
- VD18 18" Depth ..... No Charge
  - VD19 19" Depth ..... No Charge
  - VD20 20" Depth ..... \$329.00

### VAN ARM TYPE (PAIR)

- VAH Full Length, Adj. Height ..... \$230.00 E0973
- VFH Full Length, Fixed Height ..... Standard
- VAHD Desk Length, Adj. Height ..... \$230.00 E0973
- VFHD Desk Length, Fixed Height ..... No Charge

### VAN MIXED ARMS

Arms are flip-up and adjustable width.

- Right Arm  Clear
- RVAHF Right Arm, Full Length, Adj. Height ..... \$115.00 E0973
  - RVFHF Right Arm, Full Length, Fixed Height ..... No Charge
  - RVAHD Right Arm, Desk Length, Adj. Height ..... \$115.00 E0973
  - RVFHD Right Arm, Desk Length, Fixed Height ..... No Charge
- Left Arm  Clear
- LVAHF Left Arm, Full Length, Adj. Height ..... \$115.00 E0973
  - LVFHF Left Arm, Full Length, Fixed Height ..... No Charge
  - LVAHD Left Arm, Desk Length, Adj. Height ..... \$115.00 E0973
  - LVFHD Left Arm, Desk Length, Fixed Height ..... No Charge

## POSITIONING STRAPS

### CHEST POSITIONING STRAPS

Clear

- 7311BK Airline Buckle Chest Pos. Strap ..... \$100.00 E0960
- 9515 Push Button Style Chest Pos. Strap ..... \$100.00 E0960
- 7321BK Hook & Loop Chest Pos. Strap ..... \$100.00 E0960

### SEAT POSITIONING STRAPS

Clear

- 1515 Push Button Style Seat Pos. Strap ..... Standard
- 1311BK Airline Buckle Seat Pos. Strap ..... \$50.00

## DRIVE WHEEL TYPE

### DRIVE WHEELS

- DW10535 10 1/2" x 3 1/2" Pneumatic Tires ..... Standard  
Must select DW10535FF drive wheels when selecting TRRO.
- DW10535FF 10 1/2" x 3 1/2" Flat Free Tires ..... No Charge

## COLORS

### SHROUD COLORS

- RD Candy Red Pearl ..... No Charge
- BL Deep Blue Pearl ..... No Charge

## BATTERY OPTIONS

### INSTALLATION OF BATTERIES

- U1INSTAL U1 Battery Install (On Chair)(1 pair) ..... \$232.00 E2365

## ELECTRONICS

### JOYSTICK MOUNTING HARDWARE

- QL Quad Link Retractable Joystick Mount SPJ+ .... \$230.00 E1028
- LEFT Left Hand Mounted Joystick ..... No Charge

## ACCESSORIES

### MISCELLANEOUS ACCESSORIES

- LESTR Lester 8 Amp Battery Charger ..... \$320.00 E2367
- 1831M6 MK6i Programmer ..... \$350.00
- VANO2 Oxygen Holder (Solid Base Van Seat) ..... \$150.00 E2208
- 1209 Tire Pump ..... \$35.00
- SEATO2 Oxygen Holder (ASBA Seat Only) ..... \$150.00 E2208  
Only one option can be used at one time.
- HFO2-ML6 HomeFill O2 Holder-ML6 Only ..... \$150.00 K0104

## FRONT RIGGING OPTIONS

### FOOTBOARD

- 90FB Footboard ..... Standard
- L90FB Omit Footboard ..... Less \$40.00

### FOOTREST/LEGREST (HEMI)

- PHW93 Swingaway Footrests 13 1/2" - 19" ..... \$165.00
- PH904A Swingaway ELR 16" - 19 3/4" ..... \$260.00 E0990
- PHAL4A Smartleg® Art. Legrests 17 1/4" - 21 1/2" ..... \$225.00 K0053

### FOOTPLATE OPTIONS

- 1350 Extra Large Aluminum Footplates ..... \$115.00
- 1450 Aluminum Footplates Deluxe ..... No Charge
- 1650 Aluminum Footplates ..... No Charge
- 1651 Composite Flip-Up Footplates ..... No Charge
- AT5543 Adjustable Angle Flip-Up Footplates ..... \$165.00

### FOOTREST/LEGREST ACCESSORIES

Clear

- LPT2 Longer Pivot Slide Tube For 1350 ..... \$83.00
- ALPT Longer Pivot And Slide Tube ..... \$84.00
- ALPT4 Longer Pivot And Slide Tube 4" ..... \$84.00
- A1337 Calf Strap ..... \$27.00 K0038
- IMPACTG Impact Guards ..... \$51.00
- HLOOP Heel Loops ..... \$45.00 E0951

**NOTE: All specifications and dimensions are approximate.**

**For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

# Invacare® Pronto® M51™ with SureStep® (Rehab Version)

| Specifications                  | ADJASBA                     | Van Seat w/Solid Base |
|---------------------------------|-----------------------------|-----------------------|
| <b>Dimensions</b>               |                             |                       |
| Lengths (with Footboard folded) | 34"                         | 34"                   |
| Lengths (with Footboard down)   | 38" - 40"                   | 38" - 40"             |
| Widths (excluding Joystick)     | 24"                         | 24"                   |
| Heights                         | 35" - 37"                   | 35" - 37"             |
| Widths                          | 16" - 20"                   | 16" - 20"             |
| Depths                          | 16" - 20"                   | 16", 18", 20"         |
| Back Ht. without Headrest       | 16" - 20"                   | 19"                   |
| Seat to Floor Heights           | 18"-20" (seat pan to floor) | 19" - 21"             |
| <b>Weight</b>                   |                             |                       |
| With Batteries                  | 192 lbs.                    | 198 lbs.              |
| Without Batteries               | 145 lbs.                    | 151 lbs.              |
| <b>Wheel and Caster Sizes</b>   |                             |                       |
| Front Casters                   | 6" x 2"                     | 6" x 2"               |
| Drive Wheels                    | 10 1/2" x 3 1/2"            | 10 1/2" x 3 1/2"      |
| Rear Casters                    | 6" x 2"                     | 6" x 2"               |
| <b>Maximum Speed</b>            | 4 1/4 mph                   | 4 1/4 mph             |
| <b>Incline Capability</b>       | 9°                          | 9°                    |
| <b>Weight Capacity</b>          | 300 lbs.                    | 300 lbs.              |
| <b>Turning Radius</b>           |                             |                       |
| Front (with footboard down)     | 19 1/2"                     | 19 1/2"               |
| Rear                            | 19 1/2"                     | 19 1/2"               |
| Ground Clearance                | 2.6"                        | 2.6"                  |
| <b>Warranties</b>               | See Owner's Manual          | See Owner's Manual    |

**Seating and Positioning**  
 For your seating needs please refer to the following:  
 Rehab Seating and Positioning Products  
 see Price List Form 07-027

Comments

**NOTE: All specifications and dimensions are approximate.**  
 For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).



**Invacare Corporation**  
 USA  
 One Invacare Way  
 Elyria, OH USA  
 44035-4190  
 800-333-6900

[www.invacare.com](http://www.invacare.com)

©2013 Invacare Corporation  
 All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.