



Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.

Invacare® Pronto® M51™ with SureStep® (Rehab Version)

Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



Yes, you can.®

TRANSACTION TYPE:

Order Quote

Date _____ Account # _____ P.O. # _____

Company Name _____ Phone # _____

Address _____ City _____ State _____ Zip Code _____ Country _____

Contact Name _____ Comments _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

PRECONFIGURED STOCK SELECTIONS

STOCK MODEL WITH SEMI-RECLINE VAN SEAT WITH SOLID SEAT PAN

- M51PRSOLID16R 16"W x 16"D, Red \$5,986.00 K0822
- M51PRSOLID16B 16"W x 16"D, Blue \$5,986.00 K0822
- M51PRSOLID18R 18"W x 18"D, Red \$5,986.00 K0822
- M51PRSOLID18B 18"W x 18"D, Blue \$5,986.00 K0822
- M51PRSOLID20R 20"W x 18"D, Red \$5,986.00 K0822
- M51PRSOLID20B 20"W x 18"D, Blue \$5,986.00 K0822

*Total MSRP (Base Model plus Pre-Configured Features)

- Save time and money with stock models available for quick shipment!
- Pronto M51 Power Base w/ Solid Seat Pan... (MSRP) \$5454.00 K0822
 - Semi-Recline Van Seat w/ Solid Seat Pan, High Back, Light Gray Vinyl (MSRP) \$70.00
- Note: All stocked M51 (and M41) seats swivel, have a center post seat mount, and are interchangeable.*
- U1 Batteries Installed (MSRP) \$232.00 E2365
 - Full Length, Adjustable Height Arms (MSRP) \$230.00 E0973
 - Flat Free Tires
 - Flip-Up Footboard
 - MK5 NX controller with SPJ+ Joystick
 - Push-Button Style Seat Positioning Strap

FRONT RIGGING OPTIONS

- FRAKH Front Rigging Attachment Kit \$55.00
FRAKH required for use with front riggings on stock models.
- PHW93/34P Swing Away Footrests \$165.00
- PH904A/U67/34P Swing-Away Elevating Legrests \$260.00 E0990

ABSOLUTE™ CUSHION

General Purpose Cushion

- EC66 Absolute Cushion 16x16 \$104.00 E2601
- EC88 Absolute Cushion 18x18 \$104.00 E2601
- EC08 Absolute Cushion 20x18 \$104.00 E2601

COMFORT-MATE EXTRA CUSHION

Skin Protection Cushion

- CMEX66 Comfort Mate Extra 16x16 \$208.00 E2603
- CMEX88 Comfort Mate Extra 18x18 \$208.00 E2603
- CMEX08 Comfort Mate Extra 20x18 \$208.00 E2603

SWINGAWAY JOYSTICK KITS FOR STOCKED MODELS

- QLRIGHT Quad Link Joystick Mount Kits, Right Side ... \$230.00 E1028
- QLLEFT Quad Link Joystick Mount Kits, Left Side \$230.00 E1028

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Pronto® M51™ with SureStep® (Rehab Version)

FRAME SELECTION

POWER CHAIR TYPE-CONFIGURABLE

- M51PR M51 with SureStep with MK₅™ NX w/SPJ+ \$5,454.00 K0822
- M50SEAT Pronto M51 Seat Only \$740.00

TRANSPORT TIE DOWN

- TRRO Transport Ready Option \$300.00
Not available with SEMIRB seat.
- TRBKTS Wheelchair Transport Brackets \$250.00

Transport Tie Down

TRRO includes four factory-installed transport brackets and a wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNAWC Vol 1, Section 19.5.3 Frontal Impact Test requirements for wheelchairs with a 168 lbs. crash dummy, which corresponds to a person with a weight of 114 lbs. to 209 lbs.

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

*Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

SEAT SELECTIONS

SEAT SELECTION

Note: Both ADJASBA and SEMIRB seat styles on model M51PR have a 4-post seat mount for adjustable seat-to-floor height, adjustable seat slope, and can be easily retrofitted with a power tilt kit.

- ADJASBA Seat with Adj. Seat Tilt and Back Angle Standard
- SEMIRB Van Seat - Semi Recline with Solid Base \$100.00

ADJUSTABLE ASBA SEAT SELECTIONS

ADJUSTABLE SEAT WIDTH RANGE

- 1620W 16" to 20" Seat Width Range No Charge

SEAT WIDTH

- W16 16" Seat Width No Charge
- W17 17" Seat Width No Charge
- W18 18" Seat Width No Charge
- W19 19" Seat Width No Charge
- W20 20" Seat Width No Charge

ADJUSTABLE SEAT DEPTH RANGE

- 1619D 16" to 19" Seat Depth Range No Charge
- 1920D 19" to 20" Seat Depth Range No Charge

SEAT DEPTH

- D16 16" Seat Depth Standard
- D17 17" Seat Depth No Charge
- D18 18" Seat Depth No Charge
- D19 19" Seat Depth No Charge
- D20 20" Seat Depth \$529.00

UPHOLSTERY COLORS

- U240 Black Nylon Upholstery No Charge
- L8601 Omit Back Upholstery Less \$10.00
- RBP Reinforced Back Upholstery No Charge

SEAT BACK TYPE

- SECT Sectional Rear Uprights No Charge

SEAT BACK HEIGHT

- BH16 16" Back Height No Charge
Not available with TRRO.
- BH17 17" Back Height No Charge
Not available with TRRO.
- BH18 18" Back Height No Charge
- BH19 19" Back Height No Charge
- BH20 20" Back Height No Charge

ARM TYPE (PAIR)

Arm Types are flip-up and detachable.

- 28F Flip Up Desk Length Fixed Ht Conv No Charge
- 29F Flip Up Desk Length Adj. Ht Conv \$230.00 E0973
- 18F Flip Up Full Length Fixed Ht Conv No Charge
- 19F Flip Up Full Length Adj. Ht Conv \$230.00 E0973

MIXED ARMS (EACH)

The below arm selections are flip-up and detachable.

- 19FR Right Arm, Full Length, Adj. Height \$115.00 E0973
- 29FR Right Arm, Desk Length, Adj. Height \$115.00 E0973
- 19FL Left Arm, Full Length, Adj. Height \$115.00 E0973
- 29FL Left Arm, Desk Length, Adj. Height \$115.00 E0973

ARM ACCESSORIES

- Clear
- L8401 Omit Arms Less \$20.00
- L8405 Omit Arm Pads Less \$10.00

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Pronto® M51™ with SureStep® (Rehab Version)

VAN SEAT WITH SOLID SEAT PAN

VAN UPHOLSTERY COLORS

- GV Gray Vinyl Upholstery No Charge

VAN SEAT WIDTH AND DEPTH

- VW16 16" Wide Van Seat No Charge
- VD16 16" Depth No Charge
 - VD17 17" Depth No Charge
 - VD18 18" Depth No Charge
- VW18 18" Wide Van Seat Standard
- VD16 16" Depth Standard
 - VD17 17" Depth No Charge
 - VD18 18" Depth No Charge
- VW20 20" Wide Van Seat \$329.00
- VD18 18" Depth No Charge
 - VD19 19" Depth No Charge
 - VD20 20" Depth \$329.00

VAN ARM TYPE (PAIR)

- VAH Full Length, Adj. Height \$230.00 E0973
- VFH Full Length, Fixed Height Standard
- VAHD Desk Length, Adj. Height \$230.00 E0973
- VFHD Desk Length, Fixed Height No Charge

VAN MIXED ARMS

Arms are flip-up and adjustable width.

- Right Arm Clear
- RVAHF Right Arm, Full Length, Adj. Height \$115.00 E0973
- RVFHF Right Arm, Full Length, Fixed Height No Charge
- RVAHD Right Arm, Desk Length, Adj. Height \$115.00 E0973
- RVFHD Right Arm, Desk Length, Fixed Height No Charge
- Left Arm Clear
- LVAHF Left Arm, Full Length, Adj. Height \$115.00 E0973
- LVFHF Left Arm, Full Length, Fixed Height No Charge
- LVAHD Left Arm, Desk Length, Adj. Height \$115.00 E0973
- LVFHD Left Arm, Desk Length, Fixed Height No Charge

POSITIONING STRAPS

CHEST POSITIONING STRAPS

Clear

- 7311BK Airline Buckle Chest Pos. Strap \$100.00 E0960
- 9515 Push Button Style Chest Pos. Strap \$100.00 E0960
- 7321BK Hook & Loop Chest Pos. Strap \$100.00 E0960

SEAT POSITIONING STRAPS

Clear

- 1515 Push Button Style Seat Pos. Strap Standard
- 1311BK Airline Buckle Seat Pos. Strap \$50.00

DRIVE WHEEL TYPE

DRIVE WHEELS

- DW10535 10 1/2" x 3 1/2" Pneumatic Tires Standard
- Must select DW10535FF drive wheels when selecting TRRO.**
- DW10535FF 10 1/2" x 3 1/2" Flat Free Tires No Charge

COLORS

SHROUD COLORS

- RD Candy Red Pearl No Charge
- BL Deep Blue Pearl No Charge

BATTERY OPTIONS

INSTALLATION OF BATTERIES

- U1INSTAL U1 Battery Install (On Chair)(1 pair) \$232.00 E2365

ELECTRONICS

JOYSTICK MOUNTING HARDWARE

- QL Quad Link Retractable Joystick Mount SPJ+ \$230.00 E1028
- LEFT Left Hand Mounted Joystick No Charge

ACCESSORIES

MISCELLANEOUS ACCESSORIES

- LESTR Lester 8 Amp Battery Charger \$320.00 E2367
- 1831M6 MK6i Programmer \$350.00
- VANO2 Oxygen Holder (Solid Base Van Seat) \$150.00 E2208
- 1209 Tire Pump \$35.00
- SEATO2 Oxygen Holder (ASBA Seat Only) \$150.00 E2208
- Only one option can be used at one time.**
- HFO2-ML6 HomeFill O2 Holder-ML6 Only \$150.00 K0104

FRONT RIGGING OPTIONS

FOOTBOARD

- 90FB Footboard Standard
- L90FB Omit Footboard Less \$40.00

FOOTREST/LEGREST (HEMI)

- PHW93 Swingaway Footrests 13 1/2" - 19" \$165.00
- PH904A Swingaway ELR 16" - 19 3/4" \$260.00 E0990
- PHAL4A Smartleg® Art. Legrests 17 1/4" - 21 1/2" \$225.00 K0053

FOOTPLATE OPTIONS

- 1350 Extra Large Aluminum Footplates \$115.00
- 1450 Aluminum Footplates Deluxe No Charge
- 1650 Aluminum Footplates No Charge
- 1651 Composite Flip-Up Footplates No Charge
- AT5543 Adjustable Angle Flip-Up Footplates \$165.00

FOOTREST/LEGREST ACCESSORIES

Clear

- LPT2 Longer Pivot Slide Tube For 1350 \$83.00
- ALPT Longer Pivot And Slide Tube \$84.00
- ALPT4 Longer Pivot And Slide Tube 4" \$84.00
- A1337 Calf Strap \$27.00 K0038
- IMPACTG Impact Guards \$51.00
- 1600BK Heel Loops with Ankle Straps \$56.00 E0951

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Pronto® M51™ with SureStep® (Rehab Version)

Specifications	ADJASBA	Van Seat w/Solid Base
Dimensions		
Lengths (with Footboard folded)	34"	34"
Lengths (with Footboard down)	38" - 40"	38" - 40"
Widths (excluding Joystick)	24"	24"
Heights	35" - 37"	35" - 37"
Widths	16" - 20"	16" - 20"
Depths	16" - 20"	16", 18", 20"
Back Ht. without Headrest	16" - 20"	19"
Seat to Floor Heights	18"-20" (seat pan to floor)	19" - 21"
Weight		
With Batteries	192 lbs.	198 lbs.
Without Batteries	145 lbs.	151 lbs.
Wheel and Caster Sizes		
Front Casters	6" x 2"	6" x 2"
Drive Wheels	10 1/2" x 3 1/2"	10 1/2" x 3 1/2"
Rear Casters	6" x 2"	6" x 2"
Maximum Speed	4 1/4 mph	4 1/4 mph
Incline Capability	9°	9°
Weight Capacity	300 lbs.	300 lbs.
Turning Radius		
Front (with footboard down)	19 1/2"	19 1/2"
Rear	19 1/2"	19 1/2"
Ground Clearance	2.6"	2.6"
Warranties	See Owner's Manual	See Owner's Manual

Seating and Positioning
 For your seating needs please refer to the following:
 Rehab Seating and Positioning Products
 see Price List Form 07-027

Comments

NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.



Invacare Corporation
 USA
 One Invacare Way
 Elyria, OH USA
 44035-4190
 800-333-6900

www.invacare.com

©2013 Invacare Corporation
 All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.