



Yes, you can.®

Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.

Invacare® Pronto® M61™ with SureStep® and Standard Power Elevated Seat (Rehab Version)

Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



Yes, you can.®

TRANSACTION TYPE:

- Order Quote

Date _____ Account # _____ P.O. # _____
 Company Name _____ Phone # _____
 Address _____ City _____ State _____ Zip Code _____ Country _____
 Contact Name _____ Comments _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

FRAME SELECTION

POWER CHAIR TYPE

- M61R M61 w/SureStep with MK660 ACC controller
and SPJ+ Joystick and Actuator \$5454.00 K0830
 M61SEAT Pronto M61 Seat Only..... \$740.00

ELEVATING SEAT

- ELEV Standard 5" Power Elevating Seat \$1200.00
**ELEV option is required on M61. If elevating seat is not desired,
 order Power Chair model M51.**

SEAT SELECTION

- ASBA Seat w/ Adj. Seat Tilt & Back Angle No Charge
 M5SRB Van Seat Semi Recline with Solid Base \$70.00

ASBA SEAT

SEAT WIDTH

- W16 16" Seat Width No Charge
 W17 17" Seat Width No Charge
 W18 18" Seat Width No Charge
 W19 19" Seat Width No Charge
 W20 20" Seat Width No Charge

SEAT DEPTH

- D16 16" Seat Depth Standard
 D17 17" Seat Depth No Charge
 D18 18" Seat Depth No Charge
 D19 19" Seat Depth No Charge
 D20 20" Seat Depth \$529.00

UPHOLSTERY

- U240 Black Nylon Upholstery No Charge

BACK SELECTIONS

SEAT BACK TYPE

- SECT Sectional Rear Uprights No Charge

SEAT BACK HEIGHT

- BH16 16" Back Height No Charge
 BH17 17" Back Height No Charge
 BH18 18" Back Height No Charge
 BH19 19" Back Height No Charge
 BH20 20" Back Height No Charge

ARM SELECTIONS

ARM TYPE (PAIR)

Arms are flip-up and detachable.

- 28F Flip Up Desk Length Fixed Ht Conv Standard
 29F Flip Up Desk Length Adj. Ht Conv \$230.00 E0973
 18F Flip Up Full Length Fixed Ht Conv No Charge
 19F Flip Up Full Length Adj. Ht Conv \$230.00 E0973

VAN SEAT - WITH SOLID BASE

VAN UPHOLSTERY COLORS

- GV Gray Vinyl Upholstery No Charge

**NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.**

Invacare® Pronto® M61™ with SureStep® and Standard Power Elevated Seat (Rehab Version)

VAN SEAT - WITH SOLID BASE continued

VAN SEAT WIDTH AND DEPTH

- VW16 16" Wide Van Seat No Charge
 - VD16 16" Depth No Charge
 - VD17 17" Depth No Charge
 - VD18 18" Depth No Charge
- VW18 18" Wide Van Seat Standard
 - VD16 16" Depth Standard
 - VD17 17" Depth No Charge
 - VD18 18" Depth No Charge
- VW20 20" Wide Van Seat \$329.00
 - VD18 18" Depth No Charge
 - VD19 19" Depth No Charge
 - VD20 20" Depth \$329.00

VAN ARM TYPE (PAIR)

Arms are flip-up and detachable.

- VAH Desk Length Adj. Height \$230.00 E0973
- VFH Full Length, Fixed Height Standard
- VAHD Desk Length, Adj Height \$230.00 E0973
- VFHD Desk Length, Fixed Height No Charge

POSITIONING STRAPS

SEAT POSITIONING STRAPS

Clear

- 1515 Push Button Style Seat Pos. Strap Standard
- 1311BK Airline Buckle Seat Pos. Strap \$50.00

CHEST POSITIONING STRAPS

Clear

- 7311BK Airline Buckle Chest Pos. Strap \$100.00 E0960
- 9515 Push Button Style Chest Pos. Strap \$100.00 E0960
- 7321BK Hook & Loop Chest Pos. Strap \$100.00 E0960

DRIVE WHEEL TYPE

DRIVE WHEELS

- DW10535 10 1/2" x 3 1/2" Pneumatic Tires Standard
- DW10535FF 10 1/2" x 3 1/2" Flat Free Tires No Charge

COLORS

SHROUD COLORS

- RD Candy Red Pearl No Charge
- BL Deep Blue Pearl No Charge

BATTERY OPTIONS

INSTALLATION OF BATTERIES

- U1INSTAL U1 Battery Install on Chair (1 pair) \$232.00 E2365

FRONT RIGGING OPTIONS

FOOTBOARD

- 90FB Footboard Standard
- L90FB Omit Footboard Less \$40.00

FOOTREST/LEGREST (HEMI)

- PHW93 Swingaway Footrests 13 1/2" - 19" \$165.00
- PH904A Swingaway ELR 16" - 19 3/4" \$260.00 E0990
- PHAL4A Smartleg® Art. Legrests 17 1/4" - 21 1/2" \$225.00 K0053

FOOTPLATE OPTIONS

- 1350 Extra Large Aluminum Footplates \$115.00
- 1450 Aluminum Footplates Deluxe No Charge
- 1650 Aluminum Footplates No Charge
- 1651 Composite Flip Up Footplates No Charge
- AT5543 Adjustable Angle Flip Up Footplates \$165.00

FOOTREST/LEGREST ACCESSORIES

Clear

- LPT2 Longer Pivot Slide Tube For 1350 \$83.00
- ALPT Longer Pivot And Slide Tube \$84.00
- ALPT4 Longer Pivot And Slide Tube 4" \$84.00
- A1337 Calf Strap \$27.00 K0038
- IMPACTG Impact Guards \$51.00
- HLOOP Heel Loops \$45.00 E0951

MISCELLANEOUS ACCESSORIES

JOYSTICK MOUNTING HARDWARE

- QL Quad Link Retractable Joystick Mount SPJ+ \$230.00 E1028
- LEFT Left Hand Mounted Joystick No Charge

MISCELLANEOUS ACCESSORIES

- SEATO2 Oxygen Holder (ASBA Seat Only) \$150.00 E2208
Only one option can be used at a time.
- 1209 Tire Pump \$35.00
- VANO2 Oxygen Holder (Solid Base Van Seat) \$150.00 E2208
- P725 Crutch and Cane Holder \$55.00 E2207
Only one option can be used at a time.
- P770 Walker Holder \$99.00
Only one option can be used at a time.
- 1831M6 MK6i Programmer \$350.00
- LESTR Lester 8 Amp Battery Charger \$320.00 E2367

Seating and Positioning

For your seating needs please refer to the following:
Rehab Seating and Positioning Products
see Price List Form 07-027

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Pronto® M61™ with SureStep® and Standard Power Elevated Seat (Rehab Version)

Specifications	ASBA	Van Seat w/Solid Base
Dimensions		
Lengths (with Footboard folded)	34"	34"
Lengths (with Footboard down)	38" - 40"	38" - 40"
Widths (excluding Joystick)	24"	24"
Heights	35" - 37"	35" - 37"
Widths	17" - 21" (seat pan to floor)	19" - 23"
Depths	16" - 20"	16", 18", 20"
Back Ht w/o Headrest	16" - 20"	19"
Seat to Floor Heights	17" - 19"	17" - 19"
Weight		
With Batteries	192 lbs.	198 lbs.
Without Batteries	145 lbs.	151 lbs.
Wheel and Caster Sizes		
Front Casters	6" x 2"	6" x 2"
Drive Wheels	10.5" x 3 1/2"	10 1/2" x 3 1/2"
Rear Casters	6" x 2"	6" x 2"
Maximum Speed	4 1/4 mph	4 1/4 mph
Incline Capability	9 degrees	9 degrees
Weight Capacity	300 lbs.	300 lbs.
Turning Radius		
Front (w/footboard down)	19 1/2"	19 1/2"
Rear	19 1/2"	19 1/2"
Ground Clearance	2.6"	2.6"
Warranties	See Owner's Manual	See Owner's Manual

Comments

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.



Invacare Corporation

www.invacare.com

USA

One Invacare Way
 Elyria, OH USA
 44035-4190
 800-333-6900

©2013 Invacare Corporation
 All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.