



# Bariatric Bed



Yes, you can.®

- 1000 lb. client weight capacity
- Split-spring assembly for easy delivery and set-up





**Yes, you can.®**

The Invacare® BAR1000 Bariatric Bed is a heavy-duty full-electric bed frame designed for bariatric individuals, capable of supporting up to 1,000 pounds. Using the hand pendant, the client or caregiver is able to change the positioning of the head and foot sections, as well as the bed height. And, setup within the home is easy, even though this is a substantial bed.



BAR1000 with BARMATT1000 mattress.



Hand pendant

**Invacare Corporation**  
www.invacare.com

**USA**

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

**Canada**

570 Matheson Blvd. E.  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

Invacare, the Blue Medallion Design and Yes, you can. are registered trademarks and the Bariatric Solutions design and IVC are trademarks of Invacare Corporation. Specifications are subject to change without notification.

© 2005  
Invacare Corporation

Form No. 05-121

10/05



*"Invacare says 'Yes, you can.' and I couldn't have said it better myself."*

## Features

- 1,000 lb. client weight capacity
- Heavy-duty frame design ensures added strength for client support
- Split-spring design for delivery and setup
- Two locking casters on head section and two steering lock casters on foot section
- Two high-low motors, one on head section and one on foot section
- Phillips screwdriver and 9/16" open-end wrench required for assembly
- Pendant allows for dynamic client positioning
- Quiet, smooth operation
- Emergency crank manually functions to lower any raised section of bed
- UL73 certified
- Use with Invacare bariatric mattress (BARMATT1000) and bed rails (BAR6600)
- Optional Bed Cart (BAR5140) available for easy delivery

## Specifications

BAR1000 Bariatric bed

### Dimensions

Bed height *	min. 13.5" - max. 28"
Sleep surface length	80"
Sleep surface width	48"
Overall length	88"
Overall width	48"
Shipping dimensions	41" x 49" x 58"

### Limited warranty

All welds	lifetime
Frame/structure	5 years
Mechanical/electrical	1 year

### Weights

BAR1000	433 lbs.
Head section	216 lbs.
Foot section	216 lbs.
Shipping weight	510 lbs.

\*Measurement from floor to bed deck.

## Bed Packages

BAR1000	1,000 lb. Bariatric Bed includes 1 set of BAR6600 1/2 rails
BARPKG1000	BAR1000, BARMATT1000
BARPKG1002	BAR1000, BARMATT1000, BAR5140



BAR5140 Bed Cart shown w/BAR1000 Bed