

INVACARE® ROLLATORS



Invacare® Rollators will meet all your needs with a focus on style and affordability. These rollators feature easy-to-operate ergonomic hand brakes and easy folding for storage and transport. Invacare Rollators offer great maneuverability and stability, and are height adjustable for a custom fit. Baskets are standard on all rollators. Trays are standard on models P429/2 and 65550.

Invacare® FR300™ Rollator Features

Model # 67100
Tray and Basket included
Cross Brace
Streamlined design
Curved backrest
Curb climber

65650 / 65650R Rollator Features

Model # 65650 Blue, 65650R Red
Easy-to-operate, ergonomic hand brakes lock the rear wheels for security
Flip-up padded seat with zipper compartment for storage of valuables
Removable foam padded curved backrest enhances comfort and can be easily removed
Lightweight aluminum frame
Rollator folds easily with basket attached
Universal basket fits underneath seat or hooks easily on front seat brace

HCPSC Code: E0143 & E0156

65851R / 65851B Adult Rollator Features

Model # 65851R Red, 65851B Blue
Easy-to-operate, ergonomic hand brakes lock the rear wheels for security
Flip-up seat
Lightweight aluminum frame
Basket hooks easily on front seat brace
Removable foam-padded straight backrest enhances comfort and can be easily removed

65350GR Junior Rollator Features

Model # 65350GR
Easy-to-operate, ergonomic hand brakes lock the rear wheels for security
Flip-up seat
Lightweight aluminum frame
Basket hooks easily on front seat brace
Removable foam-padded straight backrest enhances comfort and can be easily removed

HCPSC Code: E0143 & E0156

NOTE: Rollators are not intended to be used as wheelchairs.

All references to HCPSC codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further details.

P429/2 Rollator Features

Model # P429/2

Easy-to-operate, ergonomic hand brakes
Rollator folds easily with tote bag attached
Indoor/outdoor 7½" x 2" solid rubber wheels
Basket standard, attaches easily to the frame
Tray standard, attaches easily to the basket
Tote bag standard, attaches easily to the frame

65550 Rollator Features

Model # 65550

Easy-to-operate, ergonomic hand brakes lock the rear wheels for security
Sturdy seat
Rollator folds easily with basket attached
Basket and tray standard

**HCPCS Code: E0143 Model # P429/2
E0143 & E0156 Model # 65550**

P429/2: Three-Wheel Rollator



65550: Four-Wheel Rollator

66550 Bariatric Rollator Features

Model # 66550

500 lb. weight capacity
Wide, deep frame
Flip-up padded seat
Removable padded backrest
Folds for ease of storage and transport
Large 8" wheels for improved mobility over rougher terrain

Height adjustable loop handbrakes with foam handgrips
Durable steel construction
Pouch included

HCPCS Code: E0149

1. Bariatric size



66550: Bariatric Rollator

NOTE: Rollators are not intended to be used as wheelchairs.

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further details.



Model #	Overall Dimensions	Dimension Inside Hand Grips	Seat Dimensions	Seat-to-Floor-Height	Caster Dimensions	Product Weight	Weight Capacity
67100	24.4" W 26.3" D	16"	15.5" W 8" D	23.6"	8"	17.5 lb.	285 lb.
65650, 65650R	24.25" W 28.28" D	18"	13.75" W 16.375" D	23.25"	7.5"	19 lb.	300 lb.
65851R, 65851B	28.5" W 27" D	16"	14" W 13" D	22"	6"	13.4 lb.	300 lb.
65350GR	26" W 24.5" D	13.5"	14" W 13" D	19.5"	6"	13 lb.	250 lb.
P429/2	26" W 24" D	14"	14" W 13" D	19.5"	6"	17 lb.	294 lb.
65550	23.25" W 27" D	16.5"	15.25" W 12.5" D	24"	7.5"	25 lb.	285 lb.
66550	27" W 29.5" D	22"	18" W 14" D	22"	7.5"	26 lb.	500 lb.



	Patient Height Range	Units Per Carton	Limited Warranty
67100	N/A	1	5 Years Brakes 1 Year Components*
65650, 65650R	5'5"-6'3"	1	5 Years Brakes 1 Year Components*
65851R, 65851B	5'6"-6'3"	1	5 Years Brakes 1 Year Components*
65350GR	5'-5'6"	1	5 Years Brakes 1 Year Components*
P429/2	5'3"-6'	1	5 Years Brakes 1 Year Components*
65550	5' 3"-6'2"	1	5 Years Brakes 1 Year Components*
66550	N/A	1	5 Years Brakes 1 Year Components*

*Excluding tires and wheels
NOTE: Rollators are not intended to be used as wheelchairs

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2012
Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and
®. All trademarks are
owned by or licensed
to Invacare Corporation
unless otherwise noted.
Specifications are subject
to change without
notification.
Form No. 05-139
Rev. 03/12
111259

