



Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.

Please consider the environment before printing this page with your order form.
For more information on features and benefits of electronic order forms, please visit www.invacare.com/eForms

Invacare® Formula™ CG Tilt Only for Storm Series® Bases

Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



TRANSACTION TYPE:

Order Quote

Date _____ Account # _____ P.O. # _____

Company Name _____ Phone # _____

Address _____ City _____ State _____ Zip Code _____ Country _____

Contact Name _____ Comments _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

POWERED SEATING BASE

- 3GAR-CG Arrow for Formula CG with MK6i \$9,775.00 K0884
The Arrow Base comes Standard with GB TrueTrack Motors for up to a top speed of 7 MPH.
- 3GRX-CG Ranger X for Formula CG with MK6i \$8,000.00 K0862
The Ranger X Base comes Standard with Heavy-Duty 4-Pole Motors, Monroe® Spring Suspension, and up to a top speed of 5 MPH.
- 3GTQ-CG Torque SP for Formula CG with MK6i \$6,850.00 K0877
The Torque SP Base comes Standard with High Performance 4-Pole Motors, and up to a top speed of 6 1/2 MPH.
- 3GTQ3-CG Torque 3 for Formula CG with MK6i \$5,970.00 K0856
The Torque 3 Base comes Standard with 4-Pole Motors and up to a top speed of 5.8 MPH.

Group 3 - K0862 Basic Package	
3GRX-CG Ranger X for Formula CG with MK6i	
400 lbs. Weight Capacity	
1431-3 14" x 3" Tires with Foam Filled Inserts	No Charge
Heavy Duty 4 Pole Motors	Standard
1228-3 8" x 2 1/4" Semi Pneumatic Casters	Standard
Flip Back Adjustable Height Arm	Standard
MK690ACC Controller with Actuator Control	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard
Seat Option	
CGT Tilt (E1002)	\$4,590.00

Group 4 - K0884 Basic Package	
3GAR-CG Arrow for Formula CG with MK6i	
300 lbs. Weight Capacity	
1431-3 14" x 3" Tires with Foam Filled Inserts	No Charge
TrueTrack GB Motors	Standard
1228-3 8" x 2 1/4" Semi Pneumatic Casters	Standard
Flip Back Adjustable Height Arm	Standard
MK6TT TrueTrack Controller	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard
Seat Option	
CGT Tilt (E1002)	\$4,590.00

Group 4 - K0877 Basic Package	
3GTQ-CG Torque SP for Formula CG with MK6i	
300 lbs. Weight Capacity	
1431-3 14" x 3" Tires with Foam Filled Inserts	No Charge
Performance 4 Pole Motors	Standard
1228-3 8" x 2 1/4" Semi Pneumatic Casters	Standard
Flip Back Adjustable Height Arm	Standard
MK690ACC Controller with Actuator Control	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard
Seat Option	
CGT Tilt (E1002)	\$4,590.00

Group 4 - K0886 Basic Package	
3GAR-CG Arrow for Formula CG with MK6i	
400 lbs. Weight Capacity	
1431-3 14" x 3" Tires with Foam Filled Inserts	No Charge
TrueTrack GB Motors	Standard
1228-3 8" x 2 1/4" Semi Pneumatic Casters	Standard
Flip Back Adjustable Height Arm	Standard
MK6TT TrueTrack Controller	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard
Seat Option	
CGT Tilt (E1002)	\$4,590.00

Group 3 - K0856 Basic Package	
3GTQ3-CG Torque 3 for Formula CG with MK6i	
300 lbs. Weight Capacity	
1431-3 14" x 3" Tires with Foam Filled Inserts	No Charge
4 Pole Motors	Standard
1228-3 8" x 2 1/4" Semi Pneumatic Casters	Standard
Flip Back Adjustable Height Arm	Standard
MK690ACC Controller with Actuator Control	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard
Seat Option	
CGT Tilt (E1002)	\$4,590.00

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.**

Invacare® Formula™ CG Tilt Only for Storm Series® Bases

OMIT SEAT FRAME

- OSF Omit Seat Frame for Base Only..... Less \$600.00

TRANSPORT TIE DOWN

- TRBKTS Wheelchair Transport Brackets \$250.00
Only available with MED or TALL seat-to-floor height. Not available with vent tray.

TRANSPORT TIE-DOWN*

TRBKTS includes four factory installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport. As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

* Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

USER WEIGHT LIMIT

Clear

Arrow

- U300 User Weight <300 lbs. No Charge
With 300 lbs. weight capacity, maximum range of tilt is 55 degrees.
- U400 User Weight <400 lbs. No Charge
With 400 lbs. weight capacity, maximum range of tilt is 45 degrees.

Ranger X

- U400 User Weight <400 lbs. No Charge
With 400 lbs. weight capacity, maximum range of tilt is 45 degrees. Comes standard with Heavy Duty Monroe® Spring Suspension.

Torque SP / Torque 3

- U300 User Weight <300 lbs. No Charge
With 300 lbs. weight capacity, maximum range of tilt is 55 degrees.

DRIVE TYPE OPTION

Clear

- TRUTK TrueTrack Upgrade \$1,899.00
TrueTrack includes Gearless Brushless GB Motors and Monroe® Spring Shock Suspension. Must upgrade to MK6TT Controller. Must upgrade to an Expandable Controller System. Not compatible with SPJ+ joystick. Available option on Ranger X, Torque SP and Torque 3. Standard on Arrow. When selecting TrueTrack motors for Torque 3 the top speed of the chair increases to 7 MPH and therefore no longer performs to the Group 3 specifications.
- PERF Performance 4 Pole Motor No Charge
PERF Standard on Torque SP. Top speed of up to 6 1/2 MPH. Not available on Torque 3

VENT TRAYS

Vent Tray on chair reduces chair weight capacity by approximately 50 lbs.

- VENT Articulating Vent Tray \$1,750.00 E1030

FRAME TYPE

- SFS Short Base No Charge
Not Available on seat depths of 1922D. Must pick Long Base.
- LFS Long Base No Charge

CASTER OPTIONS

- 1259 6" x 2" Semi-Pneumatic Casters No Charge
- 1228-3 8" x 2 1/4" Semi-Pneu Casters Lt Grey Standard
- 1251-3 8" x 2" Pneumatic Casters No Charge
- 1220-3 9" x 2 3/4" Pneumatic Casters No Charge
Not available with short base. Not available with TRBKTS.
- 1465-3 9" x 2 3/4" Foam Filled Casters No Charge
Not available with short base. Not available with TRBKTS.

CASTER SUSPENSION

Clear

The below selections are not available with 9" casters.

- 9213-6 Shock Fork For 6" Casters \$119.00 E1016
- 9213-8 Shock Fork For 8" Casters \$119.00 E1016
Not available with short base. Not available with TRBKTS.

WHEEL OPTIONS

- 1430-3 14" x 3" Pneumatic Tires Standard
- 1431-3 14" x 3" Tires with Foam Filled Inserts No Charge
- 1440 14" x 4" Pneumatic Tires \$200.00
14" x 4" wheels increase the overall width of the base by 2 1/2", does not include Storm Series design wheel or fender, not available with TRUTK TrueTrack Package.
- 1441 14" x 4" Tires with Foam Filled Inserts \$320.00
14" x 4" wheels increase the overall width of the base by 2 1/2", does not include Storm Series design wheel or fender, not available with TRUTK TrueTrack Package.

FREE WHEEL HUB

Will add 3/4" overall width. Not available with TrueTrack, does not include Storm Series design wheel or fender.

- FHUB Free Wheel Hub \$125.00

FRAME FINISH

CONTEMPORARY

- 24P Wet Black No Charge
- 60P Silver Metallic No Charge
- 105P Midnight Blue No Charge
- 115P Black Prism No Charge

TOO HOT

- 30P Sunny Yellow No Charge
- 61P Electric Red No Charge
- 119P Bubblegum Pink No Charge
- 120P Grape Madness No Charge

WAY COOL

- 62P Electric Blue No Charge
- 127P Grasshopper Green No Charge
- 116P Sunkissed Orange No Charge
- 117P Tangerine Red No Charge
- 29P White Snow No Charge

Discontinued colors 71P Silver Vein, 104P Electric Teal, 121P Lolly Pop Blue, 122P Cosmic Blue, and 125P Emerald Green can be order through specials number C-070109031327

**NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.**

Invacare® Formula™ CG Tilt Only for Storm Series® Bases

SEAT

SYSTEM TYPE

- Clear
- CGT CG Tilt Only System \$4,590.00 E1002
- Clear
- CNVACT Conventional Actuator..... Standard
- SMTACT Smart Actuator \$500.00
Must select MK690 Controller when selecting SMTACT (Smart Actuator), does not apply when selecting MK6TT. SMTACT is Standard when selecting 3GAR-CG.

SEAT-TO-FLOOR HEIGHT

- LOW 17 3/4" at 0° No Charge
Not available with Transport Tie Down option.
- MED 18 3/4" at 0° No Charge
- TALL 19 3/4" at 0° No Charge

SEAT SELECTION

- CTB Contoura™ Seating \$625.00 E2620 E2621
- OMCU Omit Contoura Seat Cushion Less \$20.00
- CTBMSH Contoura™ seating with mesh upholstery \$635.00 E2620 E2621
- OMCU Omit Contoura Seat Cushion Less \$20.00
- BPO Back Post Only, Round Canes Only No Charge
Push handles included.
- CURV Curved No Charge
Push handles included.
- CNB Conventional Black No Charge
Push handles included.
For 12" - 15" seat widths and depths must select Conventional Black or BPO. User weight is under 150 lbs. These sizes are not available with Contoura Seating Systems. Conventional back not available on widths of 16" - 20" and 20" - 24" or depths of 16" - 19" or 19" - 22".

ADJUSTABLE SEAT WIDTH RANGE

- Clear
- 1216W 12" to 16" Seat Width (Adjustable) No Charge
1216W and 1218D not available with 3GRX-CG. For 16" seat width on the 3GRX-CG, must select 16" to 20" width.
- 1620W 16" to 20" Seat Width (Adjustable) No Charge
- 2024W 20" to 24" Seat Width (Adjustable) No Charge

SEAT WIDTH SETTING

- Clear
- W12 12" wide No Charge
- W13 13" wide No Charge
- W14 14" wide No Charge
- W15 15" wide No Charge
- W16 16" wide No Charge
- W17 17" wide No Charge
- W18 18" wide No Charge
- W19 19" wide No Charge
- W20 20" wide No Charge
- W21 21" wide No Charge
- W22 22" wide No Charge
- W23 23" wide No Charge
Available only with 400 lbs. package.
- W24 24" wide No Charge
Available only with 400 lbs. package.

ADJUSTABLE SEAT DEPTH RANGE

- Clear
- 1218D 12" to 18" Deep (Adjustable) No Charge
1216W and 1218D not available with 3GRX-CG. For 16" seat width on the 3GRX-CG, must select 16" - 20" width.
- 1619D 16" to 19" Seat Depth Range No Charge
- 1922D 19" to 22" Seat Depth Range No Charge

ADJUSTABLE SEAT DEPTH SETTING

- Clear
- D12 Seat Depth 12" (Short Base) No Charge
- D13 Seat Depth 13" (Short Base) No Charge
- D14 Seat Depth 14" (Short Base) No Charge
- D15 Seat Depth 15" (Short Base) No Charge
- D16 Seat Depth 16" No Charge
- D17 Seat Depth 17" No Charge
- D18 Seat Depth 18" No Charge
- D19 Seat Depth 19" No Charge
- D20 Seat Depth 20" No Charge
- D21 Seat Depth 21" No Charge
- D22 Seat Depth 22" No Charge

SEAT TILT

- Indicates static tilt of seat. Impacts seat-to-floor height. Every 5° raises the seat-to-floor height.**
- 0SA 0° Standard
- 5SA 5° No Charge

BACK HEIGHT

- Clear
- TSBH16 Back Height 16" No Charge
- TSBH17 Back Height 17" No Charge
- TSBH18 Back Height 18" No Charge
- TSBH19 Back Height 19" No Charge
- TSBH20 Back Height 20" No Charge
- TSBH21 Back Height 21" No Charge
- TSBH22 Back Height 22" No Charge
- TSBH23 Back Height 23" No Charge
- TSBH24 Back Height 24" No Charge
- TSBH25 Back Height 25" No Charge
- TSBH26 Back Height 26" No Charge

CONTOURA BACK HEIGHT

- Clear
- TSBH22 Back Height 22" No Charge
- TSBH25 Back Height 25" No Charge
- TSBH28 Back Height 28" No Charge

**NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.**

SEAT Continued

ADJUSTABLE ARM TYPE (PAIR) Clear

- 19F Full Length Adj. Hgt No Charge
- 29F Desk Length Adj. Hgt No Charge

MIXED ADJUSTABLE ARM TYPE (EACH)

Right Arm Clear

- 19FR Full Length Adj. Hgt. Right No Charge
- 29FR Desk Length Adj. Hgt. Right No Charge

Left Arm Clear

- 19FL Full Length Adj. Hgt. Left No Charge
- 29FL Desk Length Adj. Hgt. Left No Charge

CANTILEVER ARM TYPE (PAIR) Clear

Cantilever Arms only available with 300 lbs. rated systems. Not available with BPO, Otto Bock, or 12" - 15" seat width. Stealth Lateral part number's TWB-CPR and TWB-CPRL are the only lateral support compatible with the current Cantilever Arm design. Must choose WF armpad type with the Stealth Laterals. See chart on last page for arm height adjustment range.

- CLVF Full Length Arm \$350.00
- CLVD Desk Length Arm \$350.00

CANTILEVER MIXED ARM TYPE (EACH)

Right Arm Clear

- CLVFR Full Length Arm Right \$175.00
- CLVDR Desk Length Arm Right \$175.00

Left Arm Clear

- CLVFL Full Length Arm Left \$175.00
- CLVDL Desk Length Arm Left \$175.00

ARMPAD OPTIONS (CANTILEVER ARMS) Clear

- WF Waterfall Armpads Standard
- RCL Reclining Armpads No Charge

SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Positioning Strap \$50.00
Not available with 12" - 15" seat widths.
- 1515 Push Button Style Seat Positioning Strap Standard

CHEST POSITIONING STRAPS Clear

- 7311BK Airline Buckle Chest Positioning Strap \$100.00 E0960
Not available with 12" - 15" seat widths.
- 9515 Push Button Style Chest Positioning Strap \$100.00 E0960
- 7321BK Hook and Loop Chest Positioning Strap \$100.00 E0960

POSITIONING ACCESSORIES

HEADREST TYPE (OPTIONAL) Clear

Hardware is packaged with equipment and cannot be ordered separately.

- TAR127700 Curved Headrest \$210.00 E0955
TAR27700HW Removable & Adjustable Hardware \$210.00 E1028
- TAR128000 Two-Step Headrest \$210.00 E0955
TAR128000HW Removable & Adjustable Hardware \$210.00 E1028
- T0022400 Large Headrest \$210.00 E0955
T0022400HW Removable & Adjustable Hardware \$210.00 E1028
- T0022600 Small Headrest \$210.00 E0955
T0022600HW Removable & Adjustable Hardware \$210.00 E1028

STEALTH HEADRESTS Clear

Hardware is packaged with equipment and cannot be ordered separately.

- CP175 14" Multi Axis Headrest \$225.00 E0955
CP175HW Removable & Adjustable Hardware \$225.00 E1028
- CP275 10" Multi Axis Headrest \$225.00 E0955
CP275HW Removable & Adjustable Hardware \$225.00 E1028
- CP180 14" TWB Headrest \$225.00 E0955
CP180HW Removable & Adjustable Hardware \$225.00 E1028
- CP280 10" TWB Headrest \$225.00 E0955
CP280HW Removable & Adjustable Hardware \$225.00 E1028

OTTO BOCK™ ARM PADS (PAIR) Clear

Not compatible with quadlink retractable joystick mount. If selecting Otto Bock hand pad, a large, medium, or small forearm pad must be selected.

- OB20400 One Piece Channel Armpads \$230.00 E2209
- OB20500 Large Forearm Pads \$235.00
- OB20600 Medium Forearm Pads \$235.00
- OB20700 Small Forearm Pads \$104.00

OTTO BOCK HAND PADS (PAIR) Clear

Not compatible with quadlink retractable joystick mount. If selecting Otto Bock hand pad, a large, medium, or small forearm pad must be selected.

- OB23000 Large Flat Hand Pads \$100.00
- OB23300 Medium Flat Hand Pads \$100.00

ADDUCTORS Clear

Hardware is packaged with equipment and cannot be ordered separately.

- ADD Adductors Small \$200.00 E0956
ADDHW Hardware for ADD \$400.00 E1028
- ADDL Adductors Large \$200.00 E0956
ADDHWL Hardware for ADDL \$400.00 E1028

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.**

POSITIONING ACCESSORIES Continued

LATERALS (Pair) Clear

Conventional Back / Curved Back / BPO

Hardware is packaged with equipment and cannot be ordered separately.

- TCSLM Clamp-On Swingaway Laterals Medium \$210.00 E0956
TCSLMHW Removable and Adjustable Hardware \$210.00 E1028
- TCSLL Clamp-On Swingaway Laterals Large \$210.00 E0956
TCSLLHW Removable and Adjustable Hardware \$210.00 E1028

Contoura Back Clear

For use with Contoura seating only. Hardware is packaged with equipment and cannot be ordered separately.

- CSLM Swingaway Contoura Lat. Med \$210.00 E0956
CSLMHW Removable and Adjustable Hardware \$210.00 E1028
- CSLL Swingaway Contoura Lat. Lg \$210.00 E0956
CSLLHW Removable and Adjustable Hardware \$210.00 E1028

STEALTH LATERALS (Pair) Clear

Conventional Back / Curved Back / BPO

When using Stealth laterals the arm will not swing past the lateral support. Will need to remove arm for transfer. Hardware is packaged with equipment and cannot be ordered separately.

- TWB-RPR Adj. Swing Round Lat. Med \$320.00 E0956
RPRHW Removable and Adjustable Hardware \$320.00 E1028
- TWB-RPRL Adj. Swing Round Lat. Lg \$320.00 E0956
RPRLHW Removable and Adjustable Hardware \$320.00 E1028

Contoura Back Clear

When using Stealth laterals the arm will not swing past the lateral support. Will need to remove arm for transfer. For use with Contoura seating only. Hardware is packaged with equipment and cannot be ordered separately.

- TWB-CPR Adj. Swing Contoura Lat. Med \$320.00 E0956
CPRHW Removable & Adjustable Hardware \$320.00 E1028
- TWB-CPRL Adj. Swing Contoura Lat. Lg \$320.00 E0956
CPRLHW Removable & Adjustable Hardware \$320.00 E1028

ELECTRONICS

MK6i CONTROLLERS Clear

Ranger X / Torque SP / Torque 3

- MK690ACC 90 AMP Controller with Actuator Contol Standard
- MK690 90AMP Controller No Charge
For use with Smart Actuators only.
- MK6TT TrueTrack Controller No Charge
- GTRAC G-Trac Module \$995.00 K0108
G-Trac provides Improved Tracking Capabilities and Enhanced Traction. Does not require Special Motors. Not available with SPJ+ joysticks or MK6TT.

Arrow

- MK6TT TrueTrack Controller Standard

STANDARD DRIVER CONTROL

Torque SP / Torque 3

- SPJAP SPJ+ with Power Seating Switch Standard

ACTUATOR CONTROL THROUGH SPJ+

- SPJAPM SPJ+ with Joystick Actuator Control \$1,300.00 E2310
Optional for Torque SP and Torque 3. Allows actuator operation through Driver Control. Includes port for attendant switch.

EXPANDABLE CONTROLLER SYSTEM Clear

- PWH Harness req. for Expandable Electronics \$350.00 E2313
Required to interface joystick or MK6i display with expandable controllers.
- CMPJM6U MPJ+ Multiple Drive Joystick with Color Display .. \$750.00 E2377
- PSFM6 Personal Joystick, Inductive on Front \$1,130.00 E2377
TOD-SS: Toggle On/Off Drive Select, Speed Potentiometer Standard.
- PSRM6 Personal Joystick, Inductive on Rear \$1,130.00 E2377
TOD-SS: Toggle On/Off Drive Select, Speed Potentiometer Standard. Clear
- PBOD Push Button On/Off with PB Drive Select No Charge
PSRM6 Only.
- PBSS Push Button On/Off with Speed Select No Charge
Drive select through mode switch, or mode switch and joystick.
PSRM6 Only. Clear
- OR On/Off Mounted Right No Charge
PSRM6 Only.
- OL On/Off Mounted Left No Charge
PSRM6 Only.
- ASLRDYM6 MK6 Display Only (No Driver Control) \$750.00 E2377
Display only, No Driver Control. Allows choice from order form.
- MCD Mini Color Display \$750.00 E2377
No port for remote On/Off switch. Stereo port only.

FOR YOUR ASL DRIVE SYSTEM NEEDS
PLEASE REFER TO THE
ASL SINGLE INVOICE PRICE LIST FORM # 02-097.

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS Clear

- 1500M6 Prop RIM Head Control with MK6i Display .. \$2,875.00 E2327 Clear
- SMHD Small Headrest No Charge
- LGHD Large Headrest No Charge
- 1558MM6 Compact Joystick with Display \$1,130.00 E2373

ALTERNATIVE DIGITAL CONTROL SYSTEMS Clear

Must also use either a Multiple Drive Joystick or MK6i Display.

- SNPM6 Sip-N-Puff Digital Interface Box \$1,500.00 E2325
- 5018M6 Wafer Board and Digital Interface Box \$1,550.00 E2322
- 5020M6 Mini Tash Joystick and Digital Interface Box \$1,750.00 E2321

SIP-N-PUFF KIT

- PKG32666 Therafin® Sip-N-Puff Breath Tube Kit \$380.00 E2326

ADDITIONAL CONTROL CHOICES

- 1558M6 Compact Joystick - Less Std Mount \$770.00
- 5020 Mini Tash Joystick \$405.00
Requires SNPM6 Sip-N-Puff Interface Box.
- PACM6 Prop Attendant Control with quick release mount \$850.00
- 1552M Digital Attendant Control \$850.00
Requires SNPM6 Sip-N-Puff Interface Box.

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.**

ELECTRONICS Continued

JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retract. Joystick Mnt CMPJ+ \$279.00 E1028
Includes ARM20 Standard.
- QLM Quad Link Retract. Joystick Mount SPJ+ \$279.00 E1028
Includes ARM20 Standard.
- ARM250 Stealth Height Adjustable Joystick Mount \$210.00 E1028
- LEFT Left Handed Mounted Joystick No Charge
- RIMHW RIM Hardware For 1500 \$219.00 E1028
- PKG32669 Complete Bib Assembly (Tash/1558)..... \$265.00
- GATMPJ6 Gatlin Midline for CMPJ \$600.00 E1028
- GAT1812 Gatlin Midline for 1812 \$600.00 E1028
- RCM Rear Cane Mount for Compact Joystick \$395.00

ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer w/Pro Memory Card \$475.00
- 1813SD Pro Memory Card (USBReady) \$90.00
- 1560 "T" Handle Flexible Joystick Ext. \$77.00 E2323
- 1561 Straight Handle Flexible Joystick Ext. \$77.00 E2323
- 1826 Chin Cup \$48.00 E2324
- A24VPS 24V MK6 Auxiliary Power Source \$200.00
- AUX12M6 Auxiliary Module for outputs 1 & 2 \$795.00
- AUX34M6 Auxiliary Module for outputs 3 & 4 \$795.00
- MPY Monoport "Y" Cable for CMPJ+ or Display \$15.00
- 2PB Dual Push Button \$147.00
Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options.
- 2SW Dual Toggle \$172.00
Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options.
- PBMT Push Button Mounting Kit \$252.00
- RFM Wireless (RF) Mouse Only Module \$750.00
RFM includes Transmitter, USB Radio Frequency Mouse Dongle (Receiver), and Stereo Switch with hook and loop mount for mouse clicks. Not compatible with SPJ+.

POWERED SEATING ACCESSORIES

- Clear
- SANODE One Actuator Operation thru Driver Cntrl..... \$1,300.00 E2310
Required for Tilt Operation through Driver Control - Single Actuator Systems Only. Not available with SPJ+ Joystick. Does not come with Egg switch.
- 4WSB Four Way Switch Box No Charge
 Clear
- FWT 4 Way Toggle No Charge
- QPB 4 Push Buttons No Charge
- S4WSB Multiple Actuator Interface Box \$2,600.00 E2311
 Clear
- FWT 4 Way Toggle No Charge
- QPB 4 Push Buttons No Charge



Systems

For systems with two or more actuators, NOT operating through driver control, choose 4WSB with either 4-way toggle or 4-push button.

For systems with 2 or more actuators OPERATING THROUGH DRIVER CONTROL, choose S4WSB.
Uses 4-way toggle or 4-push button as attendant control.

ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Reset-Black \$156.00
- ASL304 Light Wobble Switch \$180.00
- ASL202J Fiber Optic Switch \$732.00
Not available with SPJ+ Joystick. Must add MK6 Auxiliary Power Source.
- ASL208J Adjustable Proximity Switch \$660.00
Not available with SPJ+ Joystick. Must add MK6 Auxiliary Power Source.
- ASL504M6 Remote Emergency Stop Switch \$1,185.00
Not available with SPJ+ Joystick. Must add MK6 Auxiliary Power Source.
- ASL611 Adjustable Switch Mounting Bracket \$252.00

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Formula™ CG Tilt Only for Storm Series® Bases

FRONT RIGGINGS

HANGER BRACKET TYPE

- Clear
- ASB Center Pivot Style Bracket No Charge
- PSB Pin Style Bracket No Charge

CENTER PIVOT STYLE (PAIR)

- Clear
- 60HD 60° Swingaway Footrests 13 1/2"-17 1/2"..... No Charge
- 70NHD 70° Swingaway Footrests 13 1/4"-17" No Charge
- 70STAPER 70° Tapered Footrests 13"-17" No Charge
Not available with 19" or 20" wide and 9" casters. Standard with impact guards on the footrests.
- AT5544 Elevating Legrests \$255.00 E0990
A manual elevating legrest.

PIN STYLE (PAIR)

- Clear
- PW93 Swingaway Footrests 13 1/2" - 19"..... No Charge
- P904A Swingaway Elev. Legrests 16" - 19 3/4" \$255.00 E0990
Not available with 9" casters. Standard with impact guards on the footrests.
- PAL4A Smartleg Art. Legrests 17 1/4" - 21 1/2" \$320.00 K0053
Not available with 9" casters. Standard with impact guards on the footrests.

POWER ELEVATING LEGRESTS - ARTICULATING

- ELRPW Power ELR with Articulation (pair) \$1,895.00 E1010
Not available with SPJ+ Joystick. Must pick 4WSB or S4WSB to operate power legs. Independent powered elevating articulating legrests with adjustable angle footplates. Extension range is 13"-20". Sold as a pair.

FOOTPLATE OPTIONS (PAIR)

- Clear
- 1651 Flip-Up Composite Footplates No Charge
- 1350 Extra Large Aluminum Footplates No Charge
- AT5543 Adjustable Angle Flip-Up Footplates \$165.00 K0040

FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap \$27.00 E0038
- IMPACTG Impact Guards \$51.00
- 1600BK Heel Loops with Ankle Straps \$51.00 E0951
- LPT2 Longer Pivot Slide Tube For 1350 \$83.00
- ALPT Longer Pivot And Slide Tube \$84.00
- ALPT4 Longer Pivot And Slide Tube 4" \$84.00

MISCELLANEOUS

BATTERY BOXES

- 22NF 22NF Battery Boxes No Charge
- GP24 Group 24 Battery Boxes No Charge
Option on Torque SP and Torque 3. Must select TRBKTS.

BATTERIES

- Clear
- Batteries will be installed by Invacare Corporation at no additional charge. Chair will be shipped complete with batteries.**
- 22NFBATTERY 22NF Battery (On Chair)..... \$300.00 E2361
Not available with Arrow, Ranger X or True Track systems or chairs with TRBKTS
- 24BATTERY 24 Gel Battery (On Chair) \$410.00 E2363

BATTERY CHARGER OPTION

- Clear
- H0CHARGER 110 Volt Battery Charger No Charge
- LESTR Lester 8 AMP Battery Charger \$320.00 E2367
- L8900 Omit Battery Charger Less \$80.00

MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks \$45.00
- 9205-3 Suspension Shocks \$400.00 E1016
Standard with Ranger, Arrow, or TrueTrack Systems.
- 9205H-3 Suspension Shocks (Heavy Duty) \$400.00 E1018
Not available with 19" or 20" wide and 9" casters. Standard with impact guards on the footrests.
- FDR Fender \$320.00
Not available with 19" or 20" wide and 9" casters. Standard with impact guards on the footrests.
- EXTP Extended Anti-Tippers \$250.00

Cantilever Arm Height Adjustment			
Back Type	Height	Adjustment	Adjustment with laterals
CURV	20" or greater	11" to 16"	11" to 14"
CURV	19"	11" to 15"	11" to 13"
CURV	18"	11" to 14"	11" to 12"
CURV	17"	11" to 13"	no adjustment
CURV	16"	11" to 12"	not available
Contoura	25" to 28"	12 1/2" to 14 1/2"	no adjustment

Contoura Seat Sizes							
	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X
W24	*	*	X	X	X	X	X
*Not available with Contoura Seat Cushion							
Back Heights	22", 25" or 28"						

Curved	
Seat Widths	16" - 24"
Seat Depths	16" - 22"
Back Heights	16" - 26"

Conventional	
Seat Widths	12" - 15"
Seat Depths	12" - 18" (Conventional)
Back Heights	16" - 20" for 12" - 15" (Seat Widths Only)

Back Height Chart	
Back Type	Back Height
CTB - Contoura	22", 25" and 28"
BPO - Back Posts Only	20"-24"
CURV - Curved	16"-26"
CNB-Conventional Back	16"-20" on seat widths of 12"-15"
Note: Back Heights of 21"-26" not available with seat widths of 12"-15"	

**NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.**

