



Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.

Please consider the environment before printing this page with your order form.
For more information on features and benefits of electronic order forms, please visit www.invacare.com/eForms

Invacare® Formula™ CG Recline Only for Storm Series® Bases

Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



TRANSACTION TYPE:

Order Quote

Date _____ Account # _____ P.O. # _____

Company Name _____ Phone # _____

Address _____ City _____ State _____ Zip Code _____ Country _____

Contact Name _____ Comments _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

POWERED SEATING BASE

- 3GAR-CG Arrow for Formula CG with MK6i \$9,775.00 K0884
The Arrow Base comes Standard with GB TrueTrack Motors for up to a top speed of 7 MPH.
- 3GRX-CG Ranger X for Formula CG with MK6i \$8,000.00 K0862
The Ranger X Base comes Standard with Heavy-Duty 4-Pole Motors, Monroe® Spring Suspension, and up to a top speed of 5 MPH.
- 3GTQ-CG Torque SP for Formula CG with MK6i \$6,850.00 K0877
The Torque SP Base comes Standard with High Performance 4-Pole Motors, and up to a top speed of 6 1/2 MPH.
- 3GTQ3-CG Torque 3 for Formula CG with MK6i \$5,970.00 K0856
The Torque 3 Base comes Standard with 4-Pole Motors and up to a top speed of 5.8 MPH.

Group 4 - K0884 Basic Package	
3GAR-CG Arrow for Formula CG with MK6i	
300 lbs. Weight Capacity	
1431-3 14" x 3" Tires with Foam Filled Inserts	No Charge
TrueTrack Gearless Brushless™ GB Motors	Standard
1228-3 8" x 2 1/4" Semi Pneumatic Casters	Standard
Flip Back Adjustable Height Arm	Standard
MK6TT TrueTrack Controller	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard
Seat Option CGR Recline (E1004)	\$5,355.00

Group 4 - K0886 Basic Package	
3GAR-CG Arrow for Formula CG with MK6i	
400 lbs. Weight Capacity	
1431-3 14" x 3" Tires with Foam Filled Inserts	No Charge
TrueTrack Gearless Brushless GB Motors	Standard
1228-3 8" x 2 1/4" Semi Pneumatic Casters	Standard
Flip Back Adjustable Height Arm	Standard
MK690ACC Controller with Actuator Control	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard
Seat Option CGR Recline (E1004)	\$5,355.00

Group 3 - K0862 Basic Package	
3GRX-CG Ranger X for Formula CG with MK6i	
300 lbs. Weight Capacity	
1431-3 14" x 3" Tires with Foam Filled Inserts	No Charge
Heavy-Duty 4 Pole Motors	Standard
1228-3 8" x 2 1/4" Semi Pneumatic Casters	Standard
Flip Back Adjustable Height Arm	Standard
MK690ACC Controller with Actuator Control	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard
Seat Option CGR Recline (E1004)	\$5,355.00

Group 4 - K0877 Basic Package	
3GTQ-CG Torque SP for Formula CG with MK6i	
300 lbs. Weight Capacity	
1431-3 14" x 3" Tires with Foam Filled Inserts	No Charge
Performance 4 Pole Motors	Standard
1228-3 8" x 2 1/4" Semi Pneumatic Casters	Standard
Flip Back Adjustable Height Arm	Standard
MK690ACC Controller with Actuator Control	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard
Seat Option CGR Recline (E1004)	\$5,355.00

Group 3 - K0856 Basic Package	
3GTQ3-CG Torque 3 for Formula CG with MK6i	
300 lbs. Weight Capacity	
1431-3 14" x 3" Tires with Foam Filled Inserts	No Charge
4 Pole Motors	Standard
1228-3 8" x 2 1/4" Semi Pneumatic Casters	Standard
Flip Back Adjustable Height Arm	Standard
MK690ACC Controller with Actuator Control	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard
Seat Option CGR Recline (E1004)	\$5,355.00

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.**

Invacare® Formula™ CG Recline for Storm Series® Bases

OMIT SEAT FRAME

- OSF Omit Seat Frame for Base Only..... Less \$600.00

TRANSPORT TIE DOWN

- TRBKTS Wheelchair Transport Brackets \$250.00
Only available with MED or TALL seat-to-floor height. Not available with vent tray.

TRANSPORT TIE-DOWN*

TRRO includes four factory installed wheelchair transport brackets. TRRO has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport. As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

* Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

USER WEIGHT LIMIT

Clear

Arrow

- U300 User Weight <300 lbs. No Charge
 U400 User Weight <400 lbs. No Charge
Comes standard with Heavy Duty Monroe® Spring Suspension.

Ranger X

- U400 User Weight <400 lbs. No Charge
Comes standard with Heavy Duty Monroe® Spring Suspension.

Torque SP / Torque 3

- U300 User Weight <300 lbs. No Charge

DRIVE TYPE OPTION

Clear

- TRUTK TrueTrack Upgrade \$1,899.00
TrueTrack includes Gearless Brushless GB Motors and Monroe® Spring Shock Suspension. Must upgrade to MK6TT Controller. Must upgrade to an Expandable Controller System. Not compatible with SPJ+ joystick. Available option on Ranger X, Torque SP and Torque 3. Standard on Arrow. When selecting TrueTrack motors for Torque 3 the top speed of the chair increases to 7 MPH and therefore no longer performs to the Group 3 specifications.
- PERF Performance 4 Pole Motor No Charge
PERF Standard on Torque SP. Top speed of up to 6 1/2 MPH. Not available on Torque 3

VENT TRAY

Vent tray on chair reduces chair weight capacity by approximately 50 lbs. Not available with TRBKTS.

- VENT Articulating Vent Tray \$1,750.00 E1030

FRAME TYPE

- SFS Short Base Standard
Not available with seat depths 19" - 22". Must pick long base.
- LFS Long Base No Charge

CASTER OPTIONS

- 1259 6" x 2" Semi-Pneumatic Casters No Charge
 1228-3 8" x 2 1/4" Semi-Pneu Casters Lt Grey Standard
 1251-3 8" x 2" Pneumatic Casters No Charge
 1220-3 9" x 2 3/4" Pneumatic Casters No Charge
Not available with short base. Not available with TRBKTS.
- 1465-3 9" x 2 3/4" Foam Filled Casters No Charge
Not available with short base. Not available with TRBKTS.

CASTER SUSPENSION

Clear

Not available with 9" casters.

- 9213-6 Shock Fork For 6" Casters \$119.00 E1016
 9213-8 Shock Fork For 8" Casters \$119.00 E1016
Not available with short base.

WHEEL OPTION

- 1430-3 14" x 3" Pneumatic Tires Standard
 1431-3 14" x 3" Tires with Foam Filled Inserts No Charge
 1440 14" x 4" Pneumatic Tires \$200.00
14" x 4" wheels increase the overall width of the base by 2 1/2". Does not include Storm Series design wheel or fender. Not available with (TRUTK) TrueTrack package.
- 1441 14" x 4" Tires with Foam Filled Inserts \$320.00
14" x 4" wheels increase the overall width of the base by 2 1/2". Does not include Storm Series design wheel or fender. Not available with (TRUTK) TrueTrack package.

FREE WHEEL HUB

Will add 3/4" overall width. Not available with True Track. Does not include Storm Series design wheel or fender.

- FHUB Free Wheel Hub \$125.00

FRAME FINISH

CONTEMPORARY

- 24P Wet Black No Charge
 60P Silver Metallic No Charge
 105P Midnight Blue No Charge
 115P Black Prism No Charge

TOO HOT

- 30P Sunny Yellow No Charge
 61P Electric Red No Charge
 119P Bubblegum Pink No Charge
 120P Grape Madness No Charge

WAY COOL

- 62P Electric Blue No Charge
 127P Grasshopper Green No Charge
 116P Sunkissed Orange No Charge
 117P Tangerine Red No Charge
 29P White Snow No Charge

Discontinued colors 71P Silver Vein, 104P Electric Teal, 121P Lolly Pop Blue, 122P Cosmic Blue, and 125P Emerald Green can be order through specials number C-070109031327

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.**

Invacare® Formula™ CG Recline for Storm Series® Bases

SEAT

SYSTEM TYPE

- CGR CG Recline System Only \$5,355.00 E1004
- CNVACT Conventional Actuator Standard
- SMTACT Smart Actuator \$500.00

Must select MK690 controller when selecting SMATACT (Smart Actuator). Does not apply when selecting MK6TT. SMTACT is standard when selecting 3GAR-CG.

SEAT-TO-FLOOR HEIGHT

- LOW 17 3/4" at 0° No Charge
Not available with Transport Tie Down option.
- MED 18 3/4" at 0° No Charge
- TALL 19 3/4" at 0° No Charge

SEAT SELECTION

- CNB Conventional Back No Charge
- CTB Contoura™ Seating \$625.00 E2620 E2621
- OMCU Omit Contoura Seat Cushion Less \$20.00
- CTBMSH Contoura™ with mesh upholstery \$635.00 E2620 E2621
- OMCU Omit Contoura Seat Cushion Less \$20.00

ADJUSTABLE SEAT WIDTH RANGE

- Clear
- 1620W 16" to 20" Seat Width Range No Charge
- 2024W 20" to 24" Seat Width Range No Charge

SEAT WIDTH SETTING

- Clear
- W16 16" Seat Width No Charge
- W17 17" Seat Width (Not available with Contoura) No Charge
- W18 18" Seat Width No Charge
- W19 19" Seat Width (Not available with Contoura) No Charge
- W20 20" Seat Width No Charge
- W21 21" Seat Width (Not available with Contoura) No Charge
- W22 22" Seat Width No Charge
- W23 23" Seat Width (Not available with Contoura) No Charge
Available only with 400 lbs. package.
- W24 24" Seat Width No Charge
Available only with 400 lbs. package.

ADJUSTABLE SEAT DEPTH RANGE

- Clear
- 1619D 16" to 19" Seat Depth Range No Charge
- 1922D 19" to 22" Seat Depth Range No Charge

SEAT DEPTH SETTING

- Clear
- D16 16" Seat Depth No Charge
- D17 17" Seat Depth No Charge
- D18 18" Seat Depth No Charge
- D19 19" Seat Depth No Charge
- D20 20" Seat Depth No Charge
- D21 21" Seat Depth No Charge
- D22 22" Seat Depth No Charge

SEAT TILT

Indicates static tilt of seat. Impacts seat-to-floor height, every 5 degrees raises the seat-to-floor height.

- 0SA 0° Standard
- 5SA 5° No Charge

CONVENTIONAL BACK HEIGHT

- Clear
- TSBH20 20" Back Height No Charge
- TSBH21 21" Back Height No Charge
- TSBH22 22" Back Height No Charge
- TSBH23 23" Back Height No Charge
- TSBH24 24" Back Height No Charge
- TSBH25 25" Back Height No Charge
- TSBH26 26" Back Height No Charge

BACK HEIGHT CONTOURA

- Clear
- TSBH25 25" Back Height No Charge
- TSBH28 28" Back Height No Charge

ADJUSTABLE ARM TYPE (EACH)

- Clear
- Right Arm**
- RA19FR Full Length Adj. Hgt. Right No Charge
- RA29FR Desk Length Adj. Hgt. Right No Charge
- Left Arm**
- Clear
- RA19FL Full Length Adj. Hgt. Left No Charge
- RA29FL Desk Length Adj. Hgt. Left No Charge

SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Positioning Strap \$50.00
- 1515 Push Button Style Seat Positioning Strap Standard

CHEST POSITIONING STRAPS

- Clear
- 7311BK Airline Buckle Chest Positioning Strap \$100.00 E0960
- 9515 Auto Style Chest Positioning Strap \$100.00 E0960
- 7321BK Hook and Loop Chest Positioning Strap \$100.00 E0960

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.**

POSITIONING ACCESSORIES

HEADREST TYPE (OPTIONAL) Clear

Hardware is packaged with equipment and cannot be ordered separately.

- TAR127700 Curved Headrest \$210.00 E0955
TAR27700HW Removable & Adjustable Hardware \$210.00 E1028
- TAR128000 Two-Step Headrest \$210.00 E0955
TAR128000HW Removable & Adjustable Hardware \$210.00 E1028
- T0022400 Large Headrest \$210.00 E0955
T0022400HW Removable & Adjustable Hardware \$210.00 E1028
- T0022600 Small Headrest \$210.00 E0955
T0022600HW Removable & Adjustable Hardware \$210.00 E1028

STEALTH HEADRESTS Clear

Hardware is packaged with equipment and cannot be ordered separately.

- CP175 14" Multi Axis Headrest \$225.00 E0955
CP175HW Removable & Adjustable Hardware \$225.00 E1028
- CP275 10" Multi Axis Headrest \$225.00 E0955
CP275HW Removable & Adjustable Hardware \$225.00 E1028
- CP180 14" TWB Headrest \$225.00 E0955
CP180HW Removable & Adjustable Hardware \$225.00 E1028
- CP280 10" TWB Headrest \$225.00 E0955
CP280HW Removable & Adjustable Hardware \$225.00 E1028

OTTO BOCK™ ARM PADS (PAIR) Clear

Not compatible with quadlink retractable joystick mount. If selecting Otto Bock hand pad, a large, medium, or small forearm pad must be selected.

- OB20400 One Piece Channel Armpads \$230.00 E2209
- OB20500 Large Forearm Pads \$235.00
- OB20600 Medium Forearm Pads \$235.00
- OB20700 Small Forearm Pads \$104.00

OTTO BOCK HAND PADS (PAIR) Clear

Not compatible with quadlink retractable joystick mount. If selecting Otto Bock hand pad, a large, medium, or small forearm pad must be selected.

- OB23000 Large Flat Hand Pads \$100.00
- OB23300 Medium Flat Hand Pads \$100.00

LATERALS (Pair)

CONVENTIONAL BACK Clear

For use with channel canes only. Hardware is packaged with equipment and cannot be ordered separately.

- TMSLM Clamp-On Swingaway Laterals Medium \$210.00 E0956
TMSLMHW Removable and Adjustable Hardware \$210.00 E1028
- TMSLL Clamp-On Swingaway Laterals Large \$210.00 E0956
TMSLLHW Removable and Adjustable Hardware \$210.00 E1028

CONTOURA BACK Clear

Hardware is packaged with equipment and cannot be ordered separately.

- CSLM Swingaway Contoura Lat. Med \$210.00 E0956
CSLMHW Removable and Adjustable Hardware \$210.00 E1028
- CSLL Swingaway Contoura Lat. Lg \$210.00 E0956
CSLLHW Removable and Adjustable Hardware \$210.00 E1028

STEALTH LATERALS (Pair)

CONVENTIONAL BACK Clear

Hardware is packaged with equipment and cannot be ordered separately.

- TWB-TPR Adj. Swing Contoura Lat. Med \$320.00 E0956
TPRHW Removable and Adjustable Hardware \$320.00 E1028
- TWB-TPRL Adj. Swing Contoura Lat. Large \$320.00 E0956
TPRLHW Removable and Adjustable Hardware \$320.00 E1028

COUNTOURA BACK Clear

Hardware is packaged with equipment and cannot be ordered separately.

- TWB-CPR Adj. Swing Contoura Lat. Med \$320.00 E0956
CPRHW Removable & Adjustable Hardware \$320.00 E1028
- TWB-CPRL Adj. Swing Contoura Lat. Lg \$320.00 E0956
CPRLHW Removable & Adjustable Hardware \$320.00 E1028

ELECTRONICS

MK6i CONTROLLERS Clear

Ranger X / Torque SP / Torque 3

- MK690ACC 90 AMP Controller with Actuator Control Standard
- MK690 90AMP Controller No Charge
For use with Smart Actuators only.
- MK6TT TrueTrack Controller No Charge
- GTRAC G-Trac Module \$995.00 K0108
G-Trac provides Improved Tracking Capabilities and Enhanced Traction. Does not require Special Motors. Not available with SPJ+ joysticks or MK6TT.

Arrow

- MK6TT TrueTrack Controller Standard

STANDARD DRIVER CONTROL

Torque SP / Torque 3

- SPJAP SPJ+ with Power Seating Switch Standard

ACTUATOR CONTROL THROUGH SPJ+

- SPJAPM SPJ+ with Joystick Actuator Control \$1,300.00 E2310
Optional for Torque SP and Torque 3. Allows actuator operation through Driver Control. Includes port for attendant switch.

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.**

ELECTRONICS Continued

EXPANDABLE CONTROLLER SYSTEM Clear

- PWH Harness req. for Expandable Electronics \$350.00 E2313
Required to interface joystick or MK6i display with expandable controllers.
- CMPJM6U MPJ+ Multiple Drive Joystick with **Color Display** .. \$750.00 E2377
- PSFM6 Personal Joystick, Inductive on Front \$1,130.00 E2377
TOD-SS: Toggle On/Off Drive Select, Speed Potentiometer Standard.
- PSRM6 Personal Joystick, Inductive on Rear \$1,130.00 E2377
TOD-SS: Toggle On/Off Drive Select, Speed Potentiometer Standard. Clear
- PBOD Push Button On/Off with PB Drive Select No Charge
PSRM6 Only.
- PBSS Push Button On/Off with Speed Select No Charge
Drive select through mode switch, or mode switch and joystick.
PSRM6 Only. Clear
- OR On/Off Mounted Right No Charge
PSRM6 Only.
- OL On/Off Mounted Left No Charge
PSRM6 Only.
- ASLRDYM6 MK6 Display Only (No Driver Control) \$750.00 E2377
Display only, No Driver Control. Allows choice from order form.
- MCD Mini Color Display \$750.00 E2377
No port for remote On/Off switch. Stereo port only.

FOR YOUR ASL DRIVE SYSTEM NEEDS
PLEASE REFER TO THE
ASL SINGLE INVOICE PRICE LIST FORM # 02-097.

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS Clear

- 1500M6 Prop RIM Head Control with MK6i Display .. \$2,875.00 E2327 Clear
- SMHD Small Headrest No Charge
- LGHD Large Headrest No Charge
- 1558MM6 Compact Joystick with Display \$1,130.00 E2373

ALTERNATIVE DIGITAL CONTROL SYSTEMS Clear

- Must also use either a Multiple Drive Joystick or MK6i display.**
- SNPM6 Sip-N-Puff Digital Interface Box \$1,500.00 E2325
 - 5018M6 Wafer Board and Digital Interface Box \$1,550.00 E2322
 - 5020M6 Mini Tash Joystick and Digital Interface Box \$1,750.00 E2321

SIP-N-PUFF KIT

- PKG32666 Therafin® Sip-N-Puff Breath Tube Kit \$380.00 E2326

ADDITIONAL CONTROL CHOICES

- 1558M6 Compact Joystick - Less Std Mount \$770.00
- 5020 Mini Tash Joystick \$405.00
Requires SNPM6 Sip-N-Puff Interface Box.
- PACM6 Prop Attendant Control with quick release mount \$850.00
- 1552M Digital Attendant Control \$850.00
Requires SNPM6 Sip-N-Puff Interface Box.

JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retract. Joystick Mnt CMPJ \$279.00 E1028
Includes ARM250 standard.
- QLA Quad Link Retract. Joystick Mnt SPJ+ \$279.00 E1028
Includes ARM250 standard.
- ARM250 Stealth Height Adj Joystick \$210.00 E1028
- LEFT Left Handed Mounted Joystick No Charge
- RIMHW RIM Hardware For 1500 \$219.00 E1028
- PKG32669 Complete Bib Assembly (Tash/1558) \$265.00
- GATMPJM6 Gatlin Midline for CMPJ \$600.00 E1028
- GAT1812 Gatlin Midline for 1812 \$600.00 E1028
- RCM Rear Cane Mount for Compact Joystick \$395.00

ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer with Pro Memory Card \$475.00
- 1813SD Pro Memory Card (USBReady) \$90.00
- 1560 "T" Handle Flexible Joystick Ext. \$77.00 E2323
- 1561 Straight Handle Flexible Joystick Ext. \$77.00 E2323
- 1826 Chin Cup \$48.00 E2324
- A24VPS 24V MK6 Auxiliary Power Source \$200.00
- AUX12M6 Auxiliary Module for outputs 1 and 2 \$795.00
- AUX34M6 Auxiliary Module for outputs 3 and 4 \$795.00
- 2PB Dual Push Button \$147.00
Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options.
- 2SW Dual Toggle \$172.00
Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options.
- PBMT Push Button Mounting Kit \$252.00
- MPY Monoport "Y" Cable for CMPJ+ or Display \$15.00
- RFM Wireless (RF) Mouse Only Module \$750.00
RFM includes Transmitter, USB Radio Frequency Mouse Dongle (Receiver), and Stereo Switch with hook and loop mount for mouse clicks. Not compatible with SPJ+.

POWERED SEATING ACCESSORIES Clear

- SANODE One Actuator Operation thru Driver Cntrl..... \$1,300.00 E2310
Required for Tilt Operation through Driver Control - Single Actuator Systems Only. Not available with SPJ+ Joystick. Does not come with Egg switch.
- 4WSB Four Way Switch Box No Charge Clear
- FWT 4 Way Toggle No Charge
- QPB 4 Push Buttons No Charge
- S4WSB Multiple Actuator Interface Box \$2,600.00 E2311 Clear
- FWT 4 Way Toggle No Charge
- QPB 4 Push Buttons No Charge



Systems

For systems with two or more actuators, NOT operating through driver control, choose 4WSB with either 4-way toggle or 4-push button.

For systems with 2 or more actuators OPERATING THROUGH DRIVER CONTROL, choose S4WSB.
Uses 4-way toggle or 4-push button as attendant control.

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Formula™ CG Recline for Storm Series® Bases

ELECTRONICS Continued

ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Switch - Black \$156.00
- ASL304A Light Wobble Switch \$180.00
- ASL202J Fiber Optic Switch \$732.00
Must add MK6 Auxiliary Power Source.
- ASL208J Adjustable Proximity Switch \$660.00
Must add MK6 Auxiliary Power Source.
- ASL504M6 Remote Emergency Stop Switch \$1,185.00
Must add MK6 Auxiliary Power Source.
- ASL611 Adjustable Switch Mounting Bracket \$252.00

MISCELLANEOUS

BATTERY BOXES

- 22NF Group 22 Battery Boxes No Charge
- GP24 Group 24 Battery Boxes No Charge
Option on Torque SP and Torque 3. Must select TRBKTS.

BATTERIES

- Clear
- Batteries will be installed by Invacare Corporation. Chair will be shipped complete with batteries.**
- 22NFBATTERY 22NF Battery (On Chair) \$300.00 E2361
Not available with Arrow GB, Ranger X, or TrueTrack systems or chairs with TRBKTS.
- 24BATTERY 24 Gel Battery (On Chair) \$410.00 E2363

BATTERY CHARGER OPTION

- Clear
- 110CHARGER 110 Volt Battery Charger No Charge
- LESTR Lester 8 AMP Battery Charger \$320.00 E2367
- L8900 Omit Battery Charger Less \$80.00

MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks \$45.00
- 9205-3 Suspension Shocks \$400.00 E1016
Standard with Ranger, Arrow GB, or TrueTrack Systems.
- 9205H-3 Suspension Shocks (Heavy Duty) \$400.00 E1018
Standard with Ranger, Arrow GB, or TrueTrack Systems.
- FDR Fender \$320.00
Standard with Ranger, Arrow GB, or TrueTrack Systems.
- EXTP Extended Anti-Tippers \$250.00
- ADD Adductors Small \$200.00 E0956
- ADDHW Hardware for ADD \$400.00 E0956
Hardware is package with equipment and cannot be ordered separately.
- ADD L Adductors Large \$200.00 E1028
- ADDHWL Hardware for ADD L \$400.00 E1028
Hardware is package with equipment and cannot be ordered separately.

FRONT RIGGINGS

HANGER BRACKET TYPE

- Clear
- ASB Center Pivot Style Bracket No Charge
- PSB Pin Style Bracket No Charge

CENTER PIVOT

- Clear
- 60HD 60° Swingaway Footrests 13 1/2" - 17 1/2" No Charge
- 70NHD 70° Swingaway Footrests 13 1/4" - 17" Standard
- 70STAPER 70° Tapered Footrests 13" - 17" No Charge
Not available with 19" or 20" wide and 9" casters. Standard with impact guards on the footrests.
- AT5544 Elevating Legrests \$255.00 E0990
Manual elevating legrest.

PIN STYLE (PAIR)

- Clear
- PW93 Swingaway Footrests 13 1/2" - 19" No Charge
- P904A Swingaway Elev. Legrests 16" - 19 3/4" \$255.00 E0990
Manual elevating legrests.
- PAL4A Smartleg Art. Legrests 17 1/4" - 21 1/2" \$320.00 K0053
Manual elevating legrests.

POWER ELEVATING LEGRESTS - ARTICULATING

- Clear
- ELRPW Power ELR with Articulation 20" Ext. \$1,895.00 E1010
Not available with SPJ+ Joystick. Must pick 4WSB or S4WSB to operate the power legs. Independent powered elevating legrests including elevating and articulating legrests with angle adjustable footplates. Not available with the 16" or 17" seat depths. Extension range 13" - 20".

FOOTPLATE OPTIONS (PAIR)

- Clear
- 1651 Composite Footplates No Charge
- 1350 Extra Large Aluminum Footplates No Charge
- AT5543 Adjustable Angle Flip-Up Footplates \$165.00 K0040

FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap \$27.00 E0038
- IMPACTG Impact Guards \$51.00
- 1600BK Heel loops with ankle straps \$56.00 E0951
- LPT2 Longer Pivot Slide Tube for 1350 \$83.00
- ALPT Longer Pivot and Slide Tube \$84.00
- ALPT4 Longer Pivot and Slide Tube 4" \$84.00

Contoura Seat Sizes

	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X
W24	*	*	X	X	X	X	X

*Not available with Contoura Seat Cushion

Back Heights 25" or 29.5"

Conventional Seat Sizes

Seat Widths	16" - 24"
Seat Depths	16" - 22"
Back Heights	20" - 26"

Battery Boxes

	Standard Battery Tray*
Arrow	Group 24
Ranger X	Group 24
Torque SP/Torque 3	22 NF; Group 24 (Optional)
*Sealed Gel Cell Batteries Recommended	

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.**

Invacare® Formula™ CG Recline for Storm Series® Bases

Storm for Formula CG Recline Specifications				
Max Speed MPH	Arrow	Ranger X	Torque SP	Torque 3
TrueTrack	7	7 (option)	7 (option)	7 (option)
4-Pole	N/A	N/A	N/A	5.8
4-Pole Perf	N/A	N/A	6 1/2	N/A
4-Pole HD	N/A	5	N/A	N/A
Seat-to-Floor at 0°				
Low	17 3/4"	17 3/4"	17 3/4"	17 3/4"
Med	18 3/4"	18 3/4"	18 3/4"	18 3/4"
High	19 3/4"	19 3/4"	19 3/4"	19 3/4"
Weight Limits				
<300 lbs.			X	X
<400 lbs.		X		
<300 lbs. TrueTrack	X		X	X
<400 lbs. TrueTrack	X	X		

Memory Card - BASIC

- Standard with ALL Expandable Control Systems (CMPJ+, PSR+, PSF+, MK6i Display).
- Used to back up / restore programmed settings for only one chair.
- May be given to user or kept with provider's files for safe keeping.
- Does not contain advanced diagnostics / help screens / file structure.
- Not compatible with SPJ+ Joysticks.

Memory Card - PROFESSIONAL

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains advanced diagnostics, help screens, software upgrades, file storage / retrieval
- Used to back up multiple chairs and programming settings.
- Intended for use by qualified providers only.
- Not compatible with SPJ+ Joysticks.

After Market Backrests

The following companies also have hardware to make their backrests work with our new canes:

- AEL - planar seating system
- Freedom Designs® - planar seating system
- Canyon - planar seating system

Batteries (not included)

22NF batteries with terminal configuration (positive on the left and negative on the right) MUST be used. 22NF batteries that have the reverse terminal configuration MUST not be used. Terminals MUST have a cross hole for proper battery connection.

GP24 batteries with terminal configuration (negative on the left and positive on the right) MUST be used. GP24 batteries that have the reverse terminal configuration MUST not be used. Terminals MUST have a cross hole for proper battery connection. See Owner's Manual, part number 1114809. These recommendations MUST be followed otherwise injury and damage may occur.

Tech Tips

Operating Powered Seating through the Driver's Control. All chairs require adding or replacing one component to allow actuator operation through the driver control.

SPJ+ with PSS:

Replace joystick with SPJ+ with ACC (Actuator Control) - E2310.

Single Actuator Systems with CMPJ+, PSR+, PSF+ MK6i Display:

Add a SANODE (Single Actuator Node) to the system - E2310.

Multiple Actuator Systems:

Replace the standard 4 Way Switch box with a S4WSB (Multiple Actuator Interface Box) - E2311.

Formula CG Recline Features and Specifications

- Weight capacity of 300 lbs or 400 lbs (with heavy-duty motors).
- Recline range of 170°.
- Mechanical Shear reduction of 3".
- Conventional or Contoura style seating.
- Adjustable seat frame with Conventional style.
- Seat-To-Floor height of 17 3/4" at 0° tilt.

Comments

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare Corporation

www.invacare.com

USA

One Invacare Way
 Elyria, OH USA
 44035-4190
 800-333-6900

©2011 Invacare Corporation

All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.

