



Yes, you can.®

ValueCare Bed

- *Single Motor Design*
- *Compatible with Invacare® Universal Bed Ends*





Yes, you can.®

The Invacare® ValueCare Semi-Electric bed utilizes a single motor design making the bed lightweight and easy to assemble. It features the popular universal bed ends with the Invacare® deck design and a 9V battery backup for low cost protection from power outages. The advanced pendant design allows the head and knee sections to move either simultaneously or independently at the push of a button.



Features

Model no. VC5310
 Bed Packages: VCPKG1VC: VC5310,5085, 6629
 VCPKG2VC: VC5310,5085, 6630

- Single motor design makes the bed lightweight and easy to assemble
- 9V battery backup provides convenient low cost protection in case of power outages
- Advanced pendant design allows head and knee sections to operate simultaneously or independently
- Universal bed ends are high impact and resist impacts and scratches
- Semi-electric operation:
 - Electric head and knee gatch
 - Manual hi/lo operation
- 450 lb. weight capacity



9V battery.



Single motor.



Hand pendant.

Specifications

Model no.	Bed deck height*	Overall height	Sleep surface	Weight capacity	Limited warranty
VC5310	15" - 23"	36"W x 88"L	36"W x 80"L	450 lb.**	2 years materials, workmanship and electrical

*Measured from floor to deck.

**Inclusive of all accessories.

Invacare Corporation
 www.invacare.com

USA
 One Invacare Way
 Elyria, Ohio
 44035-4190
 (800) 333-6900

Canada
 570 Matheson Blvd. E.
 Unit 8
 Mississauga, Ontario
 L4Z 4G4 Canada
 (800) 668-5324

Invacare, the Blue Medallion Design and Yes, you can. are registered trademarks of Invacare Corporation. Specifications are subject to change without notification.

© 2006
 Invacare Corporation

Form No. 06-024

Invacare® ValueCare Bed



"Invacare says 'Yes, you can,' and I couldn't have said it better myself."