

INVACARE®
LYNX™
SCOOTER



The **Invacare® Lynx™ L-3 and L-4 Micro-Portable Scooters** allow you to accomplish daily activities with ease. You'll be amazed at how simple they are to assemble and disassemble for transport. Their compact size and seamless functionality allow them to easily transition from the tight turns of the home to the wide-open spaces of the outdoors. With many of the same attributes on both models, choose between two stylish scooters packed with great features at a competitive price.

The three-wheel **Invacare® Lynx™ L-3 Scooter** takes portability to a new level. The L-3 model conveniently disassembles into four compact pieces for easy transport. Pick the destination, stow in the trunk and go. With an estimated seven miles of range, a plug-in battery charger and flat-free tires, consumers can plan their schedules around what they want to do, not around the limitations of their scooter.

Features

- Easy 4-piece disassembly
- Estimated 7 mile travel range
- Connector-free, drop-in battery pack
- On or off-board battery charging
- Foam-filled, flat-free tires
- Comfortable 360° rotating chair
- 3-position adjustable tiller
- Ample leg room and foot clearance

HCPSC Code: K0800



1. Disassembles for transport
2. 360° rotating swivel chair
3. On or off-board battery charging
4. Blue
5. Red



Model #	Overall Dimensions	Seat Dimensions	Arm Dimensions	Width Between Arms	Wheels	Ground Clearance	Step Height
L-3R, L-3B	39.8" L 22" W	18" W 16" H 15" D	11.6" L 9" H	15"-23"	8"	1.5"	4"
	Turning Radius	Incline	Total Weight	Weight Capacity	Max Speed	Battery Range	Limited Warranty
	39.4"	8°	98 lb.	300 lb.	5 mph	7 Miles	1 Year

Do not operate on roads, streets or highways.


All references to HCPSC codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

The **Invacare® Lynx™ L-4 Scooter** offers the features of the L-3 model with the performance and bold look of a four-wheel scooter. The stable four-wheel base of the L-4 model handles varying surfaces with ease. Customize the L-4 model to fit the consumer by simply adjusting the seat height, armrest width and tiller angle. Get comfortable, get set and go! With a top speed of 5 mph, you can feel confident the Lynx L-4 scooter will get you to your destination on time.

Features

- Easy 4-piece disassembly
- Estimated 7 mile travel range
- Connector-free, drop-in battery pack
- On or off-board battery charging
- Foam-filled, flat-free tires
- Comfortable 360° rotating chair
- 3-position adjustable tiller

HCPCS Code: K0800



1. Disassembles for transport
2. 360° rotating swivel chair
3. On or off-board battery charging
4. Blue
5. Red



Model #	Overall Dimensions	Seat Dimensions	Arm Dimensions	Width Between Arms	Wheels	Ground Clearance	Step Height
L-4R, L-4B	39.8"L 22" W	18"W 16"H 15"D	11.6" L 9" H	15"-23"	8"	1.5"	4"
	Turning Radius	Incline	Total Weight	Weight Capacity	Max Speed	Battery Range	Limited Warranty
	51.2"	8°	107 lb.	300 lb.	5 mph	7 Miles	1 Year

Do not operate on roads, streets or highways.



All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

The compact **Invacare® Lynx™ L-3X Scooter** has the power, comfort and travel range of a large scooter without the unnecessary bulk. With a 300 lb. weight capacity, sliding seat and ample leg room, the L-3X model can fit a wide variety of users. You no longer have to stock several models in your showroom when the Lynx scooter will do the trick!

Features

- 360° rotating and sliding seat
- 3-position adjustable tiller
- Ample leg room and foot clearance
- Easy disassembly for travel
- Comfortable seat with soft contoured foam
- Foam-filled, flat-free tires

HCPCS Code: K0800



1. Disassembles for transport
2. Sliding seat
3. Blue
4. Red



Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2012

Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and
®. All trademarks are
owned by or licensed
to Invacare Corporation
unless otherwise noted.
Specifications are subject
to change without
notification.
Form No. 06-047
Rev. 09/12
120354



Yes, you can.®

Model #	Overall Dimensions	Seat Dimensions	Arm Dimensions	Width Between Arms	Wheels	Ground Clearance	Step Height
L-3XR, L-3XB	44" L 22" W	18" W 13.5" H 16" D	11.6" L 7" H	15"-24"	8" Front 9" Rear	3"	5"
	Turning Radius	Incline	Total Weight	Weight Capacity	Max Speed	Battery Range	Limited Warranty
	44.7"	8°	158 lb.	300 lb.	5 mph	16-20 Miles	1 Year

Do not operate on roads, streets or highways.

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.