



Yes, you can.®

Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.



Yes, you can.®

TRANSACTION TYPE:

Order Quote

Date _____ Account # _____ P.O. # _____

Company Name _____ Phone # _____

Address _____ City _____ State _____ Zip Code _____ Country _____

Contact Name _____ Comments _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

STOCK MODELS

GROUP 2 - K0823 - CAPTAIN'S SEAT

All seat depths are adjustable and can be increased up to 2".
All seats have a center-post seat mount and can swivel.

- M41SRR Red, 18"W x 18"D Office Style Seat \$5,939.00 K0823
- M41SRB Blue, 18"W x 18"D Office Style Seat \$5,939.00 K0823
- M41SR16R Red, 16"W x 16"D Office Style Seat \$5,939.00 K0823
- M41SR16B Blue, 16"W x 16"D Office Style Seat \$5,939.00 K0823
- M41SR20R Red, 20"W x 18"D Office Style Seat \$5,939.00 K0823
- M41SR20B Blue, 20"W x 18"D Office Style Seat \$5,939.00 K0823

***Total MSRP (Base Model Plus Pre-Configured Features):**

- Pronto M41 power base with Semi-Recline, High-Back Office-Style Seat(MSRP) \$5,477.00 K0823
 - 10" Flat Free Tires
 - Shark NX 40 Amp Controller with SPJ+ Joystick
 - Push-Button Style Seat Positioning Strap
 - Fold-Down Foot Platform
 - 3 Amp Off-Board Charger
- U1 AGM Batteries Installed on Chair(MSRP) \$232.00 E2365
- Full Length, Adjustable Height Arms(MSRP) \$230.00 E0973

GROUP 2 - K0822 - CAPTAIN'S SEAT WITH SOLID PAN

All seat depths are adjustable and can be increased up to 2".
All seats have a center-post seat mount and can swivel.

- M41RSOLIDR Red, 18"W x 18"D Seat \$5,939.00 K0822
- M41RSOLIDB Blue, 18"W x 18"D Seat \$5,939.00 K0822
- M41RSOLID16R Red, 16"W x 16"D Seat \$5,939.00 K0822
- M41RSOLID16B Blue, 16"W x 16"D Seat \$5,939.00 K0822
- M41RSOLID20R Red, 20"W x 18"D Seat \$5,939.00 K0822
- M41RSOLID20B Blue, 20"W x 18"D Seat \$5,939.00 K0822

***Total MSRP (Base Model Plus Pre-Configured Features):**

- Pronto M41 power base with Semi-Recline Captain's Seat with Solid Seat Pan(MSRP) \$5,477.00 K0823
 - 10" Flat Free Tires
 - Shark NX 40 Amp Controller with SPJ+ Joystick
 - Push-Button Style Seat Positioning Strap
 - Fold-Down Foot Platform
 - 3 Amp Off-Board Charger
- U1 AGM Batteries Installed on Chair(MSRP) \$232.00 E2365
- Full Length, Adjustable Height Arms(MSRP) \$230.00 E0973

GROUP 2 - K0823 - FOLD-DOWN SEAT

All fold-down seats have a center-post seat mount and swivel.

- M41FDR Red, 18"W x 17"D Dark Gray Seat \$5,676.00 K0823
- M41FDB Blue, 18"W x 17"D Dark Gray Seat \$5,676.00 K0823

***Total MSRP (Base Model Plus Pre-Configured Features):**

- Pronto M41 Power Base with medium back , Fold-Down seat(MSRP) \$5,440.00 K0823
 - Full Length, Fixed Height Arms
 - 10" Flat Free Tires
 - Shark NX 40 Amp Controller with SPJ+ Joystick
 - Push-Button Style Seat Positioning Strap
 - Fold-Down Foot Platform
 - 3 Amp Off-Board Charger
- U1 AGM Batteries Installed on Chair(MSRP) \$232.00 E2365

ABSOLUTE™ CUSHION

General Purpose Cushion

- EC66 Absolute Cushion 16x16 \$104.00 E2601
- EC88 Absolute Cushion 18x18 \$104.00 E2601
- EC08 Absolute Cushion 20x18 \$104.00 E2601

COMFORT-MATE EXTRA CUSHION

Skin Protection Cushion

- CMEX66 Comfort Mate Extra 16x16 \$208.00 E2603
- CMEX88 Comfort Mate Extra 18x18 \$208.00 E2603
- CMEX08 Comfort Mate Extra 20x18 \$208.00 E2603

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Pronto® M41

SEMI-RECLINE OFFICE STYLE SEAT ONLY

- M5SRDG 18"W x 18"D Office Style Seat \$740.00
- M5SRDG16 16"W x 16"D Office Style Seat \$740.00
- M5SRDG20 20"W x 20"D Office Style Seat \$740.00

OPTIONAL STOCKED FRONT RIGGINGS

The below options are only available on Semi-Recline.

- FRAKH Front Rigging Attachment Kit \$55.00
- PHW93/34P Swing Away Footrests \$165.00
Must order FRAKH with this option.
- PH904A Swing Away Elevating Legrests \$260.00 E0990
/U67/34P **Must order FRAKH with this option.**

OPTIONS AND ACCESSORIES

- QLRIGHT Quad Link Joystick Mount Kit, Right Side \$230.00 E1028
- QLLEFT Quad Link Joystick Mount Kit, Left Side \$230.00 E1028
- VAN02 Oxygen Holder \$150.00 E2208
Only available on Semi-Recline.
- P725 Crutch and Cane Holder \$99.00 E2207
Only available on Semi-Recline.
- P770 Walker Holder \$99.00
Only available on Semi-Recline.
- P780 Large Rear Basket \$65.00
Only available on Semi-Recline.
- P710 Safety Flag Kit \$30.00
Only available on Semi-Recline.

Pronto® M41		
Dimensions	English	Metric
Length	33.7"	85.6 cm
Width (excluding joystick)	24"	61 cm
Height (back at 17")	33"	83.8 cm
Seat-to-Floor	19" - 24"	48.3 - 60.9 cm
Seat Width	16", 18", 20"	40.6, 45.8, 51 cm
Seat Depths	16" - 20"	40.6 - 50.8 cm
Back Height	16" - 23 1/2"	47.6 cm
Weight Capacity	300 lbs.	136 kg
Wheel & Caster Sizes		
Rear & Front Casters	6" x 2"	15.2 x 5.1 cm
Drive Wheels (Flat Free)	10 1/2" x 3 1/2"	26.7 x 8.9 cm
Maximum Speed	5 mph	8 km/h
Incline Capability	7 degrees grade	7 degrees grade
Turning Radius		
Front	17 3/4"	45.07 cm
Rear	19 1/2"	49.53 cm
Weight		
With U1 Batteries (18" w)	175 lbs.	79.4 kg
W/O Batteries	124 lbs.	56.2 kg
Ground Clearance	2"	5.1 cm
Warranties		
Base Frame	60 months	60 months
Seat Frame	36 months	36 months
Motors & Electronics	12 months	12 months

Comments

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.



Invacare Corporation
USA
One Invacare Way
Elyria, OH USA
44035-4190
800-333-6900

www.invacare.com

©2013 Invacare Corporation
All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.