

INVACARE®  
**BAR750**  
BARIATRIC BED

---



---

The **Invacare® BAR750 Bariatric Bed** features a heavy-duty, full-electric bed frame capable of supporting up to 750 lb, and comes standard with a battery backup with the patient in mind, this bed expands from 39" wide to a full 48" and from 80" in length to 88". Four ultra-quiet DC actuators easily reposition the head and foot sections, as well as the bed height.

---

## Features

Heavy-duty frame design ensures added strength for client support

Split-spring design eases delivery and setup

Two locking casters on head section and two steering lock casters on foot section

Two high-low motors, one on head section and one on foot section

Pendant is designed to allow for dynamic client positioning

Quiet, smooth operation


Integrated battery backup, bed ends and half rails standard

Choose from three mattress width options

BARMATT Mattress options meet fire safety regulations 16CFR1633 and 16CFR1632

675 lb. patient weight capacity

## HCPCS Code: E0302



1. Heavy-duty frame
2. Split spring design
3. Locking casters



Model #	Description
BAR750	Bariatric Bed, with 1 pair half rails

### Bed Packages

BARPKG750-1633	BAR750/BARMATT39
BARPKG750-1-1633	BAR750/BARMATT42
BARPKG750-2-1633	BAR750/BARMATTEXP

### Accessories

BARMATT39	Bariatric Mattress, 700 lb. capacity	80" L x 39" W x 6.5" D
BARMATT42	Bariatric Mattress, 700 lb. capacity	80" L x 42" W x 6.5" D
BARMATT48	Bariatric Mattress, 700 lb. capacity	80" L x 48" W x 6.5" D
BARMATTEXP	Bariatric Mattress, 700 lb. capacity	80" x 39" - 48" W x 6" D (expandable)
BAR750BOLSTER	Bed extension from 80" to 88"	8" x 39" (expandable)

**Invacare Corporation**  
www.invacare.com

### USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

### Canada

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324



Model #	Dimensions	Sleep Surface	Weight Capacity	Shipping Weight	Carton Dimensions	Regulation Compliance	Limited Warranty
BAR750	88"-96" L 15"-29" H 51" W	80-88" L 39", 42", 48" W	Safe Working Load: 440 lb. 750 lb.* Max. Patient Load: 675 lb.		56" L 44" H 41" W	N/A	5 Years Frame/Structure 3 Years Mechanical/ Electrical

©2018  
Invacare Corporation.  
All rights reserved.  
Trademarks are  
identified by the  
symbols ™ and ®.  
All trademarks are  
owned by or licensed  
to Invacare Corporation  
unless otherwise  
noted. Specifications  
are subject to change  
without notification.  
Form No. 07-010  
Rev. 11/18  
181135

\*Includes patient, mattress and all accessories.

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

