



Price List and Order Form

Effective Dec. 1, 2008

For ease of ordering, contact customer Service Toll Free at:

800-333-6900 Fax: 800-678-4682

To order Literature: 800-828-6282

www.invacare.com

Company Name _____ Account # _____

Phone _____ P.O. # _____ Date _____

Ship to Name & Address _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

FRAME SELECTION

CHAIR TYPE

- 3GTQ3 Torque 3 Base And Seat w/MK6i™ (K0848) .. 5,970.00 ¹
- 3GSEAT Seat Only 755.00 ²
-

Group 3 – K0848 Basic Package

- 3G Torque 3 with MK6i
- 300 lb Weight Capacity
- 1431-3 14" X 3" Tires w/Foam Filled Inserts No Charge
- 4 Pole Motors.....Standard
- 1228-3 8" X 2.25" Semi Pneumatic CastersStandard
- 28F Desk Length Fixed Height Arm.....Standard
- 70NHD 70° Swingaway FootrestsStandard
- 8 Amp Off Board Battery ChargerStandard
- 1515 Push Button Style Seat Positioning StrapStandard

TRANSPORT TIE-DOWN ³

- TRRO Transport Ready Option 300.00
- TRBKTS Wheelchair Transport Brackets 250.00
-

* TRANSPORT TIE-DOWN

TRRO includes four factory-installed transport brackets and a wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESN A W C Vol 1, Section 19 Frontal Impact Test requirements for wheelchairs with a 168 lb. crash dummy, which corresponds to a person with a weight of 114 lbs. to 209 lbs.

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with W C19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

* Transport Tie-Down option is not retrofittable to existing models and is not field serviceable.

TRRO FEATURES FOR TORQUE 3

Features Required for TRRO

- Torque 3 Base & Seat
- 300 lb. user weight limit
- Short or Long Base Frame
- Adjustable Rehab Seat with Adjustable Seat Tilt & Back Angle
- Any Seat Width
- Any Seat Depth
- Low Seat-to-Floor Height
- Sectional or Straight Rear Uprights
- 18" - 24" ASBA Back Height
- 90° or 95° Back Angle
- Any Arm Type
- 8" x 2.25" Semi Pneumatic Casters
- Suspension Shocks
- 14" x 3" Tires with Foam Filled Inserts
- Rear Wheel Standard Location
- Any Color Frame Finish
- Black Nylon or Black Vinyl Upholstery
- Center Pivot Style Bracket
- 60° Swingaway, 70° Swingaway or 70° Tapered Footrests
- Group 24 Battery Box

Optional Features for TRRO

- TrueTrack Upgrade
- Fender Option
- Flip-Up Composite Footplates, Aluminum Footplates-Deluxe, Adjustable Angle Flip Up Footplates
- Calf Strap
- Impact Guards
- All PSR / PSF Joysticks
- Left Hand Mounted Joystick
- "T" Handle Flexible Joystick Ext.
- Straight Handle Flexible Joystick Ext.
- MK6i Remote Programmer
- SPJ™+ or MPJ™+ Joysticks
- Reinforced Back Upholstery
- Group 24 Batteries

Please note that all TRRO chairs are set at a 10° seat angle

FOOTNOTES

¹ Torque 3 comes standard with 300 lb. user weight capacity, 5.8 mph top speed, 4-pole motor package, MK6i SPJ+ joystick, MK690 Controller.

² When ordering a replacement seat and changing depth to 20" or deeper, order new arms.

³ Must select Group 24 battery boxes for TRRO or TRBKTS.

NOTE: All specifications and dimensions are approximate.

For the most current pricing information go to www.invacare.com.

Invacare® 3G Torque™ 3 Accessories and Options

SEAT SELECTION

USER WEIGHT LIMITS

- U300 User Weight <300 lb No Charge

FRAME TYPE

- SFS Short Base..... **Standard**
 LFS Long Base..... No Charge

ADJUSTABLE SEAT WIDTH RANGE

- 1216W 12" to 16" Seat Width Range No Charge
 1620W 16" to 20" Seat Width Range No Charge
 2022W 20" to 22" Seat Width Range No Charge

FRAME WIDTH ⁴

- W12 12" Wide 110.00 ⁵
 W13 13" Wide 110.00 ⁵
 W14 14" Wide 110.00 ⁵
 W15 15" Wide No Charge ⁵
 W16 16" Wide No Charge
 W17 17" Wide No Charge
 W18 18" Wide No Charge
 W19 19" Wide No Charge
 W20 20" Wide No Charge
 W21 21" Wide 329.00
 W22 22" Wide 329.00

TRUETRACK UPGRADE

- TRUTK TrueTrack Upgrade 1,899.00 ⁶

SEAT SELECTION

- ADJASBA Adj. Rehab Seat w/ Seat Tilt & Back Angle... **Standard**
 ADJRECL Manual Recline for Adjustable ASBA (E1226) .. 744.00
 CTS Contoura Seating 395.00 ⁷
 OMCU Omit Contoura Seat Cushion..... Less 20.00

ADJUSTABLE SEAT DEPTH RANGE

- 1218D 12" to 18" Seat Depth Range No Charge ⁸
 1619D 16" to 19" Seat Depth Range No Charge
 1922D 19" to 22" Seat Depth Range (Long Base) ... No Charge

SEAT DEPTH

- D12 Seat Depth 12" (Short Base)..... 110.00 ⁵
 D13 Seat Depth 13" (Short Base)..... 110.00 ⁵
 D14 Seat Depth 14" (Short Base)..... 110.00 ⁵
 D15 Seat Depth 15" (Short Base)..... No Charge ⁵
 D16 Seat Depth 16" (Short Base)..... No Charge
 D17 Seat Depth 17" (Short Base)..... No Charge
 D18 Seat Depth 18" (Short Base)..... No Charge
 D19 Seat Depth 19" (Long Base)..... No Charge
 D20 Seat Depth 20" (Long Base)..... No Charge
 D21 Seat Depth 21" (Long Base)..... 529.00
 D22 Seat Depth 22" (Long Base)..... 529.00

ADJASBA SEAT-TO-FLOOR HEIGHT ⁹

- LOW Seat To Floor – Low 17.75" No Charge
 MED Seat To Floor – Med 20" No Charge
 TALL Seat To Floor – Tall 21.25" No Charge

FOOTNOTES

⁴ Contoura™ Seating available in even widths 16" - 22".
⁵ Jr. Seat widths/depths require a straight rear upright and a user weight limit of 150 lbs. Manual Reclining back not available. Additional charges apply to widths/depths less than 16". Jr. seat widths have Jr. arms with an adjustment range of 6"-10" inches in height.
⁶ TrueTrack driving technology, which enables driver to track straighter regardless of the terrain, includes Monroe Spring Suspension, <300 lb. User Weight Limit, and 7 mph maximum speed. **Must upgrade to an Expandable Controller System. Not compatible with SPJ+ joystick. When selecting TrueTrack motors, the top speed increases to 7 mph and therefore no longer performs to the Group 3 specifications.**
⁷ Selected seat depth on Contoura assumes 1" of back cushion compression. For example: When choosing an 18" depth the measurement from the front of the back cushion to the front of the seat pan equals 17". The Contoura back cushion is 2" thick. Please see page 5 for seat size offerings.
⁸ Must select 12"-18" Depth Range (1218D) when ordering the 12"-16" Seat Width Range (1216W). 12" & 13" seat widths not available with 20"-24" back heights.
⁹ Based on 5° tilt and seat depth range of 16-19": add .25" for seat depth range of 19-22"

BACK SELECTION

ADJASBA BACK TYPE

- SECT Sectional Rear Uprights..... No Charge
 SRUP Straight Rear Uprights No Charge
 9720P No Push Handles (Straight Canes) 82.00
 L8661 Omit Rear Uprights..... Less 20.00

RECLINER BACK HEIGHT

- BHR18 Back Height 18" No Charge
 BHR20 Back Height 20" No Charge
 BHR22 Back Height 22" No Charge
 BHR24 Back Height 24" No Charge
 BHR26 Back Height 26" No Charge

BACK HEIGHT

ADJASBA

- BH12 Back Height 12" No Charge
 BH13 Back Height 13" No Charge
 BH14 Back Height 14" No Charge
 BH15 Back Height 15" No Charge
 BH16 Back Height 16" No Charge
 BH17 Back Height 17" No Charge
 BH18 Back Height 18" No Charge
 BH19 Back Height 19" No Charge
 BH20 Back Height 20" No Charge
 BH21 Back Height 21" No Charge
 BH22 Back Height 22" No Charge
 BH23 Back Height 23" No Charge
 BH24 Back Height 24" No Charge

CONTOURA™

- TSBH25 Back Height 25" No Charge
 TSBH29.5 Back Height 29.5" No Charge

BACK ANGLE DEGREE

- 85BA 85° No Charge
 90BA 90° No Charge
 95BA 95° No Charge
 100BA 100° No Charge
 105BA 105° No Charge

BACK ACCESSORIES

- L8601 Omit Back Upholstery No Charge
 RBP Reinforced Back Upholstery No Charge

ARM SELECTION

ARM TYPE ¹⁰

- 18F Full Length Fixed Hgt. Conv. Flip Up No Charge
 28F Desk Length Fixed Hgt. Conv. Flip Up **Standard**
 19F Full Length Adj. Hgt. Conv. Flip Up (E0973) 230.00
 29F Desk Length Adj. Hgt. Conv. Flip Up (E0973) ... 230.00

CANTILEVER ARM TYPE (PAIR) ¹¹

- CLVF Full Length Arm \$350.00
 CLVD Desk Length Arm \$350.00

CANTILEVER MIXED ARM TYPE (EACH)

- CLVFR Full Length Arm Right \$175.00
 CLVDR Desk Length Arm Right \$175.00
 CLVFL Full Length Arm Left \$175.00
 CLVDL Desk Length Arm Left \$175.00

ARMPAD OPTIONS (CANTILEVER ARMS)

- WF Waterfall Armpads **Standard**
 RCL Reclining Armpads No Charge

FOOTNOTES

¹⁰ All arms are flip up and detachable. Seat widths 12-15" not available with fixed height arms must use adjustable height arm style only.
¹¹ Cantilever Arms only available with 300 lb. rated systems. Not available with BPO, Otto Bock or 12"-15" seat width.

Invacare® 3G Torque™ 3 Accessories and Options

MK6I ELECTRONICS

MK6i CONTROLLERS

- MK690 90 Amp Controller..... **Standard**
- MK6TT TrueTrack Controller..... No Charge ¹²
- GTRAC G-Trac Module (K0108)..... 995.00 ¹³

STANDARD DRIVER CONTROL

- SPJM6 SPJ+ Single Drive Joystick **Standard**

EXPANDABLE CONTROLLER SYSTEMS

- PWH Harness req. for Expandable Electronics (E2313) 350.00 ¹⁴
- MPJM6U MPJ+ Multiple Drive Joystick (E2377) 750.00 ¹⁵
- PSFM6 Personal Joystick, Inductive on Front 1,130.00 ¹⁵
- PSRM6 Personal Joystick, Inductive on Rear 1,130.00 ¹⁵
 - PBOD Push Button On/Off w/PB Drive Select No Charge ¹⁵
 - PBSS Push Button On/Off w/Speed Select No Charge ¹⁵
 - OR On/Off Mounted Right..... No Charge ¹⁵
 - OL On/Off Mounted Left..... No Charge ¹⁵
- ASLRDYM6 MK6 Display Only(No Driver Control) (E2377) 750.00 ¹⁷

**FOR YOUR ASL DRIVE SYSTEM NEEDS
PLEASE REFER TO THE FOLLOWING:
ASL SINGLE INVOICE PRICE LIST FORM # 02-097.**

ALTERNATE PROPORTIONAL CONTROL SYSTEMS

- 1500M6 Prop RIM Head Control w/Display (E2327) 2,875.00
- SMHD Small Headrest..... No Charge
- LGHD Large Headrest..... No Charge
- 1558MM6 Compact Joystick w/Display (E2373)..... 1,130.00

ALTERNATIVE DIGITAL CONTROL SYSTEMS

- SNPM6 Sip-N-Puff Digital Interface Box (E2325)..... 1,500.00 ¹⁸
- 5018M6 Wafer Board & Digital Interface Box (E2322). 1,550.00 ¹⁸
- 5020M6 Mini Tash Jystk & Dig Interface Box (E2321). 1,750.00 ¹⁸

SIP N PUFF KITS

- PKG32666 Therafin® Sip-N-Puff Breath Tube Kit (E2326). 380.00

ADDITIONAL CONTROL CHOICES

- 1558M6 Compact Joystick - Less std Mount 770.00
- 5018 Wafer Board 310.00 ¹⁹
- 5020 Mini Tash Joystick 405.00 ¹⁹
- PACM6 Prop. Attendant Control w/quick release mount.. 850.00
- 1552M Digital Attendant Control..... 850.00 ¹⁹

JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retract. Joystick Mnt MPJ+ (E1028) 279.00 ²⁰
- QLA Quad Link Retract. Joystick Mnt SPJ+ (E1028).. 279.00 ²⁰
- ARM250 Stealth Height Adjustable Joystick Mount..... 210.00
- LEFT Left Handed Mounted Joystick No Charge
- RIMHW RIM Hardware For 1500 (E1028)..... 219.00
- PKG32669 Complete Bib Assembly (Tash/1558)..... 265.00
- GATMPJ6 Gatlin Midline for MPJ+ (E1028)..... 600.00
- GAT1812 Gatlin Midline for 1812 (E1028) 600.00
- RCM Rear Cane Mount for Compact Joystick 395.00
-

FOOTNOTES

¹² **Must upgrade to an Expandable Controller System. Not compatible with SPJ+ joysticks.**

¹³ G-Trac provides Improved Tracking Capabilities and Enhanced Traction. Does not require Special Motors. Not available w/ SPJ+ joysticks or MK6TT

¹⁴ Harness is required to interface joystick or MK6i display with expandable controllers. Option is located under Electronic Accessories in configurator.

¹⁵ TOD-SS: Toggle On/Off Drive Select, Speed Potentiometer Standard.

¹⁶ Drive select through mode switch, or mode switch and Joystick command.

¹⁷ Display only, allows choice from order form.

¹⁸ Must also use either a Multiple Drive Joystick or MK6 Display.

¹⁹ Requires SNPM6 Sip-N-Puff Interface Box.

²⁰ QLA and QLAM6 include ARM250 standard.

ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer w/Pro Memory Card 475.00
- 1813SD Pro Memory Card (USB Ready) 90.00
- 1560 "T" Handle Flexible Joystick Ext. (E2323)..... 77.00
- 1561 Straight Handle Flexible Joystick Ext. (E2323)..... 77.00
- 1826 Chin Cup (E2324)..... 48.00
- A24VPS 24V MK6 Auxillary Power Source 200.00
- AUX12M6 Auxillary Module For 1 & 2..... 795.00
- AUX34M6 Auxillary Module For 3 & 4..... 795.00
- MPY Monoport "Y" Cable for MPJ+ or Display 15.00
- 2PB Dual Push Button..... 147.00 ²¹
- 2SW Dual Toggle 172.00 ²¹
- PBMT Push Button Mounting Kit..... 252.00

ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Switch Black 65.00
- ASL304 Light Wobble Switch w/611 Mount..... 320.00
- ASL202J Fiber Optic Switch..... 610.00 ²²
- ASL208J Adjustable Proximity Switch..... 550.00 ²²
- ASL504M6 Remote Emergency Stop Switch 1,185.00 ²²

FRAME FINISH AND UPHOLSTERY COLOR

FRAME FINISH

CONTEMPORARY

- 24P Wet Black No Charge
- 60P Silver Metallic No Charge
- 71P Silver Vein No Charge
- 105P Midnight Blue..... No Charge
- 115P Black Prism..... No Charge

TOO HOT

- 30P Sunny Yellow No Charge
- 61P Electric Red..... No Charge
- 119P Bubblegum Pink No Charge
- 120P Grape Madness No Charge

WAY COOL

- 62P Electric Blue No Charge
- 104P Electric Teal No Charge
- 121P Lolly Pop Blue No Charge
- 122P Cosmic Blue..... No Charge
- 125P Emerald Green No Charge
- 127P Grasshopper Green No Charge

UPHOLSTERY TYPE

- U240 Black Nylon Upholstery **Standard**
- U67 Black Vinyl..... No Charge

FOOTNOTES

²¹ Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options.

²² Not available with SPJ+ Joystick. Must add MK6 Auxillary Power Source.

Invacare® 3G Torque™ 3 Accessories and Options

CASTER AND WHEEL SELECTION

CASTER OPTIONS

- 1259 6" X 2" Semi-Pneumatic Casters No Charge
- 1228-3 8" X 2 1/4" Semi Pneu. Casters Lt Grey **Standard**
- 1251-3 8" X 2" Pneumatic Casters No Charge
- 1220-3 9" X 2 3/4" Pneumatic Casters No Charge ²³
- 1465-3 9" X 2 3/4" Foam Filled Casters (K0097) No Charge ²⁴

CASTER SUSPENSION ²⁴

- 9213-6 Shock Fork For 6" Casters (E1016) 119.00
- 9213-8 Shock Fork For 8" Casters (E1016) 119.00 ²⁵

WHEEL OPTION

- 1430-3 14" X 3" Pneumatic Tires **Standard**
- 1431-3 14" X 3" Tires w/Foam Filled Inserts..... No Charge ^{26,27}
- 1440 14" X 4" Pneumatic Tires 200.00 ^{26,27}
- 1441 14" X 4" Tires w/Foam Filled Inserts..... 320.00 ^{26,27}

FREE WHEEL HUB ²⁸

- FHUB Free Wheel Hub..... 125.00

REAR WHEEL LOCATION ²⁹

- RWF2 Rear Wheel Forward 2" No Charge
- RWF1 Rear Wheel Forward 1" No Charge
- RWS Rear Wheel Standard..... **Standard**

FENDER OPTION ³⁰

- FDR Fender 320.00

POSITIONING STRAPS

SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Positioning Strap.....50.00 ³¹
- 1515 Auto Style Buckle Seat Positioning Strap..... **Standard**

CHEST POSITIONING STRAPS

- 7311BK Airline Buckle Chest Positioning Strap (E0960) .. 100.00 ³¹
- 9515 Auto Style Chest Positioning Strap (E0960)..... 100.00
- 7321BK Hook & Loop Chest Positioning Strap (E0960) .. 100.00

MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks45.00 ³²
- SEATO2 Oxygen Holder (ASBA Seat Only) (E2208)..... 150.00
- 9205-3 Suspension Shocks (E1016) 400.00 ³³
- 9205H-3 Suspension Shocks Heavy Duty..... 400.00

BATTERY SELECTION

BATTERY BOXES

- 22NF 22NF Battery No Charge
- GP24 Group 24 Battery Boxes No Charge ³⁴

INSTALLATION OF BATTERIES ³⁵

- 22NFBATTERY 22NF Battery (on chair) (E2361) 300.00
- 24BATTERY 24 Gel Battery (on chair) (E2363) 410.00

BATTERY CHARGER OPTION

- 110CHARGER 110 Volt Battery Charger..... No Charge
- LESTR Lester 8 AMP Battery Charger (E2367) 320.00
- L8900 Omit Battery Charger Less 80.00

FOOTNOTES

- ²³ Not available with short base. Not available with Transport Tiedown options.
- ²⁴ Not available with 9" casters.
- ²⁵ Not available with TRRO or TRBKTS. Contact Customer Service for available special.
- ²⁶ 14" x 4" wheels increase overall width of the base by 2 1/2". Not available with TrueTrack package (TRUTK).
- ²⁷ Does not include Storm Series® mag wheel.
- ²⁸ Free wheel hub increases overall width 3/4". Not needed or available with TrueTrack package (TRUTK).
- ²⁹ This location changes the center of gravity and performance. Forward wheel position will make turning and climbing easier, but more prone to "weeely".
- ³⁰ Standard with TrueTrack option.
- ³¹ 1311BK not available with 12" - 15" width.
- ³² Not available with pneumatic tires.
- ³³ Not available with 9" casters.
- ³⁴ Required for transport tie down option.
- ³⁵ Batteries will be installed by Invacare Corporation at no additional charge. Chair will be shipped complete w/batteries.

FRONT RIGGING OPTIONS

HANGER BRACKET TYPE

- ASB Center Pivot Style Bracket No Charge
- PSB Pin Style Bracket No Charge
- NHB No Hanger Bracket No Charge

FRONT RIGGING

CENTER PIVOT

- 60HD 60° Swingaway Footrests 13.5"-17.5" No Charge
- 70NHD 70° Swingaway Footrests 13.25"-17" **Standard**
- 70STAPER 70° Tapered Footrests 13"-17" No Charge
- 90SWING 70° Swingaway Footrests 3"-10" No Charge
- AT5544 Elevating Legrests (E0990) 255.00

PIN STYLE

- PAS4A Smartleg® Art. Legrests 17 1/4"-21 1/2" (K0053)320.00
- PW93 Swingaway Footrests. 13 1/2" - 19" No Charge
- P904A Swingaway Elev Legrests 16"-19 3/4" (E0990) ... 303.00
- PAL4A Smartleg Art. Legrests 17 1/4" - 21 1/2" (K0053) 320.00 ³⁶

POWER LEGS

- ELRPW Power Elev Legrests w/Articulation (E1010) ... 1,895.00 ³⁷

POWER LEGS CONTROL MODULE ³⁸

- 4WSB Four Way Switch Box No Charge
- S4WSB Multiple Actuator Interface Box (E2311) 2,600.00
 - FWT 4 Way Toggle No Charge
 - QPB 4 Push Buttons..... No Charge

For all systems with two or more actuators you must choose 4WSB at no charge or the S4WSB at an upcharge.

S4WSB Replaces the Four Way Switch Box and allows multiple actuator operation through the driver control.

FWT or QPB is required when choosing the 4WSB or S4WSB.

FOOTPLATE OPTIONS

- 1651 Flip-Up Composite Footplates No Charge
- 1350 Extra Large Aluminum Footplates No Charge
- 1450 Aluminum Footplates-Deluxe No Charge
- S1263 One Piece Wooden Footboard No Charge
- AT5543 Adjustable Angle Flip Up Footplates (K0040) 165.00

FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap (K0038) 27.00
- IMPACTG Impact Guards..... 51.00
- A9342 "H" Heel Strap (K0039) 51.00
- LPT2 Longer Pivot Slide Tube For 1350 83.00
- ALPT Longer Pivot And Slide Tube 84.00
- ALPT3 3" Extension For 90° Footrest 31.00
- ALPT4 Longer Pivot And Slide Tube 4" 84.00

SEATING AND POSITIONING

FOR YOUR SEATING NEEDS
PLEASE REFER TO THE FOLLOWING

REHAB SEATING & POSITIONING PRODUCTS
See Price List Form 07-027

FOOTNOTES

- ³⁶ For seat widths greater than or equal to 16".
- ³⁷ Not available with SPJ+ Joystick. Must pick 4WSB or S4WSB to operate power legs. Independent powered legrests include elevating and articulating legrests with angle adjustable footplates. Extension range of 13"-20". Sold as pair.
- ³⁸ These options located in the Configurator under the MK6i Electronics Folder. Category Header is labeled Power Seating Accessories. Must select when ordering ELRPW.

Invacare® 3G Torque™ 3 Accessories and Options

Chair Features and Options	
•	Independent rear wheel suspension with new rubber element design
•	5.8 mph standard speed range
•	Programmable performance to meet individuals needs
•	17 3/4" Seat-to-Floor Standard, 20" and 21.25" Optional
•	Seat location system allows for balanced weight distribution
•	Continuous tilt seat 0° - 10°
•	Dynamic anti tippers
•	Rigid seat pan
•	Convenient rear battery access
•	22NF batteries (standard) In a single 22NF battery box Single battery connector
•	Arms are flip back/removable
•	Adjustable depth rigging attachment 2" (1" increments)
•	Automatic dual mode battery charger (Lester Brand)
•	Compact size: 24 3/4" overall base width w/14" x 3" wheels
•	Weight limit <300 lb.

Weight Limit 300 Pounds	
Invacare recommends all wheelchairs be fitted and programmed to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.	
NOTE: All weight limits are inclusive of add-on equipment and seat systems.	

Footrest Extension Range			
	AT5043	AT5543	1651
60HD	N/A	12 1/4" – 16 1/2"	13 1/2" – 17 1/2"
SWING914	N/A	7 3/4" – 12 1/2"	9 1/2" – 14 1/2"
70MFX	6 3/4" – 12 1/2"	4 1/2" – 10"	N/A
70NHD	N/A	10 3/4" – 15 3/4"	11 3/4" – 17"
70STAPER	N/A	12 1/2" – 16 1/2"	13 1/4" – 17 1/4"
90SWING	4" – 9"	5" – 8"	N/A
Measured from Seat Pan to Back of Footplate			

Torque 3 Specifications		
* Overall Width	25" 27.5"	w/14"x3" wheel w/14"x4" wheel
**Overall Length	28.5" 31.5"	Short Base Long Base
***Seat-to-Floor	17.5" Standard 19.75" Optional 21" Optional	Low Medium Tall
Battery	22NF deep cycle; gel/sealed (Not supplied)	
Shipping Weight 266 lb. (18" X 18" seat w/Group 24 batteries)		
*Measurement based on wheel base. Seat width is width picked plus 7.5" for the armrest and joystick.		
**Measured from headtube to anti-tipper.		
***Based on 5° tilt and 16" seat depth; add 1/4" for additional 2" of seat depth		
Note: All specifications and dimensions are approximate.		

Manual Recline Specifications			
Seat Widths	14"-22" (1" increments)		
Seat Depths	16"-22" (1" increments) (Adjustable 16"-19", Fixed 20"-22")		
Back Heights	24" Standard, 18 1/2", 20", 22", 24", 26"		
Weight Limits	300 lb. with Standard Cylinders		
Headrest	10" Height		
Recline Range	90° to 170° (Approximate)		
Overall Back Height (w/o extension)	Low Seat	Med. Seat	High Seat
Standard Base	41 1/2"	43 1/2"	45 1/4"
(Measured with 5° Tilt and 10" for Extension. Field adjustable for 10° of total pre-tilt.)			

Contoura Seat Sizes							
	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X
W24	*	*	X	X	X	X	X
*Not available with Contoura Seat Cushion.							
Back Heights	25" or 29.5"						

Memory Card - BASIC

- Standard with ALL Expandable Control Systems (MPJ+, PSR+, PSF+, MK6i Display).
- Used to Back Up/Restore Programmed Settings for Only One Chair.
- May be given to user or kept with provider's files for safe keeping.
- Does not Contain Advanced Diagnostics/Help Screens/File Structure.
- Not compatible with SPJ+ Joysticks.

Memory Card - PROFESSIONAL

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains Advanced Diagnostics, Help Screens, Software Updates, File Storage/Retrieval.
- Used to back up multiple chairs and programming settings.
- Intended for use by qualified providers only.
- Not compatible with SPJ+ Joysticks

Invacare Corporation

www.invacare.com

USA

All rights reserved. Trademarks are identified by the symbols™ and ©. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

One Invacare Way

Elyria, OH USA

44035-4190

800-333-6900

Monroe is a registered trademark of Tenneco Automotive, Inc.

Therafin is a registered trademark of W/R Medical Electronics Co.

Products designated by ASL are Manufactured by Adaptive Switch Labs, Inc.

Specifications and prices are subject to change without notification.

©2008 Invacare Corporation

