



Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.

Invacare® Formula™ CG Tilt Only for TDX® SR Base

Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



TRANSACTION TYPE:

Order Quote

Date _____ Account # _____ P.O. # _____

Company Name _____ Phone # _____

Address _____ City _____ State _____ Zip Code _____ Country _____

Contact Name _____ Comments _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

POWERED SEATING BASE

- TDXSR-CG TDXSR for Formula CG with MK6i™1 \$9,950.00 K0877
TDXSR-CG comes standard with: 300 lbs. user weight capacity, 7 MPH top speed, GB TrueTrack motors, MK6TT Controller and Quiet Stability Lock. Base width is 25 1/2".
- TDXSR-CG-HD TDX SR HD for Formula CG \$10,500.00 K0879
TDXSR-CG-HD comes standard with: 400 lbs. user weight capacity, 7 MPH top speed, GB TrueTrack motors, MK6TT Controller and Quiet Stability Lock. Base width is 25 1/2".

BASE SELECTIONS

OMIT SEAT FRAME

OSF Omit Seat Frame for Base Only No Charge

TRANSPORT TIE DOWN

TRBKTS Wheelchair Transport Brackets \$250.00

Group 4 - K0879 Basic Package	
TDXSR-CG-HD TDXSR HD for Formula CG with MK6i	
400 lbs. Weight Capacity	\$400.00
1431-3 14" x 3" Tires with Foam Filled Inserts	Standard
Gearless Brushless Motors	Standard
Flip Back Adjustable Height Arm	Standard
MK6TT MK6 TrueTrack Controller	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard
Seat Option CGT Tilt (E1002)	\$4,590.00

TRANSPORT TIE-DOWN*

TRRO includes four factory installed wheelchair transport brackets. TRRO has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport. As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

* Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

Group 4 - K0877 Basic Package	
TDXSR-CG TDXSR for Formula CG with MK6i	
300 lbs. Weight Capacity	
1431-3 14" x 3" Tires with Foam Filled Inserts	No Charge
Gearless Brushless Motors	Standard
Flip Back Adjustable Height Arm	Standard
MK6TT MK6 TrueTrack Controller	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard
Seat Option CGT Tilt (E1002)	\$4,590.00

USER WEIGHT LIMITS

- U300 User Weight <300 lbs. No Charge
- U400 User Weight <400 lbs. \$400.00 K0870

WHEEL OPTIONS

1431-3 14" x 3" Tires with Foam Filled Inserts Standard

FORK OPTION

FKPKG One-Sided Fork Package Standard

**NOTE: All specifications and dimensions are approximate.
 For the most current pricing information go to www.invacare.com.**

Invacare® Formula™ CG Tilt Only for TDX® SR Base

FRAME FINISH

CONTEMPORARY

- 24P Wet Black No Charge
- 60P Silver Metallic No Charge
- 71P Silver Vein No Charge
- 105P Midnight Blue No Charge
- 115P Black Prism No Charge

TOO HOT

- 30P Sunny Yellow No Charge
- 61P Electric Red No Charge
- 119P Bubblegum Pink No Charge
- 120P Grape Madness No Charge

WAY COOL

- 62P Electric Blue No Charge
- 104P Electric Teal No Charge
- 121P Lolly Pop Blue No Charge
- 122P Cosmic Blue No Charge
- 125P Emerald Green No Charge
- 127P Grasshopper Green No Charge

SEAT

SYSTEM TYPE

- CGT CG Tilt Only System \$4,590.00 E1002
- SMTACT Smart Actuator Standard

SEAT-TO-FLOOR HEIGHT

Nominal measurement at front edge of seat pan to floor on 18" deep seat.

- LOW 17 1/2" at 0° No Charge
- MED 18 1/2" at 0° No Charge
- TALL 19 1/2" at 0° No Charge

SEAT SELECTION

- BPO Back Posts Only - Round Canes Only No Charge
- CNB Conventional Back No Charge
- CURV Crurved Back No Charge E1226
- CTB Contoura Seating \$625.00 E2620 E2621
- OMCU Omit Seat Cushion Less \$20.00
- CTBMSH Contoura seating with mesh upholstery \$635.00 E2620 E2621
- OMCU Omit Seat Cushion Less \$20.00

ADJUSTABLE SEAT WIDTH RANGE

- 1216W 12" to 16" Seat Width Range No Charge
Seat Widths 12" - 15" and Seat Depths 12" - 15" are only available with Back Post Only (BPO) or Conventional Seating. For users under 150 lbs.
- 1620W 16" to 20" Seat Width Range No Charge
- 2024W 20" to 24" Seat Width Range No Charge

SEAT WIDTH SETTING

Clear

- W12 Seat Width 12" (Short Base) No Charge
Available only with TDXSR-CG-HD.
- W13 Seat Width 13" (Short Base) No Charge
Available only with TDXSR-CG-HD.
- W14 Seat Width 14" (Short Base) No Charge
Available only with TDXSR-CG-HD.
- W15 Seat Width 15" (Short Base) No Charge
Available only with TDXSR-CG-HD.
- W16 Seat Width 16" (Short Base) No Charge
- W17 Seat Width 17" (Short Base) No Charge
- W18 Seat Width 18" (Short Base) No Charge
- W19 Seat Width 19" (Long Base) No Charge
- W20 Seat Width 20" (Long Base) No Charge
- W21 Seat Width 21" (Long Base) No Charge
- W22 Seat Width 22" (Long Base) No Charge
- W23 Seat Depth 23" (Long Base) No Charge
Available only with TDXSR-CG-HD.
- W24 Seat Width 24" (Long Base) No Charge
Available only with TDXSR-CG-HD.

ADJUSTABLE SEAT DEPTH RANGE

Clear

- 1218D 12" to 18" Seat Depth Range No Charge
Available only with TDXSR-CG-HD.
- 1619D 16" to 19" Seat Depth Range No Charge
- 1922D 19" to 22" Seat Depth Range No Charge

CONTOURA SEAT DEPTH RANGE

Clear

- 1619D 16" to 19" Seat Depth Range No Charge
- 1922D 19" to 22" Seat Depth Range No Charge

SEAT DEPTH SETTING

Clear

- D12 Seat Depth 12" (Short Base) \$110.00
Available only with TDXSR-CG-HD.
- D13 Seat Depth 13" (Short Base) \$110.00
Available only with TDXSR-CG-HD.
- D14 Seat Depth 14" (Short Base) \$110.00
Available only with TDXSR-CG-HD.
- D15 Seat Depth 15" (Short Base) No Charge
Available only with TDXSR-CG-HD.
- D16 Seat Depth 16" (Short Base) No Charge
- D17 Seat Depth 17" (Short Base) No Charge
- D18 Seat Depth 18" (Short Base) No Charge
- D19 Seat Depth 19" (Long Base) No Charge
- D20 Seat Depth 20" (Long Base) No Charge
- D21 Seat Depth 21" (Long Base) \$529.00
- D22 Seat Depth 22" (Long Base) \$529.00

SEAT TILT

- 0SA 0° Standard
- 5SA 5° No Charge

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information go to www.invacare.com.**

SEAT Continued

BACK HEIGHT

Clear

- TSBH16 16" Back Height No Charge
- TSBH17 17" Back Height No Charge
- TSBH18 18" Back Height No Charge
- TSBH19 19" Back Height No Charge
- TSBH20 20" Back Height No Charge
- TSBH21 21" Back Height No Charge
- TSBH22 22" Back Height No Charge
- TSBH23 23" Back Height No Charge
- TSBH24 24" Back Height No Charge
- TSBH25 25" Back Height No Charge
- TSBH26 26" Back Height No Charge

CONTOURA BACK HEIGHT

Clear

- TSBH22 22" Back Height No Charge
- TSBH25 25" Back Height No Charge
- TSBH28 28" Back Height No Charge

Back Height Chart	
Back Type	Back Height
CTB - Contoura	22", 25" and 28"
BPO - Back Posts Only	16"-24"
CURV - Curved	16"-26"
Conventional Back	16"-20" on seat widths of 12"-15"
Note: Back Heights of 21"-26" not available with seat widths of 12"-15"	

ADJUSTABLE ARM TYPE (PAIR)

Clear

- 19F Full Length Adj Hgt Conv. No Charge
- 29F Desk Length Adj Hgt Conv. No Charge

MIXED ADJUSTABLE ARM TYPE (EACH)

All arms are flip-up and detachable.

Right Arm

Clear

- 19FR Full Length Adj. Hgt. Right No Charge
- 29FR Desk Length Adj. Hgt. Right No Charge

Left Arm

Clear

- 19FL Full Length Adj. Hgt. Left No Charge
- 29FL Desk Length Adj. Hgt. Left No Charge

CANTILEVER ARM TYPE (PAIR)

Clear

Cantilever Arms only available with 300 lbs. rated systems. Not available with BPO, Otto Bock, or 12"-15" seat width. Stealth Laterals part number's TWB-CPR and TWB-CPRL are the only lateral support compatible with the current Cantilever Arm design. Must choose WF armpad type with Stealth Laterals.

- CLVF Full Length Arm \$350.00
- CLVD Desk Length Arm \$350.00

CANTILEVER MIXED ARM TYPE (EACH)

Right Arm

Clear

- CLVFR Full Length Arm Right \$175.00
- CLVDR Desk Length Arm Right \$175.00

Left Arm

Clear

- CLVFL Full Length Arm Left \$175.00
- CLVDL Desk Length Arm Left \$175.00

ARMPAD OPTIONS (CANTILEVER ARMS)

Clear

- WF Waterfall Armpads Standard
- RCL Reclining Armpads No Charge

SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Positioning Strap \$50.00
Not available with 12" - 15" seat widths.
- 1515 Push Button Style Seat Positioning Strap Standard

CHEST POSITIONING STRAPS

Clear

- 7311BK Airline Buckle Seat Positioning Strap \$100.00 E0960
- 9515 Push Button Style Chest Positioning Strap \$100.00 E0960
- 7321BK Hook & Loop Chest Positioning Strap \$100.00 E0960

POSITIONING ACCESSORIES

HEADREST TYPE (OPTIONAL)

Clear

Hardware is packaged with equipment and cannot be ordered separately.

- TAR127700 Curved Headrest \$210.00 E0955
TAR27700HW Removable & Adjustable Hardware \$210.00 E1028
- TAR128000 Two-Step Headrest \$210.00 E0955
TAR128000HW Removable & Adjustable Hardware \$210.00 E1028
- T0022400 Large Headrest \$210.00 E0955
T0022400HW Removable & Adjustable Hardware \$210.00 E1028
- T0022600 Small Headrest \$210.00 E0955
T0022600HW Removable & Adjustable Hardware \$210.00 E1028

STEALTH HEADRESTS

Clear

Hardware is packaged with equipment and cannot be ordered separately.

- CP175 14" Multi Axis Headrest \$225.00 E0955
CP175HW Removable & Adjustable Hardware \$225.00 E1028
- CP275 10" Multi Axis Headrest \$225.00 E0955
CP275HW Removable & Adjustable Hardware \$225.00 E1028
- CP180 14" TWB Headrest \$225.00 E0955
CP180HW Removable & Adjustable Hardware \$225.00 E1028
- CP280 10" TWB Headrest \$225.00 E0955
CP280HW Removable & Adjustable Hardware \$225.00 E1028

OTTO BOCK ARM PADS (PAIR)

Clear

Not compatible with retractable joystick mount. If selecting an Otto Bock hand pad a large, medium or small forearm pad must be selected.

- OB20400 One Piece Channel Arm pads \$230.00 E2209
- OB20500 Large Forearm Pads \$235.00
- OB20600 Medium Forearm Pads \$235.00
- OB20700 Small Forearm Pads \$235.00

OTTO BOCK HAND PADS (PAIR)

Clear

Not compatible with retractable joystick mount. If selecting an Otto Bock hand pad a large, medium or small forearm pad must be selected.

- OB23000 Large Flat Hand Pads \$100.00
- OB23300 Medium Flat Hand Pads \$100.00

ADDUCTORS

Clear

Hardware is packaged with equipment and cannot be ordered separately.

- ADD Adductors Small \$200.00 E0956
ADDHW Hardware for ADD \$400.00 E1028
- ADDL Adductors Large \$200.00 E0956
ADDHWL Hardware for ADDL \$400.00 E1028

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information go to www.invacare.com.**

Invacare® Formula™ CG Tilt Only for TDX® SR Base

LATERALS (Pair)

CONVENTIONAL BACK/BPO/CURVED BACK Clear

Hardware is packaged with equipment and cannot be ordered separately.

- TCSLM Clamp-On Swingaway Laterals Medium \$210.00 E0956
TCSLMHW Removable and Adjustable Hardware \$210.00 E1028
- TCSLL Clamp-On Swingaway Laterals Large \$210.00 E0956
TCSLLHW Removable and Adjustable Hardware \$210.00 E1028

STEALTH LATERALS FOR CONVENTIONAL BACK Clear

Hardware is packaged with equipment and cannot be ordered separately.

- TWB-RPR Adj. Swing Contoura Lat. Med \$320.00 E0956
RPRHW Removable & Adjustable Hardware \$320.00 E1028
- TWB-CRPL Adj. Swing Contoura Lat. Lg \$320.00 E0956
RPRLHW Removable & Adjustable Hardware \$320.00 E1028

LATERALS FOR CONTOURA SEATING Clear

When using Stealth laterals the arm will not swing past the lateral support. Will need to remove arm for transfer. For use with Contoura seating only. Hardware is packaged with equipment and cannot be ordered separately.

- CSLM Swingaway Contoura Lat. Med \$210.00 E0956
CSLMHW Removable and Adjustable Hardware \$210.00 E1028
- CSLL Swingaway Contoura Lat. Lg \$210.00 E0956
CSLLHW Removable and Adjustable Hardware \$210.00 E1028

STEALTH LATERALS FOR COUNTOURA Clear

When using Stealth laterals the arm will not swing past the lateral support. Will need to remove arm for transfer. For use with Contoura seating only. Hardware is packaged with equipment and cannot be ordered separately.

- TWB-CPR Adj. Swing Contoura Lat. Med \$320.00 E0956
CPRHW Removable & Adjustable Hardware \$320.00 E1028
- TWB-CPRL Adj. Swing Contoura Lat. Lg \$320.00 E0956
CPRLHW Removable & Adjustable Hardware \$320.00 E1028

ELECTRONICS

MK6i™ CONTROLLERS

- MK6TT TrueTrack Controller Standard

EXPANDABLE CONTROLLER SYSTEM Clear

- PWH Harness req. for Expandable Electronics \$350.00 E2313
Required to interface joystick or MK6i display with expandable controllers.
- CMPJM6U MPJ+ Multiple Drive Joystick with **Color Display** .. \$750.00 E2377
- PSFM6 Personal Joystick, Inductive on Front \$1,130.00 E2377
TOD-SS: Toggle On/Off Drive Select, Speed Potentiometer Standard.
- PSRM6 Personal Joystick, Inductive on Rear \$1,130.00 E2377
TOD-SS: Toggle On/Off Drive Select, Speed Potentiometer Standard.
- PBOD Push Button On/Off with PB Drive Select No Charge
PSRM6 Only.
- PBSS Push Button On/Off with Speed Select No Charge
Drive select through mode switch, or mode switch and joystick. PSRM6 Only.
- OR On/Off Mounted Right No Charge
PSRM6 Only.
- OL On/Off Mounted Left No Charge
PSRM6 Only.
- ASLRDYM6 MK6 Display Only (No Driver Control) \$750.00 E2377
Display only, No Driver Control. Allows choice from order form.

FOR YOUR ASL DRIVE SYSTEM NEEDS
PLEASE REFER TO THE
ASL SINGLE INVOICE PRICE LIST FORM # 02-097.

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS Clear

- 1500M6 Prop RIM Head Control with MK6i Display .. \$2,875.00 E2327
 SMHD Small Headrest No Charge
 LGHD Large Headrest No Charge
- 1558MM6 Compact Joystick with Display \$1,130.00 E2373

ALTERNATIVE DIGITAL CONTROL SYSTEMS Clear

Must also use either a Multiple Drive Joystick or MK6i Display.

- SNPM6 Sip-N-Puff Digital Interface Box \$1,500.00 E2325
- 5018M6 Wafer Board and Digital Interface Box \$1,550.00 E2322
- 5020M6 Mini Tash Joystick & Digital Interface Box ... \$1,750.00 E2321

SIP-N-PUFF KIT

- PKG32666 Therafin® Sip-N-Puff Breath Tube Kit \$380.00 E2326

ADDITIONAL CONTROL CHOICES

- 1558M6 Compact Joystick - Less Std Mount \$770.00
- 5018 Wafer Board \$310.00
Requires SNPM6 Sip-N-Puff Interface Box.
- 5020 Mini Tash Joystick \$405.00
Requires SNPM6 Sip N Puff Interface Box.
- PACM6 Prop Attendant Control with quick release mount \$850.00
- 1552M Digital Attendant Control \$850.00
Requires SNPM6 Sip N Puff Interface Box.

JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retract. Joystick Mnt CMPJ+ \$279.00 E1028
- ARM250 Stealth Height Adjustable Joystick Mount \$210.00 E1028
- LEFT Left Handed Mounted Joystick No Charge
- RIMHW RIM Hardware For 1500 \$219.00 E1028
- PKG32669 Complete Bib Assembly (Tash/1558)..... \$265.00
- GATMPJ6 Gatlin Midline for CMPJ+ \$600.00 E1028
- GAT1812 Gatlin Midline for 1812 \$600.00 E1028
- RCM Rear Cane Mount for Compact Joystick \$395.00

ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer w/Pro Memory Card \$475.00
- 1813SD Pro Memory Card (USBReady) \$90.00
- 1560 "T" Handle Flexible Joystick Ext. \$77.00 E2323
- 1561 Straight Handle Flexible Joystick Ext. \$77.00 E2323
- 1826 Chin Cup \$48.00 E2324
- A24VPS 24V MK6 Auxiliary Power Source \$200.00
- AUX12M6 Auxiliary Module for outputs 1 & 2 \$795.00
- AUX34M6 Auxiliary Module for outputs 3 & 4 \$795.00
- 2PB Dual Push Button \$147.00
Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options.
- 2SW Dual Toggle \$172.00
Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options.
- PBMT Push Button Mounting Kit \$252.00
- MPY Monoport "Y" Cable for CMPJ+ or Display \$15.00
- RFM Wireless (RF) Mouse Only Module \$750.00
RFM will be available in September, 2010. Allows wireless computer mouse operation through the driver control. Includes two button stereo switch (with hook/loop mount) for mouse clicks. Not compatible with SPJ+ joysticks.

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information go to www.invacare.com.**

ELECTRONICS Continued

POWERED SEATING ACCESSORIES

- Clear
- SANODE One Actuator Operation thru Driver Cntrl..... \$1,300.00 E2310
Required for Tilt Operation through Driver Control - Single Actuator Systems Only. Not available with SPJ+ Joystick. Does not come with Egg switch.
- 4WSB Four Way Switch Box No Charge
 - Clear
 - FWT 4 Way Toggle No Charge
 - QPB 4 Push Buttons No Charge
- S4WSB Multiple Actuator Interface Box \$2,600.00 E2311
 - Clear
 - FWT 4 Way Toggle No Charge
 - QPB 4 Push Buttons No Charge

Systems

For systems with two or more actuators, NOT operating through driver control, choose 4WSB with either 4-way toggle or 4-push button.

For systems with 2 or more actuators OPERATING THROUGH DRIVER CONTROL, choose S4WSB. Uses 4-way toggle or 4-push button as attendant control.



ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Switch - Black \$156.00
- ASL304A Light Wobble Switch \$180.00
- ASL202J Fiber Optic Switch \$732.00
Must add MK6 Auxiliary Power Source.
- ASL208J Adj Proximity Switch \$660.00
Must add MK6 Auxiliary Power Source.
- ASL504M6 Remote Emergency Stop Switch \$1,185.00
- ASL611 Adjustable Switch Mounting Bracket \$252.00

BATTERY TRAY TYPE

24TRY and 22TRY3 base width 25 1/2"

- 22TRY3 22NF Style Tray for 3 Batteries Standard
Required for use with ventilator tray.
- 24TRY 24 Group Style Tray No Charge

BATTERIES

- Clear
- 22NFBATT3 22NF Battery (On Chair) \$450.00
Available with Vent Option.
- 24BATTERY 24Gel Battery (On Chair) \$410.00 E2363
- 22NFTRY3BAT2 22NF Battery 2 Batteries (On Chair) \$350.00 E2361
Use when adding aftermarket vent tray.

BATTERY CHARGER OPTION

- Clear
- 110CHARGER 110 Volt Battery Charger No Charge
- LESTR Lester 8 AMP Battery Charger \$320.00 E2367
- L8900 Omit Battery Charger Less \$80.00

VENT TRAYS

- VENT Articulating Vent Tray \$1,275.00 E1030

MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks \$45.00

FRONT RIGGINGS

HANGER BRACKET TYPE

- Clear
- ASB Center Pivot Style Bracket No Charge
- PSB Pin Style Bracket No Charge

CENTER MOUNT

- Clear
- CTMT-F Center Mount Adj. Knee Angle (12 1/4" - 15") \$300.00
For a 14"-18" extension range, available through special #SM-071217095848 with Tall STF Only and 5 Degree Seat-Tilt Angle.
Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 12 1/2" - 15". Footplates are adjustable angle. Must choose either SM, MED or LRG footplate.
 - Clear
 - SM Small Footplates (5 1/4" x 7 1/2") No Charge
 - MED Medium Footplates (5 1/2" x 9 1/2") No Charge
 - LRG Large Footplates (5.90" x 11 1/2") No Charge
 - CALF Calf Pads For Center Mount \$200.00
- CTMT-FJ Jr. Center Mount Adj. Knee Angle (9 1/4" - 12") \$300.00
Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 9 1/4" - 12". Footplates are adjustable angle. Available on low seat-to-floor only. Must choose either SM, MED or LRG footplate.
 - Clear
 - SM Small Footplates (5 1/4" x 7 1/2") No Charge
 - MED Medium Footplates (5 1/2" x 9 1/2") No Charge
 - LRG Large Footplates (5.90" x 11 1/2") No Charge
 - CALF Calf Pads For Center Mount \$200.00

CENTER PIVOT STYLE (PAIR)

- Clear
- 70NHD 70° Swingaway Footrests 13 1/2" - 17 1/2" Standard
With seat widths of 21" to 24" and depth of 16", must pick footplate 1350. The AT5543 is not available.
- 70STAPER 70° Tapered Footrests 13"-17" No Charge
Standard with impact guards.
- AT5544 Manual Elevating Legrests \$255.00 E0990
Not available with 16" wide by 16" to 18" deep or 17" wide by 16" deep and low seat-to-floor.

PIN STYLE-HEMI (PAIR)

- Clear
- PHW93 70° Swingaway Footrests 13 1/2" - 19" No Charge
- PHAL4A Manual Art./ELR 17 1/4" - 21 1/2" \$225.00K0553
Not available on seat widths 14" and 15".

POWER ELEVATING LEGRESTS - ARTICULATING

- Clear
- MPCTMT Power Center Mount Foot Platform \$2,995.00 K0108
Includes a standard, single footboard.
 - IFP Individual Footplates \$275.00
- MPCTMT-L Power Center Mount Foot Platform Long \$2,995.00
Includes a standard, single footboard.
 - IFP Individual Footplates \$275.00
- ELRPW Power ELR with Articulation 20" Ext. \$1,895.00 E1010
Not available with SPJ+ Joystick. Independent powered elevating legrests including elevating and articulating legrest with angle adjustable footplates. Not available with the 16" or 17" seat depths. Extension range 13" - 20". Sold as pair. Not available with Sanode.
- ELRPW22 Power ELR with Articulation 22" Extension .. \$1,995.00
Not available with SPJ+ Joystick. Independent powered elevating legrests including elevating and articulating legrest with angle adjustable footplates. Not available with the 16" or 17" seat depths. Extension range 15" - 22". Sold as pair. Not available with Sanode.
 - IFP Individual Footplates \$275.00

NOTE: All specifications and dimensions are approximate. For the most current pricing information go to www.invacare.com.

FRONT RIGGINGS Continued

FOOTPLATE OPTIONS (PAIR)

Clear

- 1651 Flip-Up Composite Footplates No Charge
Not available with AT5544 and 17" wide by 16" deep and low seat-to-floor.
- 1350 Extra Large Aluminum Footplates No Charge
Not available with AT5544 and 16" wide by 16"-18" deep or 17" wide by 16" deep and low seat-to-floor.
- AT5543 Adjustable Angle Flip-Up Footplates \$165.00 K0040
- CTMTFB Footboard for Center Mount No Charge
One piece footboard for Center Mount (CTMT-F) Riggings.

FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap \$27.00 E0038
- IMPACTG Impact Guards \$51.00
- 1600BK Heel Loops with Ankle Straps \$56.00 E0951
- LPT2 Longer Pivot Slide Tube For 1350 \$83.00
- ALPT Longer Pivot And Slide Tube \$84.00
- ALPT4 Longer Pivot And Slide Tube 4" \$84.00

Power Center Mount KTH	
MPCTMT	13" - 17"
MPCTMT-L	15" - 19"

Seating and Positioning

For your seating needs, please refer to the following rehab seating and positioning products: Price List Form No. 07-027

Cantilever Arm Height Adjustment			
Back Type	Height	Adjustment	Adjustment with laterals
CURV	20" or greater	11" to 16"	11" to 14"
CURV	19"	11" to 15"	11" to 13"
CURV	18"	11" to 14"	11" to 12"
CURV	17"	11" to 13"	no adjustment
CURV	16"	11" to 12"	not available
Contoura	22", 25" or 28"	12 1/2" to 14 1/2"	no adjustment

Contoura Seat Sizes							
	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X
W24	*	*	X	X	X	X	X
*Not available with Contoura Seat Cushion							
Back Heights	22", 25" or 28"						

Conventional or BPO Seat Sizes	
Seat Widths	12" - 15"
Seat Depths	12" - 22" (BPO) 12" - 18" (Conventional)
Back Heights	16" - 20" for 12" - 15" (Seat Widths Only)

Curved or BPO Seat Sizes	
Seat Widths	16" - 24"
Seat Depths	16" - 22"
Back Heights	16" - 26"

Memory Card - BASIC

- Standard with ALL Expandable Control Systems (CMPJ+, PSR+, PSF+, MK6i Display).
- Used to back up / restore programmed settings for only one chair.
- May be given to user or kept with provider's files for safe keeping.
- Does not contain advanced diagnostics / help screens / file structure.
- Not compatible with SPJ+ Joysticks.

Memory Card - PROFESSIONAL

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains advanced diagnostics, help screens, software upgrades, file storage / retrieval
- Used to back up multiple chairs and programming settings.
- Intended for use by qualified providers only.
- Not compatible with SPJ+ Joysticks.

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information go to www.invacare.com.**

Invacare® Formula™ CG Tilt Only for TDx® SR Base

Formula CG Tilt Only Features and Specifications

- Weight capacity of 300 lbs. or 400 lbs. (with Heavy-Duty 4-Pole Motors).
- Curved, Conventional or Contoura style seating.
- Adjustable seat frame.

Tech Tips

Operating Powered Seating through the Driver's Control.
All chairs require adding or replacing one component to allow actuator operation through the driver control.

SPJ+ with PSS:

Replace joystick with SPJ+ with ACC (Actuator Control) - E2310.

Single Actuator Systems with CMPJ+, PSR+, PSF+ MK6i Display:

Add a SANODE (Single Actuator Node) to the system - E2310.

Multiple Actuator Systems:

Replace the standard 4 Way Switch box with a S4WSB (Multiple Actuator Interface Box) - E2311.

Comments

NOTE: All specifications and dimensions are approximate.
For the most current pricing information go to www.invacare.com.



Invacare Corporation

USA

One Invacare Way

Elyria, OH USA

44035-4190

800-333-6900

www.invacare.com

©2010 Invacare Corporation

All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.