



Yes, you can.®

## Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email\* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

\*Saving and emailing the eForms is available with Acrobat Professional software

### If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

**Please note:** The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

### If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit [www.adobe.com](http://www.adobe.com).



# Invacare® TDX® SC Rehab Seat Power Wheelchair

## Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



Yes, you can.®

### TRANSACTION TYPE:

☐ Order ☐ Quote

Date \_\_\_\_\_ Account # \_\_\_\_\_ P.O. # \_\_\_\_\_

Company Name \_\_\_\_\_ Phone # \_\_\_\_\_

Address \_\_\_\_\_ City \_\_\_\_\_ State \_\_\_\_\_ Zip Code \_\_\_\_\_ Country \_\_\_\_\_

Contact Name \_\_\_\_\_ Comments \_\_\_\_\_

### Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

### CHAIR TYPE

☐ TDXSC TDX SC with MK6i™ ..... \$5,720.00 K0848

**TDXSC comes standard with: 250 lbs. user weight capacity, 4.8 MPH top speed, 2-pole motor package and MK6i SPJ+ Joystick, MK660ACC Controller, and Quiet Stability Lock. Base width is 24".**

☐ TDXSCSEAT TDXSC Seat Only ..... \$755.00

#### Group 3 - K0848 Basic Package

TDX SC with MK6i and Adjustable Rehab Seat	
250 lbs. Weight Capacity	
DW1053FF 10 1/2" x 3 1/2" Flat Free Tires	Standard
2 Pole Motors	Standard
Flip Back Fixed Height Arm	Standard
MK660ACC 60 Amp Controller	Standard
SPJ+ Single Drive Joystick	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard

#### Transport Tie Down\*

TRRO includes four factory-installed transport brackets and a wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNA WC Vol 1. Section 19 Frontal Impact Test requirements for wheelchairs with a 168 lbs. crash dummy, which corresponds to a person with a weight of 114 lbs. to 209 lbs.

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

\* Transport Tie-Down option is not retrofittable to existing models and is not field serviceable.

#### TRRO Features for TDX SC

##### Required Features for TRRO

- TDXSC
- 70° Swingaway, 70° Tapered (70NHD or 70STAPER) or Manual Center Mount Front Riggings
- Any Arm Style
- Adjustable Rehab Seat with Adjustable Seat Tilt & Back Angle
- Sectional or Straight Rear Uprights or No Push Handles
- 10 1/2" x 3 1/2" Tires with Foam Filled Inserts
- 90° or 95° Back Angle
- 18" - 24" ASBA Back Height
- Any Seat Width
- Any Seat Depth
- 250 lb user weight limit
- Low Seat-to-Floor Height
- Any Color Frame Finish

##### Optional Features for TRRO

- Adjustable Angle Flip Up Footplates, Flip-Up Composite Footplates or Aluminum Footplates-Deluxe
- Wheel Locks
- Calf Strap
- Heel Loops with Ankle Straps
- Impact Guards
- One Sided Fork Package
- SPJ+, CMPJ+ Joystick
- All PSR+ / PSF+ Joysticks
- Proportional Attendant Control
- "T" Handle Flex. Joystick Ext.
- Left Hand Mounted Joystick
- Straight Handle Flexible Joystick Ext.
- Easy Remote Programmer
- Reinforced Back Upholstery
- U1 Batteries

Please note that all TRRO chairs are set at 5° seat angle

### BASE SELECTIONS

#### TRANSPORT TIE-DOWN

☐ Clear

Not available with elevating seat.

- ☐ TRRO Transport Ready Option ..... \$300.00
- ☐ TRBKTS Wheelchair Transport Brackets ..... \$250.00

#### USER WEIGHT LIMIT

- ☐ UW250 User Weight <250 lbs. .... Standard

#### WHEEL OPTION (PAIR)

- ☐ DW1053FF 10 1/2" x 3 1/2" Flat Free Tires ..... Standard

#### FORK OPTION

- ☐ FKPKG One-Sided Fork Package ..... \$500.00

**NOTE: All specifications and dimensions are approximate.**  
**For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**



# Invacare® TDX® SC Rehab Seat Power Wheelchair

## FRAME FINISH

- ☐ 24P Wet Black ..... No Charge  
☐ 61P Electric Red ..... No Charge  
☐ 62P Electric Blue ..... No Charge  
☐ 60P Silver Metallic ..... No Charge  
☐ 115P Black Prism ..... No Charge  
☐ 120P Grape Madness ..... No Charge  
☐ 127P Grasshopper Green ..... No Charge

## SEAT

### SYSTEM TYPE

- ☐ CGE 5" Elevating Seat ..... \$1,895.00

### SEAT-TO-FLOOR HEIGHT

Due to variation in options the listed seat-to-floor heights have a tolerance of +/- 0.5"

- ☐ LOW 16.5" at 0° ..... Standard  
☐ MED 17.5" at 0° ..... No Charge  
☐ TALL 18.5" at 0° ..... No Charge

### SEAT SELECTION

- ☐ ADJASBA Adj Rehab Seat with Adj Seat Tilt and Back Angle ... Standard  
☐ ADJRECL Manual Recline for Adj ASBA ..... \$744.00 E1226  
☐ CTB Contoura™ Seating ..... \$625.00 E2621 E2620  
☐ CTBMSH Contoura™ seating with mesh upholstery ..... \$635.00 E2621 E2620  
☐ OMCU Omit Contoura Seat Cushion ..... Less \$20.00  
☐ CURV Curved ..... No Charge  
☐ OMCU Omit Contoura Seat Cushion ..... Less \$20.00

### ADJUSTABLE SEAT WIDTH RANGE

- ☐ 1620W 16" to 20" Seat Width Range ..... No Charge

### ADJASBA SEAT WIDTH

☐ Clear

Contoura available in even Seat Widths from 16" - 20".

- ☐ W16 16" Seat Width ..... No Charge  
☐ W17 17" Seat Width ..... No Charge  
☐ W18 18" Seat Width ..... No Charge  
☐ W19 19" Seat Width ..... No Charge  
☐ W20 20" Seat Width ..... No Charge

### ADJUSTABLE SEAT DEPTH RANGE

☐ Clear

Contoura available in Seat Depth from 16" - 20".

- ☐ 1619D 16" to 19" Seat Depth Range ..... No Charge  
☐ 1920D 19" to 20" Seat Depth Range ..... No Charge

### ADJASBA SEAT DEPTH

☐ Clear

- ☐ D16 16" Seat Depth ..... No Charge  
☐ D17 17" Seat Depth ..... No Charge  
☐ D18 18" Seat Depth ..... No Charge  
☐ D19 19" Seat Depth ..... No Charge  
☐ D20 20" Seat Depth ..... No Charge

### FIXED ARM TYPE (PAIR)

☐ Clear

- ☐ 18F Full Length Fix. Hgt. Conv.Arms ..... No Charge  
☐ 28F Desk Length Fix. Hgt. Conv. Arms ..... Standard

### ADJUSTABLE ARM TYPE (PAIR)

☐ Clear

- ☐ 19F Full Length Adj Hgt Conv. .... \$230.00 E0973  
☐ 29F Desk Length Adj Hgt Conv. .... \$230.00 E0973

## MIXED ADJUSTABLE ARM TYPE (EACH)

### Right Arm

☐ Clear

- ☐ 19FR Full Length Adj. Hgt. Right ..... \$115.00 E0973  
☐ 29FR Desk Length Adj. Hgt. Right ..... \$115.00 E0973

### Left Arm

☐ Clear

- ☐ 19FL Full Length Adj. Hgt. Left ..... \$115.00 E0973  
☐ 29FL Desk Length Adj. Hgt. Left ..... \$115.00 E0973

## CANTILEVER ARM TYPE (PAIR)

☐ Clear

Cantilever Arms only available with 300 lbs. rated systems. Not available with BPO, Otto Bock, or 12" - 15" seat width. Must use hard back. See chart on last page for arm height adjustment range.

- ☐ CLVF Full Length Arm ..... \$350.00  
☐ CLVD Desk Length Arm ..... \$350.00

## CANTILEVER MIXED ARM TYPE (EACH)

### Right Arm

☐ Clear

- ☐ CLVFR Full Length Arm Right ..... \$175.00  
☐ CLVDR Desk Length Arm Right ..... \$175.00

### Left Arm

☐ Clear

- ☐ CLVFL Full Length Arm Left ..... \$175.00  
☐ CLVDL Desk Length Arm Left ..... \$175.00

## ARMPAD OPTIONS (CANTILEVER ARMS)

☐ Clear

- ☐ WF Waterfall Armpads ..... Standard  
☐ RCL Reclining Armpads ..... No Charge

## SEAT TILT

5° seat-to-floor height not available with Elevating Seat and not available on LOW seat-to-floor height on non-elevating seat. Add 1/2" to Med or Tall STF at 5° seat tilt.

- ☐ 0SA 0° Seat Tilt ..... No Charge  
☐ 5SA 5° Seat Tilt ..... No Charge

## SEAT BACK TYPE

- ☐ SECT Sectional Rear Uprights ..... Standard  
☐ L8661 Omit Rear Uprights ..... Less \$20.00  
☐ SRUP Straight Rear Uprights ..... No Charge  
☐ L8661 Omit Rear Uprights ..... Less \$20.00  
☐ 9720P No Push Handles (Straight Canes) ..... \$82.00

## SEAT BACK ANGLE DEGREE

- ☐ 90BA 90° ..... Standard  
☐ 95BA 95° ..... No Charge

**NOTE: All specifications and dimensions are approximate.**  
**For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**



# Invacare® TDX® SC Rehab Seat Power Wheelchair

## BACK HEIGHT

### BACK HEIGHT ADJASBA

☐ Clear

- ☐ BH16 Back Height 16" ..... No Charge
- ☐ BH17 Back Height 17" ..... No Charge
- ☐ BH18 Back Height 18" ..... No Charge
- ☐ BH19 Back Height 19" ..... No Charge
- ☐ BH20 Back Height 20" ..... No Charge
- ☐ BH21 Back Height 21" ..... No Charge
- ☐ BH22 Back Height 22" ..... No Charge
- ☐ BH23 Back Height 23" ..... No Charge
- ☐ BH24 Back Height 24" ..... No Charge

### CONTOURA BACK HEIGHT

☐ Clear

- ☐ TSBH22 Back Height 22" ..... No Charge
- ☐ TSBH25 Back Height 25" ..... No Charge
- ☐ TSBH28 Back Height 28" ..... No Charge

### BACK HEIGHT RECLINER

☐ Clear

The above selections are used only with manual reclining seat selection.

- ☐ BHR18 Back Height 18" ..... No Charge
- ☐ BHR20 Back Height 20" ..... No Charge
- ☐ BHR22 Back Height 22" ..... No Charge
- ☐ BHR24 Back Height 24" ..... No Charge
- ☐ BHR26 Back Height 26" ..... No Charge

### SEAT POSITIONING STRAPS

- ☐ 1311BK Airline Buckle Seat Positioning Strap ..... \$50.00
- ☐ 1515 Push Button Style Seat Positioning Strap ..... Standard

### CHEST POSITIONING STRAPS

☐ Clear

- ☐ 7311BK Airline Buckle Seat Positioning Strap ..... \$100.00 E0960  
Not available with 12" - 15" seat widths.
- ☐ 9515 Push Button Style Chest Positioning Strap ..... \$100.00 E0960
- ☐ 7321BK Hook and Loop Chest Positioning Strap ..... \$100.00 E0960

## ELECTRONICS

### MK6I CONTROLLERS

- ☐ MK660ACC 60 Amp Controller ..... Standard
- ☐ GTRAC G-Trac Module ..... \$995.00 K0108  
G-Trac provides improved tracking capabilities and enhanced traction. Does not require special motors. Not compatible with SPJ+ joysticks.

### STANDARD DRIVER CONTROLS

☐ Clear

- ☐ SPJPLUS SPJ+ Single Drive Joystick ..... Standard
- ☐ SPJAP SPJ+ for Powered Seating ..... No Charge  
SPJAP for use with Elevating Seat only. Egg switch with 611 mount will come standard to operate elevate feature.

### ACTUATOR CONTROL THROUGH SPJ+

- ☐ SPJAPM SPJ+ with Joystick Actuator Control ..... \$1,300.00 E2310  
SPJAM for use with Elevating Seat only. Allows actuator operation through the drive control.

### EXPANDABLE CONTROLLER SYSTEM

☐ Clear

- ☐ PWH Harness req. for Expandable Electronics ..... \$350.00 E2313  
Required to interface joystick or MK6i display with expandable controllers.
- ☐ CMPJM6U MPJ+ Multiple Drive Joystick with Color Display .. \$750.00 E2377
- ☐ ASLRDYM6 MK6 Display Only (No Driver Control) ..... \$750.00 E2377  
Display only, No Driver Control. Allows choice from order form.
- ☐ MCD Mini Color Display ..... \$750.00 E2377  
No port for remote On/Off switch. Stereo port only.

FOR YOUR ASL DRIVE SYSTEM NEEDS  
PLEASE REFER TO THE  
ASL SINGLE INVOICE PRICE LIST FORM # 02-097.

### ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS

☐ Clear

- ☐ 1500M6 Prop RIM Head Control with MK6i Display .. \$2,875.00 E2327
- ☐ SMHD Small Headrest ..... No Charge
- ☐ LGHD Large Headrest ..... No Charge
- ☐ 1558MM6 Compact Joystick with Display ..... \$1,130.00 E2373

### ALTERNATIVE DIGITAL CONTROL SYSTEMS

☐ Clear

Must also use either a Multiple Drive Joystick or MK6i Display.

- ☐ SNPM6 Sip-N-Puff Digital Interface Box ..... \$1,500.00 E2325
- ☐ 5020M6 Mini Tash Joystick and Digital Interface Box ..... \$1,750.00 E2321

**NOTE: All specifications and dimensions are approximate.**  
**For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**



# Invacare® TDX® SC Rehab Seat Power Wheelchair

## ELECTRONICS Continued

### SIP-N-PUFF KIT

- ☐ PKG32666 Therafin® Sip-N-Puff Breath Tube Kit ..... \$380.00 E2326

### ADDITIONAL CONTROL CHOICES

- ☐ 1558M6 Compact Joystick - Less Std Mount ..... \$770.00
- ☐ 5020 Mini Tash Joystick ..... \$405.00  
**Requires SNPM6 Sip-N-Puff Interface Box.**
- ☐ PACM6 Prop Attendant Control with quick release mount .... \$850.00
- ☐ 1552M Digital Attendant Control ..... \$850.00  
**Requires SNPM6 Sip-N-Puff Interface Box.**

### JOYSTICK MOUNTING HARDWARE

- ☐ QLAM6 Quad Link Retract. Joystick Mnt CMPJ+ ..... \$279.00 E1028  
**Include ARM250 standard. Not available with TRRO**
- ☐ QLA Quad Link Retract. Joystick Mnt SPJ+ ..... \$279.00 E1028  
**Include ARM250 standard. Not available with TRRO**
- ☐ ARM250 Stealth Height Adj Joystick ..... \$210.00 E1028
- ☐ LEFT Left Handed Mounted Joystick ..... No Charge
- ☐ RIMHW RIM Hardware For 1500 ..... \$219.00 E1028
- ☐ PKG32669 Complete Bib Assembly (Tash/1558) ..... \$265.00
- ☐ GATMPJM6 Gatlin Midline for CMPJ+ ..... \$600.00 E1028
- ☐ GAT1812 Gatlin Midline for 1812 ..... \$600.00 E1028
- ☐ RCM Rear Cane Mount for Compact Joystick ..... \$395.00

### ELECTRONIC ACCESSORIES

- ☐ 1813M6 MK6i Programmer with Pro Memory Card .... \$475.00
- ☐ 1813SD Pro Memory Card (USB Ready) ..... \$90.00
- ☐ 1560 "T" Handle Flexible Joystick Ext. .... \$77.00 E2323
- ☐ 1561 Straight Handle Flexible Joystick Ext. .... \$77.00 E2323
- ☐ 1826 Chin Cup ..... \$48.00 E2324
- ☐ MPY Monoport "Y" Cable for CMPJ+ or Display ..... \$15.00  
**For use in the mode switch port of the MK6i Display, CMPJ+, PSR+, or PSF+ joysticks to provide two programmable switch functions (Monoport 1 and Monoport 2 of the Calibrations menu). Not compatible with SPJ+ joysticks.**
- ☐ AUX12M6 Auxiliary Module for 1 and 2 ..... \$795.00
- ☐ AUX34M6 Auxiliary Module for 3 and 4 ..... \$795.00
- ☐ 2PB Dual Push Button ..... \$147.00  
**Option to replace Monoport "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options. "Y" Cable standard when ordering powered seating.**
- ☐ 2SW Dual Toggle ..... \$172.00  
**Option to replace Monoport "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options. "Y" Cable standard when ordering powered seating.**
- ☐ PBMT Push Button Mounting Kit ..... \$252.00
- ☐ RFM Wireless (RF) Mouse Only Module ..... \$750.00  
**RFM includes Transmitter, USB Radio Frequency Mouse Dongle (Receiver), and Stereo Switch with hook and loop mount for mouse clicks. Not compatible with SPJ+.**
- ☐ IR Infrared Module with Wireless RF Mouse Control ..... \$995.00  
**Capable of operating up to six remote Infra-Red controlled devices. Also allows wireless (RF) Mouse operation through the driver control. Compatible with RFM module for control of two wireless mice from one chair. Includes two button stereo switch (with hook/loop mount) for mouse clicks. Not compatible with SPJ+ Joysticks.**

### POWER LEGS CONTROL MODULE

☐ Clear

**These options are located in the Configurator under the MK6i Electronics Folder. The Category Header inside the folder is labeled Power Seating Accessories**

- ☐ 4WSB Four Way Switch Box ..... No Charge
- ☐ Clear
- ☐ FWT 4 Way Toggle ..... No Charge
- ☐ QPB 4 Push Buttons ..... No Charge
- ☐ S4WSB Multiple Actuator Interface Box ..... \$2,600.00 E2311
- ☐ Clear
- ☐ FWT 4 Way Toggle ..... No Charge
- ☐ QPB 4 Push Buttons ..... No Charge

#### Systems

For systems with 2 actuators NOT operating through driver control, choose ACM-E2 with either 4-way toggle or 4-push button.

For systems with 2 actuators OPERATING THROUGH DRIVER CONTROL, choose ACM-ES2.

Uses 4-way toggle or 4-push button as attendant control.

### POWERED SEATING ACCESSORIES

- ☐ SANODE One Act. Oper. through Driver Ctrl ..... \$1,300.00 E2310  
**Required for Elevate Operation through Driver Control - Single Actuator Systems. Not compatible with SPJ+ Joysticks. The Egg Switch with the 611 Mount does not come standard when ordering a SANODE. Use for operating through Driver Control.**

### ASL OPTIONAL ACCESSORIES

- ☐ EGSLK Egg Switch - Black ..... \$156.00
- ☐ ASL304A Light Wobble Switch ..... \$180.00
- ☐ ASL202J Fiber Optic Switch ..... \$732.00  
**Not available with SPJ+ joystick. Must add MK6 Auxiliary Power Source.**
- ☐ ASL208J Adj Proximity Switch ..... \$660.00  
**Not available with SPJ+ joystick. Must add MK6 Auxiliary Power Source.**
- ☐ ASL504M6 Remote Emergency Stop Switch ..... \$1,185.00
- ☐ ASL611 Adjustable Switch Mounting Bracket ..... \$252.00
- ☐ ASL804M6 24 Volt Power Source Adaptor..... \$390.00

**NOTE: All specifications and dimensions are approximate.  
For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**



# Invacare® TDX® SC Rehab Seat Power Wheelchair

## MISCELLANEOUS

### BACK ACCESSORIES

- ☐ Clear
- ☐ L8601 Omit Back Upholstery ..... No Charge
- ☐ RBP Reinforced Back Upholstery ..... No Charge

### MISCELLANEOUS ACCESSORIES

- ☐ 1303 Wheel Locks ..... \$45.00
- ☐ SEAT02 Oxygen Holder (ADJASBA Seat Only) ..... \$150.00 E2208

### ADDUCTORS

- ☐ Clear
- Hardware is packaged with equipment and cannot be ordered separately.
- ☐ ADD Adductors Small ..... \$200.00 E0956
- ADDHW Hardware for ADD ..... \$400.00 E1028
- ☐ ADDL Adductors Large ..... \$200.00 E0956
- ADDHWL Hardware for ADDL ..... \$400.00 E1028

### BATTERY

- ☐ UIINSTALL U1 Battery Installation (On Chair) ..... \$232.00
- Required pick on elevating seat.

### BATTERY CHARGER OPTION

- ☐ Clear
- ☐ 110CHARGER 110 Volt Battery Charger ..... No Charge
- ☐ LESTR Lester 8 AMP Battery Charger ..... \$320.00 E2367
- Upcharge available July 1, 2008
- ☐ L8900 Omit Battery Charger ..... Less \$80.00

## FRONT RIGGINGS

### HANGER BRACKET TYPE

- ☐ Clear
- ☐ ASB Center Pivot Style Bracket ..... No Charge
- ☐ PSB Pin Style Bracket ..... No Charge

### CENTER MOUNT

For a 13.25" - 17.25" extension range, available through special #SM-071217095848 with Tall STF Only and 5 Degree Seat-Tilt Angle.

- ☐ CTMT-F Center Mount Adj. Knee Angle (12.25" - 15") ..... \$300.00
- Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 12.5" - 15". Footplates are adjustable angle. Must choose either SM, MED or LRG footplate.

- ☐ Clear
- ☐ SML Small Footplates (5.25" x 7.5") ..... No Charge
- ☐ MDM Medium Footplates (5.5" x 9.5") ..... No Charge
- ☐ LRG Large Footplates (5.90" x 11.5") ..... No Charge
- ☐ CALF Calf Pads For Center Mount ..... \$200.00

- ☐ CTMT-FJ Jr. Center Mount Adj. Knee Angle (9.25" - 12") ..... \$300.00
- Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 9.25" - 12". Footplates are adjustable angle. Available on low seat-to-floor only. Must choose either SM, MED or LRG footplate.

- ☐ Clear
- ☐ SML Small Footplates (5.25" x 7.5") ..... No Charge
- ☐ MDM Medium Footplates (5.5" x 9.5") ..... No Charge
- ☐ LRG Large Footplates (5.90" x 11.5") ..... No Charge
- ☐ CALF Calf Pads For Center Mount ..... \$200.00

### CENTER PIVOT STYLE (PAIR)

- ☐ Clear
- ☐ 70NHD 70° Swingaway Footrests 13 1/2"-17 1/2" ..... Standard
- The composite footplate 1651 is standard.
- ☐ 70STAPER 70° Tapered Footrests 13"-17" ..... No Charge
- Standard with impact guards.
- ☐ AT5544 Manual Elevating Legrests ..... \$255.00 E0990
- Not available with 16" wide by 16" to 18" deep or 17" wide by 16" deep and low seat-to-floor.

### PIN STYLE-HEMI (PAIR)

- ☐ Clear
- ☐ PHW93 70° Swingaway Footrests 13 1/2" - 19" ..... No Charge
- Standard with composite footplates.
- ☐ PHAL4A Manual Art./ELR 17 1/4" - 21 1/2" ..... \$225.00K0053
- Standard with impact guards.

### POWER LEGS

- ☐ Clear
- Independent powered elevating legrests with Adjustable Angle Footplates. Not available with the 16" or 17" seat depths. Extension range 13" - 20" (ELRPW) 15" - 22" (ELRPW22). Sold as pair. To operate power legs, must select either 4 way switch box (4WSB) or the multiple actuator interface box (S4WSB). Not available with elevating seat. Not available with Sanode or SPJ+ joystick.
- ☐ ELRPW Power ELR with Articulation 20" Ext ..... \$1,895.00 E1010
- ☐ ELRPW22 Power ELR with Articulation 22" Ext ..... \$1,995.00 E1010

### FOOTBOARD

- ☐ 90FB Footboard ..... \$150.00
- 90FB is for VAN seat with SEMIRB. Not available with ADJASBA or Elevating Seat.

### FOOTPLATE OPTIONS (PAIR)

- ☐ Clear
- ☐ 1350 Extra Large Aluminum Footplates ..... No Charge
- Not available with AT5544 and 16" wide by 16" - 18" deep or 17" wide by 16" deep and low seat-to-floor.
- ☐ AT5543 Adjustable Angle Flip-Up Footplates ..... \$165.00 K0040
- ☐ CTMTFB Footboard for Center Mount ..... No Charge
- One piece footboard for Center Mount (CTMT-F) Riggings.

### FOOTREST/LEGREST ACCESSORIES

- ☐ 1337 Calf Strap ..... \$27.00 E0038
- ☐ HLOOP Heel Loops (pair) ..... \$45.00 E0951
- ☐ IMPACTG Impact Guards ..... \$51.00
- ☐ ALPT Longer Pivot And Slide Tube ..... \$84.00
- ☐ ALPT4 Longer Pivot And Slide Tube 4" ..... \$84.00

**NOTE: All specifications and dimensions are approximate.**  
**For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**



# Invacare® TDX® SC Rehab Seat Power Wheelchair

TDX SC Specifications		
Seat Width Range (in 1" increments)	Adult (Adjust)	
Seat Depth Range	16" - 20"	
Back Height Range	16" - 24"	
Back Angle Range	85° to 105° in 5" increments	
Seat-to-Floor (approx)	Rehab Seat	Elevating
Standard	16 1/2" at 0° (at 5° add 1/2")	16 3/4" at 0°
Medium	17 1/2" at 0° (at 5° add 1/2")	
Tall	18 1/2" at 0° (at 5° add 1/2")	
	At 5° add 1/2" to Med or Tall. Low STF height not available at 5°.	
Semi Recline Solid Seat	16" Low, 17" Med	
Overall Width of Base	24" (without Joystick)	
Overall Height (ADJASBA) (with 20" back height and 0° seat angle)		
Standard	36 1/2"	
Medium	37 1/2"	
Tall	38 1/2"	
Overall Height (Semi-Recline)		
with headrest	38"	
without headrest	45"	
Overall Length		
without Front Rigging	35 1/4"	
with 90° front rigging	43"	
Weight		
2 Pole Motor Rehab Seat	without U1 Batteries	178 lbs.
	with two U1 Batteries	225 lbs.
2 Pole Motor Elevate	without U1 Batteries	178 lbs.
	with U1 Batteries	225 lbs.
2 Pole Motor Semi Recline Van Seat	without U1 Batteries	183 lbs.
	with U1 Batteries	230 lbs.
Drive Axle	Non-Adjustable	
Drive Wheels/Tires (foam filled)		
Standard	10 1/2" x 3 1/2"	

Contoura Seat Sizes					
	D16	D17	D18	D19	D20
W16	X	X	X	X	X
W18	X	X	X	X	X
W20	X	X	X	X	X
Back Heights	22", 25" or 28"				

## Memory Card - BASIC

- Standard with ALL Expandable Control Systems (CMPJ+, MK6i Display).
- Used to back up / restore programmed settings for only one chair.
- May be given to user or kept with provider's files for safe keeping.
- Does not contain advanced diagnostics / help screens / file structure.
- Not compatible with SPJ+ Joysticks.

## Memory Card - PROFESSIONAL

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains advanced diagnostics, help screens, software upgrades, file storage / retrieval
- Used to back up multiple chairs and programming settings.
- Intended for use by qualified providers only.
- Not compatible with SPJ+ Joysticks.

## Weight Limit 250 Pounds (Standard)

Invacare recommends all wheelchairs be fitted and programmed to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

## Seating and Positioning

For your seating needs, please refer to the following rehab seating and positioning products: Price List Form No. 07-027

## Cantilever Arm Height Adjustment

Back Type	Height	Adjustment	Adjustment with laterals
CURV	20" or greater	11" to 16"	11" to 14"
CURV	19"	11" to 15"	11" to 13"
CURV	18"	11" to 14"	11" to 12"
CURV	17"	11" to 13"	no adjustment
CURV	16"	11" to 12"	not available
Contoura	25" to 28"	12 1/2" to 14 1/2"	no adjustment

Comments

**NOTE: All specifications and dimensions are approximate.**  
**For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**



**Invacare Corporation**  
**USA**

One Invacare Way  
 Elyria, OH USA  
 44035-4190  
 800-333-6900

[www.invacare.com](http://www.invacare.com)

©2010 Invacare Corporation

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Products designated by ASL are manufactured by Adaptive Switch Labs, Inc. Therafin is a registered trademark of W/R Medical Electronics Co. Freedom Designs is a registered trademark of Freedom Designs, Inc. Otto Bock is a trademark of Otto Bock Orthopaedische Industrien GMBH and Co. Specifications and prices to change without notification.