



Yes, you can.®

Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.

Invacare® TDX® SC Van Seat Power Wheelchair

Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



Yes, you can.®

TRANSACTION TYPE:

- Order Quote

Date _____ Account # _____ P.O. # _____

Company Name _____ Phone # _____

Address _____ City _____ State _____ Zip Code _____ Country _____

Contact Name _____ Comments _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

CHAIR TYPE

- TDXSCV TDX SCV with MK6i™ \$5,300.00 K0823
TDXSC comes standard with: 250 lbs. user weight limit, 4.8 MPH top speed, 2-pole motor package, MK6i SPJ+ Joystick, MK660ACC Controller and Quiet Stability Lock. Base width is 24". For Semi Recline with solid seat, see form #07-039 TDXSC Rehab Price List.
- TDXSCSEAT TDXSC Seat Only \$755.00

Group 2 - K0823 Basic Package	
TDX SC with MK6i and Van Seat	
250 lbs. Weight Capacity	
DW1053FF 10 1/2" x 3 1/2" Flat Free Tires	Standard
2-Pole Motors	Standard
Flip Back Fixed Height Arm	Standard
MK660ACC 60 Amp Controller	Standard
SPJ+ Single Drive Joystick	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard

WHEEL OPTION (PAIR)

- DW1053FF 10 1/2" x 3 1/2" Flat Free Tires Standard

FORK OPTION

- FKPKG One-Sided Fork Package \$500.00

FRAME FINISH

- 24P Wet Black No Charge
- 61P Electric Red No Charge
- 62P Electric Blue No Charge
- 60P Silver Metallic No Charge
- 115P Black Prism No Charge
- 120P Grape Madness No Charge
- 127P Grasshopper Green No Charge

VAN SEAT (STANDARD DIMENSIONS: 18"W x 16"D)

- VW16 16" Wide Van Seat No Charge
- VD16 16" Depth No Charge
- VD17 17" Depth No Charge
- VD18 18" Depth No Charge
- VW18 18" Wide Van Seat Standard
- VD16 16" Depth Standard
- VD17 17" Depth No Charge
- VD18 18" Depth No Charge
- VW20 20" Wide Van Seat \$329.00
- VD18 18" Seat Depth No Charge
- VD19 19" Seat Depth No Charge
- VD20 20" Seat Depth \$529.00
- VW22 22" Wide Van Seat \$329.00
- VD18 18" Seat Depth No Charge
- VD19 19" Seat Depth No Charge
- VD20 20" Seat Depth \$529.00

VAN UPHOLSTERY

- GV Grey Vinyl No Charge

SEAT STYLE

- SEMIR Semi Recline Standard

ARM TYPE (VAN)

- Clear
- Arms are flip-up and adjustable width.**
- VAH Full Length, Adjustable Height \$230.00 E0973
- VFH Full Length, Fixed Height Standard
- VAHD Desk Length, Adjustable Height \$230.00 E0973
- VFHD Desk Length, Fixed Height No Charge

MIXED VAN ARMS (EACH)

- Clear
- Right**
- RVFHF Full Length Fixed Height Right No Charge
- RVAHF Full Length Adjustable Height Right \$115.00 E0973
- RVFHD Desk Length Fixed Height Right No Charge
- RVAHD Desk Length Adjustable Height Right \$115.00 E0973
- Left**
- Clear
- LVFHF Full Length Fixed Height Left No Charge
- LVAHF Full Length Adjustable Height Left \$115.00 E0973
- LVFHD Desk Length Fixed Height Left No Charge
- LVAHD Desk Length Adjustable Height Left \$115.00 E0973

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® TDX® SC Van Seat Power Wheelchair

SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Positioning Strap \$50.00
- 1515 Push Button Style Seat Positioning Strap Standard

MK6i™ ELECTRONICS

MK6i CONTROLLERS

- MK660ACC 60 Amp Controller Standard
- GTRAC G-Trac Module \$995.00 K0108
G-Trac provides improved tracking capabilities and enhanced traction. Does not require special motors. Not compatible with SPJ+ joysticks.

STANDARD DRIVER CONTROL

- SPJPLUS SPJ+ Single Drive Joystick Standard

EXPANDABLE CONTROLLER SYSTEM Clear

- PWH Harness req. for Expandable Electronics \$350.00 E2313
Required to interface joystick or MK6i display with expandable controllers.
- CMPJM6U MPJ+ Multiple Drive Joystick with **Color Display** .. \$750.00 E2377

JOYSTICK MOUNTING HARDWARE

- QLM6 Quad Link Retract. Jystk Mnt CMPJ+ Van \$230.00 E1028
- QL Quad Link Retract. Jystk Mnt SPJ+ Van \$230.00 E1028
- ARM250 Stealth Height Adj Joystick Mount \$210.00 E1028
- LEFT Left Handed Mounted Joystick No Charge

ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer with Pro Memory Card \$475.00
- 1813SD Pro Memory Card (USBReady) \$90.00
- 1560 "T" Handle Flexible Joystick Ext. \$77.00 E2323
- 1561 Straight Handle Flexible Joystick Ext. \$77.00 E2323
- 1826 Chin Cup \$48.00 E2324
- MPY Monoport "Y" Cable for CMPJ+ or Display \$15.00
- RFM Wireless (RF) Mouse Only Module \$750.00
RFM includes Transmitter, USB Radio Frequency Mouse Dongle (Receiver), and Stereo Switch with hook and loop mount for mouse clicks. Not compatible with SPJ+.
- IR Infrared Module with Wireless RF Mouse Control \$995.00
Capable of operating up to six remote Infra-Red controlled devices. Also allows wireless (RF) Mouse operation through the driver control. Compatible with RFM module for control of two wireless mice from one chair. Includes two button stereo switch (with hook/loop mount) for mouse clicks. Not compatible with SPJ+ Joysticks.

ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Switch - Black \$156.00
- ASL611 Adjustable Switch Mounting Bracket \$252.00

MISCELLANEOUS

MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks \$45.00
- VAN02 Oxygen Holder (Van Seat Only) \$150.00 E2208

BATTERY

- U1INSTALL U1 Battery Installation (On Chair) \$232.00

BATTERY CHARGER OPTION Clear

- 110CHARGER 110 Volt Battery Charger No Charge
- LESTR Lester 8 AMP Battery Charger \$320.00 E2367
- L8900 Omit Battery Charger Less \$80.00

FRONT RIGGINGS

HANGER BRACKET TYPE Clear

- ASB Center Pivot Style Bracket No Charge
- PSB Pin Style Bracket No Charge

CENTER PIVOT STYLE (PAIR) Clear

- 70NHD 70° Swingaway Footrests 13 1/2" - 17 1/2" Standard
- 70STAPER 70° Tapered Footrests 13" - 17" No Charge
Not available with 17" wide by 16" deep and low seat-to-floor.
- AT5544 Manual Elevating Legrests \$255.00 E0990
Not available with 17" wide by 16" deep and low seat-to-floor.

PIN STYLE-HEMI (PAIR) Clear

- PHW93 70° Swingaway Footrests 13 1/2" - 19" No Charge
Standard with aluminum footplates.
- PHAL4A Manual Art./ELR 17 1/4" - 21 1/2" \$225.00 K0053
Standard with aluminum footplates.

FOOTBOARD

- 90FB Footboard \$150.00

FOOTPLATE OPTIONS (PAIR) Clear

- 1350 Extra Lg Aluminum Footplates No Charge
Not available with AT5544 and 17" wide by 16" deep and low seat-to-floor.
- AT5543 Adj Angle Flip-Up Footplates \$165.00 K0040
Not available with AT5544 and 16" wide by 16" - 18" deep or 17" wide by 16" deep.

FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap \$27.00 E0038
- IMPACTG Impact Guards \$56.00
- HLOOP Heel Loops (pair) \$45.00 E0951
- ALPT Longer Pivot And Slide Tube \$84.00
- ALPT4 Longer Pivot And Slide Tube 4" \$84.00

**NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.**

Invacare® TDX® SC Van Seat Power Wheelchair

Seating and Positioning

For your seating needs, please refer to the following rehab seating and positioning products: Price List Form No. 07-027

Weight Limit 250 Pounds (Standard)

Invacare recommends all wheelchairs be fitted and programmed to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

TDX SC Specifications

Seat Width Range (in 1" increments)	Adult (Adjust)
	16", 18", 20", 22"
Seat Depth Range	16" - 20"
Seat-to-Floor (approximate)	21" Standard
Overall Width of Base	24" (with out Joystick)
Overall Height	
Semi-Recline with out headrest	38"
Semi-Recline with headrest	45"
Overall Length	
with out front rigging	35 1/4"
with 90° front rigging	43"
Weight	
2 Pole Motor with out U1 Batteries	183 lbs.
2 Pole Motor with two (2) U1 Batteries	230 lbs.
Drive Axle	Non-Adjustable
Drive Wheels/Tires (foam filled)	
Standard	10 1/2" x 3 1/2"

Note: All weight limits are inclusive of add-on equipment and seating systems. All specifications and dimensions are approximate.

Comments

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.



Invacare Corporation
USA
 One Invacare Way
 Elyria, OH USA
 44035-4190
 800-333-6900

www.invacare.com

©2013 Invacare Corporation
 All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.