



## Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email\* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

\*Saving and emailing the eForms is available with Acrobat Professional software

### If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

**Please note:** The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

### If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit [www.adobe.com](http://www.adobe.com).

# Invacare® TDX® SC Van Seat Power Wheelchair

## Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



### TRANSACTION TYPE:

- Order  Quote

Date \_\_\_\_\_ Account # \_\_\_\_\_ P.O. # \_\_\_\_\_  
 Company Name \_\_\_\_\_ Phone # \_\_\_\_\_  
 Address \_\_\_\_\_ City \_\_\_\_\_ State \_\_\_\_\_ Zip Code \_\_\_\_\_ Country \_\_\_\_\_  
 Contact Name \_\_\_\_\_ Comments \_\_\_\_\_

### Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

### CHAIR TYPE

- TDXSCV TDX SCV with MK6i™ ..... \$5,300.00 K0823  
**TDXSC comes standard with: 250 lbs. user weight limit, 4.8 MPH top speed, 2-pole motor package, MK6i SPJ+ Joystick, MK660ACC Controller and Quiet Stability Lock. Base width is 24". For Semi Recline with solid seat, see form #07-039 TDXSC Rehab Price List.**
- TDXSCSEAT TDXSC Seat Only ..... \$755.00

Group 2 - K0823 Basic Package	
TDX SC with MK6i and Van Seat	
250 lbs. Weight Capacity	
DW1053FF 10 1/2" x 3 1/2" Flat Free Tires	Standard
2-Pole Motors	Standard
Flip Back Fixed Height Arm	Standard
MK660ACC 60 Amp Controller	Standard
SPJ+ Single Drive Joystick	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard

### WHEEL OPTION (PAIR)

- DW1053FF 10 1/2" x 3 1/2" Flat Free Tires ..... Standard

### FORK OPTION

- FKPKG One-Sided Fork Package ..... \$500.00

### FRAME FINISH

- 24P Wet Black ..... No Charge  
 61P Electric Red ..... No Charge  
 62P Electric Blue ..... No Charge  
 60P Silver Metallic ..... No Charge  
 115P Black Prism ..... No Charge  
 120P Grape Madness ..... No Charge  
 127P Grasshopper Green ..... No Charge

### VAN SEAT (STANDARD DIMENSIONS: 18"W x 16"D)

- VW16 16" Wide Van Seat ..... No Charge  
 VD16 16" Depth ..... No Charge  
 VD17 17" Depth ..... No Charge  
 VD18 18" Depth ..... No Charge
- VW18 18" Wide Van Seat ..... Standard  
 VD16 16" Depth ..... Standard  
 VD17 17" Depth ..... No Charge  
 VD18 18" Depth ..... No Charge
- VW20 20" Wide Van Seat ..... \$329.00  
 VD18 18" Seat Depth ..... No Charge  
 VD19 19" Seat Depth ..... No Charge  
 VD20 20" Seat Depth ..... \$529.00
- VW22 22" Wide Van Seat ..... \$329.00  
 VD18 18" Seat Depth ..... No Charge  
 VD19 19" Seat Depth ..... No Charge  
 VD20 20" Seat Depth ..... \$529.00

### VAN UPHOLSTERY

- GC Grey Cloth ..... No Charge  
 GV Grey Vinyl ..... No Charge  
 TV Tan Vinyl ..... No Charge

### SEAT STYLE

- SEMIR Semi Recline ..... Standard

### ARM TYPE (VAN)

Arms are flip-up and adjustable width.

- VAH Full Length, Adjustable Height ..... \$230.00 E0973  
 VFH Full Length, Fixed Height ..... Standard  
 VAHD Desk Length, Adjustable Height ..... \$230.00 E0973  
 VFHD Desk Length, Fixed Height ..... No Charge

### MIXED VAN ARMS (EACH)

- Right**  Clear
- RVFHF Full Length Fixed Height Right ..... No Charge  
 RVAHF Full Length Adjustable Height Right ..... \$115.00 E0973  
 RVFHD Desk Length Fixed Height Right ..... No Charge  
 RVAHD Desk Length Adjustable Height Right ..... \$115.00 E0973
- Left**  Clear
- LVFHF Full Length Fixed Height Left ..... No Charge  
 LVAHF Full Length Adjustable Height Left ..... \$115.00 E0973  
 LVFHD Desk Length Fixed Height Left ..... No Charge  
 LVAHD Desk Length Adjustable Height Left ..... \$115.00 E0973

**NOTE: All specifications and dimensions are approximate.  
 For the most current pricing information go to [www.invacare.com](http://www.invacare.com).**

# Invacare® TDX® SC Van Seat Power Wheelchair

## SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Positioning Strap ..... \$50.00
- 1515 Push Button Style Seat Positioning Strap ..... Standard

## MK6i™ ELECTRONICS

### MK6i CONTROLLERS

- MK660ACC 60 Amp Controller ..... Standard
- GTRAC G-Trac Module ..... \$995.00 K0108

**G-Trac provides improved tracking capabilities and enhanced traction. Does not require special motors. Not compatible with SPJ+ joysticks.**

### STANDARD DRIVER CONTROL

- SPJPLUS SPJ+ Single Drive Joystick ..... Standard

### EXPANDABLE CONTROLLER SYSTEM

Clear

- PWH Harness req. for Expandable Electronics ..... \$350.00 E2313  
**Required to interface joystick or MK6i display with expandable controllers.**

- CMPJM6U MPJ+ Multiple Drive Joystick with **Color Display** .. \$750.00 E2377

- PSFM6 Personal Joystick, Inductive on Front ..... \$1,130.00 E2377  
**TOD-SS: Toggle On/Off Drive Select, Speed Potentiometer Standard.**

- PSRM6 Personal Joystick, Inductive on Rear ..... \$1,130.00 E2377  
**TOD-SS: Toggle On/Off Drive Select, Speed Potentiometer Standard.**

Clear

- PBOD Push Button On/Off with PB Drive Select .... No Charge  
**PSRM6 Only.**

- PBSS Push Button On/Off with Speed Select ..... No Charge  
**Drive select through mode switch, or mode switch and joystick. PSRM6 Only.**

Clear

- OR On/Off Mounted Right ..... No Charge  
**PSRM6 Only.**

- OL On/Off Mounted Left ..... No Charge  
**PSRM6 Only.**

### JOYSTICK MOUNTING HARDWARE

- QLM6 Quad Link Retract. Jystk Mnt CMPJ+ Van ..... \$230.00 E1028
- QL Quad Link Retract. Jystk Mnt SPJ+ Van ..... \$230.00 E1028
- ARM250 Stealth Height Adj Joystick Mount ..... \$210.00 E1028
- LEFT Left Handed Mounted Joystick ..... No Charge

### ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer with Pro Memory Card .... \$475.00
- 1813SD Pro Memory Card (USBReady) ..... \$90.00
- 1560 "T" Handle Flexible Joystick Ext. .... \$77.00 E2323
- 1561 Straight Handle Flexible Joystick Ext. .... \$77.00 E2323
- 1826 Chin Cup ..... \$48.00 E2324
- MPY Monoport "Y" Cable for CMPJ+ or Display ..... \$15.00
- RFM Wireless (RF) Mouse Only Module ..... \$750.00

**RFM will be available in September, 2010. Allows wireless computer mouse operation through the driver control. Includes two button stereo switch (with hook/loop mount) for mouse clicks. Not compatible with SPJ+ joysticks.**

### ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Switch - Black ..... \$156.00
- ASL611 Adjustable Switch Mounting Bracket ..... \$252.00

## MISCELLANEOUS

### MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks ..... \$45.00
- VAN02 Oxygen Holder (Van Seat Only) ..... \$150.00 E2208

### BATTERY

- U1INSTALL U1 Battery Installation (On Chair) ..... \$232.00

### BATTERY CHARGER OPTION

Clear

- 110CHARGER 110 Volt Battery Charger ..... No Charge
- LESTR Lester 8 AMP Battery Charger ..... \$320.00 E2367
- L8900 Omit Battery Charger ..... Less \$80.00

## FRONT RIGGINGS

### HANGER BRACKET TYPE

Clear

- ASB Center Pivot Style Bracket ..... No Charge
- PSB Pin Style Bracket ..... No Charge

### CENTER PIVOT STYLE (PAIR)

Clear

- 70NHD 70° Swingaway Footrests 13 1/2" - 17 1/2" ..... Standard
- 70STAPER 70° Tapered Footrests 13" - 17" ..... No Charge  
**Not available with 17" wide by 16" deep and low seat-to-floor.**
- AT5544 Manual Elevating Legrests ..... \$255.00 E0990  
**Not available with 17" wide by 16" deep and low seat-to-floor.**

### PIN STYLE-HEMI (PAIR)

Clear

- PHW93 70° Swingaway Footrests 13 1/2" - 19" ..... No Charge  
**Standard with aluminum footplates.**
- PHAL4A Manual Art./ELR 17 1/4" - 21 1/2" ..... \$225.00 K0053  
**Standard with aluminum footplates.**

### FOOTBOARD

- 90FB Footboard ..... \$150.00

### FOOTPLATE OPTIONS (PAIR)

Clear

- 1350 Extra Lg Aluminum Footplates ..... No Charge  
**Not available with AT5544 and 17" wide by 16" deep and low seat-to-floor.**
- AT5543 Adj Angle Flip-Up Footplates ..... \$165.00 K0040  
**Not available with AT5544 and 16" wide by 16" - 18" deep or 17" wide by 16" deep.**

### FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap ..... \$27.00 E0038
- IMPACTG Impact Guards ..... \$56.00
- 1600BK Heel loops with ankle straps ..... \$51.00 E0951
- ALPT Longer Pivot And Slide Tube ..... \$84.00
- ALPT4 Longer Pivot And Slide Tube 4" ..... \$84.00

**NOTE: All specifications and dimensions are approximate.  
For the most current pricing information go to [www.invacare.com](http://www.invacare.com).**

# Invacare® TDX® SC Van Seat Power Wheelchair

### Seating and Positioning

For your seating needs, please refer to the following rehab seating and positioning products: Price List Form No. 07-027

### Weight Limit 250 Pounds (Standard)

Invacare recommends all wheelchairs be fitted and programmed to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

**Note:** All weight limits are inclusive of add-on equipment and seating systems. All specifications and dimensions are approximate.

### TDX SC Specifications

<b>Seat Width Range</b> (in 1" increments)	Adult (Adjust)
	16", 18", 20", 22"
<b>Seat Depth Range</b>	16" - 20"
<b>Seat-to-Floor (approximate)</b>	21" Standard
<b>Overall Width of Base</b>	24" (with out Joystick)
<b>Overall Height</b>	
Semi-Recline with out headrest	38"
Semi-Recline with headrest	45"
<b>Overall Length</b>	
with out front rigging	35 1/4"
with 90° front rigging	43"
<b>Weight</b>	
2 Pole Motor with out U1 Batteries	183 lbs.
2 Pole Motor with two (2) U1 Batteries	230 lbs.
<b>Drive Axle</b>	Non-Adjustable
<b>Drive Wheels/Tires (foam filled)</b>	
Standard	10 1/2" x 3 1/2"

Comments

**NOTE: All specifications and dimensions are approximate.**  
 For the most current pricing information go to [www.invacare.com](http://www.invacare.com).



**Invacare Corporation**  
 USA  
 One Invacare Way  
 Elyria, OH USA  
 44035-4190  
 800-333-6900

[www.invacare.com](http://www.invacare.com)

©2010 Invacare Corporation  
 All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.