



Yes, you can.®

### Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email\* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

\*Saving and emailing the eForms is available with Acrobat Professional software

#### If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

**Please note:** The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

#### If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit [www.adobe.com](http://www.adobe.com).

# Invacare® Formula™ CG Tilt Only for TDX® SC Base

## Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



Yes, you can.®

### TRANSACTION TYPE:

Order  Quote

Date \_\_\_\_\_ Account # \_\_\_\_\_ P.O. # \_\_\_\_\_

Company Name \_\_\_\_\_ Phone # \_\_\_\_\_

Address \_\_\_\_\_ City \_\_\_\_\_ State \_\_\_\_\_ Zip Code \_\_\_\_\_ Country \_\_\_\_\_

Contact Name \_\_\_\_\_ Comments \_\_\_\_\_

### Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

### POWERED SEATING BASE

TDXSC2-CG TDX SC for Formula CG with MK6i™ ..... \$5,300.00 K0835  
**TDXSC2-CG (Group 2) comes standard with: 250 lb. user weight capacity, 4.8 MPH top speed, 2-pole motor package, MK660ACC Controller and Quiet Stability Lock. Base width is 24".**

Group 2 - K0835 Basic Package	
TDXSC2-CG TDXSC for Formula CG with MK6i	
250 lbs. Weight Capacity	
DW10535FF 10 1/2" x 3 1/2" Tires with Foam Filled Inserts	Standard
2-Pole Motors	Standard
Flip Back Adjustable Height Arm	Standard
MK660ACC 60 AMP Controller	Standard
SPJ+ Joystick for Powered Seating	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard
<b>Seat Option</b> CGT Tilt (E1002)	\$4,590.00

### BASE

#### OMIT SEAT FRAME

OSF Omit Seat Frame for Base Only..... No Charge  
**Base only to send for aftermarket seating.**

#### TRANSPORT TIE DOWN

TRBKTS Wheelchair Transport Brackets ..... \$250.00

### TRANSPORT TIE-DOWN\*

TRBKTS includes four factory installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport. As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

\* Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

### USER WEIGHT LIMIT

UW250 User Weight <250 lbs. .... Standard

### WHEEL OPTION (PAIR)

DW10535FF 10 1/2" x 3 1/2" Flat Free Tires ..... Standard

### FORK OPTION

FKPKG One-Sided Fork Package ..... \$500.00

### FRAME FINISH

- 24P Wet Black ..... No Charge
- 61P Electric Red ..... No Charge
- 62P Electric Blue ..... No Charge
- 60P Silver Metallic ..... No Charge
- 115P Black Prism ..... No Charge
- 120P Grape Madness ..... No Charge
- 127P Grasshopper Green ..... No Charge

**NOTE: All specifications and dimensions are approximate.  
 For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

# Invacare® Formula™ CG Tilt Only for TDX® SC Base

## SEAT

### SYSTEM TYPE

CGT CG Tilt Only System ..... \$4,590.00 E1002

### SEAT-TO-FLOOR HEIGHT

**Due to variations in options the listed seat-to-floor heights have a tolerance of +/- 0.5"**

- LOW 17.5" at 0° ..... No Charge
- MED 18.5" at 0° ..... No Charge
- TALL 19.5" at 0° ..... No Charge

### SEAT TILT

**5° Seat-to-Floor height not available with Elevating Seat and not available on LOW Seat-to-Floor height on non elevating seat. Add 1" to Med STF and 1.25" to Tall STF.**

- 0SA 0° Seat Tilt ..... No Charge
- 5SA 5° Seat Tilt ..... No Charge

### SEAT SELECTION

- BPO Back Posts Only ..... No Charge  
**Push handles come standard with BPO and curved back.**
- CURV Curved ..... No Charge  
**Push handles come standard with BPO and curved back.**
- CTB Contoura™ Seating ..... \$625.00 E2620 E2621  
 OMCU Omit Contoura Seat Cushion ..... Less \$20.00
- CTBMSH Contoura™ seating with mesh upholstery ..... \$635.00 E2620 E2621  
 OMCU Omit Contoura Seat Cushion ..... Less \$20.00

### ADJUSTABLE SEAT WIDTH RANGE

1620W 16" to 20" Seat Width Range ..... No Charge

### SEAT WIDTH SETTING

- W16 16" Seat Width ..... No Charge
- W17 17" Seat Width (not available with Contoura) ..... No Charge
- W18 18" Seat Width ..... No Charge
- W19 19" Seat Width (not available with Contoura) ..... No Charge
- W20 20" Seat Width ..... No Charge

### ADJUSTABLE SEAT DEPTH RANGE

- 1619D 16" to 19" Seat Depth Range ..... No Charge
- 1920D 19" to 20" Seat Depth Range ..... No Charge

### SEAT DEPTH SETTING

- D16 16" Seat Depth ..... No Charge
- D17 17" Seat Depth ..... No Charge
- D18 18" Seat Depth ..... No Charge
- D19 19" Seat Depth ..... No Charge
- D20 20" Seat Depth ..... No Charge

### BACK HEIGHT

- TSBH16 16" Back Height ..... No Charge
- TSBH17 17" Back Height ..... No Charge
- TSBH18 18" Back Height ..... No Charge
- TSBH19 19" Back Height ..... No Charge
- TSBH20 20" Back Height ..... No Charge
- TSBH21 21" Back Height ..... No Charge
- TSBH22 22" Back Height ..... No Charge
- TSBH23 23" Back Height ..... No Charge
- TSBH24 24" Back Height ..... No Charge

### CONTOURA BACK HEIGHT

- Clear
- TSBH22 Back Height 22" ..... No Charge
- TSBH25 Back Height 25" ..... No Charge
- TSBH28 Back Height 28" ..... No Charge

Back Height Chart	
Back Type	Back Height
CTB - Contoura	22", 25", and 28"
BPO - Back Posts Only	16" - 24"
CURV Curved	16" - 24"

### ADJUSTABLE ARM TYPE (PAIR)

- Clear
- 19F Full Length Adj Hgt Conv. .... No Charge
- 29F Desk Length Adj Hgt Conv. .... No Charge

### FLIP UP MIXED ADJUSTABLE ARM TYPE (EACH)

**All arms are detachable.**

#### Right Arm

- Clear
- 19FR Full Length Adj. Hgt. Right ..... No Charge
- 29FR Desk Length Adj. Hgt. Right ..... No Charge

#### Left Arm

- Clear
- 19FL Full Length Adj. Hgt. Left ..... No Charge
- 29FL Desk Length Adj. Hgt. Left ..... No Charge

### CANTILEVER ARM TYPE (PAIR)

**Cantilever Arms only available with 300 lbs. rated systems. Not available with BPO, Otto Bock or 12" - 15" seat width. Stealth Lateral part number's TWB-CPR and TWB-CPRL are the only lateral support compatible with the current Cantilever Arm design. Must choose WF armpad type with Stealth Laterals. See chart on last page for arm height adjustment range.**

- CLVF Full Length Arm ..... \$350.00
- CLVD Desk Length Arm ..... \$350.00

### CANTILEVER MIXED ARM TYPE (EACH)

#### Right Arm

- Clear
- CLVFR Full Length Arm Right ..... \$175.00
- CLVDR Desk Length Arm Right ..... \$175.00

#### Left Arm

- Clear
- CLVFL Full Length Arm Left ..... \$175.00
- CLVDL Desk Length Arm Left ..... \$175.00

### ARMPAD OPTIONS (CANTILEVER ARMS)

- Clear
- WF Waterfall Armpads ..... Standard
- RCL Reclining Armpads ..... No Charge

### SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Positioning Strap ..... \$50.00
- 1515 Push Button Style Seat Positioning Strap ..... Standard

### CHEST POSITIONING STRAPS

- Clear
- 7311BK Airline Buckle Chest Positioning Strap ..... \$100.00 E0960
- 9515 Push Button Style Chest Positioning Strap ..... \$100.00 E0960
- 7321BK Hook and Loop Chest Positioning Strap ..... \$100.00 E0960

### HEADREST TYPE (OPTIONAL)

- Clear
- Hardware is packaged with equipment and cannot be ordered separately.**
- TAR127700 Curved Headrest ..... \$210.00 E0955  
TAR27700HW Removable & Adjustable Hardware ..... \$210.00 E1028
- TAR128000 Two-Step Headrest ..... \$210.00 E0955  
TAR128000HW Removable & Adjustable Hardware ..... \$210.00 E1028
- T0022400 Large Headrest ..... \$210.00 E0955  
T0022400HW Removable & Adjustable Hardware ..... \$210.00 E1028
- T0022600 Small Headrest ..... \$210.00 E0955  
T0022600HW Removable & Adjustable Hardware ..... \$210.00 E1028

**NOTE: All specifications and dimensions are approximate.  
For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

**SEAT continued**

**STEALTH HEADRESTS**

Clear

Hardware is packaged with equipment and cannot be ordered separately.

- CP175 14" Multi Axis Headrest ..... \$225.00 E0955  
CP175HW Removable & Adjustable Hardware ..... \$225.00 E1028
- CP275 10" Multi Axis Headrest ..... \$225.00 E0955  
CP275HW Removable & Adjustable Hardware ..... \$225.00 E1028
- CP180 14" TWB Headrest ..... \$225.00 E0955  
CP180HW Removable & Adjustable Hardware ..... \$225.00 E1028
- CP280 10" TWB Headrest ..... \$225.00 E0955  
CP280HW Removable & Adjustable Hardware ..... \$225.00 E1028

**LATERALS (PAIR)**

**LATERALS FOR CONVENTIONAL**

Clear

For use with all back types. Hardware is packaged with equipment and cannot be ordered separately.

- TCSLM Clamp-On Swingaway Laterals Medium ..... \$210.00 E0956  
TCSLMHW Removable and Adjustable Hardware ..... \$210.00 E1028
- TCSSL Clamp-On Swingaway Laterals Large ..... \$210.00 E0956  
TCSSLHW Removable and Adjustable Hardware ..... \$210.00 E1028

**STEALTH LATERALS**

Clear

For use with all back types. Hardware is packaged with equipment and cannot be ordered separately.

- TWB-RPR Adj. Swing Laterals Med ..... \$320.00 E0956  
RPRHW Removable and Adjustable Hardware ..... \$320.00 E1028
- TWB-RPRL Adj. Swing Laterals Lg ..... \$320.00 E0956  
RPRLHW Removable and Adjustable Hardware ..... \$320.00 E1028

**MK6i™ ELECTRONICS**

**MK6i™ CONTROLLERS**

- MK660ACC 60 Amp Controller ..... Standard
- GTRAC G-Trac Module ..... \$995.00 K0108  
**G-Trac provides improved tracking capabilities and enhanced traction. Does not require special motors. Not compatible with SPJ+ joysticks.**

**STANDARD DRIVER CONTROL**

- SPJAP SPJ+™ for Power Seating ..... Standard  
**Includes Egg Switch with 611 mount to operate single actuator**

**ACTUATOR CONTROL THROUGH SPJ+**

- SPJAPM SPJ+ with Joystick Actuator Control ..... \$1,300.00 E2310  
**Allows actuator operation through Driver Control. Attendant switch not included.**

**EXPANDABLE CONTROLLER SYSTEM**

Clear

- PWH Harness req. for Expandable Electronics ..... \$350.00 E2313  
**Required to interface joystick or MK6i display with expandable controllers.**
- CMPJM6U MPJ+ Multiple Drive Joystick with **Color Display** .. \$750.00 E2377
- ASLRDYM6 MK6 Display Only (No Driver Control) ..... \$750.00 E2377  
**Display only, No Driver Control. Allows choice from order form.**
- MCD Mini Color Display ..... \$750.00 E2377  
**No port for remote On/Off switch. Stereo port only.**

FOR YOUR ASL DRIVE SYSTEM NEEDS  
PLEASE REFER TO THE  
ASL SINGLE INVOICE PRICE LIST FORM # 02-097.

**ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS**

Clear

- 1500M6 Prop RIM Head Control with MK6i Display .. \$2,875.00 E2327
- SMHD Small Headrest ..... No Charge
- LGHD Large Headrest ..... No Charge
- 1558MM6 Compact Joystick with Display ..... \$1,130.00 E2373

**ALTERNATIVE DIGITAL CONTROL SYSTEMS**

Clear

**Must also use either a Multiple Drive Joystick or MK6i Display.**

- SNPM6 Sip-N-Puff Digital Interface Box ..... \$1,500.00 E2325
- 5020M6 Mini Tash Joystick and Digital Interface Box \$1,750.00 E2321

**SIP-N-PUFF KIT**

- PKG32666 Therafin® Sip-N-Puff Breath Tube Kit ..... \$380.00 E2326

**ADDITIONAL CONTROL CHOICES**

- 1558M6 Compact Joystick - Less Std Mount ..... \$770.00
- 5020 Mini Tash Joystick ..... \$405.00  
**Requires SNPM6 Sip-N-Puff Interface Box.**
- PACM6 Prop Attendant Control with quick release mount .... \$850.00
- 1552M Digital Attendant Control ..... \$850.00  
**Requires SNPM6 Sip-N-Puff Interface Box.**

**NOTE: All specifications and dimensions are approximate.  
For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

**ELECTRONICS Continued**

**JOYSTICK MOUNTING HARDWARE**

- QLAM6 Quad Link Retract. Joystick Mnt CMPJ+ ..... \$279.00 E1028  
*Include ARM250 standard.*
- QLA Quad Link Retract. Joystick Mnt SPJ+ ..... \$279.00 E1028  
*Include ARM250 standard.*
- ARM250 Stealth Height Adj Joystick ..... \$210.00 E1028
- LEFT Left Handed Mounted Joystick ..... No Charge
- RIMHW RIM Hardware For 1500 ..... \$219.00 E1028
- PKG32669 Complete Bib Assembly (Tash/1558) ..... \$265.00
- GATMPJM6 Gatlin Midline for CMPJ+ ..... \$600.00 E1028
- GAT1812 Gatlin Midline for 1812 ..... \$600.00 E1028
- RCM Rear Cane Mount for Compact Joystick ..... \$395.00

**ELECTRONIC ACCESSORIES**

- 1813M6 MK6i Programmer with Pro Memory Card ..... \$475.00
- 1813SD Pro Memory Card (USB Ready) ..... \$90.00
- 1560 "T" Handle Flexible Joystick Ext. .... \$77.00 E2323
- 1561 Straight Handle Flexible Joystick Ext. .... \$77.00 E2323
- 1826 Chin Cup ..... \$48.00 E2324
- AUX12M6 Auxiliary Module For outputs 1 and 2 ..... \$795.00
- AUX34M6 Auxiliary Module For outputs 3 and 4 ..... \$795.00
- 2PB Dual Push Button ..... \$147.00  
*Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options.*
- 2SW Dual Toggle ..... \$172.00  
*Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options.*
- PBMT Push Button Mounting Kit ..... \$252.00
- MPY Monoport "Y" Cable for CMPJ+ or Display ..... \$15.00
- RFM Wireless (RF) Mouse Only Module ..... \$750.00  
*RFM includes Transmitter, USB Radio Frequency Mouse Dongle (Receiver), and Stereo Switch with hook and loop mount for mouse clicks. Not compatible with SPJ+.*
- IR Infrared Module with Wireless RF Mouse Control ..... \$995.00  
*Capable of operating up to six remote Infra-Red controlled devices. Also allows wireless (RF) Mouse operation through the driver control. Compatible with RFM module for control of two wireless mice from one chair. Includes two button stereo switch (with hook/loop mount) for mouse clicks. Not compatible with SPJ+ Joysticks.*

**POWERED SEATING ACCESSORIES**

- SANODE One Actuator Operation thru Driver Cntrl..... \$1,300.00 E2310  
*Required for Tilt Operation through Driver Control - Single Actuator Systems Only. Not available with SPJ+ Joystick. Does not come with Egg switch.*
- 4WSB Four Way Switch Box ..... No Charge
  - FWT 4 Way Toggle ..... No Charge
  - QPB 4 Push Buttons ..... No Charge
- S4WSB Multiple Actuator Interface Box ..... \$2,600.00 E2311
  - FWT 4 Way Toggle ..... No Charge
  - QPB 4 Push Buttons ..... No Charge

**Systems**

For systems with 2 actuators NOT operating through driver control, choose ACM-E2 with either 4-way toggle or 4-push button.

For systems with 2 actuators OPERATING THROUGH DRIVER CONTROL, choose ACM-ES2.  
Uses 4-way toggle or 4-push button as attendant control.



**ASL OPTIONAL ACCESSORIES**

- EGSLBK Egg Switch - Black ..... \$156.00
- ASL304 Light Wobble Switch ..... \$180.00
- ASL202J Fiber Optic Switch ..... \$732.00  
*Not available with SPJ+ joystick. Must add ASL804M6 24V Power Source adaptor*
- ASL208J Adj Proximity Switch ..... \$660.00  
*Not available with SPJ+ joystick. Must add ASL804M6 24V Power Source adaptor*
- ASL504M6 Remote Emergency Stop Switch ..... \$1,185.00
- ASL611 Adjustable Switch Mounting Bracket ..... \$252.00
- ASL804M6 24 Volt Power Source Adaptor..... \$390.00

**MISCELLANEOUS**

**BATTERY**

- U1INSTALL U1 Battery Installation (On Chair) ..... \$232.00

**BATTERY CHARGER OPTION**

- Clear
- 110CHARGER 110 Volt Battery Charger ..... No Charge
- LESTR Lester 8 AMP Battery Charger ..... \$320.00 E2367
- L8900 Omit Battery Charger ..... Less \$80.00

**VENT TRAY**

- COMP Compact Vent Tray ..... \$1,540.00 E1030  
*For ventilators with self-contained power supply and back up.*

**MISCELLANEOUS ACCESSORIES**

- 1303 Wheel Locks ..... \$45.00

**ADDUCTORS**

- Clear
- Hardware is packaged with equipment and cannot be ordered separately.*
- ADD Adductors Small ..... \$200.00 E0956  
ADDHW Hardware for ADD ..... \$400.00 E1028
- ADDL Adductors Large ..... \$200.00 E0956  
ADDHWL Hardware for ADDL ..... \$400.00 E1028

**FRONT RIGGINGS**

**HANGER BRACKET TYPE**

- Clear
- ASB Center Pivot Style Bracket ..... No Charge
- PSB Pin Style Bracket ..... No Charge

**CENTER MOUNT**

- Clear
- For a 13.25" - 17.25" extension range, available through special #SM-071217095848 with Tall STF Only and 5 Degree Seat-Tilt Angle.*
- CTMT-F Center Mount Adj. Knee Angle (12.25" - 15") ..... \$300.00  
*Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 12.5" - 15". Footplates are adjustable angle. Must choose either SM, MED or LRG footplate.*
- Clear
- SML Small Footplates (5.25" x 7.5") ..... No Charge
- MDM Medium Footplates (5.5" x 9.5") ..... No Charge
- LRG Large Footplates (5.90" x 11.5") ..... No Charge
- CALF Calf Pads For Center Mount ..... \$200.00
- CTMT-FJ Jr. Center Mount Adj. Knee Angle (9.25" - 12") ..... \$300.00  
*Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 9.25" - 12". Footplates are adjustable angle. Available on low seat-to-floor only. Must choose either SM, MED or LRG footplate.*
- Clear
- SML Small Footplates (5.25" x 7.5") ..... No Charge
- MDM Medium Footplates (5.5" x 9.5") ..... No Charge
- LRG Large Footplates (5.90" x 11.5") ..... No Charge
- CALF Calf Pads For Center Mount ..... \$200.00

**NOTE: All specifications and dimensions are approximate.  
For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

**FRONT RIGGING OPTIONS continued**

**CENTER PIVOT STYLE (PAIR)**  Clear

- 70NHD 70° Swingaway Footrests 13 1/2"-17 1/2" ..... Standard  
**The complete footplate 1651 is standard.**
- 70STAPER 70° Tapered Footrests 13" - 17" ..... No Charge  
**Standard with impact guards.**
- AT5544 Manual Elevating Legrests ..... \$255.00 E0990  
**Not available with 16" wide by 16" to 18" deep or 17" wide by 16" deep and low seat-to-floor.**

**PIN STYLE-HEMI (PAIR)**  Clear

- PHW93 70° Swingaway Footrests 13 1/2" - 19" ..... No Charge  
**Standard with aluminum footplates.**
- PHAL4A Manual Articulating/ELR 17 1/4" - 21 1/2" ..... \$225.00 K0053

**POWER ELEVATING LEGRESTS - ARTICULATING**

- MPCTMT Power Center Mount Foot Platform ..... \$2,995.00 K0108  
**Not available with SPJ+ Joystick or Sanode.**
- IFP Individual Footplates ..... \$275.00
- MPCTMT-L Power Center Mount Foot Platform Long ... \$2,995.00  
 IFP Individual Footplates ..... \$275.00
- ELRPW Power ELR with Articulation 20" Ext. .... \$1,895.00 E1010  
**Not available with SPJ+ Joystick. Independent powered elevating legrests including elevating and articulating legrest with angle adjustable footplates. Not available with the 16" or 17" seat depths. Extension range 13" - 20". Sold as pair. Not available with Sanode.**
- ELRPW22 Power ELR with Articulation 22" Extension ..... \$1,995.00  
**Not available with SPJ+ Joystick. Independent powered elevating legrests including elevating and articulating legrest with angle adjustable footplates. Not available with the 16" or 17" seat depths. Extension range 15" - 22". Sold as pair. Not available with Sanode.**
- IFP Individual Footplates ..... \$275.00

**FOOTPLATE OPTIONS (PAIR)**  Clear

- 1350 Extra Large Aluminum Footplates ..... No Charge  
**Not available with AT5544 and 16" wide by 16" - 18" deep or 17" wide by 16" deep and low seat-to-floor.**
- AT5543 Adjustable Angle Flip-Up Footplates ..... \$165.00 K0040
- CTMTFB Footboard for Center Mount ..... No Charge  
**One piece footboard for Center Mount (CTMT-F) Riggings only.**

**FOOTREST/LEGREST ACCESSORIES**

- 1337 Calf Strap ..... \$27.00 E0038
- HLOOP Heel Loops (pair) ..... \$45.00 E0951
- IMPACTG Impact Guards ..... \$51.00
- ALPT Longer Pivot And Slide Tube ..... \$84.00
- ALPT4 Longer Pivot And Slide Tube 4" ..... \$84.00

Power Center Mount KTH	
MPCTMT	13" - 17"
MPCTMT-L	15" - 19"

**Seating and Positioning**  
For your seating needs, please refer to the following rehab seating and positioning products: Price List Form No. 07-027

Contoura Seat Sizes					
	D16	D17	D18	D19	D20
W16	X	X	X	X	X
W18	X	X	X	X	X
W20	X	X	X	X	X
Back Heights	22", 25", or 28"				
<b>Available with OMCU option only.</b>					

**Memory Card - BASIC**

- Standard with ALL Expandable Control Systems (CMPJ+, MK6i Display).
- Used to back up / restore programmed settings for only one chair.
- May be given to user or kept with provider's files for safe keeping.
- Does not contain advanced diagnostics / help screens / file structure.
- Not compatible with SPJ+ Joysticks.

**Memory Card - PROFESSIONAL**

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains advanced diagnostics, help screens, software upgrades, file storage / retrieval
- Used to back up multiple chairs and programming settings.
- Intended for use by qualified providers only.
- Not compatible with SPJ+ Joysticks.

**NOTE: All specifications and dimensions are approximate.**  
*For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).*



# Invacare® Formula™ CG Tilt Only for TDX® SC Base

Cantilever Arm Height Adjustment			
Back Type	Height	Adjustment	Adjustment with laterals
CURV	20" or greater	11" to 16"	11" to 14"
CURV	19"	11" to 15"	11" to 13"
CURV	18"	11" to 14"	11" to 12"
CURV	17"	11" to 13"	no adjustment
CURV	16"	11" to 12"	not available
Contoura	25" to 28"	12 1/2" to 14 1/2"	no adjustment

Conventional Seat Sizes	
Seat Widths	16" - 20"
Seat Depths	16" - 20"
Back Heights	16" - 24"

**FormulaCG Tilt Only Features and Specifications**

- Weight capacity of 250 lbs.
- Curved or Contoura style seating
- Adjustable seat frame

**Tech Tips**

Operating Powered Seating through the Driver's Control.  
All chairs require adding or replacing one component to allow actuator operation through the driver control.

**SPJ+ with PSS:**  
Replace joystick with SPJ+ with ACC (Actuator Control) - E2310.

**Single Actuator Systems with CMPJ+, PSR+, PSF+ MK6i Display:**  
Add a SANODE (Single Actuator Node) to the system - E2310.

**Multiple Actuator Systems:**  
Use of the M6Ki enables actuator operation through the driver control.

**After Market Backrests**

The following companies also have hardware to make their backrests work with our new canes:

- AEL - planar seating system
- Freedom Designs® - planar seating system
- Canyon - planar seating system

Curved ABS	
Seat Widths	16" - 20"
Seat Depths	16" - 20"
Back Heights	16" - 24"

Comments

**NOTE: All specifications and dimensions are approximate.**  
*For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).*



**Invacare Corporation**  
USA  
One Invacare Way  
Elyria, OH USA  
44035-4190  
800-333-6900

[www.invacare.com](http://www.invacare.com)

©2010 Invacare Corporation  
All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Products designated by ASL are manufactured by Adaptive Switch Labs, Inc. Therafin is a registered trademark of W/R Medical Electronics Co. Freedom Designs is a registered trademark of Freedom Designs, Inc. Otto Bock is a trademark of Otto Bock Orthopaedische Industrien GMBH and Co. Specifications and prices to change without notification.