



Yes, you can.®

Price List and Order Form

Effective May 18, 2009

For ease of ordering, contact Customer Service Toll Free at:

800-333-6900 Fax: 800-678-4682

To order Literature: 800-828-6282

www.invacare.com

Company Name _____ Account # _____

Phone _____ P.O. # _____ Date _____

Ship to Name & Address _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

**Group 3 – K0861
Basic Package**

- TDXSP-MCG TDXSP for Formula CG with MK6i
- 300 lb Weight Capacity
- 1431-3 14" X 3" Tires w/Foam Filled Inserts No Charge
- 4 Pole MotorsStandard
- Flip Back Fixed Height ArmStandard
- MK690ACC Controller w/ Actuator Control.....Standard
- 70NHD 70° Swingaway FootrestsStandard
- 8 Amp Off Board Battery ChargerStandard
- 1515 Push Button Style Seat Positioning StrapStandard

Seat Option

- CGTR Tilt/ Recline/Elevate (E1007)\$9,595.00
- E2-TDX Elevate to Assist for Transfers (E2300)\$1,895.00

POWERED SEATING BASE

- TDXSP-MCG TDXSP for Formula CG w/ MK6i™ (K0861)..... 6,275.00

BASE

OMIT SEAT FRAME

- OSF Omit Seat Frame for Base Only Less 700.00 ¹

USER WEIGHT LIMIT

- U300 User Weight <300 lbs (K0856) **Standard** ²

WHEEL OPTIONS (PAIR)

- 1431-3 14" x 3" with Foam Filled Inserts No Charge
- 1430-3 14" x 3" Pneumatic Tires..... **Standard**

FORK OPTION

- FKPKG One-Sided Fork Package 500.00

FRAME FINISH

CONTEMPORARY

- 24P Wet BlackNo Charge
- 60P Silver MetallicNo Charge
- 71P Silver VeinNo Charge
- 105P Midnight Blue.....No Charge
- 115P Black Prism.....No Charge

TOO HOT

- 30P Sunny YellowNo Charge
- 61P Electric Red.....No Charge
- 119P Bubblegum PinkNo Charge
- 120P Grape MadnessNo Charge

WAY COOL

- 62P Electric BlueNo Charge
- 104P Electric Teal.....No Charge
- 121P Lolly Pop Blue.....No Charge
- 122P Cosmic Blue.....No Charge
- 125P Emerald GreenNo Charge
- 127P Grasshopper GreenNo Charge

FOOTNOTES

¹ Base only to send for aftermarket seating.

² Package includes 4 Pole Motors and 5.8 mph top speed with <300 lb. user weight limit. Package includes Quiet Stability Lock.

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information go to www.invacare.com.**

Invacare® Formula™ CG Tilt/Recline/Elevate for TDX® SP Accessories and Options

SEAT

SYSTEM TYPE

- CGTRE CG Tilt/Recline/Elevate (E1007) 9,595.00 ³
- E2-TDX Elevate to Assist for Transfers (E2300) 1,895.00
- CNVACTS Conventional Actuator..... **Standard**
- SMTACTS Smart Actuators 800.00

SEAT SELECTION

- CNB Conventional Back No Charge
- CTB Contoura™ Seating (E2620/E2621) 625.00
- OMCU Omit Contoura Seat Cushion..... Less 20.00

ADJUSTABLE SEAT WIDTH RANGE

- 1620W 16" to 20" Seat Width Range No Charge
- 2022W 20" to 22" Seat Width Range No Charge

SEAT WIDTH SETTING

- W16 16" Seat Width..... No Charge
- W17 17" Seat Width (not available w/Contoura) ... No Charge
- W18 18" Seat Width..... No Charge
- W19 19" Seat Width (not available w/Contoura) ... No Charge
- W20 20" Seat Width..... No Charge
- W21 21" Seat Width (not available w/Contoura) ... No Charge
- W22 22" Seat Width..... No Charge

ADJUSTABLE SEAT DEPTH RANGE

CONVENTIONAL

- 1619D 16" to 19" Seat Depth Range..... No Charge
- 2022D 20" to 22" Seat Depth Range..... No Charge

CONTOURA

- 1618D 16" to 18" Seat Depth Range..... No Charge
- 1921D 19" to 21" Seat Depth Range..... No Charge
- 22D 22" Seat Depth..... No Charge

SEAT DEPTH SETTING

CONVENTIONAL & CONTOURA

- D16 Seat Depth No Charge ⁴
- D17 Seat Depth No Charge ⁴
- D18 Seat Depth No Charge
- D19 Seat Depth No Charge
- D20 Seat Depth No Charge
- D21 Seat Depth No Charge ⁵
- D22 Seat Depth No Charge ⁵

SEAT TILT

- OSA 0° Seat Tilt..... Standard

ADJUSTABLE ARM TYPE (PAIR) ⁶

- RA19FR Reclining Full Length Adj. Hgt. Right..... No Charge
- RA19FL Reclining Full Length Adj. Hgt. Left..... No Charge
- RA29FR Reclining Desk Length Adj. Hgt. Right..... No Charge
- RA29FL Reclining Desk Length Adj. Hgt. Left No Charge

SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Pos. Strap 50.00
- 1515 Push Button Style Seat Pos. Strap **Standard**

CHEST POSITIONING STRAPS

- 7311BK Airline Buckle Chest Pos. Strap (E0960) 100.00
- 9515 Push Button Style Chest Pos. Strap (E0960)..... 100.00
- 7321BK Hook & Loop Chest Pos. Strap (E0960) 100.00

BACK HEIGHT

CONVENTIONAL

- TSBH20 Back Height 20" No Charge
- TSBH21 Back Height 21" No Charge
- TSBH22 Back Height 22" No Charge
- TSBH23 Back Height 23" No Charge
- TSBH24 Back Height 24" No Charge
- TSBH25 Back Height 25" No Charge
- TSBH26 Back Height 26" No Charge

CONTOURA

- TSBH25 Back Height 25"..... No Charge
- TSBH28 Back Height 28"..... No Charge

POSITIONING ACCESSORIES

HEADREST TYPE (OPTIONAL)

- TAR127700 Curved Headrest (E0955) 210.00
- TAR27700HW Removable & Adj Hardware (E1028) 210.00
- TAR128000 Two-Step Headrest (E0955) 210.00
- TAR12800HW Removable & Adj Hardware (E1028) 210.00
- T0022400 Lg Headrest (E0955)..... 210.00
- TAR22400HW Removable & Adj Hardware (E1028) 210.00
- T0022600 Small Headrest (E0955)..... 210.00
- TAR22600HW Removable & Adj Hardware (E1028) 210.00

STEALTH HEADRESTS

- CP175 14" Multi Axis Headrest (E0955)..... 225.00
- CP175HW Removable & Adj Hardware (E1028) 225.00
- CP275 10" Multi Axis Headrest (E0955)..... 225.00
- CP275HW Removable & Adj Hardware (E1028) 225.00
- CP180 14" TWB Headrest (E0955) 225.00
- CP180HW Removable & Adj Hardware (E1028) 225.00
- CP280 10" TWB Headrest (E0955) 225.00
- CP280HW Removable & Adj Hardware (E1028) 225.00

OTTO BOCK ARM PADS (PAIR) ⁷

- OB20400 One Piece Channel Armpads (E2209)..... 230.00
- OB20500 Lg Forearm Pads..... 235.00
- OB20600 Med Forearm Pads 235.00
- OB20700 Small Forearm Pads..... 235.00

OTTO BOCK HAND PADS (PAIR) ⁸

- OB23000 Lg Flat Hand Pads 100.00
- OB23300 Med Flat Hand Pads 100.00

FOOTNOTES

³ Must choose the E2-TDX option with the CGTRE. If you want the Formula CG Tilt only (CGT), see form no. 06-039 – Formula CG Tilt Only pricelist.

⁴ Not available with Contoura seat cushion on 22" wide.

⁵ Not available with Contoura seat cushion on 16" or 18" wide.

FOOTNOTES

⁶ All arms are flip-up and detachable.

⁷ Not compatible with quadlink retractable joystick mount. . If selecting an Otto Bock Hand Pad, a large, medium or small forearm pad must be ordered.

Invacare® Formula™ CG Tilt/Recline/Elevate for TDX® SP

Accessories and Options

LATERALS (PAIR)

CONVENTIONAL BACK⁸

- CSLM Swingaway Laterals Medium (E0956)..... 210.00
- CSLMHW Removable & Adj Hardware (E1028)..... 210.00
- CSLL Swingaway Laterals Large (E0956)..... 210.00
- CSLLHW Removable & Adj Hardware (E1028)..... 210.00

STEALTH LATERALS FOR CONVENTIONAL BACK⁸

- TWB-CPR Adj. Swing Channel Laterals Med (E0956)..... 320.00
- CPRHW Removable & Adj Hardware (E1028)..... 320.00
- TWB-CPRL Adj. Swing Channel Laterals Lg (E0956)..... 320.00
- CPRLHW Removable & Adj Hardware (E1028)..... 320.00

CONTOURA SEATING⁹

- TMSLM Swingaway Contoura Laterals Med (E0956)..... 210.00
- TMSLMHW Removable & Adj Hardware (E1028)..... 210.00
- TMSLL Swingaway Contoura Laterals Lg (E0956)..... 210.00
- TMSLLHW Removable & Adj Hardware (E1028)..... 210.00

STEALTH LATERALS FOR CONTOURA⁹

- TWB-TPR Adj. Swing Contoura Laterals Med (E0956)..... 320.00
- TPRHW Removable & Adj Hardware (E1028)..... 320.00
- TWB-TPRL Adj. Swing Contoura Laterals Lg (E0956)..... 320.00
- TPRLHW Removable & Adj Hardware (E1028)..... 320.00

ELECTRONICS

MK6i™ CONTROLLER

- MK690ACC Controller w/ Actuator Control **Standard**
- GTRAC G-Trac Module (K0108) 995.00¹⁰

EXPANDABLE CONTROLLER SYSTEMS (REQUIRED)

- PWH Harness req. for Expandable Electronics (E2313)350.00¹¹
- MPJM6 MPJ™+ Multiple Drive Joystick (E2377) 750.00
- PSFM6 Personal Joystick, Inductive on Front 1,130.00¹²
- PSRM6 Personal Joystick, Inductive on Rear 1,130.00¹²
- PBOD Push Button On/Off w/PB Drive Select No Charge¹³
- PBSS Push Button On/Off w/Speed Select No Charge¹³
- OR On/Off Mounted Right..... No Charge¹³
- OL On/Off Mounted Left..... No Charge¹³
- ASLRDYM6 MK6 Display Only
(No Driver Control) (E2377)..... 750.00¹⁵

**FOR YOUR ASL DRIVE SYSTEM NEEDS
PLEASE REFER TO THE FOLLOWING:
ASL SINGLE INVOICE PRICE LIST FORM # 02-097.**

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS

- 1500M6 Prop RIM Head Control w/Display (E2327) 2,875.00
- SMHD Small Headrest No Charge
- LGHD Large Headrest No Charge
- 1558MM6 Compact Joystick w/Display (E2373) 1,130.00

ALTERNATIVE DIGITAL CONTROL SYSTEMS

- SNPM6 Sip-N-Puff Dig Interface Box (E2325) 1,500.00¹⁶
- 5018M6 Wafer Board & Dig Interface Box (E2322)..... 1,550.00¹⁶
- 5020M6 Mini Tash Jstck & Dig Interface Bx (E2321).. 1,750.00¹⁶

FOOTNOTES

- ⁸ For use with channel canes only.
- ⁹ For use with Contoura seating only..
- ¹⁰ G-Trac provides Improved Tracking Capabilities and Enhanced Traction. Does not require Special Motors. Not available w/ SPJ+ joysticks.
- ¹¹ Harness is required to interface joystick or MK6i display with expandable controllers. Option is located under Electronic Accessories in configurator
- ¹² TOD-SS: Toggle On/Off Drive Select, Speed Potentiometer standard.
- ¹³ PSRM6 Only.
- ¹⁴ Drive Select through mode switch, or mode switch and Joystick Commands.
- ¹⁵ Display Only, No Driver control. Allows choice from order form.
- ¹⁶ Must also use either a Multiple Drive Joystick or MK6i Display.

SIP-N-PUFF KIT

- PKG32666 Therafin® Sip-N-Puff Breath Tube Kit (E2326) . 380.00

ADDITIONAL CONTROL CHOICES

- 1558M6 Compact Joystick - Less std Mount 770.00
- 5018 Wafer Board..... 310.00¹⁷
- 5020 Mini Tash Joystick..... 405.00¹⁷
- PACM6 Prop Attendant Control w/quick release mount ... 850.00¹⁷
- 1552M Digital Attendant Control 850.00

JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retract. Joystick Mnt MPJ+ (E1028). 279.00¹⁸
- ARM250 Stealth Height Adj Joystick Mount (E1028) 210.00
- LEFT Left Handed Mounted Joystick No Charge
- RIMHW RIM Hardware For 1500 (E1028) 219.00
- PKG32669 Complete Bib Assembly (Tash/1558) 265.00
- GATMPJM6 Gatlin Midline for MPJ (E1028) 600.00
- GAT1812 Gatlin Midline for 1812 (E1028) 600.00
- RCM Rear Cane Mount for Compact Joystick 395.00

ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer w/Pro Memory Card 475.00
- 1813SD Pro Memory Card (USB Ready) 90.00
- 1560 "T" Handle Flexible Joystick Ext. (E2323) 77.00
- 1561 Straight Handle Flexible Joystick Ext. (E2323)..... 77.00
- 1826 Chin Cup (E2324)..... 48.00
- A24VPS 24V MK6 Auxiliary Power Source 200.00
- AUX12M6 Auxiliary Module for 1 & 2 795.00
- AUX34M6 Auxiliary Module for 3 & 4 795.00
- 2PB Dual Push Button..... 147.00¹⁹
- 2SW Dual Toggle 172.00¹⁹
- PBMT Push Button Mounting Kit..... 252.00
- IR Infrared Module with Wireless Mouse Control ... 995.00²⁰

POWERED SEATING ACCESSORIES

- 4WSB Four Way Switch Box No Charge
- FWT 4 Way Toggle No Charge
- QPB 4 Push Buttons..... No Charge
- S4WSB Multiple Actuator Interface Box (E2311) 2,600.00
- FWT 4 Way Toggle No Charge
- QPB 4 Push Buttons..... No Charge



FWT

For systems with two or more actuators, NOT operating through driver control, choose 4WSB with either 4-way toggle or 4-push button.



S4WSB

For systems with 2 or more actuators OPERATING THROUGH DRIVER CONTROL, choose S4WSB. Uses 4-way toggle or 4-push button as attendant control.

ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Switch - Black 156.00
- ASL304 Light Wobble Switch w/611 Mount 432.00
- ASL202J Fiber Optic Switch..... 610.00²¹
- ASL208J Adj Proximity Switch 550.00²¹
- ASL504M6 Remote Emergency Stop Switch 1,185.00
- ASL611 Adjustable Switch Mounting Bracket 252.00

FOOTNOTES

- ¹⁷ Requires SNPM6 Sip-N-Puff Interface Box
- ¹⁸ QLAM6 include ARM250 standard
- ¹⁹ Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options
- ²⁰ Capable of controlling up to six remote Infra-Red controlled devices and also Computer Mouse operation through the driver control. Includes stereo switch (with hook/loop mount) for mouse clicks
- ²¹ Must add MK6 Auxiliary Power Source.

NOTE: All specifications and dimensions are approximate.
For the most current pricing information go to www.invacare.com.

Invacare® Formula™ CG Tilt/Recline/Elevate for TDX® SP Accessories and Options

BATTERY TRAY TYPE

- 22TRY 22NF Style Tray No Charge

BATTERIES (REQUIRED) ²²

- 22NFBATTERY 22NF Battery (On Chair) (E2361) 300.00

BATTERY CHARGER OPTION

- 110CHARGER 110 Volt Battery Charger No Charge
- LESTR Lester 8AMP Battery Charger (E2367) 320.00
- L8900 Omit Battery Charger Less 80.00

BACK ACCESSORIES

- L8601 Omit Back Upholstery No Charge
- RBP Reinforced Back Upholstery No Charge

MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks 45.00
- ADD Adductors Small (E0956) 200.00
- ADDHW Hardware for ADD (E1028) 400.00
- ADD L Adductors Large (E0956) 200.00
- ADDHWL Hardware for ADD L (E1028) 400.00

FRONT RIGGING OPTIONS

FRONT RIGGING ²³

CENTER MOUNT ²⁴

- CTMT-F Center Mount Adjustable Knee Angle 300.00 ²⁵
- CTMT-FJ Jr. Center Mount Adjustable Knee to Heel 300.00 ²⁶
- SM Small Footplates (5.25" x 7.5") No Charge
- MED Medium Footplates (5.50" x 9.5") No Charge
- LRG Large Footplates (5.90" x 11.5") No Charge
- CALF Calf Pads For Center Mount 200.00

CENTER PIVOT STYLE (PAIR) ²⁷

- 70NHD 70° Swingaway Footrests 13.5"-17.5" Standard ²⁸
- 70STAPER 70° Tapered Footrests 13"-17" No Charge ²⁹
- AT5544 Manual Elevating Legrests (E0990) 255.00 ³⁰

PIN STYLE-HEMI (PAIR)

- PHW93 70° Swingaway Footrests 13.5" - 19" No Charge
- PHAL4A Smartleg Art. Legrests 17.25" - 21.5" (K0053) . . 225.00

POWER ELEVATING LEGRESTS-ARTICULATING

- MPCTMT Power Center Mount Foot Platform (K0108) ... 2,995.00 ³¹
- IFP Individual Footplates 250.00
- ELRPW Power ELR w/Articulation 20" Ext (E1010)..... 1,895.00 ³²
- ELRPW22 Power ELR w/Articulation 22" Ext (E1010)..... 1,995.00 ³³

FOOTPLATE OPTIONS (PAIR)

- 1651 Composite Flip-Up Footplates No Charge ³⁴
- 1350 Extra Lg Aluminum Footplates No Charge ³⁵
- AT5543 Adj Angle Flip-Up Footplates (K0040) 165.00
- CTMTFB Footboard for Center Mount No Charge ³⁶

FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap (E0038) 27.00
- 1600BK Heel Loops w/Ankle Straps (E0951) 56.00
- IMPACTG Impact Guards 51.00
- LPT2 Longer Pivot Slide Tube For 1350 83.00
- ALPT Longer Pivot And Slide Tube 84.00
- ALPT4 Longer Pivot And Slide Tube 4" 84.00

SEATING & POSITIONING

**FOR YOUR SEATING NEEDS
PLEASE REFER TO THE FOLLOWING
REHAB SEATING & POSITIONING PRODUCTS
See Price List Form 07-027**

Power Center Mount KTH	
MPCTMT	13" - 17"
MPCTMT-L	15" - 19"

FOOTNOTES

- ²² Sealed Gel Batteries highly recommended for use in Invacare power wheelchairs.
- ²³ All front riggings include a standard footplate .
- ²⁴ Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring.
- ²⁵ Extension range 12.5"-15". Must choose either SM, MED or LRG footplate.
- ²⁶ Extension range 9.25"-12". Must choose either SM, MED or LRG footplate.
- ²⁷ Provides up to 1 ½" of depth adjustment in ½" increments.
- ²⁸ With seat widths of 21" to 24" and depth of 16", must pick footplate 1350. The AT5543 is not available..
- ²⁹ Standard with impact guards.
- ³⁰ Not available with 16" wide by 16" to 18" deep or 17" wide by 16" deep and low seat-to-floor.

FOOTNOTES

- ³¹ Includes a standard, single footboard.
- ³² Independent powered elevating legrests including elevating and articulating legrests with angle adjustable footplates. Not available with the 16" or 17" seat depths. Extension range 13"-20". Sold as pair.
- ³³ Independent powered elevating legrests including elevating and articulating legrests with angle adjustable footplates. Not available with the 16" or 17" seat depths. Extension range 15"-22". Sold as pair.
- ³⁴ Not available with AT5544 and 17" wide by 16" deep and low seat-to-floor.
- ³⁵ Not available with AT5544 and 16" wide by 16"-18" deep or 17" wide by 16" deep and low seat-to-floor.
- ³⁶ One piece footboard for Center Mount (CTMT-F) Riggings.

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information go to www.invacare.com.**

Invacare® Formula™ CG Tilt/Recline/Elevate for TDX® SP Accessories and Options

Battery Boxes
*Standard Battery Tray
22NF with Elevate Feature
*Sealed Gel Cell Batteries Recommended

TDXSP Tilt/Recline/Elevate Specifications	
Maximum Speed MPH	5.8
4 Pole	
Seat-To-Floor at 0°	
Low	18.5"

Batteries

22NF batteries with terminal configuration (positive on the left and negative on the right) MUST be used. 22NF batteries that have the reverse terminal configuration MUST not be used. Terminals MUST have a cross hole for proper battery connection.

MEMORY CARD - BASIC

- Standard with ALL Expandable Control Systems (MPJ+, PSR+, PSF+, MK6i Display).
- Used to Back Up/Restore Programmed Settings for Only One Chair.
- May be given to user or kept with provider's files for safe keeping.
- Does not Contain Advanced Diagnostics/Help Screens/File Structure.
- Not compatible with SPJ+ Joysticks.

MEMORY CARD - PROFESSIONAL

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains Advanced Diagnostics, Help Screens, Software Updates, File Storage/Retrieval.
- Used to back up multiple chairs and programming settings.
- Intended for use by qualified providers only.
- Not compatible with SPJ+ Joysticks.

Contoura Seat Sizes							
	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X
Available with OMCU option only.							
Back Heights	25" or 29.5"						

Conventional Seat Sizes	
Seat Widths	16"-22"
Seat Depths	16"-22"
Back Heights	20"-26"

- | TDXSP Tilt/Recline/Elevate Features and Specifications |
|--|
| <ul style="list-style-type: none"> • Standard with 300 lb weight capacity. • Tilt range from 0° to 55° • Recline range of 170°. • Elevate range of 7". • Mechanical shear reduction of 3". • Conventional or Contoura style seating. • Adjustable seat frame with Conventional style. • Seat-to-floor height of 18.5" at 0° tilt. • Speed inhibit feature of 20% of maximum speed while elevated. |

- | After Market Backrests |
|---|
| <p>The following companies also have hardware to make their backrests work with our new canes:</p> <ul style="list-style-type: none"> • AEL – planar seating system • Freedom Designs® - planar seating system • AES – planar seating system • Canyon – planar seating system |

- | Tech Tips |
|---|
| <p>Operating Powered Seating through the Driver's Control.
All chairs require adding or replacing one component to allow actuator operation through the driver control.</p> <p>SPJ+ with PSS:
Replace joystick with SPJ+ with ACC (Actuator Control) -E2310.</p> <p>Single Actuator Systems with MPJ+, PSR+, PSF+ MK6i Display:
Add a SANODE (Single Actuator Node) to the system - E2310.</p> <p>Multiple Actuator Systems:
Replace the standard 4 Way Switch box with a Multiple Actuator Interface Box - E2311.</p> |

SEATING & POSITIONING
FOR YOUR SEATING NEEDS
PLEASE REFER TO THE FOLLOWING
REHAB SEATING & POSITIONING PRODUCTS
SEE PRICE LIST FORM 07-027

Invacare Corporation www.invacare.com



USA
One Invacare Way
Elyria, OH USA
44035-4190
800-333-6900

All rights reserved. Trademarks are identified by the symbols™ and®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Products designated by ASL are Manufactured by Adaptive Switch Labs, Inc. Therafin is a registered trademark of W/R Medical Electronics Co. Otto Bock is a trademark of Otto Bock Orthopaedische Industries GMBH & Co. Freedom Designs is a registered trademark of Freedom Designs, Inc. Specifications and prices are subject to change without notification.
©2008 Invacare Corporation