



Yes, you can.®

The Invacare® IVC™ Transport Chairs are designed for safe and comfortable patient transport. The frames weigh only 24.5 lb. and incorporate a fold-down back, making these chairs easy to transport and store.



## Features

- Durable, low maintenance, powder-coated carbon steel frame is long-lasting
- Footrests are now included with the chair and match frame finish
- Black nylon upholstery is flame retardant and bacteria resistant
- Full-length arms are padded for comfort
- Carry straps on seat allow for ease of folding and lifting
- Fold-down back makes transport and storage easy
- Weighs 24.5 lb. (w/o riggings)

TRAN19FR, IVC Transport chair.

## Specifications

Model Numbers	TRAN17FR	TRAN19FR
<b>Chair Dimensions</b>		
Seat width	17"	19"
Seat depth	16"	16"
Back height	19"	19"
Seat-to-floor height	19"	19"
<b>Overall Dimensions</b>		
Overall width	19"	21"
Folded width	7.25"	7.75"
Overall depth w/o riggings	27"	27"
Overall height	39"	39"
Weight w/o riggings	24.5 lb.	25 lb.
Weight limitation	250 lb.	250 lb.

Note: All specifications are approximate

**Invacare Corporation**  
www.invacare.com

### USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

### Canada

570 Matheson Blvd. E.  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

All rights reserved.  
Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.

© 2007  
Invacare Corporation

Form No. 07-058



# Transport Chairs

Yes, you can:

