



Yes, you can.®

### Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email\* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

\*Saving and emailing the eForms is available with Acrobat Professional software

#### If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

**Please note:** The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

#### If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit [www.adobe.com](http://www.adobe.com).

# Invacare® TDX® SI Rehab Seat Power Wheelchair

## Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



Yes, you can.®

### TRANSACTION TYPE:

Order  Quote

Date \_\_\_\_\_ Account # \_\_\_\_\_ P.O. # \_\_\_\_\_

Company Name \_\_\_\_\_ Phone # \_\_\_\_\_

Address \_\_\_\_\_ City \_\_\_\_\_ State \_\_\_\_\_ Zip Code \_\_\_\_\_ Country \_\_\_\_\_

Contact Name \_\_\_\_\_ Comments \_\_\_\_\_

### Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

### CHAIR TYPE

- TDXSI TDX SI with MK6i™ ..... \$5,970.00 K0848  
**TDX SI comes standard with: 300 lbs. user weight, 5 mph top speed, 2-pole motor package an MK660.**
- TDXSISEAT TDXSI Seat Only ..... \$755.00

Group 3 - K0848	
TDX SI with MK6i and Adjustable Rehab Seat	
300 lbs. Weight Capacity	
1431-3 14" x 3" Tires with Foam Filled Inserts	Standard
2-Pole Motors	Standard
Flip Back Fixed Height Arm	Standard
MK660 60 Amp Controller	Standard
SPJPLUS SPJ™+ Single Drive Joystick	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard

### TRANSPORT TIE-DOWN

Clear

- TRRO Transport Ready Option ..... \$300.00
- TRBKTS Wheelchair Transport Brackets ..... \$250.00

#### TRRO Features for TDX SI

<h5 style="text-align: center;">Required Features for TRRO</h5> <ul style="list-style-type: none"> <li>• TDX SI</li> <li>• 70° Swingaway, 70° Tapered, or Manual Center Mount Front Riggings</li> <li>• T-arms required with 12"-15" seat depth</li> <li>• 18F,19F,28F,29F for 12"-16" widths only available with 16", 17" or 18" deep</li> <li>• Adjustable Rehab Seat with Adjustable Seat Tilt and Back Angle</li> <li>• Sectional or Straight Rear Uprights or No Push Handles</li> <li>• 14" x 3" Tires with Foam Filled Inserts</li> <li>• 90° or 95° Back Angle</li> <li>• 18" - 24" ASBA Back Height</li> <li>• Any Seat Width</li> <li>• Any Seat Depth</li> <li>• 300 lbs.</li> <li>• Low Seat-to-Floor Height</li> <li>• 22NF Batteries</li> <li>• Any Color Frame Finish</li> </ul>	<h5 style="text-align: center;">Optional Features for TRRO</h5> <ul style="list-style-type: none"> <li>• Adjustable Angle Flip Up Footplates, Flip-Up Composite Footplates or Aluminum Footplates-Deluxe</li> <li>• Wheel Locks</li> <li>• Calf Strap</li> <li>• Heel Loops with Ankle Straps</li> <li>• Impact Guards</li> <li>• One Sided Fork Package</li> <li>• SPJ+, CMPJ+ Joystick</li> <li>• All PSR+ / PSF+ Joysticks</li> <li>• Proportional Attendant Control</li> <li>• "T" Handle Flex. Joystick Ext.</li> <li>• Left Hand Mounted Joystick</li> <li>• Straight Handle Flexible Joystick Ext.</li> <li>• Easy Remote Programmer</li> <li>• Reinforced Back Upholstery</li> </ul>
---	--

#### TRANSPORT TIE-DOWN\*

TRRO includes four factory-installed transport brackets and a wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNA WC Vol 1. Section 19 Frontal Impact Test requirements for wheelchairs with a 168 lbs. crash dummy, which corresponds to a person with a weight of 114 lbs. to 209 lbs.

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

\* Transport Tie-Down option is not retrofittable to existing models and is not field serviceable.

### USER WEIGHT LIMIT

U300 User Weight <300 lbs ..... Standard K0849

### WHEEL OPTIONS (PAIR)

**Free Wheel Hub will add 3/4" overall width. Does not include TDX Series design wheel. Use specials #SM-07050811239**

1431-3 14" x 3" with Foam Filled Inserts ..... Standard

### FORK OPTION

FKPKG One-Sided Fork Package ..... \$500.00

**NOTE: All specifications and dimensions are approximate.  
For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

# Invacare® TDX® SI Rehab Seat Power Wheelchair

## FRAME FINISH

### CONTEMPORARY

- 24P Wet Black ..... No Charge
- 60P Silver Metallic ..... No Charge
- 105P Midnight Blue ..... No Charge
- 115P Black Prism ..... No Charge

### TOO HOT

- 30P Sunny Yellow ..... No Charge
- 61P Electric Red ..... No Charge
- 119P Bubblegum Pink ..... No Charge
- 120P Grape Madness ..... No Charge

### WAY COOL

- 62P Electric Blue ..... No Charge
- 127P Grasshopper Green ..... No Charge
- 116P Sunkissed Orange ..... No Charge
- 117P Tangerine Red ..... No Charge
- 29P White Snow ..... No Charge

Discontinued colors 71P Silver Vein, 104P Electric Teal, 121P Lolly Pop Blue, 122P Cosmic Blue, and 125P Emerald Green can be order through specials number C-070109031327

## SEAT

### SEAT-TO-FLOOR HEIGHT

Due to variation in options the listed seat-to-floor heights have a tolerance of +/-0.5"

- LOW 16 1/2" at 0° ..... Standard
- MED 18 1/2" at 0° ..... No Charge
- TALL 20 1/2" at 0° ..... No Charge

The long extension manual center mount front rigging to provide up to 18" KTH is available using specials #C-07025041441. The seat angle is set at 5 degrees.

### SEAT SELECTION

- ADJASBA Adj. Rehab Seat with Adj. Seat Tilt and Back Angle . Standard
  - ADJRECL Manual Recline for Adjustable ASBA ..... \$744.00 E1226
  - CURV Curved ..... No Charge E2620
  - CTB Contoura Seating ..... \$625.00 E2621
  - OMCU Omit Contoura Seat Cushion ..... Less \$20.00
  - CTBMSH Contoura seating with mesh cover ..... \$635.00 E2620 E2621
  - OMCU Omit Contoura Seat Cushion ..... Less \$20.00
  - VAN Van Style Seat ..... \$70.00
- Not available with Transport Tie Down options (TRRO). Seat-to-Floor Height 21".

### ADJUSTABLE SEAT WIDTH RANGE

- 1216W 12" to 16" Seat Width Range ..... No Charge
  - 1620W 16" to 20" Seat Width Range ..... No Charge
  - 2022W 20" to 22" Seat Width Range ..... No Charge
- Must select 12" - 18" Depth Range (1218D) when ordering the 12" - 16" Seat Width Range (1216W).

## SEAT WIDTH SETTING

Contoura only available with even Seat Widths 16" - 22".

- W12 Seat Width 12" ..... \$110.00
  - W13 Seat Width 13" ..... \$110.00
  - W14 Seat Width 14" ..... \$110.00
  - W15 Seat Width 15" ..... No Charge
  - W16 Seat Width 16" ..... No Charge
  - W17 Seat Width 17" (not available with Contoura) ..... No Charge
  - W18 Seat Width 18" ..... No Charge
  - W19 Seat Width 19" (not available with Contoura) ..... No Charge
  - W20 Seat Width 20" ..... No Charge
  - W21 Seat Width 21" (not available with Contoura) ..... \$329.00
  - W22 Seat Width 22" ..... \$329.00
  - VW16 Seat Width 16" Van ..... No Charge
  - VW18 Seat Width 18" Van ..... No Charge
  - VW20 Seat Width 20" Van ..... \$329.00
  - VW22 Seat Width 22" Van ..... \$329.00
- User weight limit is 165 lbs.
- Seat Widths and Seat Depths for Van seat Semi Recline with Solid Base.

### ADJUSTABLE SEAT DEPTH RANGE

- 1218D 12" to 18" Seat Depth Range ..... No Charge
  - 1619D 16" to 19" Seat Depth Range ..... No Charge
  - 1922D 19" to 22" Seat Depth Range ..... No Charge
- Must select 12" - 18" Depth Range (1218D) when ordering the 12" - 16" Seat Width Range (1216W)

### SEAT DEPTH

- D12 Seat Depth 12" ..... \$110.00
  - D13 Seat Depth 13" ..... \$110.00
  - D14 Seat Depth 14" ..... \$110.00
  - D15 Seat Depth 15" ..... No Charge
  - D16 Seat Depth 16" ..... No Charge
  - D17 Seat Depth 17" ..... No Charge
  - D18 Seat Depth 18" ..... No Charge
  - D19 Seat Depth 19" ..... No Charge
  - D20 Seat Depth 20" ..... No Charge
  - D21 Seat Depth 21" ..... \$529.00
  - D22 Seat Depth 22" ..... \$529.00
  - VD16 Seat Depth 16" Van ..... No Charge
  - VD17 Seat Depth 17" Van ..... No Charge
  - VD18 Seat Depth 18" Van ..... No Charge
  - VD19 Seat Depth 19" Van ..... No Charge
  - VD20 Seat Depth 20" Van ..... \$529.00
- User weight limit is 165 lbs.
- Seat widths and Seat Depths for Van seat Semi Recline with Solid base.
- Seat widths and Seat Depths for Van seat Semi Recline with Solid base.
- Seat widths and Seat Depths for Van seat Semi Recline with Solid base.
- Seat widths and Seat Depths for Van seat Semi Recline with Solid base.
- Seat widths and Seat Depths for Van seat Semi Recline with Solid base.
- Seat widths and Seat Depths for Van seat Semi Recline with Solid base.

**NOTE: All specifications and dimensions are approximate.**  
 For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).

# Invacare® TDX® SI Rehab Seat Power Wheelchair

## SEAT Continued

### VAN UPHOLSTERY

- GV Grey Vinyl ..... No Charge  
*Only available with VAN seat selection and SEMIRB seat style.*

### VAN SEAT STYLE

- SEMIR Semi-Recline ..... Standard  
*Only available with VAN seat selection and GV Grey Vinyl.*

### FIXED ARM TYPE (PAIR)

- Clear  
 18F Full Length Fixed Height Conv. .... No Charge  
 28F Desk Length Fixed Height Conv. .... Standard

### ADJUSTABLE ARM TYPE (PAIR)

- Clear  
 19F Full Length Adj Hgt Conv. .... \$230.00 E0973  
 29F Desk Length Adj Hgt Conv. .... \$230.00 E0973  
 TADM T-arms Desk Lgth Adj Hgt Arm (8" - 12") ..... \$230.00 E0973  
 TAFM T-arms Desk Lgth Adj Hgt Arm (8" - 12") ..... \$230.00 E0973  
 VAH Full Length Adjustable Flip-Up ..... \$230.00 E0973  
*Arms only available with VAN seat selection and SEMIRB seat style*  
 VFH Full Length Fixed Height ..... Standard  
*Arms only available with VAN seat selection and SEMIRB seat style*  
 VAHD Desk Length Adjustable Height ..... \$230.00 E0973  
*Arms only available with VAN seat selection and SEMIRB seat style*  
 VFHD Desk Length Fixed Height ..... No Charge  
*Arms only available with VAN seat selection and SEMIRB seat style*

### MIXED ADJUSTABLE ARM TYPE (EACH)

- Right Arm**  Clear  
 19FR Full Length Adj. Hgt. Right ..... \$115.00 E0973  
 29FR Desk Length Adj. Hgt. Right ..... \$115.00 E0973  
**Left Arm**  Clear  
 19FL Full Length Adj. Hgt. Left ..... \$115.00 E0973  
 29FL Desk Length Adj. Hgt. Left ..... \$115.00 E0973

### CANTILEVER ARM TYPE (PAIR)

- Clear  
*Cantilever Arms only available with 300 lbs. rated systems. Not available with BPO, Otto Bock, or 12" - 15" seat width. Must use hard back. See chart on last page for arm height adjustment range.*  
 CLVF Full Length Arm ..... \$350.00  
 CLVD Desk Length Arm ..... \$350.00

### CANTILEVER MIXED ARM TYPE (EACH)

- Right Arm**  Clear  
 CLVFR Full Length Arm Right ..... \$175.00  
 CLVDR Desk Length Arm Right ..... \$175.00  
**Left Arm**  Clear  
 CLVFL Full Length Arm Left ..... \$175.00  
 CLVDL Desk Length Arm Left ..... \$175.00

### ARMPAD OPTIONS (CANTILEVER ARMS)

- Clear  
 WF Waterfall Armpads ..... Standard  
 RCL Reclining Armpads ..... No Charge

### SEAT BACK TYPE

- SECT Sectional Rear Uprights ..... Standard  
 L8661 Omit Rear Uprights ..... Less \$20.00  
 SRUP Straight Rear Uprights ..... No Charge  
 L8661 Omit Rear Uprights ..... Less \$20.00  
 9720P No Push Handles (Straight Canes) ..... \$82.00

### SEAT BACK ANGLE DEGREE

- 90BA 90° Back Angle ..... Standard  
 95BA 95° Back Angle ..... No Charge

### SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Positioning Strap ..... \$50.00  
*1311BK and 7311BK Not available on 12" - 15" wide.*  
 1515 Push Button Style Seat Positioning Strap ..... Standard

### CHEST POSITIONING STRAPS

- Clear  
 7311BK Airline Buckle Seat Positioning Strap ..... \$100.00 E0960  
*1311BK and 7311BK not available on 12" - 15" wide.*  
 9515 Push Button Style Chest Positioning Strap ..... \$100.00 E0960  
 7321BK Hook and Loop Chest Positioning Strap ..... \$100.00 E0960

### BACK HEIGHT

- Clear  
*Back Heights 23" and 24" not available with 12" - 16" widths. With Junior seat widths 12" - 16", the back height is adjustable.*  
 TSBH12 12" Back Height ..... No Charge  
 TSBH13 13" Back Height ..... No Charge  
 TSBH14 14" Back Height ..... No Charge  
 TSBH15 15" Back Height ..... No Charge  
 TSBH16 16" Back Height ..... No Charge  
 TSBH17 17" Back Height ..... No Charge  
 TSBH18 18" Back Height ..... No Charge  
 TSBH19 19" Back Height ..... No Charge  
 TSBH20 20" Back Height ..... No Charge  
 TSBH21 21" Back Height ..... No Charge  
 TSBH22 22" Back Height ..... No Charge  
 TSBH23 23" Back Height ..... No Charge  
 TSBH24 24" Back Height ..... No Charge

### CONTOURA BACK HEIGHT

- Clear  
 TSBH22 22" Back Height ..... No Charge  
 TSBH25 25" Back Height ..... No Charge  
 TSBH28 28" Back Height ..... No Charge

### BACK HEIGHT RECLINER

- Clear  
 BHR18 18" Back Height ..... No Charge  
 BHR20 20" Back Height ..... No Charge  
 BHR22 22" Back Height ..... No Charge  
 BHR24 24" Back Height ..... No Charge  
 BHR26 26" Back Height ..... No Charge

## ELECTRONICS

### MK6i™ CONTROLLERS

- MK660ACC Controller with Actuator Control ..... Standard  
 GTRAC G-Trac Module ..... \$995.00 K0108  
*G-trac provides Improved Tracking Capabilities and Enhanced Traction. Does not require Special Motors. Not available with SPJ+ Joysticks.*

### STANDARD DRIVER CONTROL

- SPJPLUS SPJ+ Single Drive Joystick ..... Standard

**NOTE: All specifications and dimensions are approximate.**  
*For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).*

# Invacare® TDX® SI Rehab Seat Power Wheelchair

## ELECTRONICS Continued

### EXPANDABLE CONTROLLER SYSTEM

Clear

- PWH Harness req. for Expandable Electronics ..... \$350.00 E2313  
**Required to interface joystick or MK6i display with expandable controllers.**
- CMPJM6U MPJ+ Multiple Drive Joystick with **Color Display** .. \$750.00 E2377
- ASLRDYM6 MK6 Display Only (No Driver Control) ..... \$750.00 E2377  
**Display only, No Driver Control. Allows choice from order form.**
- MCD Mini Color Display ..... \$750.00 E2377  
**No port for remote On/Off switch. Stereo port only.**

FOR YOUR ASL DRIVE SYSTEM NEEDS  
PLEASE REFER TO THE FOLLOWING:  
ASL SINGLE INVOICE PRICE LIST FORM # 02-097.

### ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS

Clear

- 1500M6 Prop RIM Head Control with MK6i Display .. \$2,875.00 E2327
- SMHD Small Headrest ..... No Charge
- LGHD Large Headrest ..... No Charge
- 1558MM6 Compact Joystick with Display ..... \$1,130.00 E2373

### ALTERNATIVE DIGITAL CONTROL SYSTEMS

Clear

**Must also use either a Multiple Drive Joystick or MK6i Display.**

- SNPM6 Sip-N-Puff Digital Interface Box ..... \$1,500.00 E2325
- 5018M6 Wafer Board and Digital Interface Box ..... \$1,550.00 E2322
- 5020M6 Mini Tash Joystick and Digital Interface Box \$1,750.00 E2321

### SIP-N-PUFF KIT

- PKG32666 Therafin® Sip-N-Puff Breath Tube Kit ..... \$380.00 E2326

### ADDITIONAL CONTROL CHOICES

- 1558M6 Compact Joystick - Less Std Mount ..... \$770.00
- 5020 Mini Tash Joystick ..... \$405.00  
**Requires SNPM6 Sip-N-Puff Interface Box.**
- PACM6 Prop Attendant Control with quick release mount .... \$850.00
- 1552M Digital Attendant Control ..... \$850.00  
**Requires SNPM6 Sip-N-Puff Interface Box.**

## JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retract. Joystick Mnt CMPJ+ ..... \$279.00 E1028  
**Comes standard with ARM250. Not available with TRRO.**
- QLA Quad Link Retract. Joystick Mnt SPJ+ ..... \$279.00 E1028  
**Come standard with ARM250. Not available with TRRO.**
- ARM250 Stealth Height Adjustable Joystick Mount ..... \$210.00 E1028
- LEFT Left Handed Mounted Joystick ..... No Charge
- RIMHW RIM Hardware For 1500 ..... \$219.00 E1028
- PKG32669 Complete Bib Assembly (Tash/1558) ..... \$265.00
- GATMPJM6 Gatlin Midline for CMPJ+ ..... \$600.00 E1028
- GAT1812 Gatlin Midline for 1812 ..... \$600.00 E1028
- RCM Rear Cane Mount for Compact Joystick ..... \$395.00

## ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer with Pro Memory Card ..... \$475.00
- 1813SD Pro Memory Card (USB Ready) ..... \$90.00
- 1560 "T" Handle Flexible Joystick Ext. .... \$77.00 E2323
- 1561 Straight Handle Flexible Joystick Ext. .... \$77.00 E2323
- 1826 Chin Cup ..... \$48.00 E2324
- A24VPS 24V MK6 Auxiliary Power Source ..... \$200.00
- AUX12M6 Auxiliary Module for outputs 1 and 2 ..... \$795.00
- AUX34M6 Auxiliary Module for outputs 3 and 4 ..... \$795.00
- 2PB Dual Push Button ..... \$147.00
- 2SW Dual Toggle ..... \$172.00  
**Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options.**
- PBMT Push Button Mounting Kit ..... \$252.00
- MPY Monoport "Y" Cable for CMPJ+ or Display ..... \$15.00
- RFM Wireless (RF) Mouse Only Module ..... \$750.00  
**RFM includes Transmitter, USB Radio Frequency Mouse Dongle (Receiver), and Stereo Switch with hook and loop mount for mouse clicks. Not compatible with SPJ+.**
- IR Infrared Module with Wireless RF Mouse Control .... \$995.00  
**Capable of operating up to six remote Infra-Red controlled devices. Also allows wireless (RF) Mouse operation through the driver control. Compatible with RFM module for control of two wireless mice from one chair. Includes two button stereo switch (with hook/loop mount) for mouse clicks. Not compatible with SPJ+ Joysticks.**

## ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Switch - Black ..... \$156.00
- ASL304A Light Wobble Switch ..... \$180.00
- ASL202J Fiber Optic Switch ..... \$732.00  
**Not available with SPJ+ joystick. For ASL202J, ASL208J must add ASL804M6.**
- ASL208J Adj Proximity Switch ..... \$660.00  
**Capable of controlling up to six remote Infra-Red controlled devices and also Computer Mouse operation through the driver control. Includes stereo switch (with hook/loop mount) for mouse clicks.**
- ASL504M6 Remote Emergency Stop Switch ..... \$1,185.00
- ASL611 Adjustable Switch Mounting Bracket ..... \$252.00

**NOTE: All specifications and dimensions are approximate.  
For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

# Invacare® TDX® SI Rehab Seat Power Wheelchair

## MISCELLANEOUS

### BACK ACCESSORIES

- Clear
- L8601 Omit Back Upholstery ..... No Charge
- RBP Reinforced Back Upholstery ..... No Charge

### MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks ..... \$45.00
- SEAT02 Oxygen Holder (ADJASBA Seat Only) ..... \$150.00 E2208

### ADDUCTORS

Hardware is packaged with equipment and cannot be ordered separately.

- Clear
- ADD Adductors Small ..... \$200.00 E0956
- ADDHW Hardware for ADD ..... \$400.00 E1028
- ADDL Adductors Large ..... \$200.00 E0956
- ADDHWL Hardware for ADDL ..... \$400.00 E1028

### BATTERY TRAY TYPE

22TRY Base Width 24".

- 22TRY 22NF Style Tray ..... Standard

### INSTALLATION OF BATTERIES

- 22NFBATTERY 22NF Battery(On-Chair)(1 pair) ..... \$300.00 E2361

### BATTERY CHARGER OPTION

- Clear
- 110CHARGER 110 Volt Battery Charger ..... No Charge
- LESTR Lester 8 AMP Battery Charger ..... \$320.00 E2367
- L8900 Omit Battery Charger ..... Less \$80.00

### VENT TRAY

Approximate dimensions are 11" W x 4" D x 13" H. Designed for use with the latest and most compact style vent systems.

- COMP Compact Vent Tray ..... \$1,540.00 E1030

## FRONT RIGGINGS

### HANGER BRACKET TYPE

- Clear
- ASB Center Pivot Style Bracket ..... No Charge
- PSB Pin Style Bracket ..... No Charge

### CENTER MOUNT

For a 13.25" - 17.25" extension range, available through special #SM-071217095848 with Tall STF Only and 5 Degree Seat-Tilt Angle.

- CTMT-F Center Mount Adj. Knee Angle (12.25" - 15") ..... \$300.00
- Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 12.5" - 15". Footplates are adjustable angle. Must choose either SM, MED or LRG footplate.

Clear

- SML Small Footplates (5.25" x 7.5") ..... No Charge
- MDM Medium Footplates (5.5" x 9.5") ..... No Charge
- LRG Large Footplates (5.90" x 11.5") ..... No Charge
- CALF Calf Pads For Center Mount ..... \$200.00

- CTMT-FJ Jr. Center Mount Adj. Knee Angle (9.25" - 12") ..... \$300.00
- Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 9.25" - 12". Footplates are adjustable angle. Available on low seat-to-floor only. Must choose either SM, MED or LRG footplate.

Clear

- SML Small Footplates (5.25" x 7.5") ..... No Charge
- MDM Medium Footplates (5.5" x 9.5") ..... No Charge
- LRG Large Footplates (5.90" x 11.5") ..... No Charge
- CALF Calf Pads For Center Mount ..... \$200.00

### CENTER PIVOT STYLE (PAIR)

Clear

- 70NHD 70° Swingaway Footrests 13 1/2" - 17 1/2" ..... Standard
- Standard with 1651 composite footplate. With seat widths of 21" to 22" and depth of 16", must pick footplate 1350. The AT5544 is not available.

- 70STAPER 70° Tapered Footrests 13" - 17" ..... No Charge
- Standard with impact guards and adjustable angle flip up footplates (AT5543).

- AT5544 Manual Elevating Legrests ..... \$255.00 E0990
- Not available with 16" wide by 16" to 18" deep or 17" wide by 16" deep and low seat-to-floor.

### PIN STYLE-HEMI (PAIR)

Clear

- PHW93 70° Swingaway Footrests 13 1/2" - 19" ..... No Charge
- Standard with 1651 composite footplate.

- PHAL4A Manual Art./ELR 17 1/4" - 21 1/2" ..... \$225.00 K0053
- Not available with 12" or 13" width seats.

### FOOTBOARD

Footboard only available with SEMIRB seat style.

- 90FB Footboard ..... \$150.00

### FOOTPLATE OPTIONS (PAIR)

All front rigging options include a standard footplate.

- 1651 Composite Flip Up Footplates ..... No Charge
- 1350 Extra Large Aluminum Footplates ..... No Charge

Not available with AT5544 and 16" wide by 16" - 18" deep or 17" wide by 16" deep and low seat-to-floor.

- AT5543 Adjustable Angle Flip-Up Footplates ..... \$165.00 K0040

- CTMTFB Footboard for Center Mounts ..... No Charge
- One piece footboard for Center Mount (CTMT-F) Riggings.

### FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap ..... \$27.00 E0038
- IMPACTG Impact Guards ..... \$51.00
- HLOOP Heel Loops (pair) ..... \$45.00 E0951
- LPT2 Longer Pivot Slide Tube For 1350 ..... \$83.00
- ALPT Longer Pivot And Slide Tube ..... \$84.00
- ALPT4 Longer Pivot And Slide Tube 4" ..... \$84.00

**NOTE: All specifications and dimensions are approximate.**  
**For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

# Invacare® TDX® SI Rehab Seat Power Wheelchair

TDX SI Specifications		
<b>Seat Width Range (in 1" increments)</b>	Adult (Adjust)	Junior (Adjust)
ADJASBA	16" - 24"	12" - 16"
SEMIRB (VAN)	16", 18", 20", 22"	
<b>Seat Depth Range</b>		
ADJASBA	12" - 22"	12" - 18"
SEMIRB (VAN)	16" - 20"	
<b>Back Height Range</b>		
	12" - 24"	12" - 22" Adj. Height
Recliner Back Heights	18", 20", 22", 24", 26"	
Contoura Back	25" and 29.5"	
SEMIRB	18" without headrest	
SEMIRB	25" with headrest	
<b>Back Angle Range</b>	85° to 105° in 5" increments except recliner	
<b>Seat-to-Floor (approx)</b>		
Standard	16.5" at 0° (at 5° add .5")	
Medium	18.5" at 0° (at 5° add .5")	
Tall	20.5" at 0° (at 5° add .5")	
<b>Overall Width of Base</b>	24" (w/o Joystick)	
<b>Overall Height</b>	(with 20" back height and 0° seat angle)	
Standard	35 1/2"	
Medium	37 1/2"	
Tall	39 1/2"	
SEMIRB	38" without headrest	
	45" with headrest	
<b>Overall Length</b>		
without Front Rigging	35 1/4"	
<b>Weight</b>		
2 Pole Motor ADJASB	without 22NF Batteries, 166 lbs.	
	with two 22NB Batteries, 264 lbs.	
SEMIRB	without 22NF Batteries, 179 lbs.	
	with two 22NB Batteries, 277 lbs.	
<b>Drive Wheels/Tires (foam filled)</b>		
Standard	14" x 3"	

Cantilever Arm Height Adjustment			
Back Type	Height	Adjustment	Adjustment with laterals
CURV	20" or greater	11" to 16"	11" to 14"
CURV	19"	11" to 15"	11" to 13"
CURV	18"	11" to 14"	11" to 12"
CURV	17"	11" to 13"	no adjustment
CURV	16"	11" to 12"	not available
Contoura	22", 25" to 28"	12 1/2" to 14 1/2"	no adjustment

### Memory Card - BASIC

- Standard with ALL Expandable Control Systems (CMPJ+, PSR+, PSF+, MK6i Display).
- Used to Back Up/Restore Programmed Settings for Only One Chair.
- May be given to user or kept with provider's files for safe keeping.
- Does not Contain Advanced Diagnostics/Help Screens/File Structure.
- Not compatible with SPJ+ Joysticks.

### Memory Card - PROFESSIONAL

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains Advanced Diagnostics, Help Screens, Software Upgrades, File Storage/Retrieval
- Used to back up multiple chairs and programming settings.
- Intended for use by qualified providers only.
- Not compatible with SPJ+ Joysticks.

**NOTE: All specifications and dimensions are approximate.**  
**For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

# Invacare® TDX® SI Rehab Seat Power Wheelchair

## Batteries

22NF batteries with terminal configuration (positive on the left and negative on the right) MUST be used. 22NF batteries that have the reverse terminal configuration MUST not be used. Terminals MUST have a cross hole for proper battery connection.

## Contoura Seat Sizes

	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X

**\*Not available with Contoura Seat Cushion.**

Back Heights | 22", 25" or 28"

## Weight Limit 300 Pounds

Invacare recommends all wheelchairs be fitted and programmed to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

**Note:** All weight limits are inclusive of add-on equipment and seating systems. All specifications and dimensions are approximate.

## Seating and Positioning

For your seating needs, please refer to the following rehab seating and positioning products: Price List Form No. 07-027

Comments

**NOTE: All specifications and dimensions are approximate.**  
*For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).*

## Invacare Corporation

USA

One Invacare Way

Elyria, OH USA

44035-4190

800-333-6900

[www.invacare.com](http://www.invacare.com)

©2010 Invacare Corporation

All rights reserved. Trademarks are identified by the symbols <sup>TM</sup> and <sup>®</sup>. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Products designated by ASL are manufactured by Adaptive Switch Labs, Inc. Therafin is a registered trademark of W/R Medical Electronics Co. Otto Bock is a trademark of Otto Bock Orthopaedische Industrien GMBH and Co. Specifications and prices to change without notification.

