



Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.

Please consider the environment before printing this page with your order form.
For more information on features and benefits of electronic order forms, please visit www.invacare.com/eForms

Invacare® TDX® SI Van Seat Power Wheelchair

Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



TRANSACTION TYPE:

Order Quote

Date _____ Account # _____ P.O. # _____
 Company Name _____ Phone # _____
 Address _____ City _____ State _____ Zip Code _____ Country _____
 Contact Name _____ Comments _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

CHAIR TYPE

- TDXSIV TDX SI with MK6i™ \$5,970.00 K0849
TDXSIV comes standard with: 300 lb. user weight, 5 MPH top speed, 2-pole motor package and MK660 controller.
- TDXSISEAT TDXSI Seat Only \$755.00

| Group 3 - K0849 | |
|--|----------|
| TDX SI with MK6i and Van Seat | |
| 300 lbs. Weight Capacity | |
| 1431-3 14" x 3" Tires with Foam Filled Inserts | Standard |
| 2-Pole Motors | Standard |
| Flip Back Fixed Height Arm | Standard |
| MK660 60 Amp Controller with Actuator | Standard |
| SPJPLUS SPJ™+ Single Drive Joystick | Standard |
| 70NHD 70° Swingaway Footrests | Standard |
| 8 Amp Off Board Battery Charger | Standard |
| 1515 Push Button Style Seat Positioning Strap | Standard |

TRANSPORT TIE DOWN

- TRBKTS Wheelchair Transport Brackets \$250.00 K0108

TRANSPORT TIE-DOWN*

TRBKTS includes four factory installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport. As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

* Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

USER WEIGHT LIMITS

- U300 User Weight <300 lbs. Standard

WHEEL OPTIONS (PAIR)

- 1431-3 14" x 3" with Foam Filled Inserts Standard

FORK OPTION

- FKPKG One-Sided Fork Package \$500.00 E1016

FRAME FINISH

CONTEMPORARY

- 24P Wet Black No Charge
 60P Silver Metallic No Charge
 105P Midnight Blue No Charge
 115P Black Prism No Charge

TOO HOT

- 30P Sunny Yellow No Charge
 61P Electric Red No Charge
 119P Bubblegum Pink No Charge
 120P Grape Madness No Charge

WAY COOL

- 62P Electric Blue No Charge
 127P Grasshopper Green No Charge
 116P Sunkissed Orange No Charge
 117P Tangerine Red No Charge
 29P White Snow No Charge

Discontinued colors 71P Silver Vein, 104P Electric Teal, 121P Lolly Pop Blue, 122P Cosmic Blue, and 125P Emerald Green can be order through specials number C-070109031327

**NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.**

Invacare® TDX® SI Van Seat Power Wheelchair

SEAT

VAN SEAT (STANDARD DIMENSIONS: 18"W x 16"D)

- VW16 16" Wide Van Seat No Charge Clear
- VD16 16" Depth No Charge Clear
- VD17 17" Depth No Charge
- VD18 18" Depth No Charge
- VW18 18" Wide Van Seat Standard Clear
- VD16 16" Depth Standard
- VD17 17" Depth No Charge
- VD18 18" Depth No Charge
- VW20 20" Wide Van Seat \$329.00 Clear
- VD18 18" Seat Depth No Charge
- VD19 19" Seat Depth No Charge
- VD20 20" Seat Depth \$529.00
- VW22 22" Wide Van Seat \$329.00 Clear
- VD18 18" Seat Depth No Charge
- VD19 19" Seat Depth No Charge
- VD20 20" Seat Depth \$529.00

VAN UPHOLSTERY

- GV Grey Vinyl No Charge Clear

VAN ARM TYPE

- VAH Full Length, Adjustable Height \$230.00 E0973
- VFH Full Length, Fixed Height Standard
- VAHD Desk Length, Adjustable Height \$230.00 E0973
- VFHD Desk Length, Fixed Height No Charge

MIXED ARM TYPE

- Right** Clear
- RVFHF Full Length Fixed Height Right No Charge
- RVAHF Full Length Adjustable Height Right \$115.00 E0973
- RVFHD Desk Length Fixed Height Right No Charge
- RVAHD Desk Length Adjustable Height Right \$115.00 E0973
- Left** Clear
- LVFHF Full Length Fixed Height Left No Charge
- LVAHF Full Length Adjustable Height Left \$115.00 E0973
- LVFHD Desk Length Fixed Height Left No Charge
- LVAHD Desk Length Adjustable Height Left \$115.00 E0973

SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Positioning Strap \$50.00
- 1515 Push Button Style Seat Positioning Strap Standard ^{A9900} _{E0978}

ELECTRONICS

MK6I CONTROLLERS

- MK690 90 Amp Controller Standard
- GTRAC G-Trac Module \$995.00 K0108
G-Trac provides Improved Capabilities and enhanced Traction. Does not require Special Motors.

STANDARD DRIVER CONTROL

- SPJPLUS SPJ+ Single Drive Joystick Standard

EXPANDABLE CONTROLLERS SYSTEMS

- PWH Harness req. for Expandable Electronics \$350.00 E2313
Required to interface joystick or MK6i display with expandable controllers.
- MPJM6U MPJ+ Multiple Drive Joystick \$750.00 E2377

JOYSTICK MOUNTING HARDWARE

- QLM6 Quad Link Retract. Joystick Mnt CMPJ+ Van \$230.00 E1028
- QL Quad Link Retract. Joystick Mnt SPJ+ Van \$230.00 E1028
- LEFT Left Handed Mounted Joystick No Charge

ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer with Pro Memory Card \$475.00
- 1813SD Pro Memory Card (USB Ready) \$90.00
- 1560 "T" Handle Flexible Joystick Ext. \$77.00 E2323
- 1561 Straight Handle Flexible Joystick Ext. \$77.00 E2323
- A24VPS 24V MK6i Auxiliary Power Source \$200.00
- MPY Monoport Cable or CMPJ+ or Display \$15.00
- RFM Wireless (RF) Mouse Only Module \$750.00
Allows wireless computer mouse operation through the driver control. Includes two button stereo switch (with hook/loop mount) for mouse clicks. Not compatible with SPJ+ joysticks.
- IR Infrared Module with Wireless RF Mouse Control \$995.00
Capable of operating up to six remote Infra-Red controlled devices. Also allows wireless (RF) Mouse operation through the driver control. Compatible with RFM module for control of two wireless mice from one chair. Includes two button stereo switch (with hook/loop mount) for mouse clicks. Not compatible with SPJ+ Joysticks.

ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Switch - Black \$156.00
- ASL611 Adjustable Switch Mounting Bracket \$252.00

MISCELLANEOUS

BATTERY TRAY TYPE

- 22TRY 22NF Style Tray Standard
22TRY: Base Width 24".

BATTERIES

- 22NFBATTERY 22NF Battery (On Chair) \$300.00 E2361

BATTERY CHARGER OPTION

- H10CHARGER 110 Volt Battery Charger No Charge
- LESTR Lester 8 AMP Battery Charger \$320.00 E2367
- L8900 Omit Battery Charger Less \$80.00

MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks \$45.00
- VAN02 Oxygen Holder (Van Seat Only) \$150.00 E2208

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® TDX® SI Van Seat Power Wheelchair

FRONT RIGGINGS

HANGER BRACKET TYPE Clear

- ASB Center Pivot Style Bracket No Charge
- PSB Pin Style Bracket No Charge

CENTER MOUNT Clear

Center mount adjustable angle (14"-18") available with tall STF only. 13.25" - 17.25" extension range available through special #SM-071217095848

- CTMT-F Center Mount Adj. Knee Angle (12 1/4" - 15") \$300.00
Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 12 1/2" - 15". Footplates are adjustable angle. Must choose either SM, MED or LRG footplate. Clear

- SM Small Footplates (5 1/4" x 7 1/2") No Charge
- MED Medium Footplates (5 1/2" x 9 1/2") No Charge
- LRG Large Footplates (5.90" x 11 1/2") No Charge
- CALF Calf Pads For Center Mount \$200.00

- CTMT-FJ Jr. Center Mount Adj. Knee Angle (9 1/4" - 12") \$300.00
Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 9 1/4" - 12". Footplates are adjustable angle. Available on low seat-to-floor only. Must choose either SM, MED or LRG footplate. Clear

- SM Small Footplates (5 1/4" x 7 1/2") No Charge
- MED Medium Footplates (5 1/2" x 9 1/2") No Charge
- LRG Large Footplates (5.90" x 11 1/2") No Charge
- CALF Calf Pads For Center Mount \$200.00

CENTER PIVOT STYLE (PAIR) Clear

- 70NHD 70° Swingaway Footrests 13 1/2"-17 1/2" Standard
The complete footplate 1651 is standard.
- 70STAPER 70° Tapered Footrests 13"-17" No Charge
Standard with impact guards. Not available with widths 19" to 24" and 16" deep seats.
- AT5544 Manual Elevating Legrests \$255.00 E0990
Not available with 16" wide by 16" to 18" deep or 17" wide by 16" deep and low seat-to-floor.

PIN STYLE-HEMI (PAIR) Clear

- PHW93 70° Swingaway Footrests 13 1/2" - 19" No Charge
Standard with aluminum footplates.
- PHAL4A Manual Articulating/ELR 17 1/4" - 21 1/2" \$225.00 K0053

FOOTBOARD

- 90FB Footboard \$150.00

FOOTPLATE OPTIONS (PAIR) Clear

All front riggings include a standard footplate.

- 1651 Composite Footplates No Charge
- 1350 Extra Large Aluminum Footplates No Charge
Not available with AT5544 and 16" wide by 16" - 18" deep or 17" wide by 16" deep and low seat-to-floor.
- AT5543 Adjustable Angle Flip-Up Footplates \$165.00 K0040
- CTMTFB Footboard for Center Mounts No Charge

| TDX SI Specifications | |
|--|---|
| Seat Width Range (in 1" increments) | Adult (Adjust) |
| | 16" - 22" |
| Seat Depth Range | 16" - 20" |
| Back Angle Range | 85° to 105° in 5° increments |
| Seat-to-Floor (approx) | 21" Standard |
| Overall Height | |
| Semi-Recline without headrest | 38" |
| Semi-Recline with headrest | 45" |
| Overall Length | |
| 35 1/4" without front rigging | |
| Weight | |
| 2-Pole Motor | without 22NF Batteries 179 lbs. with Two 22NF Batteries 277 lbs. |
| Drive Axle | |
| Non-Adjustable | |
| Drive Wheels/Tires (foam filled) | |
| Standard | 14" x 3" |

Batteries

22NF batteries with terminal configuration (positive on the left and negative on the right) MUST be used. 22NF batteries that have the reverse terminal configuration MUST not be used. Terminals MUST have a cross hole for proper battery connection.

Weight Limit 300 Pounds

Invacare recommends all wheelchairs be fitted and programmed to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

Note: All weight limits are inclusive of add-on equipment and seating systems. All specifications and dimensions are approximate.

Seating and Positioning

For your seating needs, please refer to the following rehab seating and positioning products: Price List Form No. 07-027

Comments

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.



Invacare Corporation
USA
One Invacare Way
Elyria, OH USA
44035-4190
800-333-6900

www.invacare.com

©2013 Invacare Corporation
All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.