

# Invacare® TDX® SI Van Seat Power Wheelchair



## Price List and Order Form

Effective May 18, 2009

For ease of ordering, contact Customer Service Toll Free at:

800-333-6900 Fax: 800-678-4682

To order Literature: 800-828-6282

www.invacare.com

Company Name \_\_\_\_\_ Account # \_\_\_\_\_

Phone \_\_\_\_\_ P.O. # \_\_\_\_\_ Date \_\_\_\_\_

Ship to Name & Address \_\_\_\_\_

### Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

### CHAIR TYPE

- TDXSIV TDX SI with MK6i™ (K0849)..... 5,970.00 <sup>1</sup>
- TDXSISEAT TDXSI Seat Only..... 755.00

### Group 3 – K0849

#### Basic Package

- TDX SI with MK6i and Van Seat
- 300 lb Weight Capacity
- 1431-3 14" X 3" Tires w/Foam Filled Inserts .....Standard
- 2 Pole Motors .....Standard
- Flip Back Fixed Height Arm .....Standard
- MK660 60 Amp Controller w/Actuator .....Standard
- SPJPLUS SPJ™+ Single Drive Joystick .....Standard
- 70NHD 70° Swingaway Footrests .....Standard
- 8 Amp Off Board Battery Charger .....Standard
- 1515 Push Button Style Seat Positioning Strap .....Standard

### FORK OPTION

- FKPKG One-Sided Fork Package ..... 500.00

### FRAME FINISH

#### CONTEMPORARY

- 24P Wet Black .....No Charge
- 60P Silver Metallic .....No Charge
- 71P Silver Vein .....No Charge
- 105P Midnight Blue .....No Charge
- 115P Black Prism .....No Charge

#### TOO HOT

- 30P Sunny Yellow .....No Charge
- 61P Electric Red .....No Charge
- 119P Bubblegum Pink .....No Charge
- 120P Grape Madness .....No Charge

#### WAY COOL

- 62P Electric Blue .....No Charge
- 104P Electric Teal .....No Charge
- 121P Lolly Pop Blue .....No Charge
- 122P Cosmic Blue .....No Charge
- 125P Emerald Green .....No Charge
- 127P Grasshopper Green .....No Charge

### TRANSPORT TIE-DOWN

- TRBKTS Wheelchair Transport Brackets..... 250.00

#### \* TRANSPORT TIE-DOWN

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

\*Transport Tie-Down option is not retrofittable to existing models and is not field serviceable.

### USER WEIGHT LIMIT

- U300 User Weight <300 lbs (K0849)..... **Standard**

### WHEEL OPTION (PAIR)

- 1431-3 14" x 3" w/Foam Filled Inserts ..... **Standard**

### FOOTNOTES:

1. TDXSIV comes standard with: 300 lb. user weight 5 MPH top speed, 2-pole motor package and MK660 controller.

NOTE: All specifications and dimensions are approximate.  
For the most current pricing information go to [www.invacare.com](http://www.invacare.com).

# Invacare® TDX® SI Van Seat Power Wheelchair Accessories and Options

## SEAT

### VAN SEAT (Standard Dimensions: 18"W x 16"D)

- VW16 16" Wide Van Seat ..... No Charge
  - VD16 16" Depth ..... No Charge
  - VD17 17" Depth ..... No Charge
  - VD18 18" Depth ..... No Charge
- VW18 18" Wide Van Seat ..... **Standard**
  - VD16 16" Depth ..... **Standard**
  - VD17 17" Depth ..... No Charge
  - VD18 18" Depth ..... No Charge
- VW20 20" Wide Van Seat ..... 329.00
  - VD18 18" Depth ..... No Charge
  - VD19 19" Depth ..... No Charge
  - VD20 20" Depth ..... 529.00
- VW22 22" Wide Van Seat ..... 329.00
  - VD18 18" Depth ..... No Charge
  - VD19 19" Depth ..... No Charge
  - VD20 20" Depth ..... 529.00

### VAN UPHOLSTERY

- GC Grey Cloth ..... No Charge
- GV Grey Vinyl ..... No Charge
- TV Tan Vinyl ..... No Charge

### SEAT STYLE

- VAN Van Style Seat ..... **Standard**
- SEMIR Semi Recline ..... **Standard**

### ARM TYPE <sup>3</sup>

- VAH Full Length, Adjustable Height (E0973) ..... 230.00
- VFH Full Length, Fixed Height ..... **Standard**
- VAHD Desk Length, Adjustable Height (E0973) ..... 230.00
- VFHD Desk Length, Fixed Height ..... No Charge

### MIXED ARMS

- RVAHD Flip up Desk Length Adjustable Height Right ..... 115.00
- RVFHD Flip up Desk Length Fixed Height Right ..... No Charge
- RVAHF Flip Up Full Length Adjustable Height Right ..... 115.00
- RVFHF Flip Up Full Length Fixed Height Right ..... No Charge
- LVFHF Flip up Full Length Fixed Height Left ..... No Charge
- LVAHF Flip Up Full Length Adjustable Height Left ..... 115.00
- LVFHD Flip Up Desk Length Fixed Height Left ..... No Charge
- LVAHD Flip Up Desk Length Adjustable Height Left ..... 115.00

### SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Pos. Strap ..... 50.00
- 1515 Push Button Style Seat Pos. Strap ..... **Standard**

## ELECTRONICS

### MK6i CONTROLLERS

- MK660ACC Controller w/ Actuator Control ..... **Standard**
- GTRAC G-Trac Module (K0108) ..... 995.00<sup>4</sup>

### STANDARD DRIVER CONTROL

- SPJPLUS SPJ+ Single Drive Joystick ..... **Standard**

### FOOTNOTES

<sup>3</sup> Arms are flip up and adjustable width  
<sup>4</sup> G-Trac provides Improved Tracking Capabilities and Enhanced Traction. Does not require Special Motors. Not available w/ SPJ+ joysticks.

### EXPANDABLE CONTROLLER SYSTEMS

- PWH Harness req. for Expandable Electronics (E2313) 350.00 <sup>5</sup>
- MPJM6U MPJ+ Multiple Drive Joystick (E2377) ..... 750.00 <sup>6</sup>
- PSFM6 Personal J/S, Inductive on Front ..... 1,130.00 <sup>6</sup>
- PSRM6 Personal J/S, Inductive on Rear ..... 1,130.00 <sup>6</sup>
  - PBOD Push Button On/Off w/PB Drive Select .... No Charge <sup>7</sup>
  - PBSS Push Button On/Off w/Speed Select ..... No Charge <sup>7</sup>
  - OR On/Off Mounted Right ..... No Charge <sup>7</sup>
  - OL On/Off Mounted Left ..... No Charge <sup>7</sup>

### JOYSTICK MOUNTING HARDWARE

- QLM6 Quad Link Retract. Joystick Mnt MPJ+ Van (E1028) ..... 230.00
- QL Quad Link Retract. Joystick Mnt SPJ+ Van (E1028) ..... 230.00
- LEFT Left Handed Mounted Joystick ..... No Charge

### ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer w/Pro Memory Card ..... 475.00
- 1813SD Pro Memory Card (USB Ready) ..... 90.00
- 1560 "T" Handle Flexible Joystick Ext. (E2323) ..... 77.00
- 1561 Straight Handle Flexible Joystick Ext. (E2323) ..... 77.00
- A24VPS 24V MK6i Auxiliary Power Source ..... 200.00
- MPY Monoport Cable or MPJ+ or Display ..... 15.00
- IR Infrared Module with Wireless Mouse Control ... 995.00 <sup>8</sup>

### ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Switch – Black ..... 156.00
- ASL611 Adjustable Switch Mounting Bracket ..... 252.00

## MISCELLANEOUS

### MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks ..... 45.00
- VAN02 Oxygen Holder (Van Seat Only) (E2208) ..... 150.00

### BATTERY TRAY TYPE <sup>9</sup>

- 22TRY 22NF Style Tray ..... **Standard**

### INSTALLATION OF BATTERIES

- 22NFBATTERY 22NF Battery (On Chair) (E2361) ..... 300.00

### BATTERY CHARGER OPTION

- 110CHARGER 110 Volt Battery Charger ..... No Charge
- LESTR Lester 8 Amp Battery Charger (E2367) ..... 320.00
- L8900 Omit Battery Charger ..... Less 80.00

### FOOTNOTES

<sup>5</sup> Harness is required to interface joystick or MK6i display with expandable controllers. Option is located under Electronic Accessories in configurator  
<sup>6</sup> TOD-SS: Toggle On/Off-Drive Select, Speed Potentiometer Standard.  
<sup>7</sup> PSRM6 Only.  
<sup>8</sup> Capable of controlling up to six remote Infra-Red controlled devices and also Computer Mouse operation through the driver control. Includes stereo switch (with hook/loop mount) for mouse clicks  
<sup>9</sup> 22TRY: Base Width 24":

# Invacare® TDX® SI Van Seat Power Wheelchair Accessories and Options

## FRONT RIGGINGS

### HANGER BRACKET TYPE

- ASB Center Pivot Style Bracket ..... No Charge
- PSB Pin Style Bracket ..... No Charge

### FRONT RIGGING <sup>10</sup>

#### CENTER MOUNT

- CTMT-F Center Mount Adj. Knee Angle ..... 300.00 <sup>11</sup>
- CTMT-FJ Junior Center Mount Adj. Knee Angle ..... 300.00 <sup>12</sup>
- SM Small Footplates (5.25" x 7.5").....No Charge
- MED Medium Footplates (5.50" x 9.5").....No Charge
- LRG Large Footplates (5.90" x 11.5").....No Charge
- CALF Calf Pads For Center Mount..... 200.00

#### CENTER PIVOT STYLE (PAIR)

- 70NHD 70° Swingaway Footrests 13.5"-17.5" ..... **Standard** <sup>13</sup>
- 70STAPER 70° Tapered Footrests 13"-17" ..... No Charge <sup>14</sup>
- AT5544 Manual ELR (E0990)..... 255.00 <sup>15</sup>

#### PIN STYLE HEMI (PAIR)

- PHW93 70° Swingaway Footrests 13.5" - 19"..... No Charge
- PHAL4A Manual Articulating/ELR (K0053)..... 225.00

#### FOOTBOARD

- 90FB Footboard..... 150.00

#### FOOTPLATE OPTIONS (PAIR)

- 1651 Composite Footplates..... No Charge
- 1350 Extra Large Aluminum Footplates..... No Charge <sup>16</sup>
- AT5543 Adjustable Angle Flip-Up Footplates (K0040).... 165.00
- CTMTFB Footboard for Center Mounts..... 150.00

#### FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap (K0038) ..... 27.00
- 1600BK Heel Loops w/Ankle Straps (Pair) (E0951)..... 56.00
- IMPACTG Impact Guards ..... 51.00
- LPT2 Longer Pivot Slide Tube For 1350..... 83.00
- ALPT Longer Pivot And Slide Tube..... 84.00
- ALPT4 Longer Pivot And Slide Tube 4"..... 84.00

### FOOTNOTES

- <sup>10</sup> All front rigging options include a standard footplate.
- <sup>11</sup> Allows 60 degree to 90 degree static knee angle positioning with a drop back lever for easier transferring. Extension range 12.5" to 15". Footplates are adjustable angle. Must choose either SM, MED or LRG footplate.
- <sup>12</sup> Allows 60 degree to 90 degree static knee angle positioning with a drop back lever for easier transferring. Extension range 9.25" to 12". Footplates are adjustable angle. Must choose either SM, MED or LRG footplate.
- <sup>13</sup> With seat widths of 21" to 22" and depth of 16", must pick 1350 footplates.
- <sup>14</sup> Standard with impact guards.
- <sup>15</sup> Not available with 16" wide and 16" to 18" deep or 17" wide by 16" deep and low seat-to-floor.
- <sup>16</sup> Not available with AT5544 and 16" wide by 16" - 18" deep or 17" wide by 16" deep and low seat-to-floor.

TDX SI Specifications		
<b>Seat Width Range:</b> (in 1" increments)	<b>ADULT (Adjust)</b> 16" – 22"	
<b>Seat Depth Range:</b>	16" – 20"	
<b>Back Angle Range</b>	85° - 105° in 5° increments	
<b>Seat-to-Floor (approximate)</b>	Standard: 21"	
<b>Overall Height:</b>	Semi-Recline: 38" (w/o headrest) Semi-Recline: 45" (w/ headrest)	
<b>Overall Length:</b>	35.25" w/o Front Rigging	
<b>Weight</b>	2 Pole Motor	
	w/o 22NF Batteries	179 lbs
	w/ Two 22NF Batteries	277 lbs
<b>Drive Axle:</b>	Non-Adjustable	
<b>Drive Wheels/Tires:</b>	(Foam Filled Standard)	14" x 3"
BATTERIES (NOT INCLUDED)		
22NF batteries with terminal configuration (positive on the left and negative on the right) MUST be used. 22NF batteries that have the reverse terminal configuration MUST not be used. Terminals MUST have a cross hole for proper battery connection.		

Weight Limit 300
Invacare recommends all wheelchairs be fitted and programmed to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitation of the product.

**NOTE:**

**All weight limits are inclusive of add-on equipment and seating systems.  
All specifications and dimensions are approximate.**

Invacare Corporation      www.invacare.com

USA  
One Invacare Way  
Elyria, OH USA  
44035-4190  
800-333-6900

All rights reserved. Trademarks are identified by the symbols <sup>TM</sup> and <sup>®</sup>. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.  
Products designated by ASL are manufactured by Adaptive Switch Labs, Inc.  
Specifications and prices are subject to change without notification.  
©2008 Invacare Corporation

