



Yes, you can.®

Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.

Invacare® Formula™ CG Tilt Only for TDX™ SI Base

Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



Yes, you can.®

TRANSACTION TYPE:

Order Quote

Date _____ Account # _____ P.O. # _____

Company Name _____ Phone # _____

Address _____ City _____ State _____ Zip Code _____ Country _____

Contact Name _____ Comments _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

POWERED SEATING BASE

- TDXSI-CG Formula CG with MK6i™ \$5,970.00 K0856
TDXSI-CG comes standard with: 300 lbs. user weight limit capacity, 5 MPH top speed, 2-pole motor package and MK660ACC.

Group 3 - K0856 Basic Package	
TDXSI-CG TDXSI for Formula CG with MK6i	
300 lbs. Weight Capacity	Standard
2-Pole Motors	Standard
1431-3 14" x 3" Tires with Foam Filled Inserts	Standard
Flip Back Adjustable Height Arm	Standard
MK660ACC 60 Amp Controller	Standard
SPJAP SPJ™+ For Powered Seating	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard
Seat Option CGT Tilt (E1002)	\$4,590.00

OMIT SEAT FRAME

- OSF Omit Seat Frame for Base Only Less \$400.00
Base only.

USER WEIGHT LIMIT

- U300 User Weight <300 lbs Standard
TDXSI-CG comes standard with: 300 lb user weight capacity, 5 MPH top speed, 2-pole motor package and MK660ACC.

TRANSPORT TIE DOWN

- TRBKTS Wheelchair Transport Brackets \$250.00

WHEEL OPTIONS (PAIR)

- 1431-3 14" x 3" with Foam Filled Inserts Standard

FORK OPTION

- FKPKG One-Sided Fork Package \$500.00

TRANSPORT TIE-DOWN*

TRBKTS includes four factory installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport. As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

* Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

FRAME FINISH

CONTEMPORARY

- 24P Wet Black No Charge
- 60P Silver Metallic No Charge
- 105P Midnight Blue No Charge
- 115P Black Prism No Charge

TOO HOT

- 30P Sunny Yellow No Charge
- 61P Electric Red No Charge
- 119P Bubblegum Pink No Charge
- 120P Grape Madness No Charge

WAY COOL

- 62P Electric Blue No Charge
- 127P Grasshopper Green No Charge
- 116P Sunkissed Orange No Charge
- 117P Tangerine Red No Charge
- 29P White Snow No Charge

Discontinued colors 71P Silver Vein, 104P Electric Teal, 121P Lolly Pop Blue, 122P Cosmic Blue, and 125P Emerald Green can be order through specials number C-070109031327

SEAT

SYSTEM TYPE

- CGT CG Tilt Only System \$4,590.00 E1002

SEAT-TO-FLOOR HEIGHT

Due to variation in options the listed seat-to-floor heights have a tolerance of +/- 0.5"

- LOW 16.50" at 0° No Charge
- MED 18.50" at 0° No Charge
- TALL 20.50" at 0° No Charge

SEAT TILT

- 0SA 0° No Charge
- 5SA 5° No Charge

SEAT SELECTION

- CTB Contoura™ Seating..... \$625.00 E2620 E2621
- OMCU Omit Contoura Seat Cushion Less \$20.00
- CTBMSH Contoura seating with mesh cover \$635.00 E2620 E2621
- OMCU Omit Contoura Seat Cushion Less \$20.00
- BPO Back Posts Only/Round Canes Only No Charge
For 12" - 15" seat widths and depths; must select Conventional Back or BPO. Conventional back is not available in seat widths of 16"-20" and 20"-24" or seat depths 16"-19" or 19"-22".
- CURV Curved No Charge
- CNB Conventional Back No Charge
For 12" - 15" seat widths and depths; must select Conventional Back or BPO. Conventional back is not available in seat widths of 16"-20" and 20"-24" or seat depths 16"-19" or 19"-22".

ADJUSTABLE SEAT WIDTH RANGE

Clear

- 1216W 12" to 16" Seat Width Range No Charge
Seat widths 12"- 15" and seat depths 12" - 15" are only available with Back Post Only (BPO) or Conventional Seating. For users under 165 lbs.
- 1620W 16" to 20" Seat Width Range No Charge
- 2022W 20" to 22" Seat Width Range No Charge

SEAT

SEAT WIDTH SETTING

Clear

- W12 Seat Width 12" (Not available with Contoura) No Charge
Seat Widths 12" - 15" and Seat Depths 12" - 15" are only available with Back Post Only (BPO) or Conventional Seating. For users under 165 lbs.
- W13 Seat Width 13" (Not available with Contoura) No Charge
Seat Widths 12" - 15" and Seat Depths 12" - 15" are only available with Back Post Only (BPO) or Conventional Seating. For users under 165 lbs.
- W14 Seat Width 14" (Not available with Contoura) No Charge
Seat Widths 12" - 15" and Seat Depths 12" - 15" are only available with Back Post Only (BPO) or Conventional Seating. For users under 165 lbs.
- W15 Seat Width 15" (Not available with Contoura) No Charge
Seat Widths 12" - 15" and Seat Depths 12" - 15" are only available with Back Post Only (BPO) or Conventional Seating. For users under 165 lbs.
- W16 Seat Width 16" No Charge
- W17 Seat Width 17" (Not available with Contoura) No Charge
- W18 Seat Width 18" No Charge
- W19 Seat Width 19" (Not available with Contoura) No Charge
- W20 Seat Width 20" No Charge
- W21 Seat Width 21" (Not available with Contoura) No Charge
- W22 Seat Width 22" No Charge

ADJUSTABLE SEAT DEPTH RANGE

Clear

- 1218D 12" to 18" Seat Depth Range No Charge
Seat Widths 12" - 15" and Seat Depths 12" - 15" are only available with Back Post Only (BPO) or Conventional Seating. For users under 165 lbs.
- 1619D 16" to 19" Seat Depth Range No Charge
- 1922D 19" to 22" Seat Depth Range No Charge

CONTOURA

Clear

- 1619D 16" to 19" Seat Depth Range (Adjustable)..... No Charge
- 1922D 19" to 22" Seat Depth Range (Adjustable)..... No Charge

SEAT DEPTH SETTING

Clear

Contoura Seating available in depths 16" - 22".

- D12 Seat Depth 12" No Charge
Seat Widths 12" - 15" and Seat Depths 12" - 15" are only available with Back Post Only (BPO) or Conventional Seating. For users under 165 lbs.
- D13 Seat Depth 13" No Charge
Seat Widths 12" - 15" and Seat Depths 12" - 15" are only available with Back Post Only (BPO) or Conventional Seating. For users under 165 lbs.
- D14 Seat Depth 14" No Charge
Seat Widths 12" - 15" and Seat Depths 12" - 15" are only available with Back Post Only (BPO) or Conventional Seating. For users under 165 lbs.
- D15 Seat Depth 15" No Charge
Seat Widths 12" - 15" and Seat Depths 12" - 15" are only available with Back Post Only (BPO) or Conventional Seating. For users under 165 lbs.
- D16 Seat Depth 16" No Charge
- D17 Seat Depth 17" No Charge
- D18 Seat Depth 18" No Charge
- D19 Seat Depth 19" No Charge
- D20 Seat Depth 20" No Charge
- D21 Seat Depth 21" No Charge
- D22 Seat Depth 22" No Charge

ADJUSTABLE ARM TYPE (PAIR)

Clear

- 19F Full Length Adj. Hgt No Charge
- 29F Desk Length Adj. Hgt No Charge

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.**

SEAT Continued

FLIP UP MIXED ADJUSTABLE ARM TYPE (EACH)

All arms are detachable.

- Right Arm** Clear
 19FR Full Length Adj. Hgt. Right No Charge
 29FR Desk Length Adj. Hgt. Right No Charge
- Left Arm** Clear
 19FL Full Length Adj. Hgt. Left No Charge
 29FL Desk Length Adj. Hgt. Left No Charge

CANTILEVER ARM TYPE (PAIR) Clear

Cantilever Arms only available with 300 lb. rated systems. Not available with BPO, Otto Bock, or 12" - 15" seat width. Stealth Lateral part number's TWB-CPR and TWB-CPRL are the only lateral support compatible with the current Cantilever Arm design. Must choose WF armpad type with Stealth Laterals. See chart on last page for arm height adjustment range.

- CLVF Full Length Arm \$350.00
 CLVD Desk Length Arm \$350.00

CANTILEVER MIXED ARM TYPE (EACH)

- Right Arm** Clear
 CLVFR Full Length Arm Right \$175.00
 CLVDR Desk Length Arm Right \$175.00
- Left Arm** Clear
 CLVFL Full Length Arm Left \$175.00
 CLVDL Desk Length Arm Left \$175.00

ARMPAD OPTIONS (CANTILEVER ARMS) Clear

- WF Waterfall Armpads Standard
 RCL Reclining Armpads No Charge

SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Positioning Strap \$50.00
 Not available with 12" - 15" seat widths.
 1515 Push Button Style Seat Positioning Strap Standard

CHEST POSITIONING STRAPS Clear

- 7311BK Airline Buckle Seat Positioning Strap \$100.00 E0960
 Not available with 12" - 15" seat widths.
 9515 Push Button Style Chest Positioning Strap \$100.00 E0960
 7321BK Hook and Loop Chest Positioning Strap \$100.00 E0960

BACK HEIGHT Clear

- TSBH16 Back Height 16" No Charge
 TSBH17 Back Height 17" No Charge
 TSBH18 Back Height 18" No Charge
 TSBH19 Back Height 19" No Charge
 TSBH20 Back Height 20" No Charge
 TSBH21 Back Height 21" No Charge
 TSBH22 Back Height 22" No Charge
 TSBH23 Back Height 23" No Charge
 TSBH24 Back Height 24" No Charge
 TSBH25 Back Height 25" No Charge
 TSBH26 Back Height 26" No Charge

CONTOURA BACK HEIGHT Clear

Not available with seat widths 12" - 15".

- TSBH22 Back Height 22" No Charge
 TSBH25 Back Height 25" No Charge
 TSBH28 Back Height 28" No Charge

Back Height Chart	
Back Type	Back Height
CTB - Contoura	22", 25" and 28"
BPO - Back Posts Only	16"-24"
CURV - Curved Back	16"-26"
Conventional Back	16"-20" on seat widths of 12"-15"
Note: Back Heights of 21"-26" not available with seat widths of 12"-15"	

HEADREST TYPE (OPTIONAL) Clear

Hardware is packaged with equipment and cannot be ordered separately.

- TAR127700 Curved Headrest \$210.00 E0955
 TAR27700HW Removable & Adjustable Hardware \$210.00 E1028
- TAR128000 Two-Step Headrest \$210.00 E0955
 TAR128000HW Removable & Adjustable Hardware \$210.00 E1028
- T0022400 Large Headrest \$210.00 E0955
 T0022400HW Removable & Adjustable Hardware \$210.00 E1028
- T0022600 Small Headrest \$210.00 E0955
 T0022600HW Removable & Adjustable Hardware \$210.00 E1028

STEALTH HEADRESTS Clear

Hardware is packaged with equipment and cannot be ordered separately.

- CP175 14" Multi Axis Headrest \$225.00 E0955
 CP175HW Removable & Adjustable Hardware \$225.00 E1028
- CP275 10" Multi Axis Headrest \$225.00 E0955
 CP275HW Removable & Adjustable Hardware \$225.00 E1028
- CP180 14" TWB Headrest \$225.00 E0955
 CP180HW Removable & Adjustable Hardware \$225.00 E1028
- CP280 10" TWB Headrest \$225.00 E0955
 CP280HW Removable & Adjustable Hardware \$225.00 E1028

OTTO BOCK ARM PADS (PAIR) Clear

Must select Otto Bock Hand Pad when selecting Otto Bock Forearm pads.

- OB20400 One Piece Channel Arm pads \$230.00 E2209
 OB20500 Large Forearm Pads \$235.00
 OB20600 Medium Forearm Pads \$235.00
 OB20700 Small Forearm Pads \$235.00

OTTO BOCK HAND PADS (PAIR) Clear

- OB23000 Large Flat Hand Pads \$100.00
 OB23300 Medium Flat Hand Pads \$100.00

LATERALS (PAIR)

CONVENTIONAL BACK/BPO/CURVED BACK Clear

Hardware is packaged with equipment and cannot be ordered separately.

- TCSLM Clamp-On Swingaway Laterals Medium \$210.00 E0956
 TCSLMHW Removable and Adjustable Hardware \$210.00 E1028
- TCSSL Clamp-On Swingaway Laterals Large \$210.00 E0956
 TCSSLHW Removable and Adjustable Hardware \$210.00 E1028

**NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.**

Invacare® Formula™ CG Tilt Only for TDX™ SI Base

STEALTH LATERALS

Clear

Hardware is packaged with equipment and cannot be ordered separately.

- TWB-RPR Adj. Swing Round Lat. Med \$320.00 E0956
RPRHW Removable and Adjustable Hardware \$320.00 E1028
- TWB-RPRL Adj. Swing Round Lat. Lg \$320.00 E0956
RPRLHW Removable and Adjustable Hardware \$320.00 E1028

LATERALS FOR CONTOURA SEATING

Clear

Hardware is packaged with equipment and cannot be ordered separately.

- CSLM Swingaway Contoura Lat. Med \$210.00 E0956
CSLMHW Removable and Adjustable Hardware \$210.00 E1028
- CSSL Swingaway Contoura Lat. Lg \$210.00 E0956
CSSLHW Removable and Adjustable Hardware \$210.00 E1028

STEALTH LATERALS FOR COUNTOURA

Clear

Hardware is packaged with equipment and cannot be ordered separately.

- TWB-CPR Adj. Swing Contoura Lat. Med \$320.00 E0956
CPRHW Removable & Adjustable Hardware \$320.00 E1028
- TWB-CPRL Adj. Swing Contoura Lat. Lg \$320.00 E0956
CPRLHW Removable & Adjustable Hardware \$320.00 E1028

ELECTRONICS

MK6i CONTROLLER

- MK660ACC Controller with Actuator Control Standard
- GTRAC G-Trac Module \$995.00 K0108
G-Trac provides improved tracking capabilities and enhanced traction. Does not require special motors. Not available with SPJ+ Joysticks.

STANDARD DRIVER CONTROL

- SPJAP SPJ+ for Power Seating Standard
Includes Egg Switch with 611 mount to operate single actuator

ACTUATOR CONTROL THROUGH SPJ+

- SPJAPM SPJ+ with Joystick Actuator Control \$1,300.00 E2310
Optional for Torque SP and Torque 3. Allows actuator operation through Driver Control. Includes port for attendant switch.

EXPANDABLE CONTROLLER SYSTEM

Clear

- PWH Harness req. for Expandable Electronics \$350.00 E2313
Required to interface joystick or MK6i display with expandable controllers.
- CMPJM6U MPJ+ Multiple Drive Joystick with Color Display .. \$750.00 E2377
- ASLRDYM6 MK6 Display Only (No Driver Control) \$750.00 E2377
Display only, No Driver Control. Allows choice from order form.
- MCD Mini Color Display \$750.00 E2377
No port for remote On/Off switch. Stereo port only.

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS

Clear

- 1500M6 Prop RIM Head Control with MK6i Display .. \$2,875.00 E2327
 Clear
- SMHD Small Headrest No Charge
- LGHD Large Headrest No Charge
- 1558MM6 Compact Joystick with Display \$1,130.00 E2373

ALTERNATIVE DIGITAL CONTROL SYSTEMS

Clear

Must also use either a Multiple Drive Joystick or MK6i Display.

- SNPM6 Sip-N-Puff Digital Interface Box \$1,500.00 E2325
- 5020M6 Mini Tash Joystick and Digital Interface Box \$1,750.00 E2321

SIP-N-PUFF KIT

- PKG32666 Therafin® Sip-N-Puff Breath Tube Kit \$380.00 E2326

ADDITIONAL CONTROL CHOICES

- 1558M6 Compact Joystick - Less Std Mount \$770.00
- 5020 Mini Tash Joystick \$405.00
Requires SNPM6 Sip-N-Puff Interface Box.
- PACM6 Prop Attendant Control with quick release mount \$850.00
- 1552M Digital Attendant Control \$850.00
Requires SNPM6 Sip-N-Puff Interface Box.

JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retract. Joystick Mnt CMPJ+ \$279.00 E1028
Include ARM250 standard.
- QLA Quad Link Retract. Joystick Mnt SPJ+ \$279.00 E1028
Include ARM250 standard.
- ARM250 Stealth Height Adj Joystick \$210.00 E1028
- LEFT Left Handed Mounted Joystick No Charge
- RIMHW RIM Hardware For 1500 \$219.00 E1028
- PKG32669 Complete Bib Assembly (Tash/1558) \$265.00
- GATMPJM6 Gatlin Midline for CMPJ+ \$600.00 E1028
- GAT1812 Gatlin Midline for 1812 \$600.00 E1028
- RCM Rear Cane Mount for Compact Joystick \$395.00

ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer with Pro Memory Card \$475.00
- 1813SD Pro Memory Card (USBReady) \$90.00
- 1560 "T" Handle Flexible Joystick Ext. \$77.00 E2323
- 1561 Straight Handle Flexible Joystick Ext. \$77.00 E2323
- 1826 Chin Cup \$48.00 E2324
- A24VPS 24V MK6 Auxiliary Power Source \$200.00
- AUX12M6 Auxiliary Module for outputs 1 and 2 \$795.00
- AUX34M6 Auxiliary Module for outputs 3 and 4 \$795.00
- 2PB Dual Push Button \$147.00
Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options.
- 2SW Dual Toggle \$172.00
Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options.
- PBMT Push Button Mounting Kit \$252.00
- MPY Monoport "Y" Cable for CMPJ+ or Display \$15.00
- RFM Wireless (RF) Mouse Only Module \$750.00
RFM includes Transmitter, USB Radio Frequency Mouse Dongle (Receiver), and Stereo Switch with hook and loop mount for mouse clicks. Not compatible with SPJ+.
- IR Infrared Module with Wireless RF Mouse Control \$995.00
Capable of operating up to six remote Infra-Red controlled devices. Also allows wireless (RF) Mouse operation through the driver control. Compatible with RFM module for control of two wireless mice from one chair. Includes two button stereo switch (with hook/loop mount) for mouse clicks. Not compatible with SPJ+ Joysticks.

FOR YOUR ASL DRIVE SYSTEM NEEDS
PLEASE REFER TO THE
ASL SINGLE INVOICE PRICE LIST FORM # 02-097.

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

POWERED SEATING ACCESSORIES

- SANODE One Act. Oper. through Driver Ctrl \$1,300.00 E2310
Required for Tilt Operation through Driver Control - Single Actuator Systems only. Not available with SPJ+ Joysticks. The Egg Switch with 611 Mount does not come standard when ordering a SANODE.

ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Switch-Black \$156.00
- ASL304A Light Wobble Switch \$180.00
- ASL202J Fiber Optic Switch \$732.00
Not available with SPJ+ Joystick. Must add MK6 Auxiliary Power Source.
- ASL208J Adjustable Proximity Switch \$660.00
Not available with SPJ+ Joystick. Must add MK6 Auxiliary Power Source.
- ASL504M6 Remote Emergency Stop Switch \$1,185.00
- ASL611 Adjustable Switch Mounting Bracket \$252.00

MISCELLANEOUS

BATTERY TRAY TYPE

- 22TRY 22NF Style Tray Standard
22TRY: Base Width 24".

BATTERIES

- 22NFBATTERY 22NF Battery (On Chair) \$300.00 E2361

BATTERY CHARGER OPTION

- Clear
- 110CHARGER 110 Volt Battery Charger No Charge
- LESTR Lester 8 AMP Battery Charger \$320.00 E2367
- L8900 Omit Battery Charger Less \$80.00

VENT TRAY

- COMP Compact Vent Tray \$1,540.00 E1030
Approximate dimensions are 11" W x 4" D x 13" H. Designed for use with the latest and most compact style vent systems.

MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks \$45.00

ADDUCTORS

- Clear
- Hardware is packaged with equipment and cannot be ordered separately.**
- ADD Adductors Small \$200.00 E0956
- ADDHW Hardware for ADD \$400.00 E1028
- ADDL Adductors Large \$200.00 E0956
- ADDHWL Hardware for ADDL \$400.00 E1028

FRONT RIGGINGS

HANGER BRACKET TYPE

- Clear
- ASB Center Pivot Style Bracket No Charge
- PSB Pin Style Bracket No Charge
- GHB Genius Hanger Brackets No Charge
- NHB No Hanger Brackets Less \$30.00

CENTER MOUNT

- Clear
- For a 13.25" - 17.25" extension range, available through special #SM-071217095848 with Tall STF Only and 5 Degree Seat-Tilt Angle.**
- CTMT-F Center Mount Adj. Knee Angle (12.25" - 15") \$300.00
Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 12.5" - 15". Footplates are adjustable angle. Must choose either SM, MED or LRG footplate.
- Clear
- SML Small Footplates (5.25" x 7.5") No Charge
- MDM Medium Footplates (5.5" x 9.5") No Charge
- LRG Large Footplates (5.90" x 11.5") No Charge
- CALF Calf Pads For Center Mount \$200.00

- CTMT-FJ Jr. Center Mount Adj. Knee Angle (9.25" - 12") \$300.00
Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 9.25" - 12". Footplates are adjustable angle. Available on low seat-to-floor only. Must choose either SM, MED or LRG footplate.
- Clear

- SML Small Footplates (5.25" x 7.5") No Charge
- MDM Medium Footplates (5.5" x 9.5") No Charge
- LRG Large Footplates (5.90" x 11.5") No Charge
- CALF Calf Pads For Center Mount \$200.00

CENTER PIVOT STYLE (PAIR)

- Clear
- All front riggings include a standard footplate.**
- 70NHD 70° Swingaway Footrests 13 1/2" - 17 1/2" Standard
- 70STAPER 70° Tapered Footrests 13" - 17" No Charge
Standard with impact guards.
- AT5544 Manual Elevating Legrests \$255.00 E0990

PIN STYLE-HEMI (PAIR)

- Clear
- All front riggings include a standard footplate.**
- PHW93 70° Swingaway Footrests 13 1/2" - 19" No Charge
Standard with composite footplates.
- PHAL4A Manual Art./ELR 17 1/4" - 21 1/2" \$225.00 K0053
Not available on 12" or 13" seat widths.

FOOTPLATE OPTIONS (PAIR)

- Clear
- All front riggings include a standard footplate.**
- 1651 Flip-Up Composite Footplates No Charge
- 1350 Extra Large Aluminum Footplates No Charge
- AT5543 Adjustable Angle Flip-Up Footplates \$165.00 K0040
- CTMTFB Footboard for Center Mounts No Charge
One piece footboard for Center Mount (CTMT-F) Riggings.

FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap \$27.00 E0038
- IMPACTG Impact Guards \$51.00
- HLOOP Heel Loops (pair) \$45.00 E0951
- LPT2 Longer Pivot Slide Tube For 1350 \$83.00
- ALPT Longer Pivot And Slide Tube \$84.00
- ALPT4 Longer Pivot And Slide Tube 4" \$84.00

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Formula™ CG Tilt Only for TDX™ SI Base

Batteries

22NF batteries with terminal configuration (positive on the left and negative on the right) MUST be used. 22NF batteries that have the reverse terminal configuration MUST not be used. Terminals MUST have a cross hole for proper battery connection.

Formula CG Tilt Only Features and Specifications

- Weight capacity of 300 lbs.
- Curved or Contoura style seating.
- Adjustable seat frame.

Contoura Seat Sizes

	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X
W24	*	*	X	X	X	X	X

*Not available with Contoura Seat Cushion

Back Heights 22", 25" or 28"

Conventional or BPO Seat Sizes

Seat Widths	12" - 15"
Seat Depths	12" - 22" (BPO) 12" - 18" (Conventional)
Back Heights	16" to 20" for 12" - 15"

Curved or BPO Seat Sizes

Seat Widths	16" - 24"
Seat Depths	16" - 22"
Back Heights	16" - 26"

Seating and Positioning

For your seating needs, please refer to the following rehab seating and positioning products: Price List Form No. 07-027

Tech Tips

Operating Powered Seating through the Driver's Control. All chairs require adding or replacing one component to allow actuator operation through the driver control.

SPJ+ with PSS:

Replace joystick with SPJ+ with ACC (Actuator Control) -E2310.

Single Actuator Systems with CMPJ+, PSR+, PSF+ MK6i Display:

Add a SANODE (Single Actuator Node) to the system - E2310.

Memory Card - BASIC

- Standard with ALL Expandable Control Systems (CMPJ+, MK6i Display).
- Used to back up / restore programmed settings for only one chair.
- May be given to user or kept with provider's files for safe keeping.
- Does not contain advanced diagnostics / help screens / file structure.
- Not compatible with SPJ+ Joysticks.

Memory Card - PROFESSIONAL

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains advanced diagnostics, help screens, software upgrades, file storage / retrieval
- Used to back up multiple chairs and programming settings.
- Intended for use by qualified providers only.
- Not compatible with SPJ+ Joysticks.

Cantilever Arm Height Adjustment

Back Type	Height	Adjustment	Adjustment with laterals
CURV	20" or greater	11" to 16"	11" to 14"
CURV	19"	11" to 15"	11" to 13"
CURV	18"	11" to 14"	11" to 12"
CURV	17"	11" to 13"	no adjustment
CURV	16"	11" to 12"	not available
Contoura	22", 25" to 28"	12 1/2" to 14 1/2"	no adjustment

Comments

NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.

Invacare Corporation

www.invacare.com

USA

©2010 Invacare Corporation

One Invacare Way

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Products designated by ASL are manufactured by Adaptive Switch Labs, Inc. Therafin is a registered trademark of W/R Medical Electronics Co. Freedom Designs is a registered trademark of Freedom Designs, Inc. Otto Bock is a trademark of Otto Bock Orthopaedische Industrien GMBH and Co. Specifications and prices to change without notification.

Elyria, OH USA

44035-4190

800-333-6900

