



Yes, you can.®

### Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email\* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

\*Saving and emailing the eForms is available with Acrobat Professional software

#### If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

**Please note:** The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

#### If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit [www.adobe.com](http://www.adobe.com).

**NOTE: All specifications and dimensions are approximate.**  
*For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).*

# Invacare® Formula™ PTO Plus Retrofit Kit for Pronto M71™ and M91™

## Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



Yes, you can.®

### TRANSACTION TYPE:

- Order  Quote

Date \_\_\_\_\_ Account # \_\_\_\_\_ P.O. # \_\_\_\_\_

Company Name \_\_\_\_\_ Phone # \_\_\_\_\_

Address \_\_\_\_\_ City \_\_\_\_\_ State \_\_\_\_\_ Zip Code \_\_\_\_\_ Country \_\_\_\_\_

Contact Name \_\_\_\_\_ Comments \_\_\_\_\_

### Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.  
**Note:** Existing base information (serial #) will need to be supplied to Customer Service to ensure correct configuration. Chairs with RH electronics and prior to 2004 not compatible with Retrofit Kit.

### SEAT TYPE

- M91SEAT Pronto M91 Seat Only ..... \$740.00  
 M71SEAT Pronto M71 Seat Only ..... \$740.00

### POWERED SEATING

#### POWER SEATING SYSTEM

- PTO-PLK Power Tilt Only Plus Kit ..... \$4,370.00 E1002

### FRAME SELECTION

#### SEAT SELECTION

- ADJASBA Adjustable Width and Back Angle ..... No Charge E2620 E2621  
 CTB Contoura™ Seating ..... \$625.00 E2620 E2621  
 OMCU Omit Contoura Seat Cushion .... Less \$20.00  
 CTBMSH Contoura™ seating with mesh upholstery..... \$635.00 E2620 E2621  
 OMCU Omit Contoura Seat Cushion .... Less \$20.00

#### USER WEIGHT LIMITS

The Power Tilt Only (PTO) is a weight-shifting tilt system with a weight limit of up to 275 lbs. or 375 lbs. (only for the M91). Van or Manual Reclining back seat styles are not available on the PTO.

- U275 User Weight under 275 lbs. .... No Charge  
 U375 User Weight under 375 lbs. .... No Charge

#### SEAT WIDTH RANGE

Clear

Contoura Seating available in even widths 16" - 22" and depths between 16" - 22" in 1" increments.

- 1620W 16" to 20" Seat Width Range (Adjustable) ..... No Charge  
 2022W 20" to 22" Seat Width Range (Adjustable) ..... No Charge

#### SEAT WIDTH SETTING

Clear

- W16 16" Seat Width ..... No Charge  
 W17 17" Seat Width (not available with Contoura) ..... No Charge  
 W18 18" Seat Width ..... No Charge  
 W19 19" Seat Width (not available with Contoura) ..... No Charge  
 W20 20" Seat Width ..... No Charge  
 W21 21" Seat Width (not available with Contoura) ..... No Charge  
 Not available with M71SEAT  
 W22 22" Seat Width ..... No Charge  
 Not available with M71SEAT  
 W23 23" Seat Width ..... No Charge  
 Not available with M71SEAT  
 W24 24" Seat Width ..... No Charge  
 Not available with M71SEAT

#### SEAT DEPTH RANGE

Clear

- 1619D 16" to 19" Seat Depth Range (Adjustable) ..... No Charge  
 1920D 19" to 20" Seat Depth Range (Adjustable) ..... No Charge  
 1922D 19" to 22" Seat Depth Range (Adjustable) ..... No Charge

### SEAT DEPTH SETTING

Clear

- D16 16" Seat Depth ..... No Charge  
 D17 17" Seat Depth (not available with Contoura) ..... No Charge  
 D18 18" Seat Depth ..... No Charge  
 D19 19" Seat Depth (not available with Contoura) ..... No Charge  
 D20 20" Seat Depth ..... No Charge  
 D21 21" Seat Depth (not available with Contoura) ..... No Charge  
 Not available with M71SEAT  
 D22 22" Seat Depth ..... No Charge  
 Not available with M71SEAT

### UPHOLSTERY COLORS

For M91SEAT only.

- U240 Black Nylon Upholstery ..... No Charge

### ELECTRONICS FOR PTO-PLK

Clear

- MK6 to MK6 MK6i on existing chair ..... No Charge  
 MK5 to MK6 MK5 on existing chair ..... No Charge

### EXISTING MK6i™ DRIVER CONTROL

Clear

- CMPJM6 CMPJ+ for Multiple Drive Joystick ..... No Charge  
 SPJAP SPJ+ with Power Seating Switch ..... No Charge  
 SPJAPM SPJ+ with Joystick Actuator Control ..... No Charge  
 Need new Drive Control ..... No Charge

### NEW CONTROLLER

Clear

- MK690ACC Controller with Actuator Control (M91) ..... No Charge  
 MK660ACC Controller with Actuator Control (M71) ..... No Charge

### NEW MK6i™ DRIVER CONTROL

Clear

- CMPJM6 CMPJ+ for Multiple Drive Joystick ..... \$750.00 E2377  
 Egg switch with 611 mount comes with joystick top operate tilt function.  
 SPJAP SPJ+ with Power Seating Switch ..... No Charge  
 Egg switch with 611 mount comes with joystick top operate tilt function.  
 SPJAPM SPJ+ with Joystick Actuator Control ..... \$1,300.00 E2310  
 Allows Actuator operation through the driver control. Egg Switch with 611 mount will not come standard with this joystick.  
 Need new Drive Control ..... No Charge  
 See Expandable Controller systems on Page 2.

### POWERED SEATING ACCESSORIES

- SANODE Single Actuator. Operation through Driver Control ..... \$1,300.00 E2310  
 Allows actuator operation through driver control. Not compatible with SPJAP or SPJAPM joysticks. Single actuator system only.  
 Egg Switch with 611 mount will not come with the SANODE.

**NOTE: All specifications and dimensions are approximate.**  
 For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).

**Invacare® Formula™ PTO Plus Retrofit Kits for M71™ and M91™**  
**TILT SELECTIONS**

**TILT BACK OPTION**  Clear

- CURV Curved Back with Headrest ..... No Charge
- POST Posts Only with Headrest ..... No Charge

**LATERALS**  Clear

Hardware is packaged with equipment and cannot be ordered separately.

- TCSLM Clamp-On Swingaway Laterals Medium ..... \$210.00 E0956  
TCSLMHW Removable and Adjustable Hardware ..... \$210.00 E1028
- TCSLL Clamp-On Swingaway Laterals Large ..... \$210.00 E0956  
TCSLLHW Removable and Adjustable Hardware ..... \$210.00 E1028

**BACKS**

**BACK TYPE**

- SECT Sectional Rear Uprights ..... Standard

**BACK HEIGHT**

- BH16 16" Back Height ..... No Charge
- BH17 17" Back Height ..... No Charge
- BH18 18" Back Height ..... No Charge
- BH19 19" Back Height ..... No Charge
- BH20 20" Back Height ..... No Charge
- BH21 21" Back Height ..... No Charge
- BH22 22" Back Height ..... No Charge
- BH23 23" Back Height ..... No Charge
- BH24 24" Back Height ..... No Charge
- TSBH22 22" Back Height ..... No Charge
- TSBH25 25" Back Height ..... No Charge  
25" and 28" only available with CTB (Contoura Seating).
- TSBH28 28" Back Height ..... No Charge  
25" and 28" only available with CTB (Contoura Seating).

**TILT RANGE**

- TR055 0° - 55° Tilt Range ..... No Charge
- TR560 5° - 60° Tilt Range ..... No Charge

**HEADREST TYPE (OPTIONAL)**  Clear

Hardware is packaged with equipment and cannot be ordered separately.

- TAR127700 Curved Headrest ..... \$210.00 E0955  
TAR27700HW Removable & Adjustable Hardware ..... \$210.00 E1028
- TAR128000 Two-Step Headrest ..... \$210.00 E0955  
TAR128000HW Removable & Adjustable Hardware ..... \$210.00 E1028
- T0022400 Large Headrest ..... \$210.00 E0955  
T0022400HW Removable & Adjustable Hardware ..... \$210.00 E1028
- T0022600 Small Headrest ..... \$210.00 E0955  
T0022600HW Removable & Adjustable Hardware ..... \$210.00 E1028

**ARMS**

**ARM TYPE**  Clear

- 18F Desk Length Fixed Height Flip-Up ..... No Charge
- 28F Full Length Fixed Height Flip-Up ..... No Charge
- 19F Full Length Adjustable Height Flip-Up ..... No Charge
- 29F Desk Length Adjustable Height Flip-Up ..... No Charge

**MIXED ARM TYPE**

**Right**  Clear

- R28F Desk Length Fixed Height Flip-Up..... No Charge
- R18F Full Length Fixed Height Flip-Up..... No Charge
- R29F Desk Length Adjustable Height Flip-Up ..... No Charge
- R19F Full Length Adjustable Height Flip-Up ..... No Charge

**Left**  Clear

- L28F Desk Length Fixed Height Flip-Up..... No Charge
- L18F Full Length Fixed Height Flip-Up..... No Charge
- L29F Desk Length Adjustable Height Flip-Up ..... No Charge
- L19F Full Length Adjustable Height Flip-Up ..... No Charge

**NOTE: All specifications and dimensions are approximate.**  
**For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

**POSITIONING STRAPS**

**SEAT POSITIONING STRAPS**

- 1311BK Airline Buckle Seat Positioning Strap ..... \$50.00
- 1515 Push Button Style Seat Positioning Strap ..... Standard

**CHEST POSITIONING STRAPS**  Clear

- 7311BK Airline Buckle Chest Positioning Strap ..... \$100.00 E0960
- 9515 Push Button Style Chest Positioning Strap ..... \$100.00 E0960
- 7321BK Hook and Loop Chest Positioning Strap ..... \$100.00 E0960

**EXPANDABLE CONTROLLER SYSTEM**  Clear

- PWH Harness req. for Expandable Electronics ..... \$350.00 E2313  
**Required to interface joystick or MK6i display with expandable controllers.**
- CMPJM6U MPJ+ Multiple Drive Joystick with Color Display .. \$750.00 E2377
- PSFM6 Personal Joystick, Inductive on Front ..... \$1,130.00 E2377  
**TOD-SS: Toggle On/Off Drive Select, Speed Potentiometer Standard.**
- PSRM6 Personal Joystick, Inductive on Rear ..... \$1,130.00 E2377  
**TOD-SS: Toggle On/Off Drive Select, Speed Potentiometer Standard.**
- PBOD Push Button On/Off with PB Drive Select .... No Charge  
**PSRM6 Only.**
- PBSS Push Button On/Off with Speed Select ..... No Charge  
**Drive select through mode switch, or mode switch and joystick. PSRM6 Only.**
- OR On/Off Mounted Right ..... No Charge  
**PSRM6 Only.**
- OL On/Off Mounted Left ..... No Charge  
**PSRM6 Only.**

- ASLRDYM6 MK6 Display Only (No Driver Control) ..... \$750.00 E2377  
**Display only, No Driver Control. Allows choice from order form.**
- MCD Mini Color Display ..... \$750.00 E2377  
**No port for remote On/Off switch. Stereo port only.**

FOR YOUR ASL DRIVE SYSTEM NEEDS  
PLEASE REFER TO THE  
ASL SINGLE INVOICE PRICE LIST FORM # 02-097.

**ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS**  Clear

- 1500M6 Prop RIM Head Control with MK6i Display .. \$2,875.00 E2327
- SMHD Small Headrest ..... No Charge
- LGHD Large Headrest ..... No Charge
- 1558MM6 Compact Joystick with Display ..... \$1,130.00 E2373

**ALTERNATIVE DIGITAL CONTROL SYSTEMS**  Clear

**Must also use either a Multiple Drive Joystick or MK6i Display.**

- SNPM6 Sip-N-Puff Digital Interface Box ..... \$1,500.00 E2325
- 5020M6 Mini Tash Joystick and Digital Interface Box \$1,750.00 E2321

**SIP-N-PUFF KIT**

- PKG32666 Therafin® Sip-N-Puff Breath Tube Kit ..... \$380.00 E2326

**ADDITIONAL CONTROL CHOICES**

- 1558M6 Compact Joystick - Less Std Mount ..... \$770.00
- 5020 Mini Tash Joystick ..... \$405.00  
**Requires SNPM6 Sip-N-Puff Interface Box.**
- PACM6 Prop Attendant Control with quick release mount .... \$850.00
- 1552M Digital Attendant Control ..... \$850.00  
**Requires SNPM6 Sip-N-Puff Interface Box.**

# Invacare® Formula™ PTO Plus Retrofit Kits for M71™ and M91™

## JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retractable Joystick Mnt CMPJ+ ..... \$279.00 E1028  
**Includes ARM250 standard.**
- QLA Quad Link Retractable Joystick Mnt SPJ+ ..... \$279.00 E1028  
**Includes ARM250 standard.**
- ARM250 Stealth Height Adjustable Joystick Mount ..... \$210.00 E1028
- LEFT Left Handed Mounted Joystick ..... No Charge
- RIMHW RIM Hardware For 1500 ..... \$219.00 E1028
- PKG32669 Complete Bib Assembly Tash/1558 ..... \$265.00
- GATMPJ6M Gatlin Midline for CMPJ+ ..... \$600.00 E1028
- GAT1812 Gatlin Midline for 1812 ..... \$600.00 E1028
- RCM Rear Cane Mount for Compact Joystick ..... \$395.00

## ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer with Pro Memory Card ..... \$475.00
- 1813SD Pro Memory Card (USB Ready) ..... \$90.00
- 1560 "T" Handle Flexible Joystick Ext. .... \$77.00 E2323
- 1561 Straight Handle Flexible Joystick Ext. .... \$77.00 E2323
- 1826 Chin Cup ..... \$48.00 E2324
- A24VPS 24V MK6 Auxiliary Power Source ..... \$200.00
- AUX12M6 Auxiliary Module for outputs 1 and 2 ..... \$795.00
- AUX34M6 Auxiliary Module for outputs 3 and 4 ..... \$795.00
- RFM Wireless (RF) Mouse Only Module ..... \$750.00  
**RFM includes Transmitter, USB Radio Frequency Mouse Dongle (Receiver), and Stereo Switch with hook and loop mount for mouse clicks. Not compatible with SPJ+.**
- IR Infrared Module with Wireless RF Mouse Control ..... \$995.00  
**Capable of operating up to six remote Infra-Red controlled devices. Also allows wireless (RF) Mouse operation through the driver control. Compatible with RFM module for control of two wireless mice from one chair. Includes two button stereo switch (with hook/loop mount) for mouse clicks. Not compatible with SPJ+ Joysticks.**

## ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Reset - Black ..... \$156.00
- ASL304A Light Wobble Switch ..... \$180.00
- ASL202J Fiber Optic Switch ..... \$732.00  
**Not available with SPJ+ Joystick. Must add MK6 Auxiliary Power Source.**
- ASL208J Adjustable Proximity Switch ..... \$660.00  
**Not available with SPJ+ Joystick. Must add MK6 Auxiliary Power Source.**
- ASL504M6 Remote Emergency Stop Switch ..... \$1185.00  
**Not available with SPJ+ Joystick.**
- ASL611 Adjustable Switch Mounting Bracket ..... \$252.00

## FRONT RIGGINGS

### PIN STYLE (PAIR)

Clear

- PHW93 Swingaway Footrests 13 1/2" - 19" ..... No Charge
- PH904A Swingaway Elevating Legrests 16" - 19 3/4" .... \$255.00
- PHAL4A Smartleg Articulating Legrests 17 1/4" - 21 1/2" ... \$320.00 K0053

### CENTER MOUNT

- CTMT Center Mount Adjustable ..... \$220.00
- CTMTJ Center Mount Adjustable Junior ..... \$220.00

## FRONT RIGGING SEAT DEPTH

Clear

- FR1218D 12" - 18" Seat Depth Range ..... No Charge
- FR1619D 16" - 19" Seat Depth Range ..... No Charge
- FR1922D 19" - 22" Seat Depth Range ..... No Charge

## FOOTPLATE OPTIONS (PAIR)

Clear

- CTMTFB Footboard for Center Mount ..... No Charge
- 1651 Flip-Up Composite Footplates ..... No Charge
- 1350 Extra Large Aluminum Footplates ..... No Charge  
**Not available on a 16" wide frame.**
- 1450 Aluminum Footplates-Deluxe ..... No Charge
- 1650 Aluminum Footplates ..... No Charge
- AT5543 Adjustable Angle Flip Up Footplates ..... \$165.00

## FOOTREST/LEGREST ACCESSORIES

- A1337 Calf Strap ..... \$27.00 E0038
- LPT2 Longer Pivot Slide Tube For 1350 ..... \$83.00
- IMPACTG Impact Guards ..... \$51.00
- ALPT Longer Pivot and Slide Tube ..... \$84.00
- ALPT4 Longer Pivot and Slide Tube 4" ..... \$84.00

### Seating and Positioning

For your seating needs, please refer to the following rehab seating and positioning products: Price List Form No. 07-027

### Formula PTO Plus for Pronto M71 and M91 Specifications

PTO Seat Dimensions	M71	M91
Widths	16" - 20"	16" - 24"
Depths	16" - 20"	16" - 22"
Back Heights (without headrest)	16" - 24"	16" - 24"
<b>Incline Capacity (with PTO)</b>	9 degrees	9 degrees
<b>Weight Capacity (with PTO)</b>	275 lbs.	275 lbs.
<b>Turning Radius (PW93)</b>	27 1/2"	27 1/2"
<b>Seat-to-Floor at 0°</b>	18"	18"

### Tech Tips

Operating Powered Seating through the Driver's Control. All chairs require adding or replacing one component to allow actuator operation through the driver control.

#### SPJ+ with PSS:

Replace joystick with SPJ+ with ACC (Actuator Control) - E2310.

#### Single Actuator Systems with CMPJ+, PSR+, PSF+ MK6i Display:

Add a SANODE (Single Actuator Node) to the system - E2310.

#### Multiple Actuator Systems:

Replace the standard 4 Way Switch box with a S4WSB (Multiple Actuator Interface Box) - E2311.

**NOTE: All specifications and dimensions are approximate.**

**For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

**Invacare Corporation**

[www.invacare.com](http://www.invacare.com)

**USA**

One Invacare Way  
Elyria, OH USA  
44035-4190  
800-333-6900

©2010 Invacare Corporation

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Products designated by ASL are manufactured by Adaptive Switch Labs, Inc. Therafin is a registered trademark of W/R Medical Electronics Co. Freedom Designs is a registered trademark of Freedom Designs, Inc. Otto Bock is a trademark of Otto Bock Orthopaedische Industrien GMBH and Co. Specifications and prices to change without notification.

